

Spectroscopic Ellipsometry of Compound Semiconductors: $Al_xGa_{1-x}N$ / GaN Hetero-Structures

Roland Seitz - Application Scientist - Thin Film Division

Introduction

Group III-Nitrides and their alloys are the most promising materials for short-wavelength optoelectronic devices such as LEDs, injection lasers, photodetectors, full colour displays and electronic devices like HFETs, HEMTs, etc. For their design and optimization a detailed knowledge of both the layer thickness and the optical properties are essential. Spectroscopic Ellipsometry is a non-destructive optical characterization method that allows determination of these required material parameters.

Materials

A typical AlGaIn / GaN heterostructure as used for LEDs and transistors is shown in figure 1.

AlGaIn 0.2-1 μm
GaN 1-2 μm
Sapphire substrate

Figure 1: Typical AlGaIn / GaN heterostructure

The GaN and AlGaIn films were deposited by MOCVD (metal-organic chemical vapour deposition) on sapphire substrates.

The following samples were analyzed:

Sample	Structure	Al Content in AlGaIn %
1	GaN / Sapphire	-
2	GaN / Sapphire	-
3	AlGaIn / GaN / Sapphire	7
4	AlGaIn / GaN / Sapphire	16
5	AlGaIn / GaN / Sapphire	25
6	AlGaIn / GaN / Sapphire	5
7	AlGaIn / GaN / Sapphire	9

Results

The work was performed using the HORIBA Jobin Yvon MM-16 spectroscopic ellipsometer which provides significant advantages in terms of speed, high resolution measurement and experi-

mental versatility. Ellipsometric measurements were performed at an angle of incidence of 70° in the spectral range 500nm-800nm. Both the thicknesses and optical properties were extracted simultaneously from the SE data analysis. When compared to conventional ellipsometer platforms, the Liquid Crystal Modulation Ellipsometer delivers exceptionally high accuracy for the ellipsometric angles (Ψ , Δ) across their full range in one measurement, without any dead spots.

Figure 2 shows the Ψ and Δ spectra of sample 1.




Figure 2: Ψ and Δ spectra of sample 1

The thickness and the dispersion of the GaN layer were determined by appropriate modelling. For this sample the result is the following:

1.5 nm Overlayer
2402.4 nm GaN
Sapphire substrate

Figure 3 shows the Ψ and Δ spectra of sample 6.




Figure 3: Ψ and Δ spectra of sample 6

The thickness and the dispersion of both the AlGaIn and the GaN layer were determined by appropriate modelling. For sample 6 the result is the following:

3.6 nm Overlayer
462.1 nm AlGaIn
1110.5 nm GaN
Sapphire substrate

For the modelling of the optical dispersions a classical Lorentz oscillator dispersion formula was used:

$$\epsilon(\omega) = 1 + \frac{(\epsilon_s - 1)\omega_t^2}{\omega_t^2 - \omega^2 + i\Gamma_0\omega} = \epsilon_1(\omega) + \epsilon_2(\omega)$$

where $E = h\omega$ is the photon energy.

The relation of ϵ_1 and ϵ_2 with n and k is: $\epsilon_1 = n^2 - k^2$ and $\epsilon_2 = 2nk$

The following table summarizes the results found for the samples in the wavelength range 500-800nm:

Sample	thickness GaN/nm	thickness AlGaIn/nm	n at 633nm	ϵ_s	ω_t	Γ_0
1	2402	0	2.361	5.19	6.76	0.1
2	2466	0	2.353	5.16	6.78	0.1
3	1283	332	2.331	5.07	6.85	0.2
4	1180	254	2.299	5.00	7.56	0.3
5	1128	401	2.292	4.93	7.11	0.5
6	1111	462	2.340	5.09	6.68	0.2
7	1124	602	2.336	5.08	6.75	0.2

The refractive index as a function of wavelength is shown in figure 4 for different Al concentrations.




Figure 4: Refractive index for various Al concentrations in AlGaIn

From these data a calibration curve can be set-up that allows the determination of the Al content in the AlGaIn layer by evaluating the optical dispersion of the material (figure 5):




Figure 5: Calibration curve for Al concentration in AlGaIn

The Al concentration can be calculated by the following formula:

$$Al - content(\%) = 182\epsilon_s^2 - 1943\epsilon_s + 5180$$

Conclusion

Liquid Crystal Modulation Spectroscopic Ellipsometry is an excellent technique for the highly accurate characterization of the compound semiconductor heterostructure AlGaIn / GaIn.

Using the MM-16 spectroscopic ellipsometer it is a straightforward procedure to determine the film thickness and optical dispersions of the complete structure even where the film is several microns thick.

The detailed knowledge of the optical parameters of AlGaIn alloys is crucial for example for the design of opto-electronic devices.

Furthermore, from the optical parameters a calibration curve could be constructed to provide a rapid and efficient determination of the Al content in the AlGaIn layers. Thus Spectroscopic Ellipsometry also proves a non-destructive technique for AlGaIn alloy composition determination.

This method can be equally applied to other compound semiconductors such as SiGe, II-VI semiconductors or classical III-V semiconductors.