

#### PS415 Standard Measurement on Partica LA-960V2 with PowderJet

#### Introduction

Poly-disperse glass bead standards were developed as a better test of complete system performance for laser diffraction analyzers, compared to monodisperse polystyrene latex dispersions that are not representative of most materials tested on these instruments. In this report PS415 was tested on the LA-960V2 with Powderjet. PS415 is a 10  $\mu\text{m}$  to 100  $\mu\text{m}$ , 1-gram single shot polydisperse glass bead standard produced by Whitehouse Scientific. Repeatability for all size classes was better than 2.5%, showing great repeatability.

#### Instrument Setup

##### • Instrument Configuration

- o Cell: Dry Cell
- o Nozzle: Medium

##### • Sample Information

- o Sample Name: PS415
- o Material: Glass beads
- o Source: Whitehouse Scientific
- o Refractive Index: 1.510 – 0.00i
- o Iteration mode: Manual
- o Convergence Factor: 15
- o Distribution base: Volume

##### • Advanced Measurement Conditions

- o Data acquisition times (Sample): 50000
- o Data acquisition times (Blank): 5000
- o Auto OFF after measurement: Check all items
- o Auto ON with auto blank: Check Vacuum and Air
- o T% for Sampling: Yes
  - Max T%: 99%
  - Min T%: 95
- o Start Trigger: No
- o Stop Trigger: Yes
  - T%: 99%
  - Behavior: Stop After Waiting
- o Feeder Speed: Auto
- o Air Pressure: 0.3 MPa



Figure 1. LA-960V2 with Dry Powder Feeder attached.

#### Test Procedure

1. Insert MEDIUM dispersion nozzle into dry cell.
2. Activate Vacuum and press Alignment. Inspect detector/channel baseline to determine cleanliness of dry cell glass.
3. Spread the whole bottle of standard along the chute.
4. Click Auto-Measurement.
5. Wait for measurement to complete then brush out feeder chute to prepare for next measurement.

#### Results

Verify that the D10, D50, and D90 are all within acceptable ranges.




Figure 2: An overlay from the LA-960V2 showing six runs of the PS415.

Data Name	D10 (µm)	D50 (µm)	D90 (µm)
PS415 RUN 01	33.401	53.211	85.626
PS415 RUN 02	32.633	51.268	78.694
PS415 RUN 03	32.577	52.485	83.213
PS415 RUN 04	32.907	52.171	81.618
PS415 RUN 05	32.349	51.447	81.060
PS415 RUN 06	33.097	51.804	80.737
<b>Average</b>	32.827	52.064	81.825
<b>Std. Dev.</b>	0.384	0.719	2.366
<b>CV (%)</b>	1.170	1.382	2.891
<b>Certified Diameter</b>	33.2	53.9	85.8
<b>Uncertainty</b>	3.0	4.2	6.0
<b>Upper Limit</b>	29.92	49.40	78.42
<b>Lower Limit</b>	36.48	58.40	93.18

Table 1: Summary of six runs from the LA-960. Upper and lower limits were determined using the published PS415 uncertainty and the ISO 13320:2020 tolerance criteria, 1.5% at D50, 2.0% at D10, and 2.5% at D90.

NIST Traceable

# Polydisperse Particle Standard

# 10-100 $\mu$ m

Part Number: PS415  
Nominal Weight: 1.0g x 10 bottles

Notes:

- (1) Traceability – NIST primary reference calibrated stage reference graticule, ref. 821/263573-00 and National Physical Laboratory (Teddington, UK) stage reference graticule, ref. 08A038/970127/106-66 have been used during the primary calibration of both sieving and image analysis methods.
- (2) Whitehouse Scientific Ltd operates to the principles of ISO17025 with all business operations covered by a comprehensive Quality system and supporting documentation.
- (3) This certificate is only valid if a complete single-shot bottle is used in the analysis.
- (4) Whitehouse Scientific Ltd does not accept responsibility for losses, financial of otherwise which may occur as a result of the interpretation or use of the information contained within this certificate.

## Certificate of Analysis

Whitehurst Road  
Waverham  
Dorset DT99 7PB  
UK  
+44 (0)149 330266

NIST TRACEABLE

### POLYDISPERSE PARTICLE STANDARD SODA-LIME GLASS MICROSPHERES


Bottles Numbers:  
11571-11580

#### 1. Detailed Review

(sieve analyses performed in dry state)  
(all results by weight/volume, total of 13 individual tests)

Method		Percentile						
		5	10	25	50	75	90	95
Electroformed Sieve Analysis (by weight)	Mean	28.8	34.60	44.1	55.7	70.9	87.6	97.2
	SD	0.7	0.7	0.4	0.3	0.4	0.6	0.5
Wire Sieve (by weight)	Mean	27.3	33.0	43.9	54.3	70.5	87.1	93.6
	SD	0.6	0.6	0.3	0.4	0.6	0.6	0.2
Microscopy (x2) (by volume)	Mean	25.8	31.10	40.8	50.5	63.7	80.7	89.5
	Mean	0.5	0.5	0.5	0.4	0.5	0.6	0.9


#### 2. Graphical Review



#### 2. Tabular Summary

Cumulative percent undersize	10	25	50	75	90
Final Mean Size - microns	33.2	43.3	53.9	69.1	85.8
Uncertainty - microns (95% confidence)	3.0	2.9	4.2	6.3	6.0

Issued by:



J A Miles, Laboratory Manager

