


Summary

Titanium powder is a naturally occurring mineral ranging in color in various shades of gray and black. Titanium is used in a wide range of manufactured products. Some applications where it can be found are coatings, car parts, lightweight aerospace components, and additive manufacturing (3D printing). It can also be found in medical devices and dental implants due to its high biocompatibility.

The particle size of titanium powder is important to final product surface quality, quality and robustness of coatings and mechanical strength of products. It also allows for the monitoring for elimination of large sizes which can plug nozzles in additive manufacturing.

Analytical test method dry

Sample amount: 1 gram
 Covered area: 0.5%
 Feeder Speed: 20
 Measurement Model: Area equivalent diameter
 Number of particles: 400,000
 Feeder Width: 20mm
 Objective: 0.735x




Example dry data
 Median: 35.70 µm Mean: 229.1 µm
 D(10%): 22.70 µm D(90%): 229.1 µm

Figure 1. Measurement of 10⁴Ti64 dry

Analytical test method wet

Dispersant fluid: water
 Surfactant: none
 Liquid level: medium
 Sonication: 1 minute at 50% power
 Circulation speed: 5
 Measurement Model: Area equivalent diameter
 Sample amount: 250 milligrams
 Number of particles: 1,000,000
 Objective: 1.333x



Example wet data
 Median: 38.14 µm Mean: 52.26 µm
 D(10%): 24.26 µm D(90%): 52.26 µm

Figure 2. Ti64 02 Wet