


## Analysis of films on glass by pulsed RF GDOES



Application Note  
Material Science  
GD32

Patrick Chapon, HORIBA Scientific, 16 rue du Canal, 91160 Longjumeau, France

Abstract: Some examples of elemental depth profiles of films on glass are shown coming from published data and results from the GD days illustrating the capability of the Pulsed RF GDOES instrument in this domain.

### Key words

Coatings on Glass, GDOES, Pulsed RF. Fast elemental depth profiling, Glow Discharge Optical Emission Spectroscopy

### Introduction

Coatings on glass are extensively used in multiple domains from thin films PV cells to optical coatings, IR protective layers on windows to gold coatings on the HORIBA Scientific SPRI prisms!

### Instrumentation

The GD Profiler 2 couples an advanced Pulsed RF Glow Discharge Source to a high resolution, wide spectral range Optical Emission Spectrometer. RF is of course mandatory when non conductive samples are to be measured but RF only is not sufficient as induced heat could be generated on the sample during the measurement leading to unwanted diffusion of elements. This is why pulsed RF (with patented automatic matching in pulsed mode) is crucially important for coated or treated glasses analysis (Ref 3 & 4).




Figure 1: HORIBA Scientific SPRI chips

Multiple deposition techniques are used to coat glasses (PVD, sol gel dip coating, CVD etc – Ref 1). GDOES permits to rapidly obtain the elemental depth profiles and assess the quality of the coatings and the reproducibility of the selected deposition technique (Ref 2).




Figure 2: Depth profile of a SPRI prism featuring the Au top layer on a sublayer on glass



Figure 3: GD Profiler 2 instrument

## Results

The first result is from Ref 2 in which titanium oxy-nitride films are studied as energy efficient glazing.




Figure 4: Pulsed RF GDOES profile of a TiON film on glass (from ref 2)

The second measurement has been shown at the 7th GD day (Ref 5). The sample is a Rugate filter with 46 periods of  $\text{SiO}_2/\text{Ta}_2\text{O}_5$  on glass. The depth resolution here was better than 50 nm at 9  $\mu$  sputtered depth!




Figure 5: Pulsed RF GD-OES intensity-time-profile and FE-SEM micrograph

## Conclusion

Pulsed RF GDOES is a valuable technique for fast elemental depth profile of coated and treated glasses. It is capable to measure all elements with excellent depth resolution and provide information on layers composition and thickness, diffusions or contamination at interfaces.

## References

- 1) « Thin films on glass” Back & Krause. Springer.
- 2) « Synthesis and energy modelling studies of titanium oxy-nitride films as energy efficient glazing”. Michael E.A.Warwick, Geoffrey Hyett, Ian Ridley, Fathima R.Laffir, Celia Olivero, Patrick Chapon, Russell Binions. Solar Energy Materials & Solar Cells 118(2013)149–156
- 3) “Determination of sodium migration in sol–gel deposited titania films on soda-lime glass with r.f. glow discharge optical emission spectroscopy”. B. Yuksel, E.D. Sam, O.C. Aktas, M. Urgen, A.F. Cakir. Applied Surface Science 255 (2009) 4001–4004
- 4) “Depth Profiling of Corning Gorilla Glass using RF GD-OES”. Yura Guryanov, Natalia Dawn, Jamie Vargeson, Anna Nached. USA GD day 2013.
- 5) “RF GD-OES Analyses of Optical Multilayers”. Olaf Zywitzki and Richard Belau, Fraunhofer FEP, Dresden, Germany. 7th GD day. Reims 2014.