

*Quantum-efficiency measurements using HORIBA spectrometers and detectors*

## Evaluation of Novel Photoresponsive Materials via EQE Measurements

### Introduction

Internal Quantum Efficiency (IQE) and External Quantum Efficiency (EQE) measurements are indicators of the effectiveness of a photosensitive device such as those used in telecommunications and solar cells. EQE is the ratio of the charges generated to the total amount of photons incident on the surface; a larger EQE indicates a more efficient device. IQE is the ratio of the charges generated to the number of photons that hit the surface and are absorbed by the cell. Both EQE and IQE measurements are vital in the interrogation of light-responsive devices.<sup>1-3</sup>

### Experiment

Photosensitive devices are studied using a variety of techniques, including scanning electron microscopy, electroluminescence (EL), and photoluminescence. A typical experimental setup for measuring the EQE of a device under test is pictured in Fig. 1. A white-light source (in this case, a tungsten-halogen lamp) is coupled to a HORIBA iHR550 imaging spectrometer, producing monochromatic light. The light output is then modulated using a mechanical chopper (HORIBA ACH-C or ACH-C-OPEN), and then the light strikes the sample. The voltage output from the sample is then measured using a Stanford Research Systems lock-in amplifier, which synchronizes the chopper frequency with the modulated signal. In addition, the number of photons generated at the output is measured using a calibrated photodiode to calculate a final EQE percentage (not shown).<sup>1,2</sup>

### EQE of thin-film solar cell

EQE spectra of organic photovoltaic cells were car-

ried out at the Tokyo Institute of Technology. In this experiment, an iHR550 with a 100 W halogen lamp created monochromatic light, which was passed through a chopper to the sample, a thin film of zinc octaethylporphyrin. Two measurements were made: one with an amorphous sample, and another in which the film was annealed at 473 K for one minute. As seen in Fig. 2, the annealed sample shows an overall increase across the entire range, with the biggest enhancement at ~400 nm. This enhancement is attributed to increased mobility of the charge carriers within the crystalline film.<sup>1</sup>




Fig. 1. Typical EQE experiment.

### EQE of InAs nanowires

A second study came from the University of California–San Diego in the fabrication and testing of InAs wires. Here a halogen lamp and iHR550 were used to create monochromatic light, which was subsequently chopped and focused onto the



sample. EQE results are shown in Fig. 3 for both the visible and infrared region in the sample, which was comprised of n-type InAs nanowires and a p-type Si substrate. The region of absorption below  $\sim 1200$  nm was attributed to both InAs and Si, while only the InAs nanowires contribute above  $1200$  nm.<sup>2</sup>




Fig. 2. EQE spectra of amorphous and annealed Zn(OEP) films. Reprinted with permission from reference 1. Copyright 2013 Sou Ryuzaki and Jun Onoe.




Fig. 3. EQE of n-InAs nanowires on p-Si substrate. Reprinted with permission from reference 2. Copyright 2009 American Chemical Society.

[info-sci@horiba.com](mailto:info-sci@horiba.com)  
[www.horiba.com/scientific](http://www.horiba.com/scientific)

**HORIBA**  
Scientific

USA: +1 732 494 8660  
UK: +44 (0)20 8204 8142  
Spain: +34 91 490 23 34  
Other Countries: +33 (0)1 64 54 13 00

France: +33 (0)1 64 54 13 00  
Italy: +39 0 2 5760 3050

China: +86 (0)10 8567 9966

Germany: +49 (0)89 4623 17-0  
Japan: +81 (0)3 38618231  
Brazil: +55 11 5545 1540



## EQE and IQE of PbS nanocrystals

Finally, EQE of electroluminescent samples of PbS nanocrystals, emitting in the near-infrared, were characterized at the Nano-Electronics Centre at the University of Surrey. The EL spectra were recorded using an iHR550 and Symphony InGaAs 512-element array and a Symphony  $1024 \times 512$  CCD, depending upon the wavelength region of interest. Integrated EL measurements were carried out to attain EQE values with an InGaAs integrating sphere. The EQE was measured as a function of voltage applied to the sample, and found to be  $\sim 1\%$ . The IQE was estimated from transmission measurements to be between 5% and 12%. These values indicate a highly efficient near-IR emitting device, which would be valuable in a variety of systems, including biosensing and optical circuits.<sup>3</sup>

## Conclusions

EQE measurements are an important indicator of the effectiveness of photosensitive devices. Three such studies carried out with HORIBA instrumentation are summarized herein. These studies have strong implications for a variety of applications including solar cell research and telecommunications.

## References

1. S. Ryuzaki and J. Onoe, *Nano Reviews* **4**, 21055 (2013).
2. W. Wei, X. Bao, C. Soci, Y. Ding, Z. Wang, and D. Wang, *Nano Letters* **9**, 2926 (2009).
3. K.N. Bourdakos, D.M.N.M. Dis-sanayake, T. Lutz, S.R.P. Silva, and R.J. Curry, *Appl. Phys. Lett.* **92**, 153311 (2008).