ISSN 0915-9916

# Readout

HORIBA Technical Reports

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March 1996 ■ No.12

**HORIBA**

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---

<input type="checkbox"/> Foreword	HORIBA and Medical Treatment Atushi HORIBA ——— 3
<input type="checkbox"/> Free Talking	Recent Clinical Laboratory and Equipment Manabu YOSHIMURA, Shinobu NAKANISHI, Hisao SAKAKI, Takuji SUZAKI ——— 4
<hr/>	
<input type="checkbox"/> Guest Forum	Present Status of Hematology Analysing Systems Jean-Christophe Frot ——— 12
<hr/>	
● Feature Articles	HORIBA's Clinical Testing Equipment Seiji USUI ——— 23
	The Franco-Japanese Co-Development of the Exclusive LC-220 Automated Blood Cell Analyzer Brice Prouve, Roger Le Comte, Narihiro OKU ——— 29
	Chararistic of Fully Automated Chemiluminescent Enzyme Immunoassay Equipment Luminomaster — Equipment Construction and Specifications — Takeshi KOHNO, Yoshiyuki SAWADA ——— 35
	— Reagents and Detection Systems — Takuji IKEGAMI, Yoshihiro SATOU ——— 43
<hr/>	
● New Products	Energy Dispersive X-Ray Analyzer EMAX-5770W TypeX, XR ——— 50
<hr/>	
<input type="checkbox"/> Guest Forum	Social Infrastructures and Measurement Technology Seiji TAKATA ——— 52
<hr/>	
◇ Topics	Expectations for a Laboratory Accreditation System for Assurance of the Reliability of Test Data and Mutual Acceptance Eizo TERANISHI ——— 57
<hr/>	
● Selected Topics	In-line Gas Monitor Naohito SHIMIZU, Hiroyuki EBI, Norio KADA Tetsuo SHIMIZU, Hiroyuki HARADA ——— 59
	Various Water Quality Analyzers for the Quality Control of Boiler Water at Thermoelectric Power Plants Susumu YAMAUCHI, Kazuhiko ASHIKAGA ——— 65
	Portable Gas Analyzer PG-200 Series for the Simple Measurement of Combustion Gases Toshikazu OHNISHI ——— 73
<hr/>	
◇ Patents	——— 78
◇ Technical Reports List	——— 79

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---

### Recent Clinical Laboratory and Equipment

Along with the maturing of society, the expectations we have of medicine have risen and diversified. Against this backdrop, the role that clinical laboratory plays is becoming all the more important in medical activities. On the other hand, amazing advances have been made in clinical laboratory equipment and peripheral technology. This, in turn, has stimulated a need for grasping their potential in the future.

Accordingly, we have invited three distinguished guests to discuss the latest developments in clinical laboratory and equipment. Dr. Manabu Yoshimura, a professor at the Kyoto Prefectural Medical College, is an authority on clinical examination and medicine. Dr. Shinobu Nakanishi is currently active as a practising physician at the vanguard of regional medicine. Mr. Hisao Sakaki, chief doctor at the Central Laboratory Dept. of the Akita University School of Medicine, is engaged in the building of up-to-date clinical examination systems. Takuji Suzuki, manager of development of medical equipment systems at HORIBA will act as MC.

This round-table discussion first deals with the automation of clinical equipment, then moves on to client server information systems and the standardization of equipment and reagents, and finally regional and home-based medicine.

---

## 1. Introduction

It has nearly taken two hundreds years since the recognition by VESALE of the blood as a fluid circulating through the body until the description of the two phases, one liquid and the other solid, under the microscope in the mid-eighteenth century. The solid phase is in fact the blood cells and another hundred years was necessary to classify the blood cells into populations thanks to the staining of blood films in the late nineteenth century by Giemsa, adapting the original Romanovsky stain.

Apart from the count of cells for each population, red blood cells (RBC), white blood cells (WBC), platelets (PLT), some more parameters were added; either measured like the pack cell volume (PCV) or hematocrit (Ht), hemoglobin concentration (Hb) or calculated like mean corpuscular hemoglobin concentration (MCHC), mean corpuscular hemoglobin (MCH) and mean cell volume (MCV).

If, at least in a normal patient, the RBC and PLT represent a quite homogenous population the WBC are classified in three populations as lymphocytes, monocytes and granulocytes. The latter is divided into three sub-populations as neutrophils, eosinophils and basophils. There were then five different populations of WBC described under the microscope according to their morphology and staining characteristics, the so-called white blood cell differentiation or 5 population differential.

So, until the 1950's, what we call the blood cell count (CBC) and the white blood cell differentiation were performed manually. As the CBC became one of the most common tests required to judge the status of a patient, it became obvious that the

	DEFINITION	UNITS	NORMAL RANGE ADULT
MCHC	$\frac{Hb(g/l)}{Ht(l/l)}$	g/l	32-36
MCH	$\frac{Hb(g/l)}{RBC(\times 10^{12}/l)}$	pg	27-33
MCV	$\frac{Ht(l/l) \times 1000}{RBC(\times 10^{12}/l)}$	fl	80-100

Table.1 Calculated parameters of CBC

	UNITS	NORMAL RANGE ADULT	
		MALE	FEMALE
RBC	$\times 10^{12}/l$	4.5-5.4	3.8-5.2
Hb concentration	$\times g/l$	131-165	120-150
Ht	l/l	0.40-0.50	0.35-0.45
PLT	$\times 10^9/l$	150-400	
WBC	$\times 10^9/l$	4-6	

Table.2 The basic measured parameters of CBC

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<Experience>

Born June 17th, 1949

1976-81: M.D. Resident in  
Clinical Medicine, Paris  
Hospitals, France.

1981-85: Deputy Manager at  
Clinical Chemistry Lab.  
University Hospital in  
Paris, France.  
Lecturer in Clinical  
Chemistry at School of  
Medicine in Paris,  
France. Ph.D. in  
Biochemistry.  
Research on clinical  
significance of protein  
profiles and antenatal  
diagnosis.

1985-87: Employed with Technicon,  
France.  
Sales representative for  
high end clinical chemistry  
analyzers.

1988- : Employed with ABX,  
Montpellier, France.  
Export & Marketing Direc-  
tor. Management of the  
promotion & distribution  
of ABX hematology ana-  
lyzers worldwide.

manual method was time-consuming and would no longer be applicable to large scale screening. The time for automation became mandatory.

In order to describe the present status of hematology analyzing systems, we should divide this article into three chapters: automation in blood cell counting, automation in WBC differential count and automation in sample processing and data management.

## 2. Automation in Blood Cell Counting

Automation has naturally taken as reference the manual methods to design and manufacture the system needed. Referring to blood cell counting, the first step is to know what blood cells can be compared to. The reference was particles for which we need to count whether it will be red cells, white cells or platelets. Then to differentiate between those three populations, there is a need not only to count but to size the particles as well. At present, there are two different technologies which can be used: electrical impedance and light scattering. Hemoglobin, on the other hand, is measured with a classical spectrophotometric method.

Referring now to the concentration of cells in the blood which differ in a 1/1000 ratio from RBC to WBC, there is a need for dilution. Dilution of the samples is then a basic requirement on blood cell counting automation.

Whatever the technology used, the analyzers are working on the same basis: they have at least two channels, one for RBC, PLT and Ht and another one for WBC and

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Hb. There are different ways to process the sample (see table 4) but the results are the same at the end. In one channel, a diluent is added and RBC and PLT are counted. In the second one, diluent and a lytic agent are added. The purpose of the lytic agent is to lyse the RBC, leaving the WBC intact and to provide a solution for Hb measurement. In that case, the lysing agent comprises the reagent needed for hemoglobin measurement.

Hb	Cyanmethemoglobin	in the presence of potassium cyanide, Hb is converted into cyanmethemoglobin measured at 540 nm by means of a photometer or colorimeter
Ht	Centrifugation	Ht is the proportion of a column of centrifuged blood which is occupied by the red cells

Table.3 Reference methods for Hb and Ht (ISCH & NCCLs recommendations)

## 2.1 Analyzers depending on electrical impedance

The principle refers to the very poor conductivity of blood cells. When placed in a conducting medium, cells are aspirated through an aperture across which an electrical current is applied. When the cell is going across the aperture, there is an increase in the electrical impedance. This increase is measurable and proportional to the size of the cell. Thus, the electrical pulse generated corresponds to one count and its height corresponds to the size of the cell counted. So, cells are counted and sized, giving the opportunity to differentiate the populations.

In the RBC and PLT channel, Ht is measured by adding all the pulses produced by all the cells, given then the total volume of cells. Another way is to calculate the MCV from the mean pulse and then deduct the Ht from the formula on Table 2.

In a normal sample, the effect of WBC present in the RBC channel is negligible but has to be taken into consideration when there is a marked increase like in leukemia. Also, the difference in size between PLT and RBC in a normal sample gives a big enough gap to differentiate the two populations. In pathological samples with either

	PRINCIPLE	CHARACTERISTICS
ASPIRATION NEEDLE	Aliquot of blood is transferred to dilution chambers With multi-step dilution to the counting chambers	* Low sample volume (12 $\mu$ l on LC 220) * Slow * Easy cleaning
SHEAR VALVE	The blood column is cut into aliquot  Aliquot directly flows to the counting chamber by diluent	* High sample volume (100 $\mu$ l at least) * Quick * Sensitive to blood clot, needs regular maintenance

Table.4 Sample processing device for dilution

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large PLT or small RBC, there could be an overlapping of the two populations. Several electronic flags can now alert the user of such overlapping.

The WBC channel is also used to measure the Hb through a spectrophotometer, which is most of the time placed behind the counting chamber. In a normal sample, RBC and PLT are lysed, which then gives no interference with the WBC count. In case of presence of NRBC, the remaining nucleus are counted and interfere with the WBC count. Electronic thresholds are able to fix the limits and to give flag in such conditions.

## 2.2 Analyzers depending on light scattering

The principle refers to both absorbance and light scattering of a particle when passing through a light beam. The light source may be white or laser light, or both. Absorbance is mostly used to count the cells as indicating that a cell is passing through the light beam. Scattered light is used to size the cell. Depending on the design of the optical bench used, there is an optimum angle to measure the size of the cell.

All the instruments using light scattering technology have a chamber design including sheath flow in order to avoid coincidence. Otherwise, counting and sizing the RBC, PLT and WBC gives comparable results as with impedance technology. The measurements of Ht and Hb remains the same.

Comparing the two technologies, the impedance needs no adjustment but is subject to more interferences. To the contrary, there is less interferences with light scattering but the adjustment is critical to get quality results and the heat of a laser when

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
used is still not convenient. It is foreseen that the use of laser diode will greatly improved this aspect.

### 2.3 Automated blood cell counting and cell distribution

As far as the analyzers are sizing the cell, it is then possible to classify within the same population (i.e. RBC) cells having the same volume. It is then possible to visualize cell volume distribution through a histogram. On the X axis, volume is represented, on the Y axis, relative number of cells for a defined volume (see Fig.1).

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The shape of the histogram is then a good indicator of the homogeneity of the population. It has then been possible to define new parameters for RBC and PLT, so-called red blood cell distribution width and platelet distribution width respectively.

### 3. Automation in WBC Differential Count

As already mentioned, manual WBC differential counting is time-consuming and inaccurate. The inaccuracy comes mainly from the small number of cells (100) counted under the microscope and the quality of the blood film and its stain.

Automated counters overturned this inconvenience, being able to count quickly a large number of cells, they are thus far more accurate. In addition to the count, the WBC differentiation is based upon the morphological characteristic of the cells. Size is one of those characteristics and we have seen above, that size between PLT, RBC, WBC at least in a normal sample, differs greatly enough to separate these three populations.

The problem becomes more complex with the WBC differentiation. In a normal sample, the lymphocytes are the smallest, the monocytes the biggest and the granulocytes are in between in terms of size. But granulocytes, considering neutrophils, eosinophils and basophils, are almost the same size. There is thus a need to define another characteristic being measurable to differentiate. Furthermore, in pathological samples, you may find large lymphocytes or cells which normally should not be present in the blood stream and which would have a size comparable to the normal size.

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All those conditions have led to the point that different or added technologies are now used in WBC differential automation. Two different ways were explored in the 70's: image of pattern recognition analyzers and flow cytochemistry instruments. During the 80's, improvement in the analysis of the signals given by cell sizing and light scattering gave the way to another generation of instruments. At present, those three approaches are still in use and even sometimes, combined in an attempt to give even more accurate results.

### 3.1 Pattern recognition instruments

Progress in computerization has eventually put that technology back into fashion. From a blood film, for which spreading and staining are standardized by a robotic slide-maker, the WBC differential is performed under a camera-equipped microscope. The image is then compared to the one which is memorized in a computer. Up to 250 characteristics of each possible cell to be found in the blood stream are memorized and there is the possibility to virtually identify each kind of WBC and count them.

Now, with the memory capacities and the speed of computer, up to 120 blood films with very accurate results are processed with these analyzers. But there is still the need for a blood film, which means either labor costs or investments in a blood film-maker robot. Nevertheless, for large laboratories, it is eventually a good additional option to the other instruments.

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### 3.2 Flow cytochemistry instruments

Using the right chemical, advantage is taken of a specific fixation of this chemical to a population of cells. WBC differential count is part of the full CBC given by the same instruments by adding one or more channels. The sample is then processed only once, blood film is not required and CBC and WBC differential count are given simultaneously. In order to comply with the requested productivity, the chemical used should react quickly with the cell and give a measurable signal. Flow cytochemistry is mainly used to differentiate eosinophils and neutrophils from the other three populations. Combining this technology with light scatter or volume, it is possible to build-up a matrix which shows a good separation of the different WBC populations (Fig.1).

### 3.3 Cell-sizing & light scattering instruments

The WBC react differently to lysis, depending on the strength of the lysing agent and as well to the pH of the solution. These characteristics are used to lyse some of the population and then count and size the remaining ones. The same occurs with light scattering and light absorbance which is partly due to, not only the size and shape of the cell, but also to its internal structure. Analyzing the scattered light at different angles and combining the different measured signals lead to the differentiation of the WBC. Again, in combining the two technologies, it is possible to build-up a matrix as for the flow cytochemistry instruments.

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

## 4. Automation in Sample Processing & Data Management

(A) COBAS ARGOS using rotating trays as an auto sampler

(B) COBAS VEGA using rack system

Fig. 2 Two examples of fully automated high-through put CBC+5DIFF analyzers.

Up to now, we did describe automation in counting the cells. We have pointed out that the actual instrument can perform up to 120 CBC and WBC differential simultaneously per hour. It is then obvious that sample processing and data management would become the most important labor cost of a hematology workstation.

Automation of these two items is then mandatory but has not developed the same way. The problems to be solved are in the forthcoming order: bring the sample to the instrument, present the sample to the instrument, identify the sample, aspirate the sample, recall the results, transfer the right results of the sample to the prescriber. Except for the first step, all others are now fully automated on large instruments.

### (1) Bring the sample to the instrument

Human resources are still the most common way. Manufacturers start to develop conveyor lines with computer-driven rack-systems which bring the samples automatically to the various work stations where the samples are transferred to the instruments. Two obstacles reserve such kind of organization to very large laboratories: it is space-consuming and the cost is still expensive, still more expensive than the earning in personnel. But at least, the concept will be valid when these two obstacles are erased.

### (2) Present the sample to the instrument

On small instruments, samples are presented by hand and one by one under the aspirating needle. Medium or large size instruments are now fitted with autosamplers. Tubes are positioned on a tray or racks, held by clips, the tray or racks being placed on the instruments which process the samples automatically. The instrument does not need someone to feed it and has more than one hour's autonomy. They are called "walkaway" systems.

Some instruments include an automatic rerun capabilities if the results are outside normal range, giving even more "walk away".

### (3) Identify the sample

It is always possible to enter the ID number or the name manually through the front panel keyboard. On most of the instruments, this is now done with a barcode reader with identification of the tray or rack number, the position of the tube in such and if labelled with barcode, the ID number of the sample.

### (4) Aspirating the sample

For evident safety reasons, all new generations of instruments are proposed with a closed tube sampling device, whether they are fitted or not with an autosampler. A needle pierces the rubber cap of the tube and with the help of a level detector, goes down to aspirate the required quantity of blood.

#### (5) Transfer of the results

All instruments are now fitted with a computer output allowing transmission of data to a main frame. In such case, sample identification is mandatory. The most sophisticated instruments have on board, or beside, a PC and work on two-way communications, allowing random access and selectivity. The work list is transferred to the instrument with the ID number and the requested parameters for each sample. When the barcoded sample is recognized by the instrument, it is processed and the results directly sent to the main frame. In such case, no intervention has been necessary.

### 5. Conclusion

Like in other fields in laboratory equipment, hematology has greatly taken advantage of automation. Robotic, electronic, computerization have given all the necessary tools to give efficiency, productivity and reliability to the labs. The tendency is now to consolidate the different available technologies and to combine several instruments to enhance productivity.

Major improvement in cell detection and differentiation are not foreseen in the near future, partly because the cost constraints in healthcare programs do not allow the development of expensive technology. Progress in automation will then rely more on robotic and computerization for the next years.

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## HORIBA's Clinical Testing Equipment

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**Abstract**

The increasingly higher standard of medical treatment today is making clinical testing all the more indispensable. HORIBA is applying the analysis technology it has acquired in other fields of industry to the development and commercialization of a variety of clinical testing equipment. Of the various clinical testing-related products made by HORIBA, this paper describes the measurement principles, applications and features of the electrolytic analyzer SERA series, blood cell counter LC series, fully automated chemiluminescent enzyme immunoassay equipment LEIA-2000 and the artery temperature detection type thermometer IT-500M.

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# The Franco-Japanese Co-Development of the Exclusive LC-220 Automated Blood Cell Analyzer

Brice Prouve\*, Roger Le Comte\*, Narihiro OKU

## <Abstract>

Horiba and ABX designed and co-developed the LC-220 Haematology Analyzer based on a simple basic operation for doctor's offices with 5 to 10 samples per day. The LC-220 can analyze the four parameters ; WBC, RBC, and Hct using electrical impedance method, Hgb using spectrophotometry. The concept of the system requires minimum knowledge of operation and very few maintenance.

## 1. History

Tests of Haematology are often used worldwide by laboratory technologists. Haematology refers to the study of human blood cells in term of quantity and quality. Quantity of cells is given as a concentration or number of cells per liter of whole blood. Quality of cells refers to their shape, size, and aspect.

Alteration of those two characteristics will diagnose many kinds of diseases or pathologies detected with the help of the microscope: The sampling is diluted and placed on a platform (Neubauer Haemocytometer) of ruled areas composed of 9 squares of 3 millimeters of equal size. Counting the cells is long, difficult and requires some specific skills and training. This method is still used as reference or confirmation of the pathology detected with an Automated Haematology Analyzer. Most of this alterations have been detected with the development of the Automated

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\* ABX S.A.

Fig.1 Automated blood cell analyzer LC-220

Haematology Blood Cell Analyzer since the fifties. The first cell counters were semiautomated: the technologists had to per-dilute the sample before aspiration into the analyzer inducing extra time. Progressively, the technology was adapted to run the sample directly on whole blood right after the collection. The LC-220 is one fully automated system including automatic dilution and automatic rinsing with simple basic technology and few parts.

Horiba and ABX designed and co-developed the LC-220 Haematology Analyzer (Fig.1) based on a simple basic operation for doctor's office with 5 to 10 samples per day. The concept of the system requires minimum knowledge of operation and very few maintenance. The specifications of LC-220 and summarized in the Table 1.

Parameters:	WBC,RBC,HGB,HCT.
Measuring Principles:	Aperture-Impedance for WBC,RBC,and HCT Size of aperture:80 micrometers Cyanmethemoglobin for HGB at 555 nm LED
Sample Volume:	10 microliters Dilution ratio for WBC at 1/240 Dilution ratio for RBC at 1/30,000
Throughout	30 samples per hour
Reproducibility:	WBC: ≤ 3.0% at $10 \times 10^9/\mu\text{l}$ RBC : ≤ 2.0% at $5 \times 10^9/\mu\text{l}$ Hgt : ≤ 1.5% at g/l Hct : ≤ 2.0% at 45%
Linearity:	WBC: 0.5 to $80.0 \times 10^9/\mu\text{l}$ : <±0.3 or±2.5% RBC : 0.5 to $7.5 \times 10^9/\mu\text{l}$ : <±0.07 or±2.0% Hgb : 2.0 to 25.0 g/l : <±0.3 or±2.0% Hct : 10.0 to 59.9 % : <±2.0 or±3.0% 60.0 to 70.0 % : <±5.0 or±5.0%
Carryover:	WBC: <2% RBC : <0.5% Hgb : <2% Hct : <2%
Display:	LED:Sample No. (1-99) Results (WBC,RBC,Hgb,Hct) States (Measuring,Ready,Operating Error)
Memory:	Last result only
Operating Temperature:	1.5-30 degrees Celcius
Humidity:	25-80%
Power Supply:	100V-50/60 Hz
Option:	External printer or RS232C

Table .1 LC-220 performance specifications

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## 2. How do we determine the four parameters: WBC, RBC, Hgb, and Hct ?

Blood sample is collected from the vein of the arm into a tube and mixed with an anticoagulant (EDTA) to prevent coagulation. The normal concentration of the Red Blood Cell (RBC) or Erythrocytes is about 5,000,000 per  $\text{mm}^3$ . Normal White Blood Cells (WBC) or Leukocytes' Concentration is about 7,000 per  $\text{mm}^3$ . It will be then necessary to dilute this concentration to count the cells.

The reagent for dilution is composed of a filtered isotonic solution similar to a physiological serum. The LC-220 performs two dilutions successively. The final ratio for the first dilution is 1/240 including 0.4 ml of lysing agent.

This lysing agent is a solution derived from the Saponine and has two purposes: The potassium cyanide will free the hemoglobin from the RBC's cytoplasmic membrane to form the chromogenous cyanmethemoglobin compound. This compound is then measured by spectrophotometry through the optical part of the unique chamber with a wave length of 555 nm. The optical density is proportional to the concentration of hemoglobin.

This first dilution will then be used to count the WBC's and measure the Hemoglobin. The final ration for the second dilution is 1/30,000.

The second dilution will be used to count the RBC's and the Hematocrit or Packed Red Cell Volume (PCV).

## 3. Counting principle

The LC-220 presents the exclusive technical performance of counting the white and red blood cells in one unique chamber.

The counting principle is based on an impedance variation generated by the passage of cells (non conductive particle) through a calibrated micro-aperture of 80 micrometers. The sample is suspended in an electrolytic diluent (current conductor). The conductivity of the diluent differs considerably from the conductivity of the cells. This dilution is aspirated through the calibrated micro-aperture. Two electrodes are placed on each side of the aperture. Electric current passes through the electrode continuously. When the cell passes one by one through the aperture, it displaces its own volume of electrolyte. Electric resistance (or impedance) between the two electrodes increases proportionally with the cell's volume.

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
Fig.3 The overview of the needle system


Fig.4 The overview of the bloc pump


Fig.5 The overview of the air cylinder

While the number of pulses indicates red or white blood cell counts, the amplitude of the electrical pulse produced depends on the cell's volume. Each pulse is then sorted according to its size (Fig.2).

The hematocrit is measured by adding up all the pulse heights. A mathematical process is then applied to the sum obtained to compensate for coincidence passages in the aperture.


Fig.2 The impedance principle

We will prevent clogging of the 80 micrometers calibrated aperture with a specific reagent cleaner used during the automatic rinsing cycles. The cycle of the LC-220 will inject cleaner at every end of cycles to avoid pollution and clogging of the pneumatic circuit.

#### 4. Exclusive Specifications of the LC-220 Haematology Analyzer

The LC-220 uses only three microprocessor-controlled step by step motor to achieve the blood sampling, dilution and pressure or vacuum functions, all required for blood cell counting.

The needle carriage mechanism is to carry the blood sample (Fig.3). The first step by step motor is driving the carriage from the sampling position to the chamber with a brief stop for the inside rinsing of the needle.

The bloc pump (Fig.4) is composed of three pumps: the dilution pump, the lyse pump, and the blood sample aspiration pump. The LC-220 uses the second step by step motor to activate the three pumps together. This will guarantee precision and reliability.

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

A piston or syringe will replace the usual compressor to input pressure or vacuum into the pneumatic circuit. This pump will be activated up (vacuum) and down (pressure) by the third step by step motor triggered by a micro-processor. Additionally, the air cylinder (Fig.5) will drain the counting chamber, the vacuum in the counting head, and the backflush and bubbling. Reliability, precision, silent light and ease of use reflects from the choice of this air syringe.


The quality control (QC) data of LC-220 are shown in Fig.6 and Fig.7.


(A) Linearity range R&D systems full range Hgb

(B) Linearity range R&D systems full range WBC

(C) Linearity range R&D systems full range RBC


(D) Linearity range R&D systems full range Hct

Fig.6 Linearity


These highly sophisticated, unique technologies have yielded a reduction in both size and number of parts, and have assured high precision and reliability. The LC-220 is 300 (W) × 410 (H) × 360 (D) in overall size and is as light as 13kg (excluding the weight of reagent). It is designed for quiet operations.

The LC-220 is so constructed that measurement, rinsing, washing and other basic operations can be performed by use of three switches on the front panel. This analyzer requires no special training for effective handling.


The LC-220 technology has been realized by ABX and Horiba through their Franco-Japanese joint development efforts.


(A) Haemoglobin

(B) White blood cells


(C) Red blood cells

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Fig.7 Correlation between LC-220 and ARGOS

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
Brice PROUVE

ABX S. A.  
Marketing  
Business Development Manager

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### Charastaristics of Fully Automated Chemiluminescent Enzyme Immunoassay Equipment Luminomaster —Equipment Construction and Specifications—

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#### Abstract

Various problems are met in operating automatic analysis equipment in immunoassay: centralization of analysis equipment for specific applications, the reduction of assay result report time, saving operator labor, and the assurance of data reliability.

HORIBA has developed the fully automated chemiluminescent enzyme immunoassay equipment Luminomaster™ jointly with Sankyo Co., Ltd. The Luminomaster features a reduced immuno reaction time and higher sensitivity achieved by use of a chemiluminescent reaction, a wider concentration detection range, a random access function, automatic dilution and auto re-assay functions, to mention but a few. Also, it is expected that connection to an automatic sample transfer system will reduce labor in large hospitals, reduce assay time, and help prevent contamination from samples.

This paper mainly describes a configuration of the Luminomaster™ and its features.

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図13 秋田大学付属病院設置図  
Akita university affiliated hospital installation drawing

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
Characteristics of Fully Automated Chemiluminescent Enzyme Immunoassay  
Equipment Luminomaster —Reagents and Detection Systems—

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## Abstract

The automation of immunoassay is proceeding at a rapid pace, and various automatic analysis equipment is being developed and put on the market. Sankyo Co., Ltd. has developed jointly with HORIBA the fully automated chemiluminescent enzyme immunoassay equipment Luminomaster™ as highly sensitive immunity measuring equipment compatible with the trend towards non-RIA. This system comprises fully automated measuring equipment and special reagents. The chemiluminescent reaction is utilized in the detection system. As a result, highly sensitive detection allows the measurement time to be reduced and a wide dynamic range allows a wide range from low- to high-concentration areas to be measured.

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---

## Social Infrastructures and Measurement Technology

The concept "infrastructure", now being worth notice not only in social sciences but also in various fields of knowledge, will be effectively applied to techno-science of measurement. The water level gauge having been installed along the Nile since the age of ancient Egypt can be interpreted as a typical infrastructure for measurement. The mileposts on the highway, the balances at the square and the clocks on the tower are examples of similar infrastructures known from the old times.

A number of such infrastructures, found both in the history of techno-science and in the modern society, are discussed. Emphasis is placed on the legal metrology and the traceability system of measurement.

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## Expectations for a Laboratory Accreditation System for Assurance of the Reliability of Test Data and Mutual Acceptance

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### Abstract

HORIBA has developed the IR-100, a compact in-line gas monitor that continuously measures the concentration of silane gas ( $\text{SiH}_4$ ) and phosphine gas ( $\text{PH}_3$ ) used in semiconductor device production processes. The IR-100 uses a flashing light source and a dual-type 2-wavelength pyroelectric detector in its sensor. This allows the IR-100 to be made compact at 80 mm (W)×35 mm (D)×152 mm (H) and coupling with the gas piping system to be simply configured. The IR-100 not only improves the productivity of semiconductor processes by monitoring gas concentrations and implementing final inspection to ensure complete gas purging; it is also expected to help ensure safety in preventing improper connection to gas cylinders. This paper describes the specifications, measurement principle and main features of the IR-100.

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
図7 半導体用各種ガスの赤外吸収波長の帯域  
Infrared absorption wavelength bands of gases used in the production of semiconductor devices

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


## Abstract

The diversification of power supply is making the high-grade control of boiler water used at thermoelectric power plants all the more important. This paper describes the purpose of measurement and general control standards of pH, electric conductivity, dissolved oxygen, sodium ion, silica and hydrazine, which are indispensable in controlling the quality of boiler water, and the measurement principles of each of these components. It also introduces the water quality analyzers developed by HORIBA that are used in the control of water quality at thermoelectric power plants.

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
図1 一般的な貫流形ボイラの構成と各種の水質分析計  
Structure of a general-purpose percolation type boiler and various water quality analyzers

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### Abstract

The portable gas analyzer PG-200 series has been developed mainly for the analysis of combustion gases. The PG-200 series analyzes carbon monoxide and carbon dioxide by non-dispersive infrared absorption, nitrogen oxides by cross modulation chemiluminescence, and oxygen by zirconia method or galvanic cell method. By combining an analyzer and sampling unit, the PG-200 series is capable of measuring up to four components. Moreover, it comes in a compact, lightweight design. This paper describes the configuration of the PG-200 series and its features.

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