ISSN 0915-9916

Readout

HORIBA Technical Reports

March 2001 ■ No. 22

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Feature The Technology Alliance for Analytical Instruments

Foreword

The Challenge of a Global Business 3
 Andrew MACKINTOSH

Discussion

Business Strategy and Technology Alliance of the Measurement Instruments in the 21st Century 5
 Hidehito OBAYASHI / Kozo ISHIDA

Feature Article

Technology and Product Development of HORIBA's Laboratory Analytical Instruments for R&D and Quality Control
 — Focusing on X-ray analytical instruments — 12
 Kenichi OBORI

SEM/EDX-integrated Analysis System SEMEDX Series 15
 Toshikazu YURUGI, Sukehiro ITO, Yoshinori NUMATA, Keith SYKES

Versatile Multi-Element Analysis Using X-ray Fluorescence
 The MESA-500 Series and The MDX-1000 Series 19
 Akimichi KIRA, Neal ROBSON

Horiba Application Center Provides Solutions : the X-ray Analytical Microscope 26
 Katsumi NAGASAWA, Akihiro HIRANO, Sumiyo ISHIKAWA

Guest Forum

Measurement and Control Required in the Forefront of Ferroelectric Thin Film Research 31
 Hiroshi FUNAKUBO

Selected Article

Gas Phase Analysis of MOCVD Materials by FTIR 36
 Tsukasa SATAKE

The Evaluation of Liquid Material Vaporization and Supply System for CVD 40
 Tetsuo SHIMIZU

Triple Raman Laser Spectrometer RAMANOR T64000 46
 Shoichi OKAMOTO

Panel Discussion

What was obtained by imaging of pH? and What will be obtained? 51

New Product

Light-scanning Chemical Microscope SCHEM™ 60

Selected Article

The Soil Analysis using MS-200 Fast VOCs Analyzer 61
 Yoshihiko ARITA

TPNA-200 Total Nitrogen/Phosphorous Analyzer:
 toward the Enforcement of the Fifth Water-Quality Total-Volume Restriction 66
 Yasushi KITANO, Susumu YAMAUCHI

New Product

Ambient CD Gas Monitor FE-215 71

Column

EANET is participating in activities to diagonos and treat the earth 72
 LI HU

Patents 73

Technical Report List 74

HORIBA World-Wide Network 76

The Challenge of a Global Business

Oxford Instruments began life over 40 years ago in Oxford, England as the first commercial 'spin-out' from Oxford University. Today there are many high technology businesses founded on the basis of university research but in those days there were very few examples to follow. Right from the start we knew that our small island could not sustain a large enough market for our ambition. Hence since our earliest days, the management vision, focus and direction have been aimed at the business opportunities outside the natural geographic borders of Great Britain and we export 87% of our \$300 million annual turnover. As example of our commitment to our worldwide markets we have more than 50 employees in Japan.

As our markets have grown they also matured. For all businesses like ours, the rate of technological change and investment required to compete at the leading edge of technology are increasing each year. This means that while our global vision remains valid and we remain one of the UK's most successful exporters, we need to seek new ways of refreshing that vision and enlist new partners to help us achieve ongoing success.

As part of this process, Oxford Instruments has been through some of important internal changes in its history during two years. These have all been aimed at fulfilling our objective of improving our service to customers in a sustainable way. This way enable our business to be a stronger and more flexible partner in a range of new alliances which will be a part of our international future.

The alliance between Horiba and Oxford Instruments can be understood within this context. Even rival companies can agree on cooperation for mutual benefit and this can happen in the field of analytical instruments just as usefully as in large markets which get more space in the newspaper.

An example is an important alliance between Simens and Oxford Instruments in the field of medical instruments, which has been running successfully over 10 years.

The relationship between Horiba and Oxford Instruments has developed over a long period of time. Both businesses have allowed for a long-term, maturing process. We have found that our cultures are similar, and our staffs at any different levels have put a lot of work in to understand and improve relationship with each other.

This is very important. I myself lived in Japan in my youth, and worked at Kyoto University caring out post-doctoral research where I developed close working relationships with Japanese colleagues. Since then I have come to Japan every year and have met many business leaders. Some say that doing business with customers in Japan is more difficult than in other countries, and it is true that the means to achieve success may be different markets. But I have found that a key factor for success in Japan is to build up good relationships through regular personal communication over a long period time. For this reason I always encourage direct conversation between people, even through e-mails may often easier!

The relationship between Horiba and Oxford Instruments is growing as a result of cooperation in various ways. The collaborative activities started in the field of X-ray microscopes, developed by Horiba, where there was no direct competition between our business and Oxford Instruments helped to investigate the market for these product outside Japan. The next stage was to look at an area where we had some overlapping technology. The consequence of this was to transfer a line of X-ray fluorescence products in an OEM relationship with Horiba. More recently and based on a

Dr. Andrew MACKINTOSH

Chief Executive

Oxford Instruments plc.

growing mutual confidence and desire to expand the relationship, our companies began to cooperate in the apparently competitive field of energy dispersive x-ray analyzer.

We have gone through these stages over rapid of years and at all times we have had the commitment of senior manager in each business. In particular our Oxford Instruments Japan office, which has been successfully promoting our business in Japan for 10 years, has played an important role in helping to redefine our strategy in Japan in the context of our developing relationship with Horiba.

We have been happy to find that Horiba's management and staff, like Oxford Instruments, also have a distinguished history of success through collaboration with overseas companies. We both come from islands with a culture of looking outside the natural boundaries for new markets and opportunities. This gives us confidence to continue and develop the relationship, which will only succeed if we continue to invest time and effort in improving the mutual understanding of our countries and business cultures. We appreciate the opportunity to learn from Horiba.

There is a well-known proverb: "A frog in a well does not know of the ocean". This alliance is helping to take us outside the well by helping to develop new ways of doing business. I believe that new business models such as these in this new century will help both Horiba and Oxford Instruments achieve profitable growth and sustainable leading positions in our chosen markets.

I appreciate this opportunity to express our views.

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Taking Japanese swords as an example, American metallography scientist C. S. Smith once said, "If the Japanese people had devoted all their energies to science and the Europeans had been outstanding engineers, the history of metallography would have been quite different."

At the draw of a new century, in the midst of a period of accelerated globalization, Dr. Hidehito OBAYASHI of Hitachi, Ltd., and Dr. Kozo ISHIDA of Horiba, Ltd. have discussed on the perspectives of measurement businesses and technology alliances.

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Technology and Product Development of HORIBA's Laboratory Analytical Instruments for R&D and Quality Control

— Focusing on X-ray analytical instruments —

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Abstract

"Ultra Quick Supplier :UQS" is the key word for Horiba's cooperate activities. We recognize that our most important role as a manufacturer of measuring instruments is to supply products and service quickly to customers based on a clear grasp of the market needs of 21st century which is changing drastically. This paper introduces how is developing products and technologies for analytical instruments for R&D and quality control, focusing on X-ray instruments. It also shows some examples in which alliance is effective to achieve UQS.

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
SEM/EDX-integrated Analysis System SEMEDX Series

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Abstract

In the recent laboratory, a scanning electron microscope (SEM) and an energy dispersive X-ray analyzer (EDX) are used in combination in most cases. In the past when functionally independent SEM and EDX were combined for analysis, operations were complicated and the accumulation of know-how was required for each instrument. Horiba, Ltd., Hitachi, Ltd., and Oxford Instruments plc have joined together with their own technical specialties to make a SEM/EDX-integrated analysis system SEMEDX Series, which "promises simple operation, speedy clear SEM observation, and accurate EDX analysis to any operators." In this paper I will explain its feature functions.

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
Keith SYKES

Oxford Instruments plc.

Versatile Multi-Element Analysis Using X-ray Fluorescence The MESA-500 Series and The MDX-1000 Series

Akimichi KIRA, Neal ROBSON *

* Oxford Instruments plc.

Abstract

Recent years have witnessed rapidly growing needs for the precise elemental analysis of steel, ceramics, organisms, foods, and other materials in various fields, and considerable attention has been focused on X-ray fluorescence element analyzers. As the market for these analyzers expands, demand is also growing for new features that allow easy access by any user in addition to the excellent features of X-ray fluorescence spectroscopy, which enable non-destructive analysis with minimal or no sample preparation. Horiba, Ltd. has released the MESA-500 series of element analyzer products, focusing on the ultimate target of usability based on EDX. Oxford Instruments plc has put on the market the MDX-1000 series of element analyzers, which pack the features of both EDX and WDX in a compact instrument. This article provides an introduction to the features and applications for each of these product series.

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1 Introduction

XRFS (X-ray fluorescence spectroscopy) refers to analytical methods characterized by simple or no sample preparation and the non-destructive, high-speed qualitative and quantitative analysis of many kinds of elements in a wide range of concentrations.

XRF systems are grouped into two major types: wavelength-dispersive X-ray fluorescence (WDX) and energy-dispersive X-ray fluorescence (EDX). In WDX, X-rays are irradiated toward a given sample in a vacuum or helium-filled atmosphere. The generated characteristic X-rays are separated into each components by the monochromator, and elements in the sample are detected by a proportional counter or NaI scintillator. Analysis must be done for each individual element. With EDX, on the other hand, generated characteristic X-rays

can be guided directly to a semiconductor X-ray detector where the energy levels are separated, so many elements can be analyzed simultaneously.

Highly sensitive WDX analyzers, which also require handling expertise, are used mainly for high-precision element analysis. Compact, easy-to-handle EDX analyzers, on the other hand, are used mainly as versatile multi-element analyzers.

Horiba has put on the market various EDX-type X-ray fluorescence element analyzers under the "MESA Series" name, and has established a wide customer base in many fields. Oxford Instruments has released to the market multi-dispersive X-ray fluorescence element analyzer products (MDX Series) which features both WDX and EDX packed incorporated into one instrument. These two companies are planning to expand the range of new applications for their X-ray fluorescence element analyzers. Table 1 shows the features of EDX, WDX, and MDX.

Principle	EDX	WDX	MDX
Elements	Na-U	Be-U	C-U
Measurement mode	Simultaneous	Sequential(Scanning type) or Simultaneous(Fixed type)	Simultaneous(Fixed type)
Analysis area	φ 1-10mm	φ 10-35mm	φ 20mm
Detection limit *1			
heavy matrix *1-1	100-1000 ppm	5-50 ppm	10-100 ppm
intermediate matrix *1-2	10-100 ppm	1-10 ppm	2-20 ppm
light matrix *1-3	1-20 ppm	0.2-2 ppm	0.5-5 ppm
Utility			
Cooling water	unnecessary	required	unnecessary
Liquid nitrogen	unnecessary or required	unnecessary	unnecessary
PR gas	unnecessary	required	required

*1 Detection limit in measuring of Ti-Cu (depend on each instrument)

*1-1 Typical heavy matrix material : Steel

*1-2 Typical intermediate matrix material : Cement

*1-3 Typical light matrix material : Plastics

Table1 The features of EDX, WDX, and MDX

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2 The MESA-500 Series


2.1 Instrument Configuration

MESA-500W is an energy-dispersive X-ray fluorescence element analyzer with the capability of analyzing element type and components in solid, powder, liquid, or any other sample form, at high speed and precision and without damaging the sample. The MESA-500W also requires no standard samples.

The MESA-500W analyzer consists of an analyzer unit, a data processing unit, a vacuum pump, and a computer. Air-cooled X-ray tube suited for EDX spectroscopy, a high-purity silicon X-ray detector that requires no supply of liquid nitrogen during storage, a power supply for the X-ray tube, and a high-precision current and voltage control circuit for the power supply are compactly arranged in the analyzer unit. The data-processing unit is composed of a high-rate pulse-processing circuit that measures detector signals, and a user interface. The computer is pre-installed with a Windows® operating system and dedicated MESA-500W software, which provides easy access to spectrum analysis, automatic qualitative analysis, and quantitative analysis using fundamental parameters and calibration methods. Fig.1 shows the the MESA-500W and its system configuration.

a) MESA-500W

b) System configuration of the MESA-500W

Fig.1 Energy-dispersive X-ray fluorescence element analyzer MESA-500W

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2.2 Features and Functions

The MESA-500W analyzer incorporates various features that enable the user to perform high-precision, high-efficiency analytical work with simple operations, in addition to the fundamental features of EDX spectroscopy that are used to analyze solid, liquid, or powder samples. As a result, any operator can perform qualitative and quantitative analyses of a wide range of elements without any special preparation of unknown samples.

(1) Quantitative analysis without standard samples

The MESA-500W analyzer can perform precise quantitative analysis with no preparatory work such as calibration using standard samples and the pre-registration of a standard spectrum. It can separate overlapping peaks and the backgrounds of the X-ray spectrum and calculate the precise intensity of the fluorescent X-rays using the non-negative-constraint and non-linear least square method, which utilizes the energy response functions of the detector. The MESA-500W also has the capability to conduct the quantitative analysis of a wide range of elements.

(2) Automatic qualitative analysis

Press the automatic qualitative analysis button, and the MESA-500W analyzer will display the element symbols of elements likely to be contained in a given sample, with the symbols shown at their corresponding positions in the fluorescent X-ray spectrum. The contained elements will also be indicated in a periodic table. This feature helps operators with no special knowledge of the X-ray spectrum to quickly identify elements in the sample.

(3) Automatic switching of X-ray tube anode voltage

Spectral data measured at the high and low anode voltages of the X-ray tube is used in the calculations in the fundamental parameter method, thus enabling the precise quantitative analysis of samples which contain both low and high atomic number elements.

(4) Automatic adjustment of X-ray tube anode current

This feature automatically adjusts the anode current of the X-ray tube to the best value for a given sample. As generation efficiency of fluorescent X-rays differs greatly between metal and organic film sample (the best currents can vary by a factor of more than 100 between the two), the anode current must be set to the best value for the composition of each sample. This feature controls the anode current to set the dead time of the pulse circuit to the best value, and adjusts the intensity of the X-rays that are incident upon the detector.

(5) Automatic calibration

The MESA-500W analyzer also features automatic calibration using the maintenance sample included with the analyzer. This feature will calibrate the horizontal axis (energy) and the vertical axis (X-ray intensity) of the X-ray spectrum and will at the same time calculate the energy resolution of the detector used for spectrum analysis. It is recommended that the analyzer be calibrated once every two months.

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2.3 Application

The MESA-500W analyzer has been used in various fields and applications, such as customs inspection and quality management at semiconductor plants. Fig.2 shows various applications of the MESA series.



Fig.2 Various applications of the MESA series

The following is a description of how the MESA-500W was used to survey the remains of a copper smelting facility. In the remains of the copper mines in Taka County of Hyogo Prefecture, the MESA-500 was used to analyze the elements in more than 150 samples, including kiln pieces found in the upper earth layer (dating to the middle of the 19th century) and furnace pieces found in the lower earth layer (dating from the end of the 16th century to the first half of the 17th century). The instrument was used to analyze the composition of 5 primary elements, iron (Fe), silicon (Si), aluminum (Al), manganese (Mn) and potassium (K), as well as the secondary elements sulfur (S), arsenic (As), tin (Sn), copper (Cu) and lead (Pb), in the soil, slag, road stones, and vein stones left in the furnaces and tuyeres. It was also used to classify the smelting and to analyze the refining processes. The following findings were reported. Fig.3 shows the measured result of a clod slag using MESA-500W.¹⁾

- * Flat or clod slags, mineral powder or sauce, and clay obtained from the outer and inner walls of kilns and furnaces, can be grouped by the composition ratios of Fe, Al, and Si.
- * A high level of arsenic (As) is contained in samples obtained from the kilns. (This indicates that As was oxidized and removed in the kilns.)
- * Levels of Ca differ greatly by groups, indicating that Ca was used selectively as a slag forming reagent.
- * Different compositions of furnace-wall clay suggest that smelting processes ranging from calcination to crude copper production were used in the time period corresponding to the lower-layer, and processes ranging from calcination to the smelting of silver and copper were used in the time period corresponding to the upper-layer.

a) Clad slag

b) A screen of the measurements

Fig.3 Measured result of a clod slag using MESA-500W¹⁾

3 The MDX-1000 Series

The compact, fully integrated MDX (Multi-Dispersive X-ray fluorescence) series of spectrometers provides the capacities for both EDXRF and WDXRF analysis in a single instrument. (Fig.4) This combines the flexibility and range of EDXRF with the higher resolution (for lower atomic number elements) and speed of WDXRF. Modular construction enables a choice of configurations to meet different needs.

Fig.4 Multi-Dispersive X-ray fluorescence element analyzer MDX-1000

3.1 Instrument Composition

The core of an MDX Spectrometer is the X-ray tube, multiple detection channels, and digital pulse processor (Fig.5).


Fig.5 System configuration of the MDX-1000

Simultaneous WDXRF provides the fastest possible analysis for routine analysis of known elements, with low detection limits. The MDX-1060 can be fitted with multiple monochromators for simultaneous WDX analysis of up to 12 elements. Each detection channel, or monochromator, consists of a collimator, crystal and proportional detector selected and set for the characteristic X-rays of a particular element. For example, the monochromator for sulfur analysis consists of a collimator, germanium crystal and argon methane detector which together are set to a "Bragg angle" of 110.7 degrees (2 theta) corresponding to a (sulfur) wavelength of 0.5373 nm.

For flexibility, or analysis of unknown elements (e.g. contaminants), the MDX-1080 can be configured with a solid state Energy Dispersive channel alongside the fixed channels, which uses EDXRF analysis to collect the spectrum of up to 50 elements, allowing rapid qualitative and quantitative analysis. For example, a semi-quantitative analysis of unknowns can be made automatically from the data library of element reference values. The combined technique allows each system to operate in its optimum range, i.e. wavelength dispersive for the measurement of low atomic number elements (Mg, Al, Si, S etc), where resolution is paramount, and energy dispersive for medium/high atomic numbers, where low backgrounds are particularly important.

3.2 Features and Functions

The use of a digital pulse processor with the energy dispersive channel allows high count-rates (up to 100,000 cps), reducing analysis time and increasing sample throughput.

The compact medium-power X-ray tube, which uses a rhodium target, underpins the performance of the whole system, operating continuously at 200 watts. It is positioned below the sample to allow the analysis of liquids. (A secondary safety window can be fitted to avoid any sample leakage entering the spectrometer.) An automatic, integral temperature control maintains the spectrometer at 37°C, without the need for an external cooling system, or its associated plumbing.

Sample preparation is minimal. Liquids are poured into a sample cell, powders are pressed into pellets, and metals only require grinding or turning to a flat surface.

The flexibility of the MDX spectrometers is made extremely accessible to a wide range of users through the *XpertEase* Windows®-based user interface. All functions, including selection for single or multiple samples, analysis of unknowns, and method set-up, are requested using easily identified click-on buttons. Fig.6 shows a screen for setting up semi-quantitative analysis of unknowns.


Fig.6 A screen for setting up semi-quantitative analysis of unknowns

Some of the features available include:

- * Sample spinner for inhomogeneous sample analysis.
- * Unattended operation - using the autosampler, up to 72 preloaded samples may be analysed without operator intervention.
- * Multitasking - new data may be entered during automatic sample analysis.
- * Priority interrupt - high priority samples may be inserted for individual analysis during automatic runs.
- * Password control - different users may be given different levels of access.

3.3 Application

A typical application of an MDX-1000 spectrometer is the determination of levels of sulfur (S), nickel (Ni), and vanadium (V) in fuel oil, the thick residual material resulting from the refinery distillation of crude oil. Fuel oil is the primary component for all grades of petrol (gasoline), diesel fuels for both motor vehicles and railway locomotives, aviation fuels, and finally heavy residual oils for marine diesel engines and power stations. During combustion, nickel and vanadium can form compounds that are corrosive to metal. At trace levels in petroleum they can deactivate catalysts during processing. Levels of sulphur are regulated because of its environmental impact. Using the Oxford Instruments' MDX-1060 wavelength dispersive X-ray fluorescence spectrometer, fitted with four monochromators one used for the automatic subtraction of background signals, the rapid, simultaneous analysis of these three elements gives excellent results (Table 2).

Analyte	Conc. range	Calibration standard error	Limit of detection (3sigma)	Mid-range precision (95%confidence)
S	0.1-1.0%	0.002%	0.002%	0.007%
Ni	0-50mg/kg	0.4mg/kg	1mg/kg	<1mg/kg
V	0-50mg/kg	0.5mg/kg	1mg/kg	<1mg/kg

Table 2 Specifications of the MDX-1060

4 Conclusion

This paper covers the specifications and some applications of the MESA-500 series analyzers. The ultimate goal in the design of these analyzers is to achieve simplicity, usability, and a compact instrument body for EDX spectroscopy. We have also looked at the MDX-1000 series spectrometers, which incorporate features of both EDX and WDX spectroscopy. These two series can perform non-destructive element analysis, a feature of X-ray fluorescence element analyzers, and they each have their own specialized functions and software features as well. By entering into a business alliance, Horiba and Oxford Instruments now have even greater capacity for providing products optimized to their customers' needs. We shall continue efforts to promote communications with our customers so as to be able to meet their on-going and changing needs.

Reference

- 1) Research of the history of copper produced in Harima area, Investigation Committee of the ruins around Mount Myoken (July 4, 1999)

Akimichi KIRA

Scientific Instruments R&D Dept.
Horiba, Ltd.

Neal ROBSON

Marketing Manager-XRF
Oxford Instruments plc.

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Horiba Application Center Provides Solutions:the X-ray Analytical Microscope

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Abstract

Recently, the request from a customer to the analytical instruments maker is not only to supply the good products but also to offer the know-how for analysis. An inexperienced engineer will be able to get the precise and reliable data using these know-how. The Horiba Application Center has been endeavor to accumulate the analytical know-how for various samples using Horiba's products. In this paper we will report roles of the Application Center through the customer service of the XGT-2700W X-ray analytical microscope that is very useful for elemental analysis.

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Coffee break 1

An American's View of Science and Technology

Americans think of themselves as practical people. They hold practical results in high regard and are perhaps less knowledgeable and respectful of theoreticians. Practical results can capture the imagination; theory is abstract and more difficult to understand. As the level of education of the average American improves, this attitude is changing.

Before the mid-twentieth century, most Americans recognized the practical and economic benefits of inventions. Inventors such as Morse, Bell, Marconi, Edison, and Ford were celebrated in the popular press because their work provided benefits to the common man. The terms "scientist" and "technologist" were not yet in common use.

During the Second World War, Americans saw the results of enormous projects that combined the theoretical and the practical. Theoretical and experimental physics were combined with electrical, chemical, and mechanical engineering to create new weapons and countermeasures. Ordinary Americans began to see the power, chiefly military, of combining science with technology.

Following the war, successful organized efforts to conquer many diseases made celebrated heroes out of doctors, researchers, and medical scientists. Later, the U.S. and the Soviet Union engaged in the competition known as the "space race." When an American team walked on the moon in July 1969, Americans saw one of the peacetime benefits of science and technology.

Today, technology relies on pure and applied research conducted by scientists. Large technology-based organizations are increasingly turning to partnerships with universities and institutions to assure that vital research is conducted by the best-qualified people. Science and technology are not viewed as separate disciplines, but instead as a continuum. A "revolving door" commonly connects the academic and the commercial worlds and scientists and technologists move easily between the two.

The new century may bring an extended period of peace for the world. In a world characterized by economic competition and political cooperation, Americans hope that nations will combine their science and technology to benefit all of mankind in many ways we cannot imagine today.

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
Measurement and Control Required in the Forefront of Ferroelectric Thin Film Research

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"What is the state of the art in the research and development of the FRAM, ultimate memory device?". This simple question was the subject of a lecture delivered by Dr. Hiroshi FUNAKUBO of the Tokyo Institute of Technology.

Why is the FRAM considered the dream memory device in the next generation? Which will be preferred as its material, PZT or SBT? Which method will be used to manufacture the thin film, PVD or MOCVD? What will be the conclusive factor for manufacturing at low temperature? Furthermore, what will be the future problems in measurement and control technology which support and lead these issues?

Releasing new findings one after another with the daily advancement of the R&D for ferroelectric thin film, Dr. Funakubo talked about expectations toward measurement instrument makers and said that "a global alliance of researchers and engineers beyond the boundaries of their specialties is the key for surviving the 21st century."

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Abstract

The ferroelectric thin film is going to be put to practical use as a semiconductor memory device of the next generation. Although various kinds of material and deposition technique for this film are proposed now, the optimum deposition conditions are not yet established. Then we have measured infrared spectrum of three gases that was supplied as source materials for MOCVD using Horiba's gas analyzer, the FT-730. And we have found that some intermediate products is formed in evaporation / mixture chamber. It is expected that FTIR is a very useful tool in the research laboratory or the production process of FRAM.

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
The Evaluation of Liquid Material Vaporization and Supply System for CVD

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Abstract

Semiconductor manufacturing technologies are diverse and fast-changing, and none more so than thin film formation technology. The MOCVD is one of the most feasible thin film formation technique for ferroelectrics thin film for FRAM, capacitor film, insulating films and other high-functional thin film. Currently various substances are proposed as source materials for thin film. As these materials become more diverse, the technology in which the liquid material efficiently and stably vaporizes becomes very important. This paper reviews the various liquid material vaporization and supply methods, and presents the vaporization test results for some material using the direct injection method. And automatic tetra etoxy silan delivery system that developed by STEC Ltd. is introduced.

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Abstract

Raman spectroscopy is a very effective method for the analysis of the formation and structure of materials. As the intensity of a Raman emission is very weak, however, a special high-performance optical system is required. The Jobin Yvon S.A., with an outstanding record in optics-related technology, has developed various Raman spectrometers and put them on the world market. In this paper I will introduce the RAMANOR T64000 triple laser Raman spectrometer with three monochrometers. The RAMANOR T64000 provides selection of the best measuring mode out of three to meet the properties of a given sample. Here I will show you its applications for Raman spectrum measurements of YBCO's high-temperature superconductors and hydrogen chloride gas, stress measurements of α -alumina.

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A French View of Science and Technology

Although science is a versatile word, it is generally regarded as the knowledge of rules, governing phenomena, facts, or objects, verifiable by experimental methods. Strictly speaking, it is necessary to distinguish the several levels of science between fundamental science and applied science. The fundamentalist studies nature as a simplified abstraction, in order to discover the fundamental laws of nature. The experimentalist describes the actual behavior of nature in its complexity. An experimentalist is close to application in the sense that he often conducts research in the perspective of applications. However, applying science at a technological level, by creating commercial products, is usually the role of the engineer. Engineering is not limited to science, but includes organizing and managing production units, and thinking of the product in commercial terms.

One specificity in French view of science comes from an education system, several hundred years old, that separates and opposes two institutions: universities and "ecoles." French universities are similar to those in Japan in terms of their educational goals—a high level of scientific knowledge that is kept up-to-date by the size and the diversity of the institution. University graduates are typically researchers and professors. However, the university often does not provide commercially-oriented engineering courses. Such courses are offered by ecole institutions, in which education is very specialised and technically oriented, with an emphasis on business skills.

This division of educational roles into distinct institutions is not strictly defined, but often overlapping. It is common to see universities integrating management courses in their programs, creating their own engineering schools, or conducting research programs sponsored by industry. On the other hand, ecoles attract many students who are focused on fundamental science. This is a result of a system of selection in which the main criteria is the prestige of the school. It is common for engineering students to attend university courses during or after the completion of their ecole studies.

This traditional competition between the two institutions drastically enhances the contrast in the French imagination between what is applied, and what is fundamental science. On one hand, some university scholars readily disregard any technical aspects in science for the sake of abstraction. On the other hand, some engineering graduates no longer consider themselves to be scientists, but consider to be generalists, with purpose to take up at executive positions. The middle ground of applied science, which is experimental but not industrially-oriented, is classically the weakest part of the system.

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
What was obtained by imaging of pH? and What will be obtained?

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"New tools open up a new world." Since Danish biochemist S. P. L. Sorensen defined the index of hydrogen-ion concentrates as $\text{pH} = -\log_{10}[\text{H}^+]$ in 1909, pH has been used widely as the most fundamental chemical index not only in electrochemistry but also in metal science, life science, and all other fields. In 1998 Horiba had developed the SCHEM™ light-scanning chemical microscope to measure and display two-dimensional pH distribution as image mapping.

"What information do we get by imaging of pH?" "What applications do we have?" "What perspective will be in the future?" To answer these questions, researchers engaged in the R&D of a pH imaging sensor and researchers using the SCHEM™ have been invited to Horiba and discussed their recent findings and future perspectives.

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
The Soil Analysis using MS-200 Fast VOCs Analyzer

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Abstract

In recent years, after the revelation of trichloroethylene-contaminated groundwater and the unlawful dumping of toxic waste, analytical services, construction companies, and other related businesses are taking active steps toward cleaning polluted and contaminated air, soil, and water and the rehabilitation of weakened environmental functions. Horiba has developed the MS-200 fast VOCs (Volatile Organic Compounds) analyzer based on the time of flight mass spectrometry. With features of high speed, sensitivity, and portability, the MS-200 can measure concentrates of VOCs in a wide range of processes from pollution surveys to purification. In this paper, I will report some application data from, and discuss future assignments on, soil analysis and the MS-200.

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Yoshihiko ARITA

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Coffee break 3

An English View of Science and Technology

The terms "science" and "technology" are to a large extent identical. "Science" comes from the Latin *Scientia*, that from *Scire* - "to know." Science traditionally includes an array of disciplines which aspire to rational, objective thought. What is, and what is not science has been a long-running debate between the "hard" sciences of physics, chemistry, and biology, and pseudo-scientific disciplines such as psychology and philosophy. This debate is fuelled by a critique of the levels of scientific methodology, but a better defining feature of a science is whether it has technology. Technology is the result of the applied sciences and is derived from the Greek *Tekhnologia*.

Nowadays technology has become almost synonymous with science. However, in England subtle distinctions remain and a discernible hierarchy continues to persist. The use of the term science implies a creative skill which is of comparable status to the arts. However because English society is dominated by the arts, science is often differentiated from an art because it involves technical skills and technology. This is an erroneous distinction given the subjectivity of technical skill, but it indicates the conceptual closeness of modern science and technology.

In most western countries, technology is not only synonymous with science but also with engineering. In England the theoretical tends to enjoy a superior social status to the practical. This prejudice places the scientist above the engineer. The original medieval Latin *Engeniare* means "engine." An English dictionary definition of an engineer is of a person who makes or maintains machines. In England we have heating engineers and telephone engineers who would be termed technicians in most other countries. This confusion between the technician who merely maintains a machine, and is often not university educated, and the highly educated engineer remains, and diminishes the status of the professional engineer. This confusion extends to whether someone is called an engineer or a scientist. Accordingly, it is common in England for a person who has a doctorate in engineering to be referred to as a scientist in the press because of the superior status this title affords.

Ultimately the distinction between science and technology in England is one grounded upon academic and cultural snobbery.

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TPNA-200 Total Nitrogen/Phosphorous Analyzer:
toward the Enforcement of the Fifth Water-Quality Total-Volume Restriction

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

Abstract

Conventionally, the total volume of drainage has been regulated, having used COD as the index in order to prevent eutrophication in the lake and sea areas in Japan. Our water quality has not improved yet, however, as seen in frequently occurring red tides and water blooms. In the Fifth Water-Quality Total-Volume Restriction scheduled to take effect in the spring of 2001, nitrogen and phosphorous, two causes of eutrophication, are going to be added newly to the list of restricted items, and the need is quickly growing for an easy-to-use low-cost total nitrogen/phosphorous analyzer. To meet these needs, Horiba has already developed the TPNA-200 that can simultaneously measure total nitrogen/phosphorous concentrates. This paper will introduce the outline of the Fifth Water-Quality Total-Volume Restriction and the features and measurement results of the TPNA-200 analyzer.

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Susumu YAMAUCHI

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Coffee break 4

A Japanese View of Science and Technology

The Japanese have put trust in science and technology for many years. In the last days of the Tokugawa shogunate, we studied them ardently and applied them to practical uses. The Meiji government sent many students to Western countries to learn various branches of science and technology. The Meiji government also invited many men of distinguished talent from Western countries to encourage the advancement of science and technology. When Dr. Hideki Yukawa was the first Japanese person to be awarded the Nobel prize in physics, we were greatly encouraged at the news. Many citizens formed a deep attachment to science and technology. With the aggravation of our environment, accidents at nuclear reactors, and occurrences of quality problems at private industries, however, we have become aware of the fact that science and technology have a dark side too.

Ladies and gentlemen, think of it like this. Neither science nor technology can speak out. They will not assert themselves nor justify themselves. We human beings have used them as we want. Recently science and technology have not been received favorably, but this is because we human beings have made mistakes in dealing with them. If we neglect our duties, science and technology will strike back without fail. It is said that a good design engineer is a good communicator with science and technology. I agree, too.

Let's leave the detailed definitions of science and technology to critics and government officials. I think like this: If we human beings can converse with nature (by which I mean science and technology) to our satisfaction and such conversation reaches the level of art, we will be able to harvest a lot of fruits from nature. It is not through our words that we can understand nature. It is through our efforts to pursue the truth that we can understand nature. In this sense, there is no boundary in nature. If we understand nature (that is, science and technology) with an appreciation of different cultures and ideologies of different nations and ethnic groups, I believe we could bear a true alliance of technologies.

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KOREA ATAGO SPECTROSCOPY CO.,LTD
HORIBA INSTRUMENTS (SINGAPORE) PTE. LTD.
HORIBA LTD. BEIJING REPRESENTATIVE OFFICE