
Readout

HORIBA Technical Reports

No.40
MARCH 2013

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We hope this issue of HORIBA Readout, Featuring the application technology in analysis and measurement in HORIBA Gr.

I visited Hokkaido, a northern island, looking for pure scenery. A few Japanese red-crowned cranes flying across the unstepped snowfield seem to indicate the advent of spring.
-Photographer Hideo MATSUI-
(Member of NIKA Association of Photographers)

Name of the book

This book is named "Readout" in the hope that "the products and technology we have created and developed will be read out and so become widely known".

Application technology in analysis and measurement

Review

- 6** Engine Emission Measurement Technologies Related to Emission Regulation Trends
Kaori INOUE
- 12** Aiming at Development of New Applications and Professionalism.
Bunji HASHIMOTO

Feature Article Pittcon

- 21** Development of On-Site Measurement of Calcium in Food Samples by Use of Ion Selective Electrode
example for Using <LAQUAtwin> Ca²⁺ Ion meter
Hisashi YAMANOUCHI
- 27** Potentiometric Determination of pH Values of Low Conductivity Solutions with the Glass Combination Electrode Equipped with Ionic Liquid Salt Bridge
Manabu SHIBATA
- 30** Portable X-ray Fluorescence Analyzer MESA-50 with New Developed Optics and New Spectrum Processing
Tomoki AOYAMA
- 38** Measurement of Automobile Exhaust N₂O by Mid-IR Laser Spectroscopy
Kenji HARA, Montajir RAHMAN

Feature Article Application

- 42** Non-Destructive Evaluation of SiC wafer for Power Device
Defect detection in μm~nm scale using optical analytical technique
Nobuyuki NAKA, Tomoya SHIMIZU, Hiroki BABA, Tomoko NUMATA, Junichi AOYAMA, Shogo AWATA, Tomoaki HATAYAMA
- 49** Analytical Evaluation for Research & Development and Manufacturing process of Lithium ion battery
Jun HIROSE
- 55** Development of an Immunoaffinity Column "SmartColumn" for Mycotoxin Analysis
Mikiko UCHIGASHIMA

Selected Article

- 59** Accurate Determination of pH by Use of Ionic Liquid Salt Bridge
Manabu SHIBATA
- 64** Atomic Resolution Holography with SEM: Internal-Detector Electron Holography
Akio UESAKA, Kouichi HAYASHI, Tomohiro MATSUSHITA, Shigetoshi ARAI
- 69** Study of the Improvement of Calculation Method for Particle Size Distribution Using a Laser Diffraction Method.
Tatsuo IGUSHI

Product Introduction

- 74** Integrated Emission Measurement Instrument Platform
HORIBA ONE PLATFORM
Katsumi URATANI
- 81** Development of Transient EGR Ratio Analyzer Utilizing Heated NDIR
Tomoshi YOSHIMURA
- 86** New Drive Recorder Systems (DRT-7000, DR-2000, HsE)
Masatomo ISHIKURA, Yasuyuki SADO, Akihiro HASHIMOTO

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Engine Emission Measurement Technologies Related to Emission Regulation Trends

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Recently, drastic changes in engine/vehicle emission regulations lead to complicated requirements for emission measurement technologies. One of the emission regulation trends is introducing limit values to new components/categories, e.g. particle number (PN) in EU countries and nitrous oxide (N_2O) as a greenhouse gas in US. To make such new standards effective, measurement techniques to be adopted must be carefully investigated and validated. Test procedures and measurement conditions defined in regulations have been revised in response to increasing interest in real-world emission behavior and diffusion of new-generation vehicles such as hybrid electric vehicles. This article gives an outline of some emission regulation trends, from the view-point of relation with measurement technologies.

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
Figure 1 Schematic of Full Flow Dilution System for Particulate Mass Emission Measurement

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
Figure 2 Schematic of Recommended Particle Sampling System by PMP^[3]

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Figure 3 Example of Vehicle Test Cell with a Chassis Dynamometer

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
Figure 4 Example of Driving Cycle for Light-duty Vehicle Test(Japanese JC08 mode)^[7]


Figure 5 Example of Control Area for Off-cycle Emission Test(WNTe)^[8]

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
Figure 6 Battery State of Charge during Driving a Plug-in Hybrid Electric Vehicle

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- [1] Regulation(EC)No 715/2007 of the European Parliament and of the Council of 20 June 2007
- [2] Regulation(EC)No 595/2009 of the European Parliament and of the Council of 18 June 2009
- [3] UNECE Regulation No. 83 rev 4(April 26, 2011)
- [4] Federal Register Vol. 75, No. 88, 25324-25728(May 7, 2010)
- [5] Federal Register Vol. 76, No. 179, 57106-57513(September 15, 2011)
- [6] Federal Register Vol. 77, No. 199, 62624-63200(October 15, 2012)
- [7] 国土交通省告示第1268号(平成18年11月1日)
- [8] Global Technical Regulation No. 10(September 9, 2009)
- [9] Federal Register Vol. 70, No. 113, 34594-34626(June 14, 2005)
- [10] Commission Regulation(EU)No 582/2011 of 25 May 2011
- [11] 中央環境審議会：今後の自動車排出ガス低減対策のあり方について(第十次答申)(平成22年7月28日)

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Aiming at Development of New Applications and Professionalism.

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Everyday activity of the application center includes: development of the new characterization methods for particular instruments, demonstration for visiting customers, sample analysis and user's trainings. On the other hand, academic activity includes: scientific articles submission, patent applications and presentations at various conferences. Recently, it became very difficult to satisfy customer's needs, by following the same way. To cover the customer's demand for versatile characterization, it became necessary to introduce application approach. In this review, application center's present activity and the vision for the near future are presented, specifically emphasizing the development of application approach and professionalism of the application engineers.

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Figure 1 Journal of The Analytical Society for Analytical chemistry "Bunseki"

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
Figure 3 Dependence of voltages generated by the calcium ion meter on the hydrogen ion activity.


Figure 4 Relationship between measurement value using by calcium meter and ICP. Four samples were measured with pretreatment(◇)and without pretreatment(□).

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Feature Article

Pittcon

Development of On-Site Measurement of Calcium in Food Samples by Use of Ion Selective Electrode example for Using <LAQUAtwin> Ca²⁺ Ion meter

Hisashi YAMANOUCHI

The concentration of calcium in food samples is commonly measured by use of Inductively Coupled Plasma (ICP), Ion Chromatography (IC) and Atomic Absorption Spectrophotometry (AAS). However, these methods are not appropriate for on-site measurement at production lines and farms because these instruments are too large to carry out. In this paper, we propose new measurement method which dissociates the compound including calcium into the calcium ion adjusting pH of samples before the potentiometric measurement with the compact calcium ion meter (abbreviated as Ca²⁺-meter). The results obtained by use of Ca²⁺-meter with pretreatment were good agreement with that obtained by use of ICP. We can measure all amount of calcium in milk products by use of Ca²⁺-meter with pretreatment which generate free calcium ions from the calcium compound in milk products.

Introduction

Calcium is one of the most important elements for our body. As a major mineralization of bone and teeth, calcium constitute a few percent of the body's weight. It has an important role as a signal for many cellular processes such as nerves work and muscle movement. We can take it into our body by eating which milk products and fish have much calcium..

The concentration of calcium in food samples is commonly measured by use of Inductively Coupled

Plasma (ICP), Ion Chromatography (IC) and Atomic Absorption Spectrophotometry (AAS).^[1] However, these methods are not appropriate for on-site measurement at production lines and farms because these instruments are too large to carry out and they need for complicated pretreatment. Recently, for on-site measurement of calcium, the minimization and simplification of instrument is demanded. An ion selective electrode is one of the suitable sensors for the on-site measurement of calcium.

A calcium ion selective electrode (abbreviated as Ca²⁺-ISE) consists of a working electrode and a reference electrode. The cell voltage between two electrodes varies depending on the concentration of calcium ion in the solution. The potentiometric measurement by use of Ca²⁺-ISE is cheaper and smaller than that used for above mentioned methods. However the Ca²⁺-ISE is not sensitive to calcium binding protein and various organic acids in food because an ion selective electrode can respond only dissociated ions in aqueous solution. Accordingly, when we determine total amount of calcium in food by use of


Figure 1 Illustration of the calcium ion meter (B-751)

Ca^{2+} -ISE, the special pretreatment which makes the calcium in the compound into a free calcium ion is needed.

In this paper, we propose new measurement method which dissociates the compound including calcium into the calcium ion adjusting pH of samples before the potentiometric measurement with Ca^{2+} -ISE. By use of this method, we can easily measure the concentration of calcium and the measured values correlate well with that obtained by use of ICP. In addition, we introduce the compact calcium ion meter (abbreviated as Ca^{2+} -meter, Horiba, Ltd, LAQUAtwin B-751) which newly released in September 2012.

Experimental

Sample.

Milk products contained much calcium were prepared. Two milks and two lactic acid were kept at room temperature (22~24 °C) during the measurements.

Instrument and standard solution.

Figure 1 illustrates the structure of Ca^{2+} -meter which consists of the calcium ion sensor, the potentiometer and the display. The ion selective membrane and liquid junction of reference electrode were located on the same plane surface whose thickness was about 0.8 mm. This instrument enables us to measure small volume sample (about 300 μL).

The ion meter was calibrated with two standard solutions before measurements. The 150 ppm and 2000 ppm calcium standard solutions (3.74 mmol dm^{-3} and 49.9


Figure 2 Illustration of milk product sample after pretreatment

mmol dm^{-3} at CaCl_2) were used for the two point calibration. 0.1 mol dm^{-3} KCl was added into the standard solutions as supporting electrolyte. The intercept (zero point) and slope for the calibration curve of the cell voltage against the concentration of calcium ion were determined from two point calibrations. Compared with the Ca^{2+} -meter, the concentration of calcium in same samples was measured with ICP (Horiba, Ltd., ACTIV-M ICP).

Pretreatment of sample

Milk products contain calcium bonding casein protein and lactic acid compounds. The pH of these sample was controlled to 2 with adding strong acid in order to break these bonds, and then it was adjusted to among 4.3 to 4.6 by diluting with ion-exchange water and adding appropriate amount of Tris powder. There are two reasons why we control pH in sample. First, the interference of hydrogen ions causes the measurement errors.^[2] Second, the free calcium ion dissociated from calcium compound may recombine casein protein and lactic acid when the pH value in sample is more than 4.6.

Supernatant liquid of samples after pretreatment was measured. The details of pretreatment procedures are as follows.

Pretreatment procedures

1. Preparation of 5 mL sample.
2. Adding of 100 μL of 5 mol dm^{-3} HCl to sample.
3. Dilution by adding ion exchange water to 50 mL.
4. Additoin of 0.05 g of Tris (hydroxymethyl) aminomethane to diluting sample to control pH4.3~4.6.
5. Keeping on leaving diluting sample for a few minutes. (Figure 2) The casein protein precipitated in pH 4.3~4.6 because the isoelectric point of casein is 4.6.
6. Removing 5 mL of supernatant liquid from diluting sample.
7. Dissolving 0.0375 g of KCl in 5 mL of supernatant liquid.

The concentration of calcium in samples was calculated by multiplying experimental values obtained by use of Ca^{2+} -meter by ten because the tenfold dilution was made of the sample through the pretreatment. In the case of the measurement with ICP, the thousandfold dilution was made of sample by ion exchange water.

Results and Discussion

Figure 3 shows the difference cell voltage (ΔE_{pH}) dependence on pH with ion meter at 1 mmol dm^{-3} CaCl_2 and 0.1 mol dm^{-3} KCl. The difference between E values at

other pH CaCl₂ solution except pH 7.5, $E_{pH(X)}$ and at pH 7.5 CaCl₂ solution, $E_{pH7.5}$, $\Delta E_{pH} = E_{pH(X)} - E_{pH7.5}$, is given in Figure 3. The measurement of E at each pH was repeated with three different ion meters. Experimental values were obtained from the average of E values and the range of error bar shows the difference between maximum and minimum values for three measurements. The pH of each solution was controlled to among 2 to 12 with adding HCl or NaOH.

The ΔE_{pH} value was independent on pH in the range of pH 5~12. The ΔE_{pH} value changed due to the interference of hydrogen ion on the Ca²⁺-ISE in a solution whose pH was less than 4. The change of pH from 4 to 2 in sample decrease the cell voltage by about 60 mV which corresponds to the 1/100 dilution of the concentration of calcium ion. Therefore, pH of samples must be controlled above 4 to prevent the interference of pH on the Ca²⁺-ISE. Figure 4 shows comparison of experimental values obtained by use of Ca²⁺-meter and ICP at four milk products. Measurements by use of Ca²⁺-meter by with


Figure 3 Dependence of voltages generated by the calcium ion meter on the hydrogen ion activity.


Figure 4 Relationship between measurement value using by calcium ion meter and ICP. Four samples were measured with pretreatment (◊) and without pretreatment (◻).

pretreatment (◻) and without pretreatment (◊) were carried out to confirm the effect of pretreatment on measurement of calcium. The results obtained by use of Ca²⁺-meter with pretreatment were good agreement with that obtained by use of ICP. The correlation coefficient, r² was 0.952. Experimental values obtained by use of Ca²⁺-meter without pretreatment were less than a half of that obtained with the pretreatment.

Thus, we can measure all amount of calcium in milk products by use of Ca²⁺-meter with pretreatment which generate free calcium ions from the calcium compound in milk products.

Conclusion

By use of Ca²⁺-meter with pretreatment, calcium concentration in milk products were determined. This method enables us to measure more easily total amount of calcium in milk products on-site at production lines and farms than the expensive instrument such as ICP, IC and AAS. It is expected that this method will be widely used on-site measurement.

On-site measurement of calcium will be more required in the food industry because of the increase of interest in healthy, safety and peace of mind. We expect this method can be helpful in these demands.

References

- [1] JIS K 0400-50-20, JIS K 0400-52-30
- [2] J. RUZICKA, E. H. HANSEN, J. CHR, TJELL, 160-173

Hisashi YAMANOUCHI

Liquid and Water Quality R&D Dept.
Application R&D Center
Reserch & Development Division
HORIBA, Ltd.

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Feature Article

Pittcon

Potentiometric Determination of pH Values of Low Conductivity Solutions with the Glass Combination Electrode Equipped with Ionic Liquid Salt Bridge

Manabu SHIBATA

It is difficult to measure accurately and quickly pH of low conductivity solutions by use of pH glass combination electrodes based on a KCl salt bridge (KCISB). The reasons are the variation of the liquid junction potential (LJP) between a KCISB and a sample solution, and the change of pH in the sample solution due to the leakage of the concentrated KCl from the liquid junction. In order to solve the intrinsic problems, a new reference electrode, where an ionic liquid was used as the material of liquid junction, has proposed. The pH glass combination electrode equipped with an ionic liquid salt bridge (PUREIL electrode) enables us to measure accurately and quickly pH of less than 10 mS m^{-1} sample solutions such as pure water and boiler water. In this paper, we show the pH measurement of 10 mS m^{-1} sulfuric acid solution by use of PUREIL electrode.

Introduction

Accurate determination of pH is of fundamental importance in not only science and technology but also many facets of our life and environments. Potentiometry by use of an electrochemical cell that consists of a glass electrode and a reference electrode equipped with a salt bridge^{*1} based on a concentrated KCl solution (abbreviated as KCISB) has been used as a reliable and convenient method of pH measurements^[1,2]. However, the potentiometry with KCISB is not accurate enough for solutions of low conductivity. The reasons are the liquid junction potential (LJP) between a KCISB and a dilute solution and the change of ionic strength in sample solution due to the leakage of a concentrated KCl solution. Although the problem was pointed out in many studies^[3-16] conducted in 1970s-1990s, it remains unsolved. A new salt bridge made of an ionic liquid^{*2} (abbreviated as ILSB) was proposed as an alternative to the KCISB^[17-20]. For the working principle of salt bridge, you can refer to p. 59 in

this paper. In this paper, we report that the pH glass combination electrode equipped with ILSB (HORIBA, Ltd. PUREIL, 9600-10D) enable us to measure accurately pH of dilute solutions.

*1: Salt bridge: The third solution inserted between different two solutions in order to eliminate the liquid junction potential between the two solutions.

*2: Ionic liquid: the molten salt at the room temperature

Experimental

Reagents

The ionic liquid, tributyl (2-methoxyethyl) phosphonium bis(pentafluoroethanesulfonyl)amide (TBMOEPC₂C₂N) was obtained from Kanto Chemical Co., Inc. and used without further purification. $20 \mu\text{mol dm}^{-3}$ (10 mS m^{-1}) H₂SO₄ was prepared by diluting with MilliQ water a standardized sulfuric acid solution, which was certified to be $(5.00 \pm 0.01) \times 10^{-2} \text{ mol dm}^{-3}$ by coulometric titration with NaOH solution (Nacalai Tesque, Inc., Japan). A phosphate standard solution ($0.025 \text{ mol kg}^{-1}$ KH₂PO₄ + $0.025 \text{ mol kg}^{-1}$ Na₂HPO₄, pH

Figure 1 The picture of PUREIL electrode

= 6.865 ± 0.005 at $25\text{ }^\circ\text{C}$) were obtained from Kanto Chemical Co., Inc. A 0.05 mol kg^{-1} citrate buffer solution ($\text{pH} = 3.776$ at $25\text{ }^\circ\text{C}$) was prepared by dissolving 11.41 g of $\text{KH}_2\text{C}_6\text{H}_5\text{O}_7$ (Kanto Chemical Co., Inc. 99 %) in pure water and diluting it to $1.0000 \pm 0.0004\text{ dm}^3$.

ILSB-Type Combination Electrode (PUREIL Electrode)

Figure 1 illustrates the structure of a combination electrode which consists of a glass electrode and a reference electrode equipped with a ring of gelled ILSB. TBMOEPC₂C₂N was used as the ILSB. The structures of the cation and anion constituting TBMOEPC₂C₂N are given in Figure 2. TBMOEPC₂C₂N was gelled by dissolving 8 g of P (VdF-HFP) and 0.008 dm^3 TBMOEPC₂C₂N in 0.1 dm^3 acetone. The mixture was dried to remove acetone for one week at room temperature to obtain a membrane of 2.5 mm thickness. The ring-shaped membrane was cut out from the membrane. The ring-shaped membrane of the gelled IL was mounted with a silicone O-ring to the cylindrical body of the combination electrode^[21]. In the ILSB-type reference electrode, the inner cell was composed of a Ag/AgCl electrode in a 0.1 mol dm^{-3} KCl saturated with a TBMOEPC₂C₂N and AgCl. In the glass electrode, the inner cell was composed of a Ag/AgCl electrode in a 0.1 mol dm^{-3} KCl saturated with a AgCl and 0.04 mol dm^{-3} KH_2PO_4 + 0.16 mol dm^{-3} Na_2HPO_4 . The composition of the glass was the same as that of the glass electrode (HORIBA, Ltd. 9618).

pH measurement method

The electrochemical cell employed for the pH determination with the glass electrode and the ILSB-type reference electrode is represented as


Figure 2 Structure of TBMOEPC₂C₂N

I	II	III	IV	V	VI	VII	VIII	VIV
Ag	AgCl	0.1 mol dm ⁻³ KCl saturated with TBMOEP C ₂ C ₂ N and AgCl	gelled TBMOEP C ₂ C ₂ N	pH standards or x μmol dm ⁻³ H ₂ SO ₄	glass	0.1 mol dm ⁻³ KCl saturated with AgCl 0.04 mol dm ⁻³ KH ₂ PO ₄ 0.16 mol dm ⁻³ Na ₂ HPO ₄	AgCl	Ag (A)

The single vertical bar indicates the phase boundary, and the single dashed vertical bar indicates the liquid junction between two electrolyte solutions of different compositions.

The cell voltage, E , i.e., the potential of the right-hand-side terminal referred to that of the left in cell (A), was measured with a pH meter (Horiba, Ltd., F53, resolution: 0.1 mV) at a sampling rate of 0.3 Hz . The polypropylene container containing 20 dm^3 sample solution (V in cell A) was set in a water bath kept at $25.0 \pm 0.1\text{ }^\circ\text{C}$. Cell A was calibrated with two pH standard buffers before the measurement of H_2SO_4 .

Calculation of pH Values from Cell Voltage

When cell A is calibrated with a standard buffer, pH_s , an unknown pH value of H_2SO_4 solutions, pH_x , in V in cell A is written

$$\text{pH}_x = \text{pH}_s - \frac{[E_{\text{ex}} - E_s - (E_{j(x)} - E_{j(s)})]F}{RT \ln 10} \quad \dots \quad (1)$$

where E_{ex} and E_s are the readings of the pH electrode for the H_2SO_4 solution and the standard buffer whose pH values are pH_x and pH_s , respectively, $E_{j(x)}$ and $E_{j(s)}$ are the liquid junction potentials (LJPs) at ILSB | $x\text{ } \mu\text{mol dm}^{-3}$ H_2SO_4 and ILSB | the standard buffer solution interfaces, F is the Faraday constant, R is the gas constant, and T is the absolute temperature. Glass electrodes may exhibit the pH response smaller than the theoretical value, $RT \ln 10 / F$ volts per pH unit. The common procedure is to use the operative Nernst slope, k' , which is obtained from the two-point calibration of the electrode,

$$k' = \frac{E_{S1} - E_{S2}}{\text{pH}_{S2} - \text{pH}_{S1}}, \quad \dots \quad (2)$$

where E_{S1} and E_{S2} are the pH cell voltages for the buffers pH_{S1} and pH_{S2} , respectively. When k' is used, eq(1) is represented by

$$\text{pH}_x = \text{pH}_s - \frac{E_{\text{ex}} - E_s - (E_{j(x)} - E_{j(s)})}{k'} \quad \dots \quad (3)$$

If the ILSB works ideally, $E_{j(x)}$ is equal to $E_{j(s)}$ and eq(3) reduces to

$$\text{pH}_x = \text{pH}_s - \frac{E_{\text{ex}} - E_s}{k'} \quad (4)$$

The pH values of the sample solution were obtained from the measured E values with eq(4).

In order to make a comparison between the PUREIL electrode and the KCl type electrode, the pH of H_2SO_4 was measured by use of the KCl-type combination electrode (HORIBA, Ltd, 9681) according to the same procedure.

Result

Figure 3 shows the time courses of pH for 15 min at 10 mS m^{-1} H_2SO_4 solution. In Figure 3, ● and ■ are the pH values measured by use of PUREIL and KCl-type electrodes, respectively, and red line is the calculated pH value^[20] of 10 mS m^{-1} H_2SO_4 solution.

The response time of PUREIL electrode is shorter than that of KCl-type. In the case of PUREIL electrode, the difference between experimental and calculated pH values is within 0.03 pH unit. On the other hand, the pH value measured by use of KCl-type is higher by 0.1 pH unit than the calculated pH value. Two reasons are considered. First, in the case of KCl-type, the pH in sample solution changes with time due to the leakage of a concentrated KCl from liquid junction. Second, the LJP between a KCISB and the low conductivity solution is unstable. Thus, we can measure accurately and quickly pH of the low conductivity solutions with a PUREIL electrode since the contamination of a sample solution is less and the LJP between the ILSB and sample is stable.


Figure 3 The time course of pH at 10 mS m^{-1} ($20 \mu\text{mol dm}^{-3}$) sulfuric acid solution. ●: PUREIL electrode, ■: KCl-type electrode, red line: the calculated pH value

Conclusions

PUREIL electrodes enable us to accurately and quickly determine the pH values of the low conductivity solution whose pH values can not be accurately measured by use of a KCl-type electrode used over 100 years. In this paper,

although we measured the dilute H_2SO_4 solution, the accurate and high stable pH measurements of sample solutions such as rain water, river water, boiler water, and pure water are possible by use of PUREIL electrode. It is anticipated that the PUREIL electrode is certified as the standard electrode for the pH measurement of low conductivity solutions. Although we omitted details of the results due to space limitation, you can refer to our paper in preparation for these details^[22].

References

- [1] R. G. Bates, *Determination of pH*, Wiley, New York, 1973.
- [2] R. P. Buck, S. Rondinini, A. K. Covington, F. G. K. Baucke, C. M. A. Brett, M. F. Camoes, M. J. T. Milton, T. Mussini, R. Naumann, K. W. Pratt, P. Spitzer and G. S. Wilson, *Pure Appl. Chem.*, 2002, **74**, 2169-2200.
- [3] J. N. Galloway, B. J. Cosby and G. E. Likens, *Limnol. Oceanogr.*, 1979, **24**, 1161-1165.
- [4] S. Y. Tyree, *Atmos. Environ.*, 1981, **15**, 57-60.
- [5] G. Marinenko and W. F. Koch, *Environ. Int.*, 1984, **10**, 315-319.
- [6] A. K. Covington, P. D. Whalley and W. Davison, *Pure Appl. Chem.*, 1985, **57**, 877-886.
- [7] W. Davison and C. Woof, *Anal. Chem.*, 1985, **57**, 2567-2570.
- [8] W. Davison and M. J. Gardner, *Anal. Chim. Acta*, 1986, **182**, 17-31.
- [9] W. F. Koch, G. Marinenko and R. C. Paule, *J. Res. Natl. Bur. Stand (U.S.)*, 1986, **91**, 23-32.
- [10] R. C. Metcalf, *Analyst*, 1987, **112**, 1573-1577.
- [11] D. Midgley, *Atmos. Environ.*, 1987, **21**, 173-177.
- [12] W. Davison, A. K. Covington and P. D. Whalley, *Anal. Chim. Acta*, 1989, **223**, 441-447.
- [13] R. C. Metcalf, D. V. Peck and L. J. Arent, *Analyst*, 1990, **115**, 899-905.
- [14] R. A. Durst, W. Davison and W. F. Koch, *Pure Appl. Chem.*, 1994, **66**, 649-658.
- [15] T. Ozeki, Y. Tsubosaka, S. Nakayama, N. Ogawa and T. Kimoto, *Anal. Sci.*, 1998, **14**, 749-756.
- [16] R. Kadis and I. Leito, *Anal. Chim. Acta*, 2010, **664**, 129-135.
- [17] Kakiuchi, T.; Yoshimatsu, T. *Bull. Chem. Soc. Jpn.* **2006**, **79**, 1017-1024.
- [18] Yoshimatsu, T.; Kakiuchi, T. *Anal. Sci.* **2007**, **23**, 1049-1052.
- [19] Sakaida, H.; Kitazumi, Y.; Kakiuchi, T. *Talanta* **2010**, **83**, 663-666.
- [20] M. Shibata, H. Sakaida and T. Kakiuchi, *Anal. Chem.*, 2011, **83**, 164-168.
- [21] M. Shibata, M. Yamanuki, Y. Iwamoto, S. Nomura, H. Sakaida and T. Kakiuchi, *Anal. Sci.*, 2010, **26**, 1203-1206.
- [22] M. Shibata, M. Kato, Y. Iwamoto, S. Nomura, and T. Kakiuchi, in preparation.

Manabu SHIBATA

Liquid and Water Quality R&D Dept.
Application R&D Center
Reserch & Development Division
HORIBA, Ltd.
Dr. Eng.

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Portable X-ray Fluorescence Analyzer MESA-50 with New Developed Optics and New Spectrum Processing

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We have developed a portable XRF which has almost same performance as conventional desktop XRF. RoHS, soil and particulate matter are typical application examples for rapid and accurate analysis with portable XRF. Detection limits of several ppms have been achieved with Pb, As, Cd, Br, Se and Hg in plastics and soil. New developed optics has achieved highly accurate and rapid measurement. In addition, new spectrum processing method has reduced the effect of overlaps between peaks. Detail of new developed optics and spectrum processing method are presented.

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Measurement condition	
Sample	Cu alloy (ZBY921-925)
Voltage	50 kV
Current	Auto
Filter	5 element Filter
Process time	P2 (2 μ s)
Measurement time	100 s

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Table 3 MESA-50の仕様

Model name	MESA-50
X-ray tube	50 kV
Detector	SDD
Measurable elements	Al-U
Signal processor	Digital Pulse Processor
X-ray filter changer	4 position (including open)
X-ray bench (optional)	Atmosphere
Mass and Size	12 kg/202 (W) × 280 (L) × 178 (H)
Radiation area	7 mm, 3 mm, 1.2 mm (Automatic changer)
Sample image	CCD camera
Software	Reporting, Judge OK/NG
Correction algorithm	FPM Thickness, Material, Pb/As peak separation, Residual component estimation
Power supply	Ni-MH battery or AC adapter (100-240 V 50/60 Hz)
Battery life	6 hours (Ni-MH battery)

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
Figure 4 Total and Modal Mass Emission of N₂O

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- [1] Ballantyne, V., Howes, P., and Stephanson, L., "Nitrous Oxide Emissions from Light Duty Vehicles," SAE Technical Paper 940304, 1994, doi: 10.4271/940304.
- [2] IPCC/UNEP/OECD/IEA. Revised 1996 IPCC Guidelines for National Greenhouse Gas Inventories. Paris: Intergovernmental Panel on Climate Change, United Nations Environment Programme, Organization for Economic Co-Operation and Development, International Energy Agency(1997).
- [3] Environmental Protection Agency, "Electronic Code of Federal Regulation, Title 40, Parts 85, 86, 600, 1033, 1036, 1037, 1039, 1065, 1066, and 1068.
- [4] Environmental Protection Agency, "Electronic Code of Federal Regulation, Title 49, CFR Parts 523, 534, and 535.
- [5] Montajir, R., "Development of an Ultra-Low Concentration N₂O Analyzer Using Quantum Cascade Laser (QCL)," SAE Technical Paper 2010-01-1291, 2010, doi: 10.4271/2010-01-1291.
- [6] "Electronic Code of Federal Regulation, Title 40: Protection of Environment, Part 86: Control of Emission from New and In-Use Highway Vehicles and Engines, Subpart B, Section 86-132-00.
- [7] "Electronic Code of Federal Regulation, Title 40: Protection of Environment, Part 86: Control of Emission from New and In-Use Highway Vehicles and Engines, Subpart B, Section 86. 144-94.

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**Montajir RAHMAN**

HORIBA Instruments Incorporated,
Ann Arbor Facility
Ph. D

Feature Article

Pittcon

Measurement of Automobile Exhaust N₂O by Mid-IR Laser Spectroscopy

Kenji HARA

Montajir RAHMAN

Nitrous Oxide (N₂O) emission reduction has gained large prominence recently due to its contribution to the climate change as a greenhouse gas. The United States Environment Protection Agency (EPA) together with the United States Department of Transport (DOT) has already regulated the N₂O emissions from Light-Duty Vehicles (LDV). For LDV, N₂O measurement should be done from sample storage bags over the light-duty FTP drive cycles. N₂O emission standard of 0.10 g/bhp-hr for Heavy-Duty Engines (HDE) is also finalized. The final N₂O standard becomes effective in 2017 model year for diesel engines. Usually raw or diluted exhaust is measured for HDE emission testing. Therefore, an analyzer capable of measuring N₂O from bag and from diluted sample continuously is required to support both LDV and HDE regulations. In this study, comparative measurement of automobile exhaust N₂O in continuous dilute and in sample storage bag has been attempted. A modern vehicle meeting 2008 Japanese emission standards has been evaluated. In this letter, it is found that the laser based exhaust gas analyzer has sufficient detection capability for bag measurement and fast response for continuous dilute measurement.

INTRODUCTION

Nitrous oxide (N₂O) is a trace gas that exists naturally in the atmosphere with a level of about 300 ppb. It is the most common tropospheric nitrogen species aside from molecular nitrogen. Because it is a stable molecule, it has very long life time which has been recorded as being from 130 to 170 years.^[1] Since N₂O has an even greater propensity for absorbing infrared radiation than CO₂, there is concern that N₂O contributes significantly to global warming.^[2]

The US-EPA together with the DOT has already regulated the N₂O emissions from Light-Duty Vehicles (LDV) under a harmonized and consistent National Program. Under the national program, automobile manufacturers will be able to build a single light-duty national fleet that satisfies all requirements while ensuring that consumers still have a full range of vehicle choices. N₂O emission from light duty LDV has been limited to 0.010 g/mile.^[3] For LDV, N₂O measurement should be done from sample storage bags. Similar to the CO₂ standard approach, the

N₂O emission level of a vehicle would be a composite of the light-duty FTP cycles. The US-EPA is also finalizing a N₂O emission standard of 0.10 g/bhp-hr for heavy-duty engines (HDE).^[4] The final N₂O standard becomes effective in 2017 model year for diesel engines. Usually raw or diluted exhaust is measured for HDE emission testing. An analyzer capable of measuring N₂O from bag and from diluted sample continuously is required for supporting both the LDV and HDE regulations.

An instrument based on mid-infrared quantum cascade laser spectroscopy had been developed for measuring ultra-low level N₂O in automobile exhaust gas sampled in a sample storage bag.^[5] This instrument is a pulsed-mode high-resolution narrow band mid-infrared analyzer. It was reported that it has very low detection limit (<5 ppb), excellent linearity, and the response time (t₁₀-t₉₀) is shorter than 4 second (sample flow rate of 4 L/min). Negligible amount of interference from co-existing gases was also confirmed by using real engine exhaust.^[5] In this study, the instrument is applied for measuring N₂O emission from a modern gasoline vehicle that meets 2008 Japanese

emission standards. Correlation between N₂O measurement from sample storage bag hereafter called as “bag measurement” and continuous N₂O measurement from CVS hereafter called as “continuous dilute measurement” has been attempted.


Figure 1 Block diagram of the N₂O measurement system

SYSTEM CONFIGURATION

Figure 1 shows the block diagram of the analyzer. It mainly consists of a gas cell, a laser unit, a detector, and a vacuum pump. The gas cell is maintained at 50 °C and under high vacuum. The cell volume is about 500 ml where a total optical path length of about 30 m was achieved by multiple reflections between two mirrors in order to achieve low detection limit. The sample gas flow rate is maintained at about 4 L/min at room temperature for bag measurement in this letter.

VEHICLE TEST PROCEDURE

The test vehicle is equipped with a 1.5 L port injection gasoline engine and with dual-TWC after treatment


Figure 2 System schematic for vehicle test on chassis dynamometer with CVS

system. The tests were performed on a chassis-dynamometer under FTP test cycles under cold start and hot start condition. In the cold start test, the vehicle was preconditioned and soaked in similar way every time before the test according to the EPA recommended procedure.^[6] The hot start test is the test started immediately after completion of a cold start test without any preconditioning and soak. The schematic of the test system is shown in Figure 2. The sample storage bag was purged with nitrogen gas well in advance. In addition, sampling-to-measurement time was also kept constant in order to minimize the variation in bag environment under which transformation of NO_x components in the sample storage bag occurs.

Measurement of N₂O was attempted both from CVS for continuous dilute measurement and from sample storage bag. The dilution air was supplied from ambient through a charcoal plus HEPA filter unit. Diluted sample and dilution air was sampled in two separate bags for background correction. The full FTP driving cycle is separated into 3 phases; phase-1: Cold Transient, phase-2: Stabilized, and phase-3: Hot Transient. In between phase-2 and phase-3 there is a predetermined soak period. Therefore a total of six bags were sampled, 2 bags for each phase in the full driving cycle. Calculation of N₂O


Figure 3 Real-time N₂O emission of the FTP driving cycle

mass emission of was done according to EPA procedure.^[7] Only continuous dilution measurement was performed without any bag sampling in this case.

RESULTS AND DISCUSSIONS

Figure 3 shows the real-time N₂O emission in different phases of the FTP driving cycle under cold and hot start condition. The upper side graph shows the real-time N₂O emission in the cold transient phase (Phase-1) while the bottom side graph shows the real-time N₂O emission in the hot transient phase (Phase-3). The graph in the middle shows the real-time N₂O emission in the stabilized phase (Phase-2). There is a predetermined soak period between the stabilize phase and the hot transient phase.

Figure 4 shows the total and individual phase based mass emission of N₂O calculated from the bag measurement data. Calculation of mass emission was done according to the “40CFR Part 86 Subpart B.”^[6] It shows that during cold start condition most of the N₂O is emitted in cold transient phase of FTP driving cycle. Concentration of N₂O in diluted samples stored in Bag-1 during this phase is about 600 ppb which gives mass emission of about 4.50 mg/mile. However when the test is done under hot started condition, concentration of N₂O in diluted samples stored in Bag-1 during this phase is about 330 ppb which gives a mass emission of about 0.17 mg/mile.

Phase-2 (stabilized phase) shows no N₂O emission at all even if the vehicle is started at either cold or hot condition. Concentration of N₂O in diluted samples stored in Bag-2 during this phase is within 290 to 300 ppb while the ambient background is about 320 ppb. Therefore it gives almost zero mass emission according to the calculation procedure proposed in the “40CFR Part 86 Subpart B.”^[6] Concentration of N₂O in the Bag-2 may become significantly lower than dilution air background depending on the engine type and driving pattern. Phase-3 shows some positive N₂O mass emissions and hot start condition has almost the same N₂O mass emissions


Figure 4 Total and Modal Mass Emission of N₂O

as that of cold start condition. It shows that the total mass emission of N₂O from this test vehicle is about 0.95 mg/mile under cold start condition while it is about 0.14 mg/mile under hot start condition. The cold start emission from this test vehicle is one tenth of the EPA regulatory value.

CONCLUSION

A modern clean vehicle has been tested under FTP test cycles. A quantum cascade mid-infrared laser spectroscopic motor exhaust gas analyzer is used to measure the N₂O emission. The analyzer showed promising potential for ultra-low concentration N₂O measurement from bag as well as faster response for continuous dilute measurement. Results revealed the following:

1. Vehicle equipped with conventional port injection gasoline engine and TWC emits very low level of N₂O and concentration of N₂O in diluted samples stored in bag sometimes goes below ambient level. N₂O emission level of the test vehicle is one tenth of the regulatory value. Therefore, compliance test of such low emission vehicle requires very sensitive analyzer with low detection limit.
2. For vehicle equipped with conventional port injection gasoline engine and TWC, most of the N₂O is emitted during first 100 to 200 seconds of any driving cycle if started from cold condition. There is almost no N₂O emission when the vehicle is warmed up.
3. The mid-infrared quantum cascade laser spectroscopic analyzer used in this study has sufficient detection capability for measurement of very low concentration of N₂O from bag. The analyzer has fast response for continuous dilute measurement from CVS. There is excellent consistency between bag measurement and continuous dilution measurement by the analyzer used in this study.

References

- [1] Ballantyne, V., Howes, P., and Stephanson, L., “Nitrous Oxide Emissions from Light Duty Vehicles,” *SAE Technical Paper 940304*, 1994, doi: 10.4271/940304.
- [2] IPCC/UNEP/OECD/IEA. Revised 1996 IPCC Guidelines for National Greenhouse Gas Inventories. Paris: Intergovernmental Panel on Climate Change, United Nations Environment Program, Organization for Economic Co- Operation and Development, International Energy Agency (1997).
- [3] Environmental Protection Agency, “Electronic Code of Federal Regulation, Title 40, Parts 85, 86, 600, 1033, 1036, 1037, 1039, 1065, 1066, and 1068.
- [4] Environmental Protection Agency, “Electronic Code of Federal Regulation, Title 49, CFR Parts 523, 534, and 535.
- [5] Montajir, R., “Development of an Ultra-Low Concentration N₂O Analyzer Using Quantum Cascade Laser (QCL),” *SAE Technical Paper 2010-01-1291*, 2010, doi: 10.4271/2010-01-1291.
- [6] “Electronic Code of Federal Regulation, Title 40: Protection of Environment, Part 86: Control of Emission from New and In-Use Highway Vehicles and Engines, Subpart B, Section 86-132-00.
- [7] “Electronic Code of Federal Regulation, Title 40: Protection of Environment, Part 86: Control of Emission from New and In-Use Highway Vehicles and Engines, Subpart B, Section 86. 144-94.

Kenji HARA

Energy System Analysis R&D Dept.
Application R&D Center,
Research & Development Division
HORIBA, Ltd.
Ph. D

Montajir RAHMAN

Ann Arbor Facility,
HORIBA Instruments Incorporated,
Ph. D

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Non-Destructive Evaluation of SiC wafer for Power Device
Defect detection in μm ~ nm scale using optical analytical technique

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Power devices using Silicon Carbide (SiC) were already commercialized and stepping into mass production. But commercial SiC wafers have a problem which several types of defects still exist. Analysis techniques for a variety of defects are important because these defects make yield ratio of SiC device worse. In this work, we carried out complex analysis for sub-millimeter and nanometer scale defect and stress in epitaxial/bulk SiC wafer. Laser scattering method is able to detect defects of submicron size on a whole surface of wafer about 3 min / 4 inch wafer. Our original CL imaging system is capable of evaluating crystal defects of nanometer-size by sub-millimeter scale panchromatic image acquired in short time and Photoluminescence (PL) lifetime measurement helps us to understand carrier lifetime. Moreover, we reveal that Raman spectroscopy can evaluate stress distribution around defects detected in CL image non-destructively.

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
Figure 1 Schematic diagram of defect type in SiC wafer.

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
Figure 2 Apparatus and Optical configuration of particle/defect detection (PR-PD2).


Figure 3 Imaging CL system(DF-100)and light emitting by electron beam irradiation.

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
Figure 4 Apparatus of fluorescence life time measuring system (FluoroCube 5000U).

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Figure 5 Apparatus of Raman scattering measuring system (LabRAM HR Evolution).

Table 1 Main specification of SiC wafer (Size : 4inch).

Sample	Wafer type	Surface	Crystal axis	Dopant conc. (cm ⁻³)	Polished surface	Thickness (μm)
Sample A	4H-SiC Epitaxial film on bulk	Si	[11-20] 8° off	4 × 10 ¹⁴	One side	20
Sample B	4H-SiC bulk	Si	[11-20] 4° off	>1 × 10 ¹⁸	Both side	—
Sample C	4H-SiC bulk	Si	[11-20] 8° off	>1 × 10 ¹⁸	One side	—


Figure 6 Particle/Defect detection result of Sample B (Bulk wafer).


Figure 7 Particle/Defect detection result and CL image of Sample C (Bulk wafer).


Figure 8 CL image and CL spectrum for stacking fault of Sample A (Epitaxial wafer).

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


(a) Panchromatic CL image (b) Optical image after KOH etching
Figure 9 Defect evaluation result of Sample A (Epitaxial wafer).

(a) Panchromatic CL image (b) Optical image after KOH etching
Figure 10 Defect evaluation result of Sample B (Bulk wafer).

(a) Panchromatic CL image (b) Optical image after KOH etching
Figure 11 Defect evaluation result of Sample C (Bulk wafer).


Figure 12 Evaluation result in depth direction for stacking fault of Sample B (Bulk wafer).

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
Figure 13 Fluorescence spectrum of Sample A, B and C.


Figure 14 Lifetime curve of Sample A, B and C.


Figure 15 Raman spectrum of Sample B(Bulk wafer).

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
(a) Panchromatic CL image (b) Raman peak shift map

Figure 16 CL image and stress evaluation result of Sample B (Bulk wafer)

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Analytical Evaluation for Research & Development and Manufacturing process of Lithium ion battery

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To promote the performance or security of Lithium ion battery, it is very important to evaluate the property of its parts using analysis methods. In this paper, I introduce the useful analytical methods to solve the deterioration mechanism for Research & Development · Manufacturing process of Lithium ion battery and their applications.

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
Figure 1 Charge-discharge mechanism of lithium ion battery


Figure 2 Cross-sectional SEM and element mapping images of LiFePO₄ system positive active material.(F-K, Fe-L)
(a)SEM and normal element mapping images at 10 kV, (b)SEM and normal element mapping images at 3 kV, (c)SEM and element mapping images by TruMap at 3 kV.

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
Figure 3 Depth profile and cross-sectional element mapping image of fluorine(F-K) on LiFePO₄ system positive active material.

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
Figure 4 Analysis result of foreign material(iron powder in 35μm size)on separator using XGT
 (a)Optical image of separator, (b)Element mapping image in 50 mm square(100μm beam / measurement time: 30 min), (c)Focused element mapping image(10μm beam / measurement time: 10 min), (d)Spectrum(measurement time: 10 sec)


Figure 5 Elemental mapping image of laminated lithium ion battery negative electrode after charge-discharge cycle using XGT


Figure 6 Depth profile of LiCoO₂ system lithium ion battery negative electrode after charge-discharge cycle using GDS
(a)Depth profile from surface of negative electrode to collector electrode,
(b)Depth profile around surface of negative electrode


Figure 7 Raman spectrum changes of LiCoO₂ system positive active material with charge-discharge cycle


Figure 8 Cross-sectional Raman image and Spectrum of Li(Co,Mn,Ni)O₂ system lithium ion battery positive active material
(a)-1 Raman spectrum before charge-discharge cycle test,
(a)-2 Raman image based on A1g peak intensity before charge-discharge cycle test,
(b)-1 Raman spectrum after charge-discharge cycle test,
(b)-2 Raman image based on A1g peak intensity after charge-discharge cycle test

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
Figure 9 Surfacial Raman image and spectrum of LiCoO₂ system lithium ion battery positive electrode after charge-discharge cycle
(a) Raman image, (b)-1 Raman spectrum in red (Lithium cobalt oxide including Cobalt oxide), (b)-2 Raman spectrum in green (Carbon), (b)-3 Raman spectrum in red (Lithium cobalt oxide),

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
Figure 10 In-situ cell for Raman analysis

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Development of an Immunoaffinity Column “SmartColumn” for Mycotoxin Analysis

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“SmartColumn” products consisting in “AFLAKING” and “OCHRAKING” were developed for mycotoxin analysis. This product is an immunoaffinity column whose performance depends on the antibody used. The newly developed IAC was the best regarding organic solvent-tolerance at the time the product was launched. Both “AFLAKING” and “OCHRAKING” can be used for acetonitrile having a concentration of up to 20%.

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
Figure 2 Reactivity of the antibody with AF in ic-ELISA

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
Figure 3 Acetonitrile tolerance of MoAb2-3 reactivity with AFB1 in ic-ELISA : 1, 10, 20, 30, 40, and 50% acetonitrile. Acetonitrile percentages are the final acetonitrile concentrations in the competitive reaction mixture.

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
Figure 4 Acetonitrile tolerance of the prepared IAC : AFB1, AFB2, AFG1, and AFG2

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Table 1 Mean total TAF recoverie extracted with 90% acetonitrile from spiked food samples by means of the prepared IAC cleanup

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Table 2 Recovery of OTs in spiked food samples by the prepared IAC.

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Accurate Determination of pH by Use of Ionic Liquid Salt Bridge

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A salt bridge is a device indispensable in electrochemical measurement. Since the invention of that, a concentrated KCl solution has fulfilled the role over 100 years. The excellent performance of a concentrated KCl salt bridge (KCISB), however, does not always work satisfactory. The reason why a KCISB has been used for a long time is mainly the lack of better alternatives. It is anticipated that a newly proposed salt bridge based on moderately hydrophobic ionic liquids solves the intrinsic problems of a KCISB and bring a new dimension to not only pH measurement but also electroanalytical chemistry. In this paper, we introduce the fundamental property of an ionic liquid salt bridge and the application of it to pH measurements.

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Table 1 Liquid junction potentials between dilute solutions and a saturated KCl solution at 25 °C in millivolts.

Molarity	Equimolar Na acetate + acetic acid	Potassium hydrogen phthalate	Sodium acetate	HCl	KOH	KCl
10 ⁻²	3.20	3.49	3.23	2.85	1.92	2.78
10 ⁻³	4.15	4.06	4.21	3.97	3.22	3.93
10 ⁻⁴	5.00	4.87	5.27	4.77	4.48	5.02
10 ⁻⁵	5.80	5.78	6.29	5.69	5.75	6.10
10 ⁻⁶	6.72	6.71	7.23	6.70	6.88	7.07

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
Bis(trifluoromethanesulfonyl)amide ($C_1C_1N^-$)

Figure 1 Structure of $C_8mimC_1C_1N$


Figure 2 Schematic view of the ionic liquid-type reference electrode. 1, lead wire; 2, cap; 3, plastic tube; 4, 3.3 mol dm^{-3} KCl aqueous solution; 5, Ag/AgCl electrode; 6, O-ring; 7, gelled ionic liquid; 8, cylindrical cap


Figure 3 E of IL-type and KCl-type electrodes after 3 min for measurements of pH standard solutions. open circle, IL-type; open square, KCl-type.

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
Figure 4 Illustration of the electrochemical cell using a ILSB sandwiched by two hydrogen electrodes. 1: glass tube for introducing hydrogen gas; 2: hydrogen electrode; 3: glass tube for exhausting hydrogen gas; 4: silicon rubber stopper; 5: U-type glass cell; 6: 500 μmol dm⁻³ H₂SO₄; 7: x μmol dm⁻³ H₂SO₄; 8: TBMOEPC₂C₂N

Table 2 The experimental and calculated pH value of 20-200 μmol dm⁻³ H₂SO₄.

Molarity of H ₂ SO ₄ / μmol dm ⁻³	Molarity of H ₂ SO ₄ / μmol kg ⁻¹	m_{H^+} / μmol kg ⁻¹	γ_{H^+}	pH _{cal}	Mean pH _{ex} ± 95% confidence interval	pH _{ex} -pH _{cel}
20	20.06	40.04	0.9910	4.401	4.410 ± 0.021	0.009
50	50.15	99.85	0.9860	4.007	4.008 ± 0.018	0.001
100	100.30	198.87	0.9804	3.710	3.713 ± 0.026	0.003
150	150.45	297.12	0.9762	3.538	3.543 ± 0.020	0.005
200	200.59	394.67	0.9727	3.416	3.421 ± 0.029	0.005

Table 3 Experimental pH values obtained by use of ILSB and pH values by use of a Harned cell or KCISB at reference buffer solutions studied at 37 °C.

buffer solution	pH _{Harned}	pH _{KCISB}	Mean pH _{ex} ± 95% confidence interval	pH _{KCISB} -pH _{Harned}	pH _{ex} -pH _{Harned}
1:3.5 phosphate	7.297	7.267	7.282 ± 0.012	-0.030	-0.015
1:3 Tris	7.427	7.381	7.395 ± 0.015	-0.046	-0.032
1:2 HEPES	7.661	7.631	7.613 ± 0.008	-0.030	-0.048
1:2 TES	7.535	7.503	7.484 ± 0.017	-0.032	-0.051

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Atomic Resolution Holography with SEM: Internal-Detector Electron Holography

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Internal-Detector Electron Holography, which is one of the Atomic Resolution Holography, can be proved with Scanning Electron Microscope. With Internal-Detector Electron Holography, we can analyze the local structure around the specific atoms in the laboratory. In this paper, I explain the theory of Internal-Detector Electron Holography and the measurement results of a bulk sample (SrTiO₃) and a thin film sample (Pt).

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
Figure 1 Sample Surface irradiated by Electron

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
Figure 2 Principle of the Atomic Resolution Holography.
(a)Normal Method(b)Inverse Method.

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
Figure 3 $f(E)$ of Ti-K α for the electron energies.^[14]


Figure 4 $f(E)$ of Pt-M α for the electron energies.^[14]


Figure 5 The adequate thickness for various energies and IMFP.^[14]


Figure 6 Measurement setup.

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
Figure 7 (a) Measured Hologram (AccV=6.00 kV). (b) Simulated Hologram (6.00 keV, 25 Å).^[4, 5]


Figure 8 Reconstructed atomic image with measured holograms.^[4, 5]


Figure 9 Measured Hologram. (a) AccV=4.00 kV (b) AccV=10.00 kV.


Figure 10 Reconstructed atomic image with measured holograms.

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Study of the Improvement of Calculation Method for Particle Size Distribution Using a Laser Diffraction Method.

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The paper is a summary of thesis titled “Study of the Improvement of Calculation Method for Particle Size Distribution Using a Laser Diffraction Method”. A motivation of the study is to investigate the influence of optical model and inversion problem on results of the particle size distribution by the laser diffraction method. The final goal of the study is to improve accuracy of particle size distribution by laser diffraction method. the effect on the kinds of inversion algorithm and these parameters, the number of detectors with same optical models, and theoretical calculation of fundamental uncertainty region based on the maximum and/or the minimum size in the preparation of standard reference particles were researched in the study, respectively. As the results, the research is extremely useful for improving the accuracy of measurement with new design of optics and inversion algorithm for the laser diffraction method.

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
Figure 1 (a)Schematic block diagram of instrument for the laser diffraction method.
(b)Schematics of a Fourier optics used in laser diffraction method and array detector.

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
Figure 2 Schematic of the computer simulation process.

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
Figure 3 Comparison of simulated inverted particle size distributions. (a) Expected particle size dispersion, (b) results of the conventional method, and (c) results of the iteration method.

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
Figure 4 Relation between the square of residual between calculated and expected signal and the number of detectors by computer simulation. The given particle size distribution (PSD) is the log-normal distributions at (a) standard deviation (SD) of 0.27, and (b) SD of 0.67, respectively.


Figure 5 Relation between the number of detectors and the error ratio of calculated and expected signal with log-normal distribution of mono disperse (SD of 0.01) with bimodal PSD. The first mode of mean diameter at (a) 0.546 μm with SD of 0.01 and (b) 48.1 μm with SD of 0.01, respectively.

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
Figure 6 Measured particle size distribution of silica glasses beads with uncertainty region (solid lines). Open circles show count base distribution measured by TEM images. Filled circles show volume base distribution calculated from count base distribution. Open squares show PSD of the silica particles measured by laser diffraction method.

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- [1] K. S.Shifrin, V. F.Truchin, L. S.Turovtseva, and V. A. Gashko, *Atmos. Oceanic. Physics.*, **8**, 739(1972).
- [2] R. Xu, "Particle characterization: Light scattering methods", Kluwer academic publishers, Dordrecht(2001).
- [3] Y. Mori: *Huntai to Kogyo*, **38 B**, 35(2006)
- [4] S. Tohno, *Earozoru Kenkyu*, **12**, 281(1997).
- [5] ISO13320: 2009 International Organization of Standardization, Geneva(2009).
- [6] T. Igushi, and H. Yoshida: *Rev. Sci. Inst.*, **82**, 015111-1(2011)
- [7] T. Igushi, and H. Yoshida: *Adv. Part. Tech.*, **23**, 682(2012)
- [8] T. Igushi, and H. Yoshida: *Rev. Sci. Inst.*, **83**, 055103-1(2012)
- [9] H. Yoshida, T. Igushi, T. Yamamoto and H. Masuda: *Adv. Part. Tech.*, **22**, 43,(2010)
- [10] Proceeding of the 8th international conference on measurement and control of granular materials (MCGM2009), Shenyang, China, 110(2009)

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Integrated Emission Measurement Instrument Platform
HORIBA ONE PLATFORM

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HORIBA ONE PLATFORM is the integrated measurement instrument platform for exhaust gas measurement system, which includes exhaust gas analyzer, constant volume sampler and other instruments. This has been released in May 2012 at the same time with MEXA-ONE. MEXA-ONE is the successor of MEXA-7000 which has been recorded more than 5000 sales in the world. MEXA-ONE is one of the instrument application created on the HORIBA ONE PLATFORM. MEXA-ONE application was developed using the concept of HORIBA ONE PLATFORM and takes advantage of its functionalities and features. This integrated measurement instrument controller provides full support for daily operation such as auto start-up sequence, etc. And furthermore, HORIBA ONE PLATFORM can provide higher test cell efficiency and availability by integrating preventive maintenance functions and self-diagnostic function such as monitoring core parameter of analytical devices, etc.

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Wider Application/Flexible System

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
Figure 1 Plug in/out Operation screen


Figure 2 Multi-point Measurement(MEXA-ONE+3*OVN+QL-NX)


Figure 3 Multi-point EGR Measurement(MEXA-ONE+2*EGR)

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
Figure 4 Integration device(s) & Data Management

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
Figure 5 System Control-Auto Sequence Operation for Daily work

Figure 6 Status View for System Device(s)


Figure 7 Self Notification-Periodic inspections/adjustments

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
Figure 8 Self Notification-Replacement of consumable parts


Figure 9 Precaution-Flow rate Monitoring


Figure 10 Self Precaution-Analyzer Sensitivity Monitoring

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

Figure 11 Extension for Trouble Shooting

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
Figure 12 HORIBA ONE PLATFORM Architecture Overview

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- ① Open Gui Platform(以後OGP)
- ② Open Application Platform(以後OAP)

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Development of Transient EGR Ratio Analyzer Utilizing Heated NDIR

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Most of the recent clean diesel engines are generally equipped with exhaust gas recirculation (EGR) technology in order to meet the strict criteria of NO_x emission regulations. In the field of emission control, accurate and fast transient EGR ratio operation is becoming very critical. In this study, a transient EGR analyzer has been developed and evaluated to confirm the advantage of the transient EGR analyzer by comparing it to a conventional exhaust gas analyzer system. A faster response has been observed from the transient EGR analyzer without significant transient error during transient test cycles. Furthermore, there was a high correlation between the EGR valve behavior and the EGR ratio measured by the transient EGR analyzer. These results show the advantages of the transient EGR analyzer when measuring the EGR ratio during the transient cycle.

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
Figure 2 Schematic diagram and general installation of conventional system

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
Figure 3 Schematic diagram of new EGR analyzer


Figure 4 Analyzer response comparison




Figure 5 Comparison under step change of EGR valve

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Table 1 Response time of each line

Analyzer	New EGR Analyzer		Conventional System	
Line	Intake	Exhaust	Intake	Exhaust
T_0-T_{10}	3.4	3.4	11.1	12.5
$T_{10}-T_{90}$	2.0	2.0	2.1	2.4
T_0-T_{90}	5.4	5.4	13.2	14.9
T_0-T_{99}	8.4	7.8	24.0	25.6

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


Figure 7 EGR measurement comparison under ETC test

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New Drive Recorder Systems
(DRT-7000, DR-2000, HssE)

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Personal standards on road safety vary among individuals. Previously there were no definite road safety standards; as a result, many evaluations have been done on the basis of personal senses or common belief. As a solution to these problems, Drive Recorder^{*1} have been developed to obtain quantitative records of driving, and systems for data evaluation have accordingly been pursued. In this report we introduce the most recent models of these devices and systems.

*1: Drive Recorder: An Drive Recorder is a device to record driving data and picture images on a memory card or other storage media when triggered by an impact shock or sudden braking caused by an accident or a near-miss incident. The driving data contains records on acceleration, braking ,speed, turn signals and GPS while the picture images shows the moments before, during and after the shock.

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Reno Office

605 Spice Island Drive, #5, Sparks, NV 89431
Phone : (1)775-358-2332 Fax : (1)775-358-0434

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HORIBA INSTRUMENTS (SHANGHAI) CO., LTD

No.200, Taitao Rd, Anting Town, Jiading District, Shanghai 201814

Phone : (86)21-6952-2835 Fax : (86)21-6952-2823

HORIBA (China) Trading Co., Ltd.

Unit D, 1F, Building A, Synnex International Park, 1068 West Tianshan Road, 200335, Shanghai

Phone : (86)21-6289-6060 Fax : (86)21-6289-5553

Beijing Office

Room 1801, SK Tower, Tower 1, No. 6 Jia, Jianguomenwai Ave., Chaoyang District, Beijing 100022

Phone : (86)10-8567-9966 Fax : (86)10-8567-9066

Guangzhou Office

Room 1611/1612, Goldlion Digital Network Center, 138Tiyu Road East, Guangzhou 510620

Phone : (86)20-3878-1883 Fax : (86)20-3878-1810

Shanghai Service Center

Rm.303, No.84, Lane887, Zu-Chong-Zhi Rd., Zhangjiang Hi-tech Park, Shanghai, 201203

Phone : (86)21-5131-7150 Fax : (86)21-5131-7660

Shanghai Technical Center

No.200, Taitao Rd, Anting Town, Jiading District, Shanghai, China 201814

Phone : (86)21-6289-6060 Fax : (86)21-6289-5553

INDIA

HORIBA India Private Limited

246, OKHLA INDUSTRIAL ESTATE, PHASE 3 NEW DELHI - 110020

Phone : (91)11-4669-5001/(91)11-4646-5000

Fax : (91)11-4669-5010/(91)11-4646-5020

Bangalore Office

Kamadhenu, No.17/1 - 32, Bannerghatta Road, Audugodi, Bangalore - 560030

Phone : (91)80-22210071

Pune Office

502, 5th Floor, Purushottam Plaza, Baner Road, Baner, Pune - 411045

Phone : (91)20-40766000 Fax : (91)20-40766010

INDONESIA

HORIBA Instruments (Singapore) Pte Ltd.

Jakarta Office

Menara Bidakara 2 Unit 11-04, Jl. Jend. Gatot Subroto Kav. 71-73, Jakarta Selatan, 12870

Phone : (62)21-2906-9419/(62)21-2906-9420

Fax : (62)21-2906-9421

KOREA

HORIBA Automotive Test Systems Ltd.

Room #906, World Meridian Venture Centre I, 60-24

Gasam-Dong, Geumcheon-Gu, Seoul, 153-781

Phone : (82)2-562-7706 Fax : (82)2-562-7630

HORIBA KOREA Ltd.

202-501, Bucheon Techno Park, 192,Yakdae-Dong, Wonmi-ku, Bucheon Kyunggido Korea

Phone : (82)32-621-0100 Fax : (82)32-621-0105

Seoul Office

10, Dogok-ro 6-gil, Gangnam-gu, Seoul, 135-270

Phone : (82)2-753-7911 Fax : (82)2-756-4972

Suwon Office

D2f Hae Sung Bldg 307-9 Yuljun-Dong Changan-Ku Suwon, Kyunggi-Do

Phone : (82)31-296-7911 Fax : (82)31-296-7913

Ulsan Office

601, Kaya Bldg, 193-5, Sinjeong 3Dong, Nam-Ku, Ulsan

Phone : (82)52-275-0122 Fax : (82)52-276-0136

HORIBA STEC KOREA, Ltd.

110, Suntech-City, 513-15, Sangdaewon, Jungwon-Ku, Sungnam-City, Kyungki-Do

Phone : (82)31-777-2277 Fax : (82)31-777-2288

SINGAPORE

HORIBA Instruments (Singapore) Pte Ltd.

10, UBI CRESCENT #05-12 LOBBY B UBI TECHPARK (408564)

Phone : (65)6-745-8300 Fax : (65)6-745-8155

TAIWAN

HORIBA Taiwan, Inc.

3F., No.18, Ln. 676, Zhonghua Rd, Zhubei City, Hsinchu County 302

Phone : (886)3-656-1160 Fax : (886)3-656-8231

Tainan Office

1F., No.117, Chenggong Rd., Shanhua Township, Tainan Country 741

Phone : (886)6-583-4592 Fax : (886)6-583-2409

THAILAND

HORIBA ABX Diagnostics (Thailand) Ltd.

393, 395, 397, 399, 401, 403 Ladyah Road, Somdejchaopraya Sub-district, Klongsan District, Bangkok 10600

Phone : (66)2-861-59-95 Fax : (66)2-861-52-00

VIETNAM

HORIBA Instruments (Singapore) Pte Ltd.

Hanoi Office

Unit 10, 4 Floor, CMC tower, Day Tan Street, Dich Vong Hau Ward, Cau Giay district, Hanoi

Phone : (84)4-3795-8552 Fax : (84)4-3795-8553

Readout HORIBA Technical Reports MARCH 2013 No.40

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