
Readout

HORIBA Technical Reports

April 2015 No. **44**

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This volume introduces the many products and wealth of experience at Horiba Medical. Horiba Medical products for the blood testing market are available in over 150 countries and five continents. This volume also describes contributions made to clinical testing at regional hospitals and research institutions.

Around the time cherry petals began to fall in the Tohoku region, I visited Hirosaki Castle. Cherry petals covered the moat of the castle as if a carpet was spread. Petals that pleased the eyes of many tourists stuck to leaves of striped bamboo to make visitors please to the last moment. I was fascinated by the sight for a while.

-Photographer Hideo MATSUI-
(Member of NIKA Association of Photographers)

Name of the book

This book is named "Readout" in the hope that "the products and technology we have created and developed will be read out and so become widely known".

Contribution of Diagnostics to Total Medical Care / Healthcare

Review

- 6** ~ Together, let's build the future ~
Arnaud PRADEL

Guest Forum

- 12** The Actual Situation of the International Standardization and External Quality Assurance on Blood Cell Count
Hiroshi KONDO

Feature Article Contribution of Diagnostics Solution on Laboratory Testing

- 16** POCT-compatible Blood Glucose Monitor ANTSENSE ROSE
Introduction and Application in Our Hospital
Osamu YAMADA
- 21** HELO: Horiba Evolutive Laboratory Organization
Laurent ARAUD, Olivier POU
- 26** HORIBA contribution to athletes' healthcare at London 2012 Olympic Games
Mandy CAMPBELL, Manuela PASTORE

Topics

- 29** Specialized Hematology: An Ambition for HORIBA Medical
Jean-Michel GARCIA, Elina ALATERRE, Sébastien RAIMBAULT
Guillaume CARTRON, Guilhem REQUIRAND, Jerome MOREAUX
Caroline BRET, Jean-François SCHVED

Feature Article Contribution of Diagnostics on Infectious Diseases / Immunological Diseases

- 35** Judicious Use of Antibiotics for Pediatric Infection
—Global Strategies to Prevent the Increase of Bacterial Resistance—
Kazunobu OUCHI
- 39** Neonatal Infections
—Clinical Manifestation and Diagnosis—
Masahiro HAYAKAWA
- 43** CRP (C reactive protein)
—Basics and Clinical Application of CRP—
Masaki NAKANISHI, Toru INABA, Naohisa FUJITA

Selected Article

- 47** Evaluation methods of carbon-nanotube dispersed aqueous solution
Yasushi NAKATA, Sakiko AKAJI, Yoko SHINOZAKI

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Arnaud PRADEL

Director of Products Development
HORIBA ABX SAS

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The Actual Situation of the International
Standardization and External Quality
Assurance on Blood Cell Count

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In the hematology tests by using automated hematology analyzer, there is not internationally certified reference material. Therefore, an internationally recognized reference method enables the clarification of traceability of measurements since the use of fresh blood is necessary for the assurance of accuracy in hematology tests. The guidelines of the reference method are formulated with under international framework. The external quality assurance acts as an invaluable reference source for evaluation the reliability of measurement system. It's significant to standardize the method of implementation for external quality assurance.

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POCT-compatible Blood Glucose Monitor ANTSENSE ROSE
Introduction and Application in Our Hospital

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Our hospital introduced POCT-compatible instruments for blood glucose testing in our wards and other locations for point of care testing. Additionally, we utilized a data management system to handle both instrument management and data management. As the use of electronic chart systems increases in medical institutions, improving the quality of blood glucose testing with points of care and applying measurement data have become important issues. By combining sophisticated measurement instruments with a data management system, we succeeded in pointing the way toward one integrated solution.

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HELO: Horiba Evolutive Laboratory Organization

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Over the past 5 years, HORIBA Medical has been proposing the “HaemCell” as the high-range solution for Hematology Laboratory organization which was based on software connections only and avoiding any mechanical link. To improve patient care and adapt to a highly evolving market, HORIBA Medical decided to make the HaemCell evolve toward the HELO concept and exploit at their best the analyzer’s performances and scalability for healthcare needs. Still driven by the LEAN philosophy the new HELO concept extends the flexibility of configuration by adding the possibility to work with a track system and offers around both the work-cell and the dedicated track a multiplicity of peripheral solutions to accommodate modern requirements of the Hematology laboratory as well as related disciplines.

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Laurent ARAUD

Data Management Corporate Product Manager
HORIBA ABX SAS

Olivier POU

Junior Director
Distributor Network
HORIBA ABX SAS

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HORIBA contribution to athletes' healthcare at London 2012 Olympic Games

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The last Olympic and Paralympic Games took place in London in 2012. Spire Healthcare, Laboratory Network in charge of the healthcare during this event, built a state-of-the-art facility within the Olympic Village to give over 10,000 Olympic and 6,200 Paralympic athletes, as well as associated team officials, access to 24-hour healthcare services. Since the London 2012 Olympics were heralded as one of the most sustainable, Spire was required to reflect this ethos and approached HORIBA Medical to secure a loan and service support for an ABX Pentra 80 analyzer. About 250 samples were tested during this event. Finally, the NHS Trust concluded that the use the pathology services of the Olympic Village Polyclinic were able to satisfy demands and that a similar structure of facilities and expertise in one clinic was recommended when planning future mass participation sporting events.

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- [1] Vanhegan, Ivor S. et al. 'The London 2012 Summer Olympic Games: an analysis of usage of the Olympic Village 'Polyclinic' by competing athletes. *Br J Sports Med*, 47, 415(2013)

Mandy CAMPBELL

Product Manager
HORIBA UK Limited

Manuela PASTORE

Scientific Manager
Marketing Dept.
HORIBA ABX SAS
PhD

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Specialized Hematology: An Ambition for HORIBA Medical

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With the increase of cancer occurrences, its management, requires novel diagnostic tools improved performance and adaptation to personalized medicine. Such technological development can only be successful by the convergence of expertise between the medical field and the experience of industry in translating innovation into products. In that spirit, HORIBA Medical teamed up with Montpellier University Regional Hospital. That collaboration is realized by the presence of a HORIBA laboratory within the hospital working closely with local teams. After a successful initial project on Minimal Residual Disease, the strategy was recently extended to other hemopathies aiming at the discovery of new oncological biomarkers that could advance earlier diagnosis as well as the constitution of an original bio-bank focused on hemopathies and in particular lymphoma and myeloma.

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- [1] EuroFlow: Resetting leukemia and lymphoma immunophenotyping. Basis for companion diagnostics and personalized medicine. van Dongen, J JM et Orfao, A. 2012, Leukemia, Vol. 26, pp. 1899-1907.
- [2] Harmonization in clinical flow cytometry (FC): The HARMONEMIA project. Porwit, Anna, et al. October 2014, ICCS 2014, p. accepted.
- [3] Multiparameter Immunophenotyping by Flow Cytometry in Multiple Myeloma: Advantages for Diagnosis and Minimal Residual Disease Monitoring. G. REQUIRAND, J. MOREAUX, A. CARAUX, S. RAIMBAULT, B. KLEIN. September 2012, Technical Journal "Readout" English Edition, Vol. 39, pp. 46-51.
- [4] Prognostic value of deep sequencing method for minimal residual disease detection in multiple myeloma. Joaquin Martinez-Lopez, Juan J. Lahuerta, François Pepin, Marcos González, Santiago Barrio, Rosa Ayala, Noemí Puig, María A. Montalban, Bruno Paiva, Li Weng, Cristina Jiménez, María Sopena, Martin Moorhead, Teresa Cedena, Immaculada Rapado, María Victoria. March 2014, Blood. DOI 10.1182/blood-2014-01-550020.
- [5] Detection and Molecular Monitoring of Minimal Residual Disease in Chronic Myelogenous leukemia. Yin, Cameron C. 4, 2012, Clinical & Experimental Pathology, Vol. 2, p. e110.
- [6] Beyond morphology: minimal residual disease detection in acute myeloid leukemia. DiNardo, CD et Luger, SM. 2, 2012, Current Opinion in Hematology, Vol. 19, pp. 82-88.
- [7] High Prognostic Impact of Flow Cytometric Minimal Residual Disease Detection in Acute Myeloid Leukemia: Data From the HOVON/SAKK AML 42A Study. Terwijn, Monique, et al. 31, 2013, Journal of Clinical Oncology, Vol. 31, pp. 3889-3897.
- [8] Minimal residual disease monitoring in chronic lymphocytic leukaemia patients. A comparative analysis of flow

cytometry and ASO IgH RQ-PCR. Raponi, Sara, et al. 3, 2014, British Journal of Haematology, Vol. 166, pp. 360-368.

- [9] Immunoglobulin and T Cell Receptor Gene High-Throughput Sequencing Quantifies Minimal Residual Disease in Acute Lymphoblastic Leukemia and Predicts Post-Transplantation Relapse and Survival. Logan, Aaron C, et al. 9, 2014, Biology of Blood and Marrow Transplantation, Vol. 20, pp. 1307-1313.
- [10] Deep-sequencing approach for minimal residual disease detection in acute lymphoblastic leukemia. Faham, Malek, et al. 26, 2012, Blood, Vol. 120, pp. 5173-5180.
- [11] Impact of Minimal Residual Disease, Detected by Flow Cytometry, on Outcome of Myeloablative Hematopoietic Cell Transplantation for Acute Lymphoblastic Leukemia. Bar, Merav, et al. 2014, Leukemia Research and Treatment, p. 421723.
- [12] Acute lymphoblastic leukemia: monitoring minimal residual disease as a therapeutic principle. Brüggemann, Monika, Gökbuget, Nicola et Kneba, Michael. 1, 2012, Seminars in Oncology, Vol. 39, pp. 47-57.
- [13] Minimal disseminated disease in high-risk Burkitt's lymphoma identifies patients with different prognosis. Mussolin, Lara, et al. 13, 2011, Journal of clinical Oncology, Vol. 29, pp. 1779-1784.
- [14] Early assessment of minimal residual disease identifies patients at very high relapse risk in NPM-ALK-positive anaplastic large-cell lymphoma. Damm-Velk, Christine, et al. 3, 2014, Blood, Vol. 123, pp. 334-337.
- [15] Bone marrow minimal disseminated disease (MDD) and minimal residual disease (MRD) in childhood T-cell lymphoblastic lymphoma stage III, detected by flow cytometry (FC) and real-time quantitative polymerase chain reaction (RQ-PCR). Stark, Batia, et al. 1, 2009, Pediatric Blood Cancer, Vol. 52, pp. 20-25.
- [16] Prognostic relevance of minimal residual disease in colorectal cancer. Bork, Ulrich, et al. 30, 2014, World Journal of Gastroenterology, Vol. 20, pp. 10296-10304.
- [17] Minimal residual disease in breast cancer: in blood veritas. Siravegna, Giulia et Bardelli, Alberto. 10, 2014, Clinical Cancer Research, Vol. 20, pp. 2505-2507.
- [18] Detection of minimal residual disease in lung cancer. Chudacek, Josef, et al. 2, 2014, Biomed Pap Med Fac Univ Palacky Olomouc Czech Repub, Vol. 158, pp. 189-193.
- [19] Flow cytometry detection of minimal residual disease in multiple myeloma: Lessons learned at FDA-NCI Roundtable Symposium. Landgren, Ola, et al. American Journal of Hematology. (Accepted Article) DOI: 10.1002/ajh.23831.
- [20] Minimal residual disease monitoring: the new standard for treatment evaluation of haematological malignancies? M, Hauwel et T., Matthes. 2014, Swiss Medical Weekly, Vol. 144, p. w13907.
- [21] Next generation sequencing: new tools in immunology and hematology. Mori, Antonio, et al. 4, 2013, Blood Research, Vol. 48, pp. 242-249.
- [22] Defining minimal residual disease in acute myeloid leukemia: which platforms are ready for "Prime Time"? Grimwade, David et Freeman, Sylvie D. 05-577593, 2014, Blood, Vol. [Epub ahead of print]. PMID: 25049280.
- [23] Is Minimal Residual Disease Monitoring Clinically Relevant in Adults with Acute Myelogenous Leukemia? Carlson, Karen-Sue B et Guzman, Monica L. 2, 2013, Current Hematologic Malignancy Reports, Vol. 8, pp. 109-115.
- [24] Biology, Risk Stratification, and Therapy of Pediatric Acute Leukemias: An Update. Pui, Ching-Hon, et al. 5, 2011, Journal of Clinical Oncology, Vol. 29, pp. 551-565.

Jean-Michel GARCIA

Clinical Research, Department of Innovation
HORIBA ABX SAS
PhD

Elina ALATERRE

Clinical Research, Department of Innovation
HORIBA ABX SAS

Sébastien RAIMBAULT

Clinical Research, Department of Innovation
HORIBA ABX SAS
PhD


Prof. Guillaume CARTRON

Department of Clinical Hematology
University Hospital Center (CHRU) of Montpellier
MD, PhD

Guilhem REQUIRAND

STI, Department of Biological Hematology
University Hospital Center (CHRU) of Montpellier

Jerome MOREAUX

STI, Department of Biological Hematology,
University Hospital Center (CHRU) of Montpellier
PhD

Caroline BRET

Department of Biological Hematology
University Hospital Center (CHRU) of Montpellier
PharmD, PhD

Prof. Jean-François SCHVED

Department of Biological Hematology
University Hospital Center (CHRU) of Montpellier
MD, PhD

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Judicious Use of Antibiotics for Pediatric Infection

—Global Strategies to Prevent the Increase of Bacterial Resistance—

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Principle of antimicrobial therapy in children is to select and use the most sensitive antibiotics supposing the causative organism under special characteristics of children. Resistant organisms have been increasing these days, so it is very important for us to use proper antibiotics under considering the trends of antimicrobial susceptibility and PK/PD. We are careful to use antibiotics based on proper differential diagnosis exploiting rapid diagnose tests, CBC, CRP and so on. The most practical manner for judicious antibiotic use is to follow the newly revised guidelines for treating each patient. Under considering shortage of new antimicrobial agents in the future, we must keep judicious antibiotic use for maximum effect and long life span of antimicrobial agents.

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Neonatal Infections

—Clinical Manifestation and Diagnosis—

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The clinical characteristics of infections in neonates are different from those in children and adults. Especially, some neonates suffer from opportunistic infections and critical infections because neonates are physiologically immunocompromised. It is important that neonatologists must make an early diagnosis and begin early treatments. However, early diagnosis sometimes is difficult because clinical symptoms are non-specific in early stage. Neonatologists should suspect neonatal infections when neonates have subtle symptoms. The diagnosis should confirm by laboratory data. Serial evaluation of C-reactive protein is useful for diagnosis of neonatal infections.

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[11] Connell TG, Rele M, Cowley D, et al, "How reliable is a negative blood culture result? Volume of blood submitted for culture in routine practice in a children's hospital", *Pediatrics*, **119**, 891(2007)

[2] Polin RA, Parravicini E, Regan JA, "Bacterial sepsis and meningitis", *Avery's Disease of the Newborn*, 8th ed. Elsevier, 551-577(2005)

[3] Newman TB, Puopolo KM, Wi S, Draper D, Escobar GJ, "Interpreting complete blood counts soon after birth in newborns at risk for sepsis", *Pediatrics*, **126**, 903(2010)

[4] Guida JD, Kunig AM, Leef KH, et al, "Platelet count and sepsis in very low birth weight neonates: is there an organism-specific response?", *Pediatrics* **111**, 1411(2003)

[5] Manzoni P, Mostert M, Galletto P, et al, "Is thrombocytopenia suggestive of organism-specific response in neonatal sepsis?", *Pediatr Int*, **51**, 206(2009)

[6] Benitz WE, Han MY, Madan A, et al, "Serial serum C-reactive protein levels in the diagnosis of neonatal infection" *Pediatrics*, **102**, E141(1998)

[7] Garges HP, Moody MA, Cotten CM, et al, "Neonatal meningitis: what is the correlation among cerebrospinal fluid cultures, blood cultures, and cerebrospinal fluid parameters?", *Pediatrics*, **117**, 1094(2006)

[8] Johnson CE, Whitwell JK, Peth K, et al "Term newborns who are at risk for sepsis: are lumbar punctures necessary?", *Pediatrics*, **99**, E10(1997)

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CRP (C reactive protein)

—Basics and Clinical Application of CRP—

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Accurate identification of infected organs, an infecting organism, the necessity of an antimicrobial and the degree of severity is important for the treatment of infection, for which it is essential to understand accurately the severity of inflammation. However, even if doctors perform careful evaluations, we occasionally come across the cases in which it is difficult to understand exactly the inflammation. These situations require a simple and accurate inflammatory marker. Although many inflammatory markers have been developed, there is no biomarker useful by itself for all of diagnosis, prognostic prediction, response evaluation. CRP is considered to be useful in the treatment of infection as a biomarker that reflects the clinical course rapidly, at a low cost and in a balanced manner. Moreover, if we understand the mechanism of CRP production and its weak point, it has a great deal of potential in the daily medical practice.

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
Figure 3 Role of CRP in human body

Table 1 Main causes of false-positive and -negative results of PCT measurement. (cited from ref. [9])

<p>Common causes of false-negative and false-positive results:</p> <p>False-positive (i.e., falsely high levels in the absence of a bacterial infection):</p> <ul style="list-style-type: none"> • Newborns (physiologically) during first days of life • Acute respiratory distress syndrome acute attacks of plasmodium falciparum malaria systemic fungal • Infections (e.g., candidiasis, aspergillosis) • Severe mechanical trauma • Following surgical trauma • Administration of monoclonal or polyclonal anti-thymocyte globulin in the treatment of care • Rejection after transplantation • Chemical pneumonitis • Severe burns and heat strokes • Patients with medullary thyroid cancer, small cell cancer of the lung, carcinoid tumours with • Paraneoplastic hormone production • Inflammation associated with "cytokine storms", e.g., ILβ, in familial Mediterranean fever, therapeutic • Infusions of TNFα for melanoma <p>False-negative (i.e., falsely low levels in the presence of a bacterial infection):</p> <ul style="list-style-type: none"> • Early course of infections • Localized infections • Subacute endocarditis

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
Figure 4 Time-course changes in acute inflammatory markers. (cited from ref. [11])

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
Figure 5 Time-course changes in ESR, CRP, and WBC in children with bone marrow and joint infections. (cited from ref. [13])

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- [1] Tillett WS, Francis T, "Serological reactions in pneumonia with a non-protein somatic fraction of pneumococcus", *J Exp Med* **52**, 561 (1930)
- [2] 宮本恵宏, "循環器疾患の炎症マーカー (高感度CRPとアディポネクチン)", *生物試料分析*, **32**(2), 121 (2009)
- [3] 中村治雄, 山下毅, 本間優, "日本人の基準範囲と動脈硬化のリスク度評価", *臨床検査*, **46**(9), 951 (2002)
- [4] Volanakis JE, "Human C-reactive protein: expression, structure and function", *Mol Immunol*, **38**, 189 (2001)
- [5] Marnell L, Mold C, Du Clos TW, "C-reactive protein: ligands, receptors and role in inflammation", *Clin Immunol*, **117**, 104 (2005)
- [6] Becker KL, Nylén ES, White JC, Müller B, Snider RH Jr, "Procalcitonin and the calcitonin gene family of peptides in inflammation, infection and sepsis: a journey from calcitonin back to its precursors", *J Clin Endocrinol Metab*, **89**(4), 1512 (2004)
- [7] Linscheid P, Seboek D, Nylén ES, et al, "In vitro and in vivo calcitonin I gene expression in parenchymal cells: a novel product of human adipose tissue", *Endocrinology*, **144**(12), 5578 (2003)
- [8] Müller B, Peri G, Doni A, et al, "High circulating levels of the IL-1 type II decoy receptor in critically ill patients with sepsis: association of high decoy receptor levels with glucocorticoid administration", *J Leukoc Biol*, **72**(4), 643 (2002)
- [9] Christ-Crain M, Müller B, "Procalcitonin in bacterial infections-hype, hope more or less?", *Swiss Med Wkly*, **135**(31-32), 451 (2005)
- [10] Bouadma L, Luyt CE, Tubach F, et al, "PRORATA trial group. Use of procalcitonin to reduce patients' exposure to antibiotics in intensive care units (PRORATA trial)", *a multicentre randomised controlled trial. Lancet*, **375**(9713), 463 (2010)
- [11] Jaye DL, Waites KB, "Clinical applications of C-reactive protein in pediatrics", *The Pediatric Infectious Disease Journal*, **16**(8), 735 (1997)
- [12] Williamson MA, "Wallach's Interpretation of Diagnostic Tests", North American Edition (Interpretation of Diagnostic Tests, 9th edition) Lippincott Williams & Wilkins, USA, 2011
- [13] Pääkkönen M, Kallio MJ, Kallio PE, Peltola H, "Sensitivity of erythrocyte sedimentation rate and C-reactive protein in childhood bone and joint infections", *Clin Orthop Relat Res*, **468**(3), 861 (2010)

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Evaluation methods of carbon-nanotube dispersed aqueous solution

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Carbon nanotube is well known nano-material with many superior features as thermo stable, toughness, low gravity, unique electrical properties, and biocompatibility and applied to practical use in many fields. Especially, for electronics application such as transparent conductive film, antistatic coating, touch panel display, and FET (Field Effect Transistor), dispersion technique of CNT in solution is now becoming more important. Aggregation degree, length, semiconductor and metallic features and conductivity of CNT are important properties on its practical application. For its bulk material, elemental analysis and thermogravimetric analysis are used, and for nano-structure observation, TEM and AFM are often used. Otherwise, for dispersed solution of CNT, analysis method used with light probe or spectroscopic methods are more useful. In this report, the analysis method will be introduced through reporting the analysis results of CNT dispersed solution which was prepared by HORIBA.

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
Figure 1 Structure of SWCNT
Carbon hexagonal direction are expressed by chiral angle θ against the tube axis of longitudinal direction expressed as long dashed short dashed lines which correspond to $\theta = 0$ degree and $\theta = 30$ degree.

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
Figure 3 NIR-PL of SWCNT by Alcohol CVD method
(a) EEM map (b) Chirality distribution. A pair of numbers is chiral indices of SWCNT.
Data courtesy of Prof. S. Maruyama, the University of Tokyo.


Figure 4 Comparison of Raman spectra on the process of purification SWCNT, which was dispersed in 20% cycloamylose aqueous solution, was purified by centrifugation. F-1 and F-2 were fraction of before and after purification respectively. F-3 was residue fraction after purification.
Data courtesy of Prof. S. Kitamura, Osaka Prefecture University.

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
Figure 5 Nanoparticle Analyzer nano Partica SZ-100


Figure 6 Time variation of particle size distribution of SWCNT in 0.2% carboxymethyl cellulose aqueous solution by ultrasonication
Time values at left-side of figures indicate process time of sonication.

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
Figure 7 Particle size distribution of SWCNT/SDBS measured by DLS (Dynamic Light Scattering)method


Figure 8 D and G Raman bands of SWCNT/SDBS
He-Ne Laser (633 nm) for excitation.


Figure 9 RBM Raman bands of SWCNT/SDBS
He-Ne Laser (633 nm) for excitation.
A pair of numbers in parenthesis is chiral indices of SWCNT.


Figure 10 Absorption Spectra of SWCNT/SDBS
A pair of numbers in parenthesis is chiral indices of SWCNT

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
Figure 11 NIR-PL EEM map of SWCNT/SDBS
A pair of number in parenthesis is chiral indices of SWCNT

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
Figure 12 Peak shift of NIR PL bands of SWCNT/SDBS
NIR-PL bands were assigned as (10,2), (9,4) and (8,6) for each. Excitation wavelength was 720nm. The left figure is the initial spectrum. Every bands shifted about 8nm longer wavelength after 214 days.

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There is growing emphasis on POCT (Point of Care Testing) which brings in-situ test results to the treatment because it enables the provision of immediate and effective medical care. HORIBA, Ltd. has started the sales of a new model Antsense Duo, a small electrode type blood glucose analyzer deriving from the conventional cartridge-type analyzer. Hereunder, we describe the features, improvements and the performance evaluation.

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Kyoto Fukuchiyama Technology Center

11-1 Ecotopia Miwa Koaza, Miwa, Fukuchiyama
620-1445
Phone : (81) 773-59-2070 Fax : (81) 773-59-2074

HORIBA TECHNO SERVICE Co., Ltd.

2, Miyanohigashi, Kisshoin, Minami-ku,
Kyoto 601-8305
Phone : (81)75-325-5291 Fax : (81)75-315-9972

BRAZIL

HORIBA Instruments Brasil, Ltda.

Rua: Presbítero Plínio Alves de Souza, 645,
Loteamento Polo Multívias Barro Medeiros,
Jundiá, São Paulo 13212-181
Phone : (55)11-2923-5400 Fax : (55)11-2923-5490

TCA/HORIBA Sistema de Testes Automotivos Ltda.

Rua Goiás 191 Vila Oriental, Diadema, São Paulo,
09941-690
Phone : (55)11-4224-0200 Fax : (55)11-4227-3133

CANADA

HORIBA Canada, Inc.

Unit102, 5555 North Service Road Burlington,
Ontario, L7L 5H7
Phone : (1)905-335-0234 Fax : (1)905-331-2362

London Office

347 Consortium Court, London, Ontario, N6E 2S8
Phone : (1)519-668-6920 Fax : (1)519-668-8437

U.S.A.

HORIBA International Corporated

9755 Research Drive, Irvine, CA 92618
Phone : (1)949-250-4811 Fax : (1)949-250-0924

HORIBA Instruments Incorporation

9755 Research Drive, Irvine, CA 92618
Phone : (1)949-250-4811 Fax : (1)949-250-0924

Albany Office

Suite (04), 58 Clifton Country Road, Clifton Park, NY 12065
Phone : (1)-518-331-1371

Alvin TX Office

5318 West Farm to Market 517 Road, Alvin, TX 77511
Phone : (1)-281-482-4334 Fax : (1)-281-614-0303

Austin Office

9701 Dessau Road, Suite 605, Austin, TX 78754
Phone : (1)512-836-9560 Fax : (1)512-836-8054

Chicago Office

1725 Roosevelt Road #115, West Chicago, IL 90185
Phone : (1)630-562-2298 Fax : (1)630-562-2324

Edison Office

3880 Park Avenue, Edison, NJ 08820
Phone : (1)732-494-8660 Fax : (1)732-549-5125

Houston Office

240 Spring Hills Drive, Suite 410, Spring, TX 77386
Phone : (1)281-367-7422 Fax : (1)281-367-7423

New Hampshire Office

315 Derry Road, Suite 13 Hudson, NH 03051
Phone : (1)603-886-4167 Fax : (1)603-886-4267

Portland Office

10110 South West Nimbus Avenue, Suite B-11,
Portland, OR 97223
Phone : (1)503-624-9767 Fax : (1)503-968-3236

Reno Office

605 Spice Island Drive, #5, Sparks, NV 89431
Phone : (1)775-358-2332 Fax : (1)775-358-0434

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HORIBA INSTRUMENTS (SHANGHAI) CO., Ltd.

No.200, Taitao Road, Anting Town, Jiading District, Shanghai 201814
Phone : (86)21-6952-2835 Fax : (86)21-6952-2823

HORIBA (China) Trading Co., Ltd.

Unit D, 1 Floor, Building A, Synnex International Park, 1068 West Tianshan Road, Shanghai, 200335
Phone : (86)21-6289-6060 Fax : (86)21-6289-5553

Beijing Office

12F, Metropolis Tower, No. 2, Haidian Dong 3 street Beijing 100080
Phone : (86)10-8567-9966 Fax : (86)10-8567-9066

Guangzhou Office

Room 1611/1612, Goldlion Digital Network Center, 138 Tiyu Road East, Guangzhou 510620
Phone : (86)20-3878-1883 Fax : (86)20-3878-1810

Shanghai Service Center

Room 303, No.84, Lane887, Zu-Chong-Zhi Road, Zhangjiang Hi-tech Park, Shanghai, 201203
Phone : (86)21-5131-7150 Fax : (86)21-5131-7660

Shanghai Technical Center

No.200, Taitao Road, Anting Town, Jiading District, Shanghai, 201814
Phone : (86)21-6289-6060 Fax : (86)21-6289-5553

INDIA

HORIBA India Private Ltd.

246, Okhla Industrial Estate, Phase 3, New Delhi, 110020
Phone : (91)11-4669-5001/(91)11-4646-5000

Fax : (91)11-4669-5010/(91)11-4646-5020

Bangalore Office

Kamadhenu, No.17/1 - 32, Bannerghatta Road, Audugodi, Bangalore - 560030
Phone : (91)80- 22210071

Chennai office

The ground floor New No. 41, Old No. 15/3, Velachery Road, Little Mount, Saidapet, Chennai - 600015
Phone : +(91)-44-42077899

Pune Office

502, 5 Floor, Purushottam Plaza, Baner Road, Baner, Pune - 411045
Phone : (91)20-40766000 Fax : (91)20-40766010

Haridwar Factory

Plot No. 26, sector-7, IIE, Sidcul, Uttarakhand 249403
Phone : (91)-1334-239139

INDONESIA

PT HORIBA Indonesia

JL. Jalur Sutera Blok 20A, No. 16-17, Kel. Kunciran, Kec. Pinang Tangerang 15144
Phone : (62)21-3044-8525 Fax : (62)21-3044-8521

KOREA

HORIBA Automotive Test Systems Ltd.

Room #906, World Meridian Venture Centre 1, 254, Beotkot-Ro, Geumcheon-Gu, Seoul, 153-781
Phone : (82)2-562-7706 Fax : (82)2-562-7630

HORIBA KOREA Ltd.

202-501, Bucheon Techno Park, 388, Songnae-Daero, Wonmi-Gu, Bucheon-Si Gyunggi-Do, 420-831
Phone : (82)32-621-0100 Fax : (82)32-621-0105

Seoul Office

10, Dogok-Ro 6-Gil, Gangnam-Gu, Seoul, 135-860
Phone : (82)2-753-7911 Fax : (82)2-756-4972

Suwon Office

D2f Hae Sung Building 107-5 Yuljeon-Ro, Jangan-Gu, Suwon-Si Gyeonggi-Do, 440-825
Phone : (82)31-296-7911 Fax : (82)31-296-7913

Ulsan Office

601, Kaya Building, 93, Sinjeong 3 Dong, Nam-Gu, Ulsan
Phone : (82)52-275-0122 Fax : (82)52-276-0136

HORIBA STEC KOREA, Ltd.

110, Suntech-City, 474, Dunchon-daero, Jungwon-Gu, Seongnam, Gyeonggi-Do
Phone : (82)31-777-2277 Fax : (82)31-777-2288

SINGAPORE

HORIBA Instruments (Singapore) Pte. Ltd.

3 Changi Business Park Vista #01-01 Akzonobel House, Singapore, 486051
Phone : (65)6-745-8300 Fax : (65)6-745-8155

West Office

83 Science Park Drive #02-02A, The Curie Singapore 118258
Phone : (65)6-908-9660

TAIWAN

HORIBA Taiwan, Inc.

3 Floor, No.18, Ln. 676, Zhonghua Road, Zhubei, Hsinchu country 302
Phone : (886)3-656-1160 Fax : (886)3-656-8231

Tainan Office

1 Floor, No.117, Chenggong Road, Shanhua, Township Tainan country 741
Phone : (886)6-583-4592 Fax : (886)6-583-2409

THAILAND

HORIBA (Thailand) Ltd.

393, 395, 397, 399, 401, 403 Latyah Road, Somdetchaopraya, Klongsan, Bangkok 10600
Phone : (66)2-861-59-95 Fax : (66)2-861-52-00

VIETNAM

HORIBA Instruments (Singapore) Pte Ltd.

Hanoi Office

Unit 10, 4 Floor, CMC tower, Duy Tan Street, Dich Vong Hau Ward, Cau Giay, Hanoi
Phone : (84)4-3795-8552 Fax : (84)4-3795-8553

Readout HORIBA Technical Reports April 2015 No.44

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