# Readout

HORIBA Technical Reports

August 2018 No. **49**

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Photonic instrumentation is getting more important as biotechnology application progresses in life science. Like the lighthouse in the jet black ocean, we hope that the latest analytical instruments introduced in this issue will continue to support the reliable future of life science.

Maritime traffic in the Seto Inland Sea has been extensive for a long time. Many ships pass through even in the night. Operation of those ships in pitch dark ocean relies on the light house serving to ensure the safety of the ships. The light it emits makes me feel like looking at the guiding light of our life path.

-Photographer Hideo MATSUI-  
(Member of Nikakai Association of Photographers)

#### Name of the book

This book is named "Readout" in the hope that "the products and technology we have created and developed will be read out and so become widely known".

## Photonic Instrumentation in Life Science

### Foreword

- 4** Bringing the Benefits of Photonic Instrumentation to the Field of Life Science  
James THEPOT

### Review

- 6** Importance of and Features in Time-Resolved Fluorescence Spectrometry Traces in HORIBA's Time-Resolved Fluorescence Spectrometer  
Kiyooki HARA
- 13** Finding the Time for Fluorescence  
Its Measurement and Applications in Life Science  
David J S Birch, David McLoskey

### Guest Forum

- 21** *In situ* Raman Spectroscopy to Monitor Interfacial Electron Transfer Process: from Fundamental to Application Case Studies  
Akiyoshi KUZUME, Ph.D.

### Feature Article Photonic Instrumentation in Life Science

- 30** Life Science Analytical Tools for Real Life Problems  
Marinella SANDROS, Fran ADAR
- 36** Spectroscopic Analysis of Red Wines with A-TEEM Molecular Fingerprinting  
Sakiko AKAJI
- 44** Impact of Raman Microscopy on Biomedical Research  
Fran ADAR, Catalina DAVID, Marinella SANDROS
- 54** Tip-Enhanced Raman Spectroscopy  
Finds a Growing Number of Applications in Biology  
Marc CHAIGNEAU, Maruda SHANMUGASUNDARAM, Fran ADAR

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**James THEPOT**

President, HORIBA FRANCE SAS  
Corporate Officer, HORIBA, Ltd.

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Importance of and Features in Time-Resolved Fluorescence Spectrometry  
Traces in HORIBA's Time-Resolved Fluorescence Spectrometer

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Time-resolved fluorescence measurements reveal intermolecular interactions depending upon the environment of the fluorophore. This paper describes the basic background and some examples of how the fluorescence lifetime is more informative than steady-state fluorescence measurements. HORIBA has many fluorescence technologies and resources.

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
Figure 2 Excimer formation in de-gassed sample: Pyrene in cyclohexane (2.2 mMol). (a) Monomer Emission, (b) Excimer emission.

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$$I(t) = A_1 \exp(-t/\tau_1) + A_2 \exp(-t/\tau_2) \dots\dots\dots (7)$$

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
Figure 3 Anisotropy decay: (a)  $I_{\parallel}(t) + 2 I_{\perp}(t)$ , (b)  $I_{\parallel}(t)$ , (c)  $I_{\perp}(t)$ , (d)  $I_{\parallel}(t) - I_{\perp}(t)$ , (e) Lamp profile. Sample was provided by Dr. Tunesha Araiso, Hokkaido Univ.


Figure 4 Excitation: 470nm, SDS conc.: 5 nM (a) Residual, (b) Excitation profile (c) Obs.630 nm,  $\tau_1=2.48$  ns (41.2%),  $\tau_2=9.03$  ns (58.8%) (d) Obs.535 nm,  $\tau_1=1.41$  ns (54%),  $\tau_2=3.48$  ns (45.5%)


Figure 5 Excitation 337 nm( $N_2$  laser) Observation 660 nm. Sample (a)  $\tau_1=680$  ps(83.1%),  $\tau_2=1.26$  ns(16.9%) Sample (b)  $\tau_1=199$  ps(83.8%),  $\tau_2=1.16$  ns(16.2%) (c) Laser pulse profile. Sample: provided by Dr. Shigekazu Minagawa, Central Research Laboratory, Hitachi Ltd.

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
Figure 6 Polarized excitation and fluorescence detection.

$$I_{\parallel}(t) = e^{-t/\tau} (1 + 2r_0 e^{-t/\theta}), I_{\perp}(t) = e^{-t/\tau} (1 - r_0 e^{-t/\theta}) \dots\dots\dots (13)$$

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
Table 1 Observable and Molecular property in TRFM.

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- [1] J. R. Lakowicz "Principle of Fluorescence Spectroscopy" Third Edition, Springer, 2006.
- [2] K. Awasthi, D. Moriya, T. Nakabayashi, L. Li and N. Ohta, J. Photochem. Photobiol. B 165, 256-265, 2016.
- [3] 中林孝和, 太田信廣, 日レ医誌(JJSLSM), 441, Vol.3, No.4, 2010.
- [4] W.R.Ware, Transient Luminescence measurement, In Creation and detection of the excited state, Vol. 1A, pp-213-302. Ed. AA Lamola. Marcel Dekker, New York, 1971.
- [5] D. V. O'Conner and D. Phillips, "Time-correlated Single Photon Counting", Academic Press, 1984.
- [6] 平山鋭, 原清明訳, ナノ・ピコ秒の蛍光測定解析法~時間相関単一光子計数法~, 学会出版センター, 1988.
- [7] I. Yokoyama, Readout, No.2, 27-36,1991.
- [8] 北海道大学学長, 試料の発光寿命測定装置, 特許1395515, 1987.8.24.
- [9] 木下一彦, 御橋廣真編: 蛍光測定~生物科学への応用~, 学会出版センター, 1983.
- [10] N.Miyoshi, K.Hara, I.Yokoyama, G.Tomita and M.Fukuda, Photochem.Photobiol. Vol.47,685-688,1988.

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Finding the Time for Fluorescence  
Its Measurement and Applications in Life Science

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We summarise how developments in technology have brought fluorescence lifetime spectroscopy from being the preserve of the specialist to becoming a major tool for research across many science and engineering disciplines. We highlight the advantages which fluorescence lifetime measurements can bring, not only to underpin research, but also through application in helping to solve real-world problems. We illustrate this with recent examples in cancer and Alzheimer's research, which are aimed at improving disease understanding, diagnosis and therapeutics.

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Figure 1 Model 5000F coaxial nanosecond flashlamp (a) and DeltaDiode picosecond LED/LD (b).

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

Figure 2 Model 5000W SAFE and its 16 channel TCSPC multiplexing module.

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
Figure 3 DeltaFlex dedicated lifetime system (a) and a hybrid FluoroMax-TCSPC lifetime system (b).

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
Figure 4 Fluorescence decay of melanin. In the presence of 4 μM Cu<sup>2+</sup> ions the decay becomes a single exponential of τ ~ 6 ns at 480 nm emission, suggesting a simplified structure.


Figure 5 Relative amplitude of two decay parameters ( $p_1/p_2$ ) when a 3-exponential decay model  $i(t) = \sum p_n \exp(-t/\tau_n)$ , for  $n=3$ , is used to describes the tyrosine fluorescence decay of Aβ during the early stages of aggregation.  $F_1/F_0$  describes the fluorescence intensity of the sheet-sensing probe ThT during this time and it can be seen to provide little information on the cytotoxic oligomer formation which occurs as a precursor to Alzheimer's disease.


Figure 6 GNR hairpin sensor sprung by binding to c-DNA and the consequential increase in the average fluorescence lifetime of oligonucleotide labelled with 6-carboxyfluorescein as hybridization occurs. Excitation was at 474 nm using the NanoLED laser diode.

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- [1] *A single-photon counting fluorescence decay-time spectrometer.* D J S Birch and R E Imhof. J. Phys. E, **10**, 1044-9, 1977.
- [2] *Coaxial nanosecond flashlamp.* D J S Birch and RE Imhof, Rev. Sci. Instrum. **52**, 1206-12, 1981.
- [3] *A new sub-nanosecond LED at 280 nm: application to protein fluorescence.* C D McGuinness, K Sagoo, D McLoskey and D J S Birch. Meas. Sci. Technol. **15**, L19-22, 2004.
- [4] *Fast time-correlated single-photon counting fluorescence lifetime acquisition using a 100 MHz semiconductor excitation source.* D. McLoskey, D. Campbell, A. Allison and G. Hungerford. Meas. Sci. Technol. **22**, 067001, 2011.
- [5] *Pulse fluorometry using simultaneous acquisition of fluorescence and excitation.* D J S Birch, R E Imhof and A Dutch, Rev. Sci. Instrum. **55**, 1255-64, 1984.
- [6] *A Multiplexed single photon instrument for routine measurement of time-resolved fluorescence anisotropy.* D J S Birch, A S Holmes, J R Gilchrist, R E Imhof, S M Al-Alawi and B Nadolski, J. Phys. E. Sci. Instrum. **20**, 471-3, 1987.
- [7] *Multiplexed array fluorometry.* D J S Birch, A S Holmes, R E Imhof, B Z Nadolski and K Suhling, J. Phys. E.: Sci. Instrum. **21** 415-7, 1988.
- [8] *Multiplexed single-photon counting 1: A time-correlated fluorescence lifetime camera.* D McLoskey, D J S Birch, A Sanderson, K Suhling, E Welch and P J Hicks, Rev. Sci. Instrum. **67**, 2228-37, 1996.
- [9] *Eumelanin Kinetics and Sheet Structure.* J Sutter, T Bidláková, J Karolin and D J S Birch. App. Phys. Letts., **100**, 113701 (4 pages), 2012.
- [10] *Metal ion influence on eumelanin fluorescence and structure.* J-U Sutter and D J S Birch. Methods Appl. Fluoresc. **2**, 024005, 2014.
- [11] *Early detection of amyloid aggregation using intrinsic fluorescence.* O J Rolinski, M Amaro and D J S Birch. Biosensors and Bioelectronics. **25** (10), 2249-52, 2010.
- [12] *Radiative decay engineering: biophysical and biomedical applications.* J. R. Lakowicz, Anal. Biochem. **298**(1), 1-24, 2001.
- [13] *Hairpin DNA-functionalized gold nanorods for mRNA detection in homogenous solution.* W Guoke, Y Jun, W Jinliang, G Peng, D J S Birch and Yu Chen. J. Biomed. Opt. **21** (9), 097001, 2016.

**David J S Birch**

President  
HORIBA JOBIN YVON IBH Ltd.  
Professor of the University of Strathclyde  
Ph.D.

**David McLoskey**

Managing Director  
HORIBA JOBIN YVON IBH Ltd.  
Ph.D.

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In situ Raman Spectroscopy to Monitor Interfacial Electron Transfer Process:  
from Fundamental to Application Case Studies

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In situ Raman spectroscopy is capable of probing electron transfer processes under electrochemical condition by monitoring local structure and reactivity properties of solid/liquid interfaces. Therefore, in situ Raman spectroscopy is now recognized as one of the powerful surface characterisation techniques for the investigation of the fundamental and application studies, such as development of energy conversion and bio-mediating devices. In addition, an introduction of nanostructured surfaces increases the sensitivity of Raman signals, capable of observing spectroscopic features from molecular level. This article gives a brief history of Raman spectroscopy in the introduction, followed by an overview of the recent in situ Raman spectroscopic studies in our group, particularly focusing on the electrochemical and microbial electron transfer processes at solid/liquid interfaces under electrochemical condition.

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Life Science Analytical Tools for Real Life Problems

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When light interacts with matter (i.e. molecule/sample) many interesting phenomena occur. Light may get absorbed by your sample and depending on the sample properties you will observe light emission, reflection and scattering. HORIBA Scientific builds tools that measure these phenomena allowing our users to gain meaningful information about their samples. This information is used to solve outstanding problems in research or guide scientists to monitor product development in various fields. Life Science, in particular, is an attractive market because, over the years, it has experienced exponential growth, and continues to do so. All the divisions that engage in the scientific study of living organisms including plants, animals and human beings fall under the life science umbrella. This article will highlight the key contributions our products offer in the life science sub-divisions - pharmaceutical and health sciences - and will compare their advantages over existing technologies.

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
Figure 1 A schematic comparison between small molecule (TAXOL) and protein (antibody IMMUNOGLOBULIN) drugs.

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
Figure 2 Raman Analysis of lysozyme (200 mg/mL) in 20 mM Citrate-PBS buffer at pH 4.0 before and after the addition of 20 % ethanol. Sample was excited with a 532 nm laser and a grating of 1800 lines/mm was used on the XploRA, a 200 mm focal length instrument.<sup>[1]</sup>


Figure 3 Comparison of the Trp band at 1555  $\text{cm}^{-1}$  full width half maximum (FWHM) width (top left) and  $I_{856}/I_{830}$  tyrosyl ratios (top right) at concentrations ranging from 10 to 200 mg/mL. A schematic representation of the microenvironment around IgG molecules changes as the concentration increases (bottom (i)-(iii)).<sup>[2]</sup>


Figure 4 Monitoring cell proliferation over 4 days using ATEEM. The EEM data were analyzed using PARAFAC analysis. As a result, four spectral components were successfully extracted and assigned to tyrosine, tryptophan, NAD(P)H and Riboflavin.<sup>[3]</sup>

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
Figure 5 A schematic representation of the three different configuration of peptides (ligands) used on the microarray chip: linear (L), branched (B) and Tandem (T).<sup>[5]</sup> (Reproduced with permission from ref.5, Copyright 2017, Elsevier)

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

Figure 6 SPRi study of alpha lactalbumin (LAC) antibody binding to the microarray peptide pre-functionalized with various orientation of alpha lactalbumin IgE epitopes LAC1\_L, LAC1\_B, LAC1\_T. The SPRi association (A-C; association rates were deduced from 3 min of association) and dissociation curves (D-F, dissociation rates were deduced from 4 h of dissociation) after the injection of 8 different concentrations (0.98-500 nM) of the antibody (colored plots).<sup>[5]</sup> (Reproduced with permission from ref.5, Copyright 2017, Elsevier)

**Raman**

**SPRi**

**A-TEEM**

Figure 7

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- [1] Marinella G. Sandros and Fran Adar. Assessing biotherapeutics stability using Raman Spectroscopy. (Application note) <http://www.horiba.com/scientific/products/horiba-life-science-solutions/application-notes/>
- [2] Chikashi Ota, Shintaro Noguchi, Satoru Nagatoishi and Kouhei Tsumoto. Assessment of the Protein-Protein Interactions in a Highly Concentrated Antibody Solution by Using Raman Spectroscopy, *Pharm. Res.*, 33, 956-969(2016).
- [3] Yuichi Kitagawa, Takumi Moriyama, Daisuke Irikura, and Yasushi Nakata. Holistic Analysis of Mammalian Cell Proliferation using Fluorescence Spectroscopy. (Application note) <http://www.horiba.com/scientific/products/horiba-life-science-solutions/application-notes/>
- [4] D.W.Zabriskie and A.E.Humphrey, Estimation of Fermentation Biomass Concentration by Measuring Culture Fluorescence. *Appl. Environ. Microbiol.*, 35, 337-343(1978).
- [5] Alessandro Gori, Marina Cretich, Renzo Vanna, Laura Sola, Paola Gagni, Giulia Bruni, Marta Liprino, Furio Gramatica, Samuele Burastero, Marcella Chiari *Analytica Chimica Acta* 983, 189-197(2017).

#### **Marinella SANDROS, Ph.D.**

Business Development Manager for Life Sciences  
SPRi Product Manager  
HORIBA Instruments Inc.

#### **Fran ADAR, Ph.D.**

Worldwide Raman Applications Manager  
HORIBA Instruments Inc

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Spectroscopic Analysis of Red Wines with A-TEEM Molecular Fingerprinting

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The patented (US patent 8901513) Aqualog which facilitates simultaneous Absorbance- Transmission and fluorescence Excitation-Emission Matrix (A-TEEM) technology, provides rapid access to a wide range of parameters of significance in water treatment, drug and protein analysis as well as the food and beverage industry. We think of Aqualog as an instrument for water analysis, but find that it is adept at wine analysis and more. Aqualog acquires a complete UV-VIS spectrum including the wine industry-standard Absorbance wavelength values at 280, 420, 520 and 620 nm which are important to evaluate a wine's phenolic content, and derive characteristic Hue and Intensity values. Aqualog also reports the Transmission spectrum which can be used to determine the CIElab Tri-Coordinate color descriptions. Aqualog reports a NIST-traceable EEM fingerprint which can be evaluated using multivariate statistics such as PARAFAC (Parallel Factor Analysis) and PCA (Principal Components Analysis). Most importantly, A-TEEM fingerprints yield qualitative and quantitative composition of key flavor and color determinants in grape juice and wine that are not discernible with simple Absorbance or Transmission data analysis.

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
Figure 3 Simultaneously recorded EEMs (A and B) and Absorbance and percent Transmission spectra (C and D) for a typical Italian red wine from a freshly opened bottle Before (A and C) and After a one week exposure to air (B and D).

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
Figure 4 Comparison of the Absorbance parameters (A) and Hue and Intensity parameters (B) defined above, measured with Aqualog for a series of freshly opened red wines from various countries.

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
Figure 5 CIE 1931 x and y coordinates for each of the wine samples from Figure 4.

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Table 3 Absorbance, Hue, Intensity and CIE Lab color index parameters for the Italy wine sample shown in Figure 3 Before and After oxidation.

Absorbance	Before	$\sigma$	After	$\sigma$	$\Delta$	$\sigma$
A620	1.14E+00	1.09E-02	2.11E+00	1.17E-02	-9.78E-01	1.60E-02
A520	4.12E+00	1.50E-02	7.15E+00	1.70E-02	-3.03E+00	2.27E-02
A420	4.14E+00	1.86E-02	7.60E+00	1.14E-02	-3.46E+00	2.18E-02
A280	5.30E+01	1.38E-01	8.21E+01	1.36E-01	-2.92E+01	1.93E-01
HUE	1.01E+00	5.82E-03	1.06E+00	2.99E-03	-5.69E-02	6.54E-03
Intensity	9.39E+00	2.63E-02	1.69E+01	2.36E-02	-7.47E+00	3.53E-02
CIE Lab						
X	1.10E+02	3.00E-02	1.06E+02	3.39E-02	4.79E+00	4.53E-02
Y	1.09E+02	3.29E-02	1.03E+02	3.85E-02	5.69E+00	5.07E-02
Z	1.06E+02	5.09E-02	9.72E+01	2.14E-02	8.57E+00	5.52E-02
x	3.40E-01	2.71E-05	3.45E-01	2.25E-05	-5.50E-03	3.52E-05
y	3.35E-01	1.59E-05	3.37E-01	3.30E-05	-2.25E-03	3.66E-05
X+Y+Z	3.25E+02	1.14E-01	3.06E+02	8.89E-02	1.91E+01	1.44E-01
L*	1.03E+02	1.21E-02	1.01E+02	1.46E-02	2.12E+00	1.89E-02
a*	1.18E+01	4.74E-03	1.32E+01	7.58E-03	-1.40E+00	8.94E-03
b*	6.57E+00	1.12E-02	8.45E+00	1.80E-02	-1.88E+00	2.12E-02
C*ab	1.35E+01	9.48E-03	1.57E+01	7.02E-03	-2.17E+00	1.18E-02
h*ab	5.09E-01	5.68E-04	5.70E-01	1.16E-03	-6.11E-02	1.29E-03
S*	1.31E-01	1.07E-04	1.55E-01	6.70E-05	-2.42E-02	1.26E-04
Q*	1.33E+02	1.08E-02	1.31E+02	1.31E-02	1.91E+00	1.70E-02


Figure 6 IFE corrected excitation and emission spectral contour loadings for the five component model evaluated for all wine sample replicates before and after oxidation. The PARAFAC model (n=84) was described with an  $r^2=0.997$  and a split half validation score  $r^2=0.906$  (where the model is tested with two halves of the dataset separately).

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Table 4 Qualitative assignment of PARAFAC components

PARAFAC Component	Excitation Max, nm	Emission Max, nm	Name
C1	278	340	Caffeic Acid
C2	263	380	Flavonol Like
C3	280	300	Epicatechin
C4	315	405	Gentisic Acid
C5	445	568	Anthocyanin


Figure 7 Comparison of the five PARAFAC component scores in the Italy wine samples (n= 3 replicates per sample) before and after oxidation. Component scores are reported normalized water Raman scattering units (RSU).


Figure 8 Comparison of the three major component cluster plots for PARAFAC (A) and PCA (B) analyses of the freshly opened wine samples ( $n=3$  replicates per sample).


Figure 9 Schematic of the relationships among the simultaneous Absorbance and EEM data acquisition and analysis with respect to the significance of wine sample and component resolution pertaining to this study.

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- [1] <https://www.mordorintelligence.com/industry-reports/wine-market> [accessed 21<sup>st</sup>, October 2017]
- [2] A. L. Waterhouse et al. Understanding Wine Chemistry, ed. Wiley and Sons 2016.
- [3] M. Sternad Lemut et al., *J Sci Food Agric.*, 2013, **93**, 3670–3681.
- [4] C. Coelho et al., *Anal. Chem.*, 2015, **87**, 8132–8137.
- [5] <https://www.etslabs.com/library/44> [accessed 21<sup>st</sup>, October 2017]
- [6] D. Airado-Rodriguez et al., *Journal of Food Composition and Analysis*, 2011, **24**, 257–264.
- [7] A. Gilmore, A. Akaji, K. Castorday, “Spectroscopic Analysis of Red Wines with A-TEEM Molecular Fingerprinting”, *Readout English Edition*, 2017, **49**, 41–48.

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Impact of Raman Microscopy on Biomedical Research

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Raman spectroscopy is a vibrational technique complimentary to Infrared (IR) Absorption. They both provide information on materials via their molecular and crystalline vibrations. While IR spectroscopy has been used more widely, Raman has some distinctly different characteristics that provide advantages when making measurements on biological materials. For this article we have selected three topics where there is groundbreaking activity in Raman microscopy worldwide and where the implementation has potential to solve outstanding problems. Information on instrumentation innovations that are critical to some of this work (especially lab-on-a-chip) will be mentioned. The references, while not exhaustive, will provide the interested reader entry to the literature. The topics that have been selected are identification of bacterial micro-organisms, disease diagnosis, and Raman imaging.

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
Figure 1 Scheme proposed by European Union for identifying infectious organism and preparing the best antibiotic defense. This figure is reproduced from <http://www.hemospec.eu/index.php/technology>.

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
Figure 2 (A) (B) Raman map and the characteristic Raman spectra of native sperm (the video image of the sperm in the inset). Three distinct regions in the sperm are observed: proximal head, distal head and tail, as well as small details as vacuoles. (C) Overlay of human sperm showing the locations of spectra found to be indicative of intact (purple) and damaged (aqua) DNA. (D) Averaged spectra showing the distinct shift at 1042 cm<sup>-1</sup> (arrow) in spectra of irradiated (red) compared with that of the untreated human sperm (green). (E) Scores cluster plot derived from PCA: native sperm (red and blue, for two observers) and UVB-irradiated sperm (green).

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
Figure 3 Histopathologic HE-stained section of normal (A and C) and of adenocarcinomatous (E) colonic tissues; B, D and F display pseudo-color Raman maps based on 11-means cluster (B and D) and 8-means cluster (F) on the same areas of the adjacent unstained section. Each cluster (consisting of similar spectra) was assigned to one color and corresponded to the extracted model.

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
Figure 4 Biochemical distribution of the individual constituent (proteins, collagen, lipids, DNA, and mucus) between normal, adenocarcinoma colonic tissues produced by applying the MLS algorithm. To determine the molecular composition of tissue composition, the relative scattering intensities of the various components were normalized to add up to 100%. The spatial and biochemical information obtained from these maps can be used to identify which biochemical markers could discriminate between normal and cancer tissues

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- [1] Z. Movasaghi, S. Rehman, and I.U. Rehman, Raman Spectroscopy of Biological Tissues, *Appl. Spectrosc. Rev.* 42, 493-541(2007)
- [2] G.J. Puppels, F.F.M. Demul, C. Otto, J. Greve, M. Robertnicoud, D.J. Arndtjovin and T.M. Jovin, Studying single living cells and chromosomes by confocal Raman microspectroscopy, *Nature* 327, 301-303, 1990
- [3] K.C. Schuster, I. Reese, E. Urlaub, J.R. Gapes and B. Lensl, Multidimensional Information on the Chemical Composition of Single Bacterial Cells by Confocal Raman Microspectroscopy, *Anal. Chem.*, 72, 5529-5534(2000)
- [4] W.E. Huang, R.I. Griffiths, I.P. Thompson, M.J. Bailey and A.S. Whiteley, Raman Microscopic Analysis of Single Microbial Cells, *Anal Chem.* 76, 4452-4458(2004)
- [5] W.E. Huang, K. Stoecker, R. Griffiths, L. Newbold, H. Daims, A.S. Whiteley and M. Wagner, Raman-FISH: combining stable-isotope Raman spectroscopy and fluorescence in situ hybridization for the single cell analysis of identity and function, *Environmental Microbiology* 9(8), 1878-1889(2007)
- [6] M. Li, J. Xu, M. Romero-Gonzalez, S.A. Banwart and W.E. Huang, Single cell Raman spectroscopy for cell sorting and imaging, *Curr. Opin. Biotechnol.* 23, 56-63(2012)
- [7] Y. Song, H. Yin and W.E. Huang, Raman activated cell sorting, *Curr Opin Chem Biology* 33, 1-8,(2016)
- [8] K. Maquelin, C. Kirschner, L.P. Choo-Smith, N. van den Braak, H.P. Endtz, D. Naumann and G.J. Puppels, Identification of medically relevant microorganisms by vibrational spectroscopy, *J. Microbiol Methods* 51, 255-271(2002)
- [9] M. Krause, B. Radt, P. Rosch and J. Popp, The identification of single living bacteria by a combination of fluorescence staining techniques and Raman spectroscopy, *J. Raman Spectrosc.* 38, 369-372(2007)
- [10] M. Harz, P. Rosch, J. Popp, Vibrational Spectroscopy - A Powerful Tool for the Rapid Identification of Microbial Cells at the Single-Cell Level, *Cytometry A*, 75A,104-113(2009)
- [11] J.F.M. Almarashi, N. Kapel, T.S. Wilkinson and H.H. Telle, Raman Spectroscopy of Bacterial Species and Strains Cultivated under Reproducible Conditions, *Spectroscopy: An International Journal* 27, 5-6, 361-365(2012)
- [12] S. Stockel, J. Kirchhoff, U. Neugebauer, P. Rosch and J. Popp, The application of Raman spectroscopy for the detection and identification of microorganisms, *J. Raman Spectrosc.* 47, 89-109(2016)
- [13] J.J. Barage, M.S. Feld, and R.P. Rava, In Situ Optical Histochemistry of Human Artery Using near-Infrared Fourier Transform Raman Spectroscopy, *Proc. Natl. Acad. Sci.* 89, 8, 3473-3477(1992)
- [14] J. Motz, M. Fitzmaruice, A. Miller, S. Gandhi, A. Haka, L. Galindo, R.R. Dasari, J. Kramer, and M. Feld, In vivo Raman Spectral Pathology of Human Atherosclerosis and Vulnerable plaque, *J. Biomed. Opt.* 11, 2, 021003(2006)
- [15] C. Matthaus, S. Dochow, G. Bergner, A. Lattermann, B.F.M. Romeike, E.T. Marple, C. Krafft, B. Dietzek, B.R. Brehm and J. Popp, In vivo Characterization of Atherosclerotic Plaque Depositions by Raman-Probe Spectroscopy and in Vitro Coherent Anti-Stokes Raman Scattering Microscopic Imaging on a Rabbit Model, *Anal. Chem.* 84, 7845-7851(2012)
- [16] J.A. Timlin, A. Carden and M.D. Morris, Chemical Microstructure of Cortical Bone Probed by Raman Transects, *Appl. Spect.* 53, 11 1429-1435(1999)
- [17] G.S. Mandair and M.D. Morris, Contributions of Raman spectroscopy to the understanding of bone strength,

- BoneKEy Reports 4, Article #620(2015)
- [18] O. Akkus, A. Polyakova-Akkus, F. Adar and M. Schaffler, Aging of Microstructural Compartments in Human Compact Bone, *J. Bone and Mineral Res.* 18, 6 1012-1019(2003)
- [19] N. Stone, C. Kendall and H. Barr, Raman Spectroscopy as a Potential Tool for Early Diagnosis of Malignancies, pp. 203-230, in *Esophageal and Bladder Tissues*, in *Vibrational Spectroscopy for Medical Diagnosis*, Ed. M. Diem, P.R. Griffiths and J.M. Chalmers, John Wiley ISBN 978-0-470-01214-7(2008)
- [20] K. Kong, C. Kendall, N. Stone and I. Notingher, Raman spectroscopy for medical diagnostics – from in-vitro biofluid assays to in-vivo cancer detection, *Adv. Drug Del Reviews*, 89 121-134(2015)
- [21] A.S. Haka, Z. Volynskaya, J.A. Gardecki, J. Nazemi, J. Lyons, D. Hicks, M. Firzmaurice, R.R. Dasari, J.P. Crowe, and M.S. Feld, In vivo Margin Assessment during Partial Mastectomy Breast Surgery Using Raman Spectroscopy, *Cancer Res* 66, 6 3317-3322(2006)
- [22] A. Beljebbar, S. Dukic, N. Amharref and M. Manfait, Ex vivo and in vivo diagnosis of C6 glioblastoma development by Raman spectroscopy coupled to a microprobe, *Anal Bioanal Chem* 398, 477-487(2010)
- [23] R. Kast, G. Auner, S. Yurgelevic, B. Broadbent, A. Raghunathan, L.M. Poisson, T. Mikkelsen, M.O. Rosenblum and S.N. Kalkanis, Identification of regions of normal grey matter and white matter from pathologic glioblastoma and necrosis in frozen sections using Raman imaging, *J. Neurooncol* 125, 287-295(2015)
- [24] J.W. Chan, D.S. Taylor, S.M. Lane, T. Zwerdling, J. Tuscano and T. Huser, *Anal. Chem.* 80, 2180-2187(2008)
- [25] S. Dochow, C. Krafft, U. Neugebauer, T. Bocklitz, T. Henkel, G. Mayer, J. Albert and J. Popp, Tumor cell identification by means of Raman spectroscopy in combination with optical traps and microfluidic environments, *Lab Chip* 11, 1484-1490 (2011)
- [26] J.J. McGarvey and J. Renwick Beattie, Raman Imaging, Chapter 7, *Raman Microscopy : A Versatile Approach to Bio-Imaging*, ISBN 978-3-642-28251-5, 2012
- [27] J. Renwick Beattie, Josephine V. Glenn, Michael E. Boulton, Alan W. Stitt, John J. McGarvey, Effect of signal intensity normalization on the multivariate analysis of spectral data in complex 'real-world' datasets, *JRS*, Volume 40, Issue 4, April 2009, Pages 429–435
- [28] Clark D, Sasić S. Chemical images: technical approaches and issues. *Cytometry A*. 2006 Aug 1; 69(8): 815-824.
- [29] Lavine B, Workman J. Chemometrics. *Anal Chem.* 2008 Jun 15; 80(12): 4519-31
- [30] Gendrin C, Roggo Y, Collet C. Pharmaceutical applications of vibrational chemical imaging and chemometrics: a review. *J Pharm Biomed Anal.* 2008 Nov 4; 48(3): 533-553
- [31] Schlücker S, Schaeberle MD, Huffman SW, Levin IW. Raman microspectroscopy: a comparison of point, line, and wide-field imaging methodologies, *Anal Chem.* 2003 Aug 15; 75(16): 4312-8
- [32] Krafft C, Knetschke T, Funk RH, Salzer R. Studies on stress-induced changes at the subcellular level by Raman microspectroscopic mapping. *Anal Chem.* 2006 Jul 1; 78(13): 4424-9
- [33] Bykov S, Lednev I, Ianoul A, Mikhonin A, Munro C, Asher SA. Steady-state and transient ultraviolet resonance Raman spectrometer for the 193-270 nm spectral region. *Appl Spectrosc.* 2005 Dec; 59(12): 1541-1552
- [34] Huang Z, McWilliams A, Lui H, McLean DI, Lam S, Zeng H. Near-infrared Raman spectroscopy for optical diagnosis of lung cancer, *Int J Cancer.* 2003 Dec 20; 107(6): 1047-1052.
- [35] C. Mallidis, J. Wistuba, B. Bleisteiner, O.S. Damm, P. Groß, F. Wubbeling, C. Fallnich, M. Burger, and S. Schlatt, In situ visualization of damaged DNA in human sperm by Raman microspectroscopy, *Human Reproduction*, Vol.26, No.7 pp. 1641–1649, 2011
- [36] A. Beljebbar, O. Bouché, M.D. Diébold, P.J. Guillou, J.P. Palot, D. Eudes, M. Manfait, Identification of Raman spectroscopic markers for the characterization of normal and adenocarcinomatous colonic tissues, *Critical Reviews Oncology/Hematology* 72(2009)255–264

### Fran ADAR, Ph.D.

Worldwide Raman Applications Manager  
HORIBA Instruments Inc.

### Catalina DAVID, Ph.D.

Application Scientist  
HORIBA FRANCE SAS

### Marinella SANDROS, Ph.D.

Business Development Manager for Life Sciences  
SPRI Product Manager  
HORIBA Instruments Inc.

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Tip-Enhanced Raman Spectroscopy  
Finds a Growing Number of Applications in Biology

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Tip-enhanced Raman spectroscopy (TERS) combines the high chemical specificity of Raman scattering with the nanoscale spatial resolution of Scanning Probe Microscopy (SPM). It enables label-free detection of surface components in multi-component samples with ultra-high spatial resolution (better than 20 nm). This article describes how TERS can be applied to study biological systems ranging from pure components like amino acids to membrane receptors found in complex environments such as cell surfaces.

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
Figure 1 Reflection and transmission TERS configurations.


Figure 2 TERS image of 100 nm × 100 nm (75 × 75 pixels) of a carbon nanotube showing an optical spatial resolution down to 8 nm; total map acquisition time < 9 minutes.

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
Figure 3 (a) AFM topography of engineered DNA, (b)-(c) corresponding TERS spectral mapping of over 50 × 20 nm<sup>2</sup> showing clear differentiation of spectral regions of pattern and size consistent with the expected (b) A/T and (c) G/C homopolymeric blocks, (d) horizontally averaged spectral map from the previous TERS maps, showing a good agreement with the (e) original sequence. (Data courtesy of Dr Noah Kolodziejski, Radiation Monitoring Devices.)

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**Marc CHAIGNEAU**

AFM/Raman Product Manager  
Sales & Marketing Division  
HORIBA FRANCE SAS  
Ph.D.

**Maruda SHANMUGASUNDARAM**

Application Scientist  
HORIBA Instruments Inc.

**Fran ADAR**

Raman Principal Scientist  
HORIBA Instruments Inc.  
Ph. D.

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Spectroscopic Ellipsometry Application in Life Science

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Nataliya NABATOVA-GABAIN


Figure 1 Setup of spectroscopic ellipsometer (UVISEL)


Figure 2 Correlation between refractive index and cell proliferation rate<sup>[3]</sup> (Correlation coefficient "r" = 0.71, Significant difference "P" = 0.0004)

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
Figure 3 (a) Optical model of tooth and (b) difference of effective roughness after brushing as function of dip into saliva time..

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- [1] 平川 誠一 他, HORIBA Technical Reports, 21, 26(2000)
- [2] 特開2014-57578
- [3] Y. Murayama, *et al.*, Abstract of 2016 Hong Kong International Conference on Engineering and Applied Sciences, HKICEAS-4905
- [4] 特開2016-128780

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