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Simple and Quick Quality Confirmation Method of Aqueous Urea Solution for Marine SCR System
Evaluation of Alkalinity of Aqueous Urea Solution with LAQUAtwin pH Meter

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The quality requirements of aqueous urea solution at 40% concentration for Marine Selective Catalytic Reduction (SCR) system is defined by International Standard (ISO-18611). One of them is the alkalinity which determined by potentiometric titration of free ammonia in the urea solution with hydrochloric acid solution. However, the titration is time-consuming work, and there are no crew being familiar with the chemical experiment, and no operation place in the ship. This paper shows compact water quality meter LAQUAtwin pH can be used for the simple and quick evaluation of aqueous urea solution quality.

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
$$\text{Alkalinity}(\%) = (0.017V) / m_s \dots\dots\dots (6)$$

V : The volume of the hydrochloric acid solution (0.01 mol/L) used for the titration (ml)
 m_s : The mass of the test portion of urea solution (g)

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
Figure 1 Model for the titration curve of 40% urea solution with 0.01 mol/L HCl


Figure 2 Comparison of neutralizing efficacy of HCl and H₂SO₄ in aqueous urea solution. (a) Neutralizing titration in fresh urea solution; (b) Neutralizing titration in degraded urea solution. 0.01 mol/L hydrochloric acid (○); 0.005 mol/L sulfuric acid (△)
 The degraded urea solution was which the fresh urea solution was stored under 60°C for 2 weeks.


Figure 3 Aqueous urea solution test kit

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
Figure 4 Correlation between ISO method and test kit method
The two symbols (●) and (○) show their individual series of the first and second time examinations.
The alkalinity 0.46% is regarded as 0.5% by rounding a number^[4].
The "X" represents the disagreement between the data of two methods.

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- [1] ISO 18611-1:2014, Ships and marine technology - Marine NOx reduction agent AUS 40 - Part 1: Quality requirements
- [2] ISO 18611-2:2014, Ships and marine technology - Marine NOx reduction agent AUS 40 - Part 2: Test methods
- [3] ISO 18611-3:2014, Ships and marine technology - Marine NOx reduction agent AUS 40 - Part 3: Handling, transportation and storage
- [4] JIS Z8401:1999 数値の丸め方

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HORIBA Advanced Techno, Co., Ltd.