Product Introduction

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Duetta: A New Fluorescence and Absorbance Spectrometer Providing Unique Two-in-One Benefits, Plus High Speed and NIR Fluorescence Detection

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In 2018, the HORIBA Instruments Fluorescence Division introduced a mid-market analytical spectrometer, called Duetta™, the first such “value, mid-market” instrument from the division. For decades now HORIBA has had a global leading market share of high-end research instruments, and the introduction of this new instrument opens up new market opportunities for HORIBA. The design of this new spectrometer concept, the sales success in the market in just its first year, and the key markets and applications of Duetta are all presented.

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Figure 1 Duetta two-in-one fluorescence and absorbance spectrometer


Figure 2 A fluorescence excitation spectrum (blue) and an emission spectrum (purple).

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

Figure 3 (Left) Duetta sample compartment shows the transmission/absorbance (red) detection optical path and the simultaneous fluorescence optical path (green) at 90 degrees from the excitation beam (blue). (Right) Duetta acquires multiple types of data in a single instrument.


Figure 4 (Left) Top view diagram of a cuvette with sample solution inside, showing incident light (I_0), transmitted light (I) and the path of fluorescence detection at 90 degrees (F). Primary and secondary IFEs on the detected fluorescence intensity are shown with blue and yellow stars, respectively. (Right) The IFE is a well-known physical phenomenon. This experiment shows that the fluorescence of fluorescein at varying concentrations exhibits a non-linear response due to IFE reabsorption (red), and the corrected linear response of Duetta (blue) when IFE correction is implemented. Thus, extending the concentration range wherein the fluorescence signal is quantitative.

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
Figure 5 On the left, Image A, Kinetic spectral scans of native BSA protein with rapidly added ANS (3×10^{-6} M). Spectra taken every 100 ms after ANS addition. As ANS binds to native BSA, the BSA Trp emission decreases and the ANS emission increases as a result of Förster resonance energy transfer (FRET), due to proximity of excited Trp and ANS. On the right, Image B, Kinetic spectral scans of SDS-denatured BSA after addition of ANS. Much higher ANS concentration (4×10^{-5} M) is required to affect binding to denatured BSA; no FRET is observed due to increased distances between ANS and Trp.

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
Figure 6 Neodymium laser glass emission spectrum. Ex: 575 nm, 5 nm band pass, Em: 700-1100 nm on fluorometer with PMT (red) and Duetta (blue). Standard PMT sensitivity falls off around 850 nm, while the Duetta CCD detector has sensitivity up to 1100 nm.

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
Figure 7 Fluorescence EEMs of two different concentrations of quinine sulfate in tonic water diluted in 0.1 M perchloric acid (aq.) with and without IFE corrections applied.

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
Figure 8 Fluorescence EEM spectra of Ag_A (top left), Ag_B (top right), Ag_C (bottom left) and the crude sample (bottom right). Dotted lines/arrows indicate the relative position of observed emission maxima. All plots are normalized to an intensity of 1 and were corrected for IFEs. The quantum yields of emission for Ag A, B and C are 0.323%, 0.232% and 0.083%, respectively. Reprinted with permission from [1] H. Ramsay, D. Simon, E. Steele, A. Hebert, R. D. Oleschuk and K. G. Stamplecoskie, RSC Adv., 8, 2018, 42080, Copyright 1999 the Royal Society of Chemistry.

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
Figure 9 Absorption (dashed), excitation (dotted) and emission (solid) spectra of (Pyr₂POAc)₃Yb(thf) in various solvents (2.6×10^{-6} mol L⁻¹). Excitation collected for emission centered at 970 nm. Emission collected with excitation at 355 nm. Slit widths: 5 nm (excitation/emission). Reprinted with permission from [2] Min Deng, Nathan D. Schley, and Gael Ung, Inorg. Chem., 57,15399(2018), Copyright 2018 American Chemical Society

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Cary Joseph DAVIES

Global Product Line Manager,
Fluorescence Division,
HORIBA Instruments Inc.

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R&D Dept. HORIBA, LTD.