# Readout

HORIBA Technical Reports

October 2021 No. **55**

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
The theme of 2021 Masao Horiba Awards is analytical and measurement technologies for development and production of biopharmaceuticals with high social needs, and vaccines or therapeutic drugs of today's great interest on the pandemic of COVID-19. This research field is directly related to human health and life, and is a field that is attracting global interest and investment.

I visited a forest detached from a famous sightseeing spot searching for the color of true autumn.

When I suddenly looked up, there was another world that seems out of this world with patterns created by red and green leaves.

For a while, I was so absorbed in taking pictures that I didn't notice the passage of time.

-Photographer MATSUI Hideo-  
(Member of Nikakai Association of Photographers)

#### Name of this Journal

This Journal is named "Readout" in the hope that "the products and technology we have created and developed will be read out and so become widely known".

## 2021 Masao Horiba Awards Optical/Spectroscopic Measurement Technologies for Life Science

### Foreword

- 4** HORIBA's Contribution to the Future of Pharmaceutical Innovation and Human Health  
*NISHIKATA Kentaro*

### Review

- 6** Eligible Fields and Award Winners for 2021 Masao Horiba Awards  
*NOGUCHI Shintaro*

### Feature Articles by 2021 Masao Horiba Awards Winners

- 8** Award Details
- 11** Development of Innovative Bio-Measurement Technology by Micro-Flow Light-Induced Acceleration  
*IIDA Takuya*
- 18** Development of Ultrafast Machine Vision-Activated Cell Sorters and Its Applications  
*OTA Sadao*
- 22** Elucidation of the Mechanism of Near-Infrared Light-Induced Cell Death and Method Establishment for Measuring Therapeutic Effects  
*SATO Kazuhide*
- 28** Real-Time Monitoring and Control of Pharmaceutical Production Processes Using Spectroscopic Data  
*Sanghong KIM*

### Guest Forum by 2021 Maso Horiba Awards Judges

- 34** The Role of Particle Design Studies in Developing Pharmaceutical Dosage Forms  
*TAKEUCHI Hirofumi*
- 41** Current Status and Objectives on Technologies in Pharmaceutical Industry  
*IBUKI Rinta*
- 46** Analytical Chemistry for Biopharmaceutics  
*TSUMOTO Kouhei*
- 54** Principles of Good Practice for Laboratory Medicine  
*MAEKAWA Masato*
- 59** Translational Raman Spectroscopic Approaches for Clinical Routine  
*Jürgen POPP*

### Feature Article

- 66** Improving Speed and Safety of Biopharmaceutical Manufacturing with HORIBA Analytical Solutions  
Innovative spectroscopic (A-TEEM and Raman) and optical techniques (ViewSizer) offer faster, accurate and less costly alternatives to chromatographic quantitative analysis in biotechnology  
*Linda Kidder YARLOTT, Adam M. GILMORE, Sean TRAVERS*
- 74** Fluorescent Bioprobes for Life Science Applications  
*Florian FORMANEK*
- 82** Providing Solutions for the Life Science Field Using Spectroscopic Analyzers  
*UCHIGASHIMA Mikiko*
- 88** Introduction of "Rapica" for Sensitive and Rapid Detection of ATP at the Level of One Microbial Cell  
*KOMAKI Naoto, NAKAYAMA Hideki, FUKAO Yoshiki, NAKAI Yoko*
- 93** History of HORIBA Medical Products to Contribute to the In Vitro Diagnostic Testing  
*NOMURA Naoyuki*
- 99** Disruptive Technologies in Haematology for POCT Market  
*Damien ISEBE, Philippe PIEDCOQ, Benoit THOUY, Aurélien DAYNES, Michael BRUCKNER*

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## Development of Innovative Bio-Measurement Technology by Micro-Flow Light-Induced Acceleration

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Dr. IIDA Takuya

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In the measurement of biological samples (proteins, nucleic acids, microorganisms, etc.) on the medical treatment, pharmaceutical and public health area, the conventional method is often complicated and requires long-term operation, and there is room for improvement in sensitivity. The authors have developed a microflow-type light-induced acceleration system (LAC-SYS) to realize rapid and highly sensitive measurement, where the measurement targets and probe particles are condensed by synergistic utilization of the “optical pressure” and “light-induced convection” in a narrow space of the microchannel for the acceleration of reaction when photoresponsive materials (substrate, particles) in liquid are irradiated with light. As a result, we achieved tens to hundreds of times higher sensitivity than the conventional method, and succeeded in quantitatively evaluating a very small amount (femtogram-level) of proteins in a few minutes. This technology innovates a wide range of bio-measurement such as drug development, personalized medicine, and food inspection.

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### Introduction

In recent years, a rapid, highly sensitive, low-cost, and compact biological sample detection system has been eagerly desired for the establishment of tailor-made drug discovery for personalized medical care, telemedicine using ICT technology, preventive medical care, and smart healthcare. Specifically, the cost reduction of preventive medical care and drug discovery, which leads to the reduction of medical expenses in Japan, which is increasing due to the declining birthrate and aging population, including the diagnosis of cancer and adult diseases, has become an important issue. In addition, due to the recent rapid growth of industrial activities in emerging countries and the problem of liberalization of imports and exports, there is also a need for technology such as genetic testing of food

and beverage production areas and the rapid and highly sensitive analysis and detection of foreign substances such as allergens and pathogenic bacteria. As conventional methods, an enzyme-linked immunosorbent assay (ELISA) using an antigen-antibody reaction is known for protein testing<sup>[1]</sup>, and a polymerase chain reaction (PCR) method is often used for nucleic acid testing such as DNA and RNA<sup>[2]</sup>. However, the conventional techniques require a high degree of professional pretreatment, and a process such as long-term incubation or complicated cleaning is required, and the detection time often takes several hours or more. There are also problems that a large amount of sample is required for analysis of a dilute target, and that reagents and devices used for fluorescent staining are expensive. In addition, since proteins cannot be amplified like nucleic acids, a breakthrough for the condensation of sample was

required for trace detection.

In order to solve these problems, technologies for non-destructive and non-contact condensation of biological samples for high-throughput analysis and measurement have been desired<sup>[3-5]</sup>. It is worth remembering that the technology of optical tweezers, which uses laser light to trap biomaterials contactlessly<sup>[6]</sup>, was awarded the Nobel Prize in Physics in 2018. Although the use of infrared wavelengths in this technology allows for low-damage manipulation, it was limited to precise motion control of a small number of biomolecules and cells. On the other hand, the authors theoretically predicted the possibility of controlling “plasmonic super-radiance” in which a large number of metallic nanoparticles (NPs) are densely arranged by optical trapping under thermal fluctuations in water at room temperature<sup>[7]</sup>. The plasmonic super-radiance is a quantum-mechanical effect, where the localized surface plasmon (LSP) in each NP strongly interact with each other via the light electromagnetic field for the increase of the radiative relaxation rate (the probability that an excited electron will emit energy as light). Based on this prediction, by focusing on the broadening of the LSP spectrum and the red-shift to the infrared region due to the increase of the radiative relaxation rate in a high-density metallic NPs, the principle of a photothermal biosensor was proposed, where pg-level proteins can be detected by accumulating them in the laser focal area due to the synergistic effect of the light-induced convection generated by the photothermal effect and the optical trapping of gold NP (AuNPs) fixed beads<sup>[8]</sup>. However, since many proteins and biological substances are sensitive to heat, it is necessary to condense them while controlling the thermal effect.

In this study, we have developed the “Light-induced Acceleration System (LAC-SYS)”, which enables high-throughput bioanalysis and measurement by controlling various biochemical reactions involving proteins, double-strand formation of nucleic acids, and condensation of microbes using laser irradiation. Particularly, we have focused on the development of a novel system that integrates microfluidic channels and LAC-SYS to increase the probability of interaction between the laser and NPs or microparticles (MPs) as probes in a narrow space. Furthermore, we also promoted the elucidation of the mechanism of optical condensation without damaging the target biomaterials, and aimed to carry out the entire process from the acquisition of new principles to the development of systems that will innovate spectroscopic analysis and measurement technology in the field of life science (Figure 1).


Figure 1 Schematic images of the light-induced acceleration system (LAC-SYS); (a) Basic principle, (b) microflow-type LAC-SYS.

## Detection of trace amounts of proteins by microflow optical condensation

As a main part of this research, we developed a micro-flow LAC-SYS by constructing a system that can introduce target molecules and small particles as probes, and irradiate them with a laser using pressure-driven flow by a syringe pump in a microchannel as shown in Figure 1(b)<sup>[9][10]</sup>. For example, in order to elucidate the basic principle in our literature<sup>[9]</sup>, we investigated the phenomenon that occurs when a protein is accumulated on the surface of a bubble (light-induced bubble) by using “optical pressure” derived from dissipative force, which is a kind of light-induced force, to accumulate metal NPs, and convection generated by the photothermal effect was also used as an auxiliary effect (Figure 2(a)). Specifically, an aggregate of AuNPs fixed on the surface of plastic beads (AuNP-fixed beads) were introduced as probe particles into a narrow microchannel of 400  $\mu\text{m}$  in width and 100  $\mu\text{m}$  in height, which is the same scale as blood vessels (for example, the thickness of small veins and arteries is 100-200  $\mu\text{m}$ ). They were used as probe to detect the target protein (in this case, albumin derived from egg white), and the research was carried out based on the idea that highly efficient trace detection would be possible by increasing the probability of interaction with laser light. In the search for conditions, we conducted experiments to launch and accumulate

AuNP-fixed beads on the ceiling of microchannel by optical pressure. Theoretical calculations confirmed that high-density metallic NPs on the surface of the beads can generate strong optical pressure even when using a laser in the infrared wavelength range, which is less damaging to living organisms, and we confirmed that AuNPs can efficiently generate optical pressure and light-induced convection in the preliminary experiments before our main research. As a result, AuNP-fixed beads and proteins were transported to the surface of the light-induced bubbles using the optical pressure and the convection generated by the photothermal effect (light-induced convection). As shown in the schematic image of Figure 2(b)(i), in the high concentration region (50  $\mu\text{g}/\text{mL}$  to 500  $\mu\text{g}/\text{mL}$ ), where the amount of protein in the observation region ranged from 340 pg to 3.4 ng, a large number of bubbles were generated and remained stable for a long time (Figure 2(c)).


Figure 2 (a) Schematic image of microflow-type LAC-SYS using AuNP-fixed beads as probe particles. (b) Schematic images of the phenomena depending on the amount of proteins. (c) Optical transmission image of multi-bubble generation when a large amount of proteins are added.

On the other hand, in the low-concentration region (0.5 ng/mL to 5  $\mu\text{g}/\text{mL}$ ), where the amount of protein in the observation region is 3.4 fg to 34 pg, a single bubble disappeared in 3 to 10 minutes when the laser was turned off after 10 seconds of irradiation from bubble generation.

In this case, the larger the amount of protein, the larger the aggregates remained (Figure 3(a), (b)). The amount of introduced target protein (in the observation region) is plotted on the horizontal axis, and the cross-sectional area of the bubble observed by optical microscopy is plotted on the vertical axis in Figure 3(c), and the cross-sectional area obtained from the transmission image of the aggregate is plotted on the vertical axis in Figure 3(d). It was found that the larger the amount of protein in the range of 1 fg to 100 pg, the larger the bubble was generated and the larger aggregate remained, which positively correlated with the size of the bubble, indicating the possibility of using this method for quantitative analysis of proteins. In particular, as shown in the inset of Figure 3(d), it shows high linearity and can be used for quantitative evaluation even when the amount of protein is only a few fg to a few hundred fg in the observation region, indicating the possibility of using it for trace analysis of proteins. We have previously reported that proteins as small as 1 pg can be detected by synergistically using light-induced force and light-induced convection in a quiescent fluid<sup>[8]</sup>, and this result shows that we have succeeded in detecting proteins with two to three orders higher sensitivity. Since this is similar to the aggregate formation when bubblegum shrinks, we call this new type of biosensor as a “bubblegum sensor”. Considering that protein adsorption occurs even after laser irradiation is stopped due to pressure-driven flow, we expect that thermal damage can be avoided since the concentration process is based on cooled bubbles. Furthermore, we have succeeded in the quantitative detection of antigens and antibodies on the order of sub-pg/mL concentration by accelerating the antigen-antibody reaction through non-thermal optical condensation using optical pressure as the main driving force<sup>[10]</sup>, and have elucidated the possibility of performing specific detection in just a few minutes with one to two orders higher sensitivity than conventional methods such as ELISA, suggesting the possibility of a breakthrough in protein measurement.

### Light-induced acceleration of DNA double strand formation

DNA has a strand-like structure consisting of four types of bases (adenine (A), thymine (T), guanine (G), and cytosine (C)). Base pairs with high complementarity, such as A-T and G-C, selectively bind to each other through a molecular recognition mechanism to form a double strand (hybridization). On the other hand, bases of the same species show low complementarity, called mismatching. AuNPs with a diameter of 30 nm were used as probe particles <I> chemically modified with 12-base, 4.08 nm long single-stranded DNA with a thiol group attached to the end as shown in Figure 4 (a), and the same size AuNPs modified with the same length of single-stranded DNA with the chemical


Figure 3 Relationship between the aggregate size and protein amount after the shrinkage of light-induced bubbles (for low concentrations).

(a), (b) Optical transmission images of the aggregates after bubble contraction for each protein concentration. Dependence of (c) the cross section (CS) of the bubble and (d) the area of the aggregate obtained after bubble contraction on the mass of the protein present in the observation region (inset: magnified view and approximate curve in the low concentration region)

structure of the end reversed were used as probe particles <I>. Each of them was dispersed in phosphate buffer solution and used. Dispersions of these probe particles <I> and <II> and a dispersion of single-stranded target DNA (24 bases) as a target were dropped onto a glass substrate, mixed, and allowed to diffuse sufficiently, and we try to control hybridization with the single-stranded DNA on the surface of the probe particles by irradiation of an infrared laser focused strongly with a high-magnification objective lens.

When the laser beam is focused near the droplet surface as a non-equilibrium open system, the convection due to the photothermal effect of the AuNPs and the condensation due to the evaporation are expected to occur, and the probe DNA on the surface of AuNP and the target DNA are expected to form a double strand efficiently. Therefore, even with a dilute concentration of 100 pM (about 300 million DNA strands were contained the 5 µL solution) of complementary strand DNA, a macroscopic network-like assembly of about 0.1 mm can be formed after only 2 min of laser irradiation (Figure 4(b)). On the other hand, no large aggregates were obtained when the same concentration


Figure 4 Light-induced acceleration of DNA double strand formation.

(a) Schematic image of the principle. Optical transmission images of light-induced network formation for (b) full-match DNA and (c) mismatch DNA added to DNA on the surface of the probe particle as target DNA. (d) Dependence of the peak wavelength of the absorbance spectra for each sequence on the laser irradiation time.

of fully mismatched DNA was used (Figure 4(c)), and it was confirmed that the difference in the base sequence also affects the aggregation phenomenon induced by light irradiation. Also, we measured the change in the local extinction spectrum of the white light transmitted through the region surrounded by the dashed lines, and observed a significant broadening and red-shift in a few minutes only when the complementary DNA was added, which is similar tendency to the theoretical prediction in the literature<sup>[7]</sup>. Even in the same full-match case, the spectral shift was larger for sequences with higher complementarity, such as “sequence A” containing G-C bonds, than for “sequence B” with only A-T bonds (Figure 4(d)). Furthermore, the number of DNA in this photometric region before laser irradiation was estimated to be about 5 zmol estimated from the above concentration (zmol: ~10<sup>2</sup> molecules), indicating that the formation of macroscopic aggregates that can be fully observed by optical microscopy was achieved in only a few minutes by triggering the double-strand formation of a very small amount of DNA. On the other hand, any aggregate could not be observed in a mixture of probe particles and target DNA at the same concentration under an optical microscope simply by letting the mixture stand

for 3 hours, but macroscopic aggregates could be formed very quickly by accelerating the reaction under optical condensation. This is an important achievement that shows that specific detection of trace DNA is possible through this process.

### Principle of damage-free light-induced assembly of microbes by bubble-mimetic substrate

We have measured the number of bacteria in liquid based on the optical condensation of microbes such as bacteria using “light-induced convection” and “bubble” generated when a substrate coated with a metallic nanofilm exhibiting a finite absorption rate in the wavelength range of the irradiated light is locally heated by laser irradiation<sup>[12]</sup>, and we have also performed the optical condensation of useful bacteria with the high viability using a multi-pore substrate<sup>[13]</sup>. In these studies, there were some problems such as the reduction of the survival rate when the laser power was increased. However, we also clarified the principle of “damage-free light-induced assembly” by coating the upper part of the solid-liquid interface with a metallic nanofilm to mimic a bubble and minimize the thermal effect on the trapping site<sup>[14]</sup>.


Figure 5 Low-damage light-induced acceleration of microbial condensation by controlling the position of heat sources.

- Schematic of optical condensation using an imitation bubble,
- Schematic of optical condensation using a conventional flat photothermal substrate.
- Electron micrograph of an example of microbe.
- Fluorescence images of optically-condensed bacteria.
- Laser power dependence of the number of aggregated bacteria.
- Laser power dependence of the survival rate of assembled bacteria.

As shown in Figure 5(a), polystyrene beads with a diameter of 100 μm were chemically fixed onto a glass bottom dish as an imitation bubble, and a thin platinum (Pt) film was coated on top of the bead by ion sputtering to create a “bubble-mimetic substrate (BMS)”. In particular, *Lactobacillus casei* (*L. casei*), a type of lactic acid bacteria, was used as a target, and green fluorescent dye (SYTO 9), which stains all the bacteria including live ones, and red fluorescent dye (PI), which stains only dead ones, were used to determine the viability. The optical system used for the experiment was an inverted microscope as shown in Figure 1(b), where a laser was defocused from below onto a Pt thin film coated on the top of a 100 μm diameter bead, and optical condensation was performed in the power range of about 15-60 mW without syringe pump.

Figure 5(d) shows the main results of the light-induced assembly of fluorescently stained *L. casei*. We found that about 40,000 bacteria could be aggregated after 5 minutes of irradiation with 27 mW laser, and a high aggregation rate of the order of  $10^4$ - $10^5$  bacteria was achieved after only a few minutes of laser irradiation (Figure 5(e)). Furthermore, we found that the survival rate of higher than 95%, which is almost the entire population, can be maintained while using convection due to the photothermal effect (Figure 5(f)). In order to elucidate the mechanism of the highly efficient light-induced assembly and the high survival rate, we evaluated the convection velocity by simulations combining theories of electro-dynamics and thermo-fluid dynamics, and estimated the number of bacteria entering the region around the imitation bubble to be about 50,000, which is consistent with the order of the number of bacteria trapped in the experiment. Furthermore, it was found that the bubble-mimetic substrate causes less damage even with high laser power because the heat source and trapping site are spatially separated. On the other hand, as shown in Figure 5(b), the conventional method of generating a bubble on the surface of a flat substrate<sup>[12]</sup> inevitably causes damage to the microbes due to the high temperature near the laser irradiation point, and the survival rate decreased with higher power (Figure 5 (f)).

### Conclusion

Based on the accelerating assembly process of biological nanomaterials and probe particles introduced into a micro-flow system using “optical condensation” by the electro-magnetics and thermohydrodynamics actions of light, we have clarified that it is possible to achieve a dramatically high-throughput in terms of measurement time, sensitivity, and sample volume for the first time in the world. In particular, we clarified the highly efficient photothermal effect of densely assembled metallic NPs arising from quantum effect and utilized it for “optical condensation

of biological samples,” and clarified a new principle that enables the detection of proteins on the order of fg ( $10^{-15}$  gram) in just a few tens of seconds to a few minutes by condensing them with laser irradiation. Remarkably, we have obtained initial results on “light-induced acceleration” of antigen-antibody reaction by extending this principle<sup>[10]</sup>, and have found that it can speed up protein detection by 2-3 orders of magnitude faster (from 5 hours to 3 minutes), increase sensitivity by 1-2 orders of magnitude higher (from pg-order to fg-order), and drastically reduce sample volume (from several hundred  $\mu\text{L}$  to several hundred nL) by eliminating the multi-step process that has been the bottleneck in conventional ELISA and a significant reduction in sample volume (from several hundred  $\mu\text{L}$  to several hundred nL). Moreover, we succeeded in light-induced acceleration of DNA double strand formation on the order of zmol (hundreds of DNA molecules), which could be applied to rapid and highly sensitive genetic testing. In addition, they succeeded in developing a substrate that can be used to measure the number of bacteria in just a few minutes, which takes several days by the conventional culture method.

These results have a high potential to open up an important

core technology for the construction of a next-generation platform for spectroscopic analysis and measurement technology for biological science, which boasts overwhelming speed and sensitivity. As shown in Figure 6, our developed micro-flow LAC-SYS can contribute to the creation of a platform for spectroscopic analysis and measurement technology in the life science field, which will lead to ultra-smart applications in drug discovery, medical testing, food testing, and environmental measurement, as well as to scientific and technological innovations that transcend the fields of physics, chemistry, and biology regarding the optical control of various biological functions.

The basic structure of the system is based on a simple optical system, which makes it easy to be compact and portable, and it is expected to be developed into a “future spectroscopic analysis and measurement technology” useful for smart drug discovery, self-care medicine and telemedicine, and on-site food inspection. Furthermore, our technology will contribute to the development of low-cost and simple devices for evaluating the local effects of rare and expensive new drugs and fluorescent probes in cancer treatment and screening, and to environmental measurements such as the measurement of biological genes (DNA, RNA) and harmful particles in the environment, including advanced measurement technology for biological cells and microbes.

### Acknowledgement

The author would like to express his sincere thanks to Prof. Shiho Tokonami, the Deputy Director of the RILACS of Osaka Prefecture University, for her support through joint research, and to Dr. Mamoru Tamura, Dr. Yushi Nishimura, Dr. Yasuyuki Yamamoto, Ms. Mayu Ueda, Mr. Kota Hayashi, and other alumni and graduate students for their cooperation in carrying out the research. Also, the author would like to express his deep gratitude to Prof. Hajime Ishihara and Prof. Syoji Ito of Graduate School of Engineering Science in Osaka University, and Prof. Ikuhiko Nakase, the Assistant Director of RILACS, for their fruitful discussions. This work was partially supported by JST-Mirai Program (No. JPMJMI18GA, No. JPMJMI21G1), Grant-in-Aid for Scientific Research (A) (No. 17H00856, No. 21H04964), Grant-in-Aid for Scientific Research (B) (No. 18H03522), Grant-in-Aid for Scientific Research on Innovative Areas (No. 16H06507) from JSPS KAKENHI, AMED (No. JP20he0622017), and the Key Project of Osaka Prefecture University.

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Figure 6 Future image of innovation in medical testing and drug discovery process by micro-flow type LAC-SYS.

## References

- [1] E. Engvall, P. Perlmann, "Enzyme-linked immunosorbent assay (ELISA) quantitative assay of immunoglobulin G", *Immunochemistry*, 8, 871-874 (1971).
- [2] R. K. Saiki, S. Scharf, F. Faloona, K. B. Mullis, G. T. Horn, H. A. Erlich, N. Arnheim, "Enzymatic amplification of beta-globin genomic sequences and restriction site analysis for diagnosis of sickle cell anemia", *Science*, 230, 1350-1354 (1985).
- [3] E. A., Raiber, R. Hardisty, P. van Delft, et al. "Mapping and elucidating the function of modified bases in DNA", *Nat. Rev. Chem.* 1, 0069 (2017).
- [4] A. Altaalar, J. Munoz, A. Heck, "Next-generation proteomics: towards an integrative view of proteome dynamics", *Nat. Rev. Genet.* 14, 35 (2013).
- [5] Y. Bian, R. Zheng, F. P. Bayer, et al. "Robust, reproducible and quantitative analysis of thousands of proteomes by micro-flow LC-MS/MS". *Nat. Commun.* 11, 157 (2020).
- [6] A. Ashkin, J. M. Dziedzic, J. E. Bjorkholm, S. Chu, "Observation of a single-beam gradient force optical trap", *Opt. Lett.* 11, 288 (1986).
- [7] T. Iida, "Control of Plasmonic Superradiance in Metallic Nanoparticle Assembly by Light-Induced Force and Fluctuations", *J. Phys. Chem. Lett.*, 3, 332 (2012).
- [8] Y. Nishimura, K. Nishida, Y. Yamamoto, S. Ito, S. Tokonami, T. Iida, "Control of Submillimeter Phase Transition by Collective Photothermal Effect", *J. Phys. Chem. C*, 118, 18799 (2014).
- [9] M. Ueda, Y. Nishimura, M. Tamura, S. Ito, S. Tokonami, T. Iida, "Microflow-mediated optical assembly of nanoparticles with femto-gram protein via shrinkage of light-induced bubbles", *APL Photon.* 4, 010802 (2019).
- [10] T. Iida, S. Tokonami, I. Nakase, "Detection method for substance to be detected and detection system for substance to be detected", *PCT/JP2020/032758* (2020).
- [11] T. Iida, Y. Nishimura, M. Tamura, K. Nishida, S. Ito, S. Tokonami, "Submillimetre Network Formation by Light-induced Hybridization of Zeptomole-level DNA", *Sci. Rep.* 6, 37768 (2016).
- [12] Y. Yamamoto, E. Shimizu, Y. Nishimura, T. Iida, S. Tokonami, "Development of a rapid bacterial counting method based on photothermal assembling", *Opt. Mater. Exp.*, 6, 1280 (2016).
- [13] S. Tokonami, S. Kurita, R. Yoshikawa, K. Sakurai, T. Suehiro, Y. Yamamoto, M. Tamura, O. Karthaus, T. Iida, "Light-induced assembly of living bacteria with honeycomb substrate". *Sci. Adv.* 6, eaaz5757 (2020).
- [14] K. Hayashi, Y. Yamamoto, M. Tamura, S. Tokonami, T. Iida, "Damage-free Light-induced Assembly of Intestinal Bacteria with a Bubble-mimetic Substrate", *Commun. Biol.*, 4, 385 (2021).

Dr. IIDA Takuya

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Professor, (Concurrent) Director of RILACS  
 Department of Physical Science,  
 Graduate School of Science,  
 Osaka Prefecture University.  
 (Concurrent) Research Institute for Light-induced  
 Acceleration System (RILACS)

## Development of Ultrafast Machine Vision-Activated Cell Sorters and Its Applications

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OTA Sadao

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A cell sorter that performs image information analysis in real-time to separate a large number of cells were desired. However, there is always a tradeoff between the processible information per cell (quality) and the number of cells per time (quantity). Therefore, it has been a great challenge to enable high throughput imaging cell sorter that holds the advantages of optical microscope (high quality) and flow cytometry (high quantity simultaneously). In order to overcome this challenge, we proposed and realized a new approach named “ghost cytometry”, originating from a concept that image reconstruction is not always necessary in image analysis when performed by machines, not by humans. This method utilizes the motion of each cell in microchannels to acquire its compressed image signal by a single pixel detector, and directly applies AI-based analysis to the signal, resulting in the development of the ultrafast and accurate image-free “imaging” cell sorter. This technology is expected to be widely applied in various biotechnologies and cell-based clinical methods.

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### Introduction

Imaging and analyzing many single cells hold the potential to significantly increase our understanding of heterogeneous systems involved in various complex life systems and diseases. Many key applications in these fields require accurate and high throughput isolation of specific populations of cells according to information contained in the high content images. However, there is always a trade-off between the processible information per cell (quality) and processible number of cells per time (quantity). This

trade-off is depicted in Figure 1. For cytometry technologies, this trade-off caused critical challenges that have to be simultaneously solved: (1) ultrafast imaging technology that simultaneously meets the needs of high sensitivity, polychromaticity, high speed and continuous acquisition, and (2) ultrafast information technology that continuously performs both computational image production and analysis which is costly in terms of both time and money. As a result, there has been no high throughput imaging-activated cell sorting technology that holds advantages of both optical microscopy (high quality) and flow cytometry (high

quantity).

In a series of works, we proposed and experimentally demonstrated that directly applying machine learning methods to compressed imaging signals measured using a single pixel detector enables ultrafast, sensitive, and accurate image-free (without image production), imaging-based cell analysis as well as sorting in real time, which we call ghost cytometry<sup>[1]</sup> (GC, Figure 1). GC has enabled the fast cell sorting based on the cells' image information in both fluorescence and label-free modes, which we name machine-vision based cell sorters (ViCS). Furthermore, using this ViCS, we have developed an ultrafast pooled platform for cell image-based phenotypic screening of genes, which can be scaled up at a large scale.


Figure 1 Ghost cytometry simultaneously achieved speed and rich information by “image analysis without image reconstruction”

Moreover, based on the invention, we founded a start-up company ThinkCyte Inc., which develops these ViCS machines through the co-development with renowned companies. Further, we are now striving to realize its practical use in the fields of medical diagnosis, drug discovery, and cell therapy, and a series of joint research with several pharmaceutical companies is underway.

### Key feature of our technology 1: ghost cytometry

Our key concept was ghost cytometry (GC), which is an approach of image analysis without image production. In GC, we applied machine learning methods to compressed imaging signals measured using a single pixel detector for enabling ultrafast, image-free, “image”-based cell analysis as well as sorting in real time (Figure 2). Thanks to skipping the time-consuming image reconstruction process, the time required for image inferring was reduced down to 10 microseconds which is orders of magnitudes shorter than other methods relying on image reconstructions. With GC, we developed a series of ultrafast cell sorters which classifies cell based either fluorescence image information or label-free image information. In the case of the fluorescence mode, fluorescent signals are detected by converting image information into temporal waveforms compressively. In the case label-free modes, refractive index distribution of

cells is detected; these label-free morphological modalities obtained are similar label free microscopy images including darkfield and bright field images.

Figure 3 shows an example workflow of GC utilizing a machine learning model to classify the cells by directly analyzing the imaging waveforms. It first starts with preparing a training data set: for each cell, we simultaneously measure the GC imaging waveforms and molecular labels that reflect types, states, phenotypes, and other characteristics of cells. Using this training data set, we then develop a model based on machine learning methods. Finally this model predicts the labels directly from the GC imaging signals of unknown cells.

In the fluorescence modes, GC can classify localization of fluorescence molecules in the cell morphology, even if their total intensities are the same or if the same kind of cell lines are used. This is particularly useful when users know what image pattern to recognize and what kind of specific molecules are their targets which is often the case in drug screening or biological research. With these pre-knowledge, one can chemically or genetically label the molecule of interest to visualize the fluorescent image pattern of interest with a good signal-to-noise-ratio. Using various optical probes, the GC can classify a variety of fluorescent patterns effectively, including nuclear translocation of molecules, organelle morphologies, and protein-protein interactions. On the other hand, that label-free mode can be useful especially when users don't want to stain cells with chemical or biological labels or damage them before the downstream uses. In addition, this approach is also effective when users are not clear about what the target molecules are or what image patterns should be specifically recognized. This label-free GC has demonstrated accurate classification of cells based on their kinds, states, and functions without observing molecular labels.


Figure 2 Image-free imaging cytometry driven by machine learning: In conventional workflows of image-based cytometry, imaging technologies were to reconstruct a clear 2D/3D images for human recognition, followed by analysis based on human recognition. Our radical idea of ghost cytometry (GC) is, “apart from human recognition-based imaging cytometry, reconstructing 2D images is not necessary for machine learning-based analysis (AI)”. In GC, its compressive imaging technology works to extract cell’s imaging information as a temporal waveform of signals measured by a single pixel detector, followed by analysis without image reconstruction using machine learning methods. Thereby, GC enabled the world’s first and fastest image-free “imaging” cell sorter by skipping the most time-consuming and computer-intensive image reconstruction in the data-intensive flow cytometry.

Example workflow of GC utilizing supervised machine learning


Figure 3 Example workflow of ghost cytometry utilizing supervised machine learning methods

Key feature of our technology 2: machine vision-activated cell sorting

By combining GC’s efficient signal processing implemented in a field programmable gate array with a microfluidic device, we realized ultrafast and accurate cell sorting based on real-time analysis of “imaging” data. Figure 4 shows its first prototype and, since then, this machine vision-activated cell sorter (ViCS) has evolved to be more stable and more flexible one through development by ThinkCyte Inc. together with excellent other companies.


Figure 4 The first prototype of ultrafast machine vision-based cell sorter (ViCS), enabled by combining the optical ghost imaging signal acquisition, the real-time machine learning-based analysis implemented in FPGA, a microfluidic cell sorting technology.

Application

In the field of drug discovery, we have developed a pooled cell phenotypic screening method using ViCS. In this application, the fluorescence mode ViCS is used when the cell phenotype is specifically known (i.e. the localization of a specific molecule or a specific intracellular organelle), and the label-free mode ViCS is used when no specific explicit imaging feature is available. While high-content screening used to take long time (several months) and high cost, our method can significantly accelerate it at lower cost. We believe we can contribute to society by greatly accelerating drug discovery research and development.

In the field of cell therapy, regenerative medicine, and cell-based production, the label-free GC and ViCS can evaluate cell types and functions without using stains or antibodies, and to enrich cells of interest. This method can be thus used for automatic monitoring of cell production lines and improving cell qualities by removing unwanted cells or contaminating particles. Similarly, it can be used in the evaluation of cells that produces antibodies and drugs. While remarkably effective cell-based drugs currently seem to put pressure on medical costs, such a system may greatly contribute to its reduction.

## Conclusion

In conclusion, we proposed and realized “ghost cytometry”, an approach based on AI-assisted cell analysis of their “image” information without image production. More concretely, by directly applying AI-based analysis to a compressed “image” signals in temporal domains, we skip a time-consuming image reconstruction process and realized the ultrafast and accurate image-free “imaging” cell sorter. We expect this technology to be widely adopted in biological sciences, and industrial applications including biotechnology, drug discovery, and cell therapy and regenerative medicine.

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**OTA Sadao**

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Associate Professor  
Research Center for Advanced Science and Technology,  
The University of Tokyo  
Ph.D.

## Elucidation of the Mechanism of Near-Infrared Light-Induced Cell Death and Method Establishment for Measuring Therapeutic Effects

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### SATO Kazuhide

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
Near-infrared photoimmunotherapy (NIR-PIT) is a new cancer treatment that involves the application of a near-infrared light probe (IR700) to a specific antibody that binds to antigens on the surface of cancer cells. However, despite its impressive efficacy and innovation, the detailed mechanism of cell death has remained unclear. From the viewpoints of analytical chemistry, condensed matter chemistry, and optics, it has been clarified that NIR-PIT cell death is the world's first photochemical-induced cell death, and the cell death can be measured and quantified by IR700 near-infrared fluorescence spectroscopy at 700 nm, thus establishing the basis for a therapeutic biomarker. Due to the novelty of this cell death, this therapy was recognized as a new cancer modality and was implemented clinically in Japan ahead of the rest of the world.

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### Introduction

Near-Infrared photoimmunotherapy (NIR-PIT) is a next-generation cancer-targeting therapy in which antibodies specific to antigens on the surface of cancer cells are attached to probes that react to near-infrared light and are then locally irradiated with near-infrared light (Figure 1)<sup>[1]</sup>. This photo-targeted therapy is a futuristic, and multidisciplinary treatment method in which probes (chemistry) and antibodies (pharmacology and biology) are conjugated under optimal conditions (conjugation chemistry), and near-infrared light (optics and physics) is added (Figure 1). Irradiation with near-infrared light at 690 nm, the peak absorption of the probe, immediately expands the cancer cells and induces cell death. It is currently undergoing international Phase III trials (LUZERA-301)<sup>[2,3]</sup> and was approved by PMDA in Japan in September 2020 for the treatment of recurrent and previously treated stage IV head and neck cancer<sup>[4]</sup>. It is now being recognized as

the fifth cancer treatment after surgery, chemotherapy, radiotherapy, and cancer immunotherapy. So far, various application studies of NIR-PIT therapy have been conducted all over the world, and further development of this therapy is expected from here. Although this therapy has attracted attention for its innovation, the detailed mechanism of cell death in NIR-PIT had not been fully elucidated at the time of the introduction of Phase III trials. This was a barrier to its approval and implementation. In this review, we will provide an overview of NIR-PIT, which is attracting attention as a “fifth cancer therapy”, its cell death mechanism, and the prediction of its efficacy by near-infrared spectroscopic imaging derived from its elucidation.

**Figure 1** Overview of NIR-PIT  
NIR-PIT is a new therapeutic method in which an antibody-IR700 complex, in which IR700DX, a light-absorbing substance, is added to the antibody, binds to cancer-specific antigens that are highly expressed on the surface of target cancer cells, and is then exposed to NIR light.

## About NIR-PIT

There are four known treatment methods for cancer: (1) surgery (surgical and endoscopic resection), (2) chemotherapy (anti-cancer drugs), (3) radiotherapy, and (4) cancer immunotherapy. These existing cancer treatments cause injury not only to cancer cells but also to surrounding normal tissues and organs. In anticancer drug therapy, chemical substances that cause damage to both normal cells and cancer cells are administered, but since the dosage is based on an extremely narrow dosage design that allows normal cells to recover earlier than cancer cells, the dosage is strictly determined, and if an overdose occurs, not only the cancer but also the individual may die. In addition, surgery requires the resection of normal tissues surrounding the cancer as margins, and radiation therapy inevitably causes injury to normal cells within and through the irradiation area. Therefore, it would be ideal if there were a cancer-targeted therapy that could superspecifically damage only cancer cells in the body without damaging normal cells and tissues.

Based on the above idea, Dr. Hisataka Kobayashi and his colleagues at the National Cancer Institute (NCI) in USA announced the development of Near-Infrared Photoimmunotherapy (NIR-PIT) in 2011<sup>[1]</sup>. NIR-PIT is a photo-targeted therapy that combines near-infrared light and antibody photo-probe adducts, in which IR700, a light-absorbing probe, is covalently attached to antibodies that bind to cancer-specific antigens. NIR-PIT is super-selective because it is doubly targeted by NIR-light and antibodies (Figure 1).

NIR-light is used in TV remote controls, CD players, etc., and is a safe light that has almost no adverse effects when irradiated to the human body. IR700 is a probe that absorbs

light in the near-infrared region of 690 nm and emits fluorescence at 700 nm, and was originally developed by LICOR as a dye for microscopic observation in biology. Antibodies are biomaterials that have been actively used in medicine, biology, and engineering in recent years, and are among the existing biomolecules that can bind specifically and stably to molecules on the surface of cell membranes. Among the existing biomolecules, antibodies can bind specifically and stably to cell surface molecules<sup>[4]</sup>. As described above, NIR-PIT is a multidisciplinary therapy that combines biology, pharmacology, chemistry, optics, engineering, and physics.

## Mechanisms of cell death in NIR-PIT

It has been proposed that oxidative stress is the main mechanism of cell death using photosensitive substances<sup>[5]</sup>. This pathway is often considered to be one in which the energy difference between photoexcitation and photoemission of a photosensitive substance acts on mitochondria and other organs, resulting in cell injury by oxidative stress agents such as singlet oxygen and free radicals, leading to cell apoptosis. The detailed mechanism of cell death in NIR-PIT has been unknown for a long time, which has been a barrier to its clinical implementation. Since the cell death reaction in NIR-PIT proceeded even after the cell function was stopped at 4°C, and the inhibition of cell death was not sufficiently effective even when oxidative stress inhibitors (free radical scavengers) were added, it was inferred that the oxidative stress was not the main mechanism for NIR-PIT. Therefore, it is inferred that the mechanism of NIR-PIT is not based on oxidative stress as in conventional photodynamic therapy (PDT).<sup>[5]</sup>

NIR-PIT has a very different therapeutic mechanism compared to conventional cancer cell death. When the antibody-IR700 conjugate bound to the cell surface is irradiated with near-infrared light (690 nm), the cell surface is rapidly injured, the cell morphology changes, and then the cell bursts, leaking the cytoplasmic contents and causing cell death within a few minutes. Thus, the mechanism is necrosis (rapid destruction of cells) rather than apoptosis, which is the conventional slow death of cancer cells. Since near-infrared light irradiation itself is harmless to cells, it can only damage the target cells to which the antibody-IR700 conjugate binds. In fact, when cells with and without antibody-IR700 conjugates were co-cultured and irradiated with NIR light at the same time, cell death was observed only in the cells with antibody-IR700 conjugates, while the neighboring cells without antibody-IR700 conjugates were not injured at all. The antibody-IR700 conjugate can selectively destroy not only cells in vitro but also target cells in vivo in space and time.<sup>[6]</sup>

### Aggregation by photochemical reaction of IR700 and detection of released products by mass spectrometry


Figure 2 Aggregation by photochemical reaction of IR700 and detection of released products by mass spectrometry

- The photochemical reaction of IR700 as the starting point of NIR-PIT action: the silanol group attached to the silicon of the silicon phthalocyanine ring of IR700 is released by near-infrared light irradiation, and as a result, the hydrophilic nature of IR700 suddenly becomes hydrophobic and aggregates. As a result of aggregation, the fluorescence is lost.
- The silanol groups released from the supernatant after irradiation were detected by mass spectrometry (IT-TOF MS), proving the reaction.
- IT-TOF-MS analysis showed a peak of silanol group only at 16 J/cm<sup>2</sup> near-infrared light irradiation. By product scan (MS3) analysis, the structure was determined and proved to be silanol group.

Originally, IR700 dye (IR700 DX) was developed and marketed with the merit that its fluorescence hardly fades compared to other dyes (AlexaFluor680, etc.) during the number of scans in the fluorescence observation environment of microscope slide samples. However, the accumulation of *in vitro* and *in vivo* experiments with NIR-PIT clearly showed that the fluorescence of IR700 faded, which was a contradiction. As a hypothesis to resolve this discrepancy, we hypothesized a change in the chemical nature of IR700 due to photochemical reactions. The hypothesis is that the chemical nature (hydrophilic) of IR700 is rapidly changed to hydrophobic by the release of the silanol side chain by near-infrared light. The hydrophobized IR700 forms aggregates in the aqueous solution, and as a result of the aggregation, the fluorescence is lost from the solution in a series of photochemical reactions (Figure 2a)<sup>[7]</sup>. It was thought that the aforementioned discrepancy could be explained by the difference between the inclusion state of the slide samples using non-water-soluble inclusion material and the water-rich biological environment.

Based on these hypotheses, in order to prove this mechanism, we attempted to measure the silanol groups (molecular weight 500) above and below the silicon phthalocyanine ring of IR700, which are thought to be released upon aggregation, by mass spectrometry. IT-TOF-MS was able to detect the peak of the molecular weight only in the case of light irradiation (Figure 2b). In addition, it has been found that this liberation reaction is enhanced in the

absence of oxygen and in an environment rich in electron donors, while the presence of electron absorbers inhibits the reaction. This indicates that this reaction is independent of oxygen and free radicals, and has the advantage that the reaction proceeds more easily in tumors under hypoxic conditions.

Next, in order to prove that the above reaction also occurred in the antibody-IR700 adduct, cetuximab-IR700 (cet-IR700) was irradiated in the near-infrared to prove that silanol could be measured and detected by mass spectrometry. In addition, the solution of cet-IR700 was confirmed by the product scan method to be of IR700 origin, not antibody origin, with the structure of the silanol group (Figure 2c).

SDS-PAGE analysis was performed to study the effect of photochemical reaction of IR700 on antibody IR700 adduct and antibody IR700 adduct-antigen protein complex. First, we found that NIR light irradiation (peak 690 nm NIR light laser irradiation) caused the antibody (cetuximab) band to smear on SDS-PAGE, resulting in the loss of IR700 fluorescence. The photochemical reaction of IR700 resulted in the protein denaturation, i.e., aggregation of the added antibody proteins (Figure 3a). As shown in the figure, the denaturation increased in a near-infrared light-dependent manner, and the aggregation was almost complete at 64 J/cm<sup>2</sup>. It is also clear that the near-infrared fluorescence is lost due to aggregation. In order to investigate the effect

**Examination of protein aggregation and loss of fluorescence of cet-IR700 by near-infrared irradiation and changes in the antibody-IR700 adduct-antigen protein complex using SDS-PAGE**


Figure 3 Examination of protein aggregation and loss of fluorescence of cet-IR700 by near-infrared irradiation and changes in the antibody-IR700 adduct-antigen protein complex using SDS-PAGE  
 a. The band of cet-IR700 became a smear in a light-dependent manner and lost its fluorescence, indicating that it was aggregated.  
 b. When cet-IR700 was bound to an antigen protein (EGFR) and irradiated with near-infrared light, the band of the antigen protein also became smeared, indicating that not only cet-IR700 but also its binding partner antigen was aggregated. With aggregation, fluorescence was also lost.

of the photochemical change of cet-IR700 on the binding partner antigen protein, we conducted a near-infrared irradiation experiment with cet-IR700 bound to EGFR protein (antibody-IR700 adduct-antigen protein complex) (peak 690 nm near-infrared laser irradiation measured at 16 J/cm<sup>2</sup>). SDS-PAGE showed that not only cetuximab but also EGFR protein, the antigen, was aggregated (Figure 3b). In addition, the loss of NIR fluorescence of IR700 was observed. From the above, it was confirmed that the complex of antibody IR700 adduct and antigen was aggregated by light irradiation. This indicates that the photochemical aggregation of IR700 causes the antibody to aggregate and also causes the antigen protein, the binding

partner of the antibody, to aggregate.

Although we were able to prove the photochemical reaction by mass spectrometry above, it is necessary to capture the change visually. This mechanism can be clarified by imaging the changes in a single antibody molecule and quantitatively measuring the photochemical reaction. Therefore, we attempted to perform nano-imaging in liquid by atomic force microscopy (AFM). Panitumumab-IR700 (pan-IR700) was adhered to a mica substrate and the top of the substrate was filled with 10 mM phosphate buffer containing 50 mM MgCl<sub>2</sub> for nano-imaging in liquid. Macroscopic views of 1.0 μm and microscopic images of

**Observation and Quantification of Changes in panitumumab-IR700 (pan-IR700) under Near-Infrared Light Irradiation by Nano-Imaging in Liquid Using FM-AFM**


Figure 4 Observation and Quantification of Changes in panitumumab-IR700 (pan-IR700) under Near-Infrared Light Irradiation by Nano-Imaging in Liquid Using FM-AFM The changes of pan-IR700 before and after near-infrared light irradiation were investigated by nano-imaging in liquid using FM-AFM.  
 a. The pan-IR700 molecules depicted on the mica substrate aggregated and became larger after light irradiation, and adhered to the mica substrate.  
 b. pan-IR700 molecules depicted on the mica substrate have a Y-shape, but they aggregate and completely change their structure after light irradiation.  
 c. Results of change quantification by liquid nano-imaging using FM-AFM; the height and area significantly increased only in the group where Pan-IR700 was irradiated with NIR-light.

Relationship of IR700 fluorescence and in vitro cell death, in vivo anti-tumor effects, respectively


Figure 5 Relationship of IR700 fluorescence and in vitro cell death, in vivo anti-tumor effects, respectively  
 a. IR700-fluorescence in vitro was detected before and after NIR-PIT with flow cytometry.  
 b. Cytotoxicity in vitro was detected before and after NIR-PIT with luciferase activities.  
 c. Positive correlations between decrease of IR700 fluorescence and cytotoxicity is detected ( $r = 0.9779$ ).  
 d. in vivo NIR-PIT regimen is showed.  
 e. Mouse tumor inoculation outline and light dose.  
 f. Representative image of treated mice.  
 g. Positive correlations between decrease of IR700 fluorescence and anti-tumor effects is detected ( $r = 0.6316$ ).

41.7 nm were taken. In the macroscopic image, a single molecule of pan-IR700 attached to the mica substrate was agglomerated and expanded as an oligomer by near-infrared light irradiation (measured value: 16 J/cm<sup>2</sup>). In the micro-image, a single molecule of pan-IR700 in the shape of a Y-shape underwent structural changes upon light irradiation, and expanded, deformed, and aggregated as an oligomer (Figure 4a, b). As a control, pan-IR700 alone (no light irradiation) and antibody (pan) plus light irradiation (16 J/cm<sup>2</sup>) were tested in the same way for quantification. As a result of quantification of the height and Area values over time, only the combination of pan-IR700 and light irradiation (16 J/cm<sup>2</sup>) was found to be dominant, with the height and Area values expanding over time (Figure 4c). Therefore, nano-imaging confirmed that the antibody-IR700 agglomerates upon near-infrared light irradiation.

In order to analyze the relationship between photochemical reactions and cell death in vitro, we analyzed the relationship between the loss of IR700 fluorescence on the cell membrane before and after light irradiation and cell death. It is known that near-infrared fluorescence is lost in a dose-dependent manner, and measuring the near-infrared fluorescence of IR700 in IR700 antibody adducts was thought to reflect the degree of protein aggregation. Changes in near-infrared fluorescence were quantified by flow cytometry (Figure 5a), and cell death was quantified by luciferase activity (Figure 5b), and the respective changes were analyzed as a scatter plot. The results showed a positive correlation between the decrease in IR700 fluorescence and the increase in cell death ( $r = 0.9779$ ) (Figure 5c). This result suggests that cell death can

be measured from the amount of decrease in near-infrared fluorescence in target cells bound with IR700 antibody adducts, proving that IR700 photochemistry is strongly associated with cell death in NIR-PIT.

To analyze the relationship between the photochemical reaction and the in vivo anti-tumor effect, we analyzed the relationship between the loss of IR700 fluorescence in tumors before and after light irradiation and the anti-tumor effect (Figure 5d). As shown in Figure 5e, three tumors were created in the same individual and irradiated with 0, 50, and 100 J/cm<sup>2</sup> of near-infrared light, respectively, and the amount of change was measured. The change in IR700 fluorescence was quantified by the 700 nm channel of the PEARL imager, and the anti-tumor effect was quantified by the luciferase activity, which was analyzed in relation to the scatter plot (Figure 5f). The results showed a positive correlation between the decrease in IR700 fluorescence and the increase in antitumor effect ( $r = 0.6316$ ) (Figure 5g). This result suggests that the anti-tumor effect can be measured from the amount of decrease in IR700 fluorescence (700 nm) across the target tumor bound with IR700 antibody adducts, proving that the photochemical reaction of IR700 is strongly related to the anti-tumor effect of NIR-PIT. Since this decrease in IR700 fluorescence may predict therapeutic efficacy, optimization of real-time imaging of IR700 fluorescence is underway.

Conclusion

The above results suggest that the initiation point of cell death in NIR-PIT is the hydrophobization of IR700 by

## Mechanism of cell destruction by NIR-PIT


Figure 6 Mechanism of cell destruction by NIR-PIT

The antibody-IR700 complex binds to cancer-specific antigens on the cell surface. When the cells are irradiated with near-infrared light, the silanol group is released from IR700 and rapid hydrophobization occurs, resulting in the formation of aggregates including antibody and antigen proteins. This causes physical damage to the cell membrane. The osmotic disparity causes water to flow into the cell from outside the cell, causing the cell to rupture. Thus, the cell contents are disseminated to the surrounding area as necrosis to produce active cell destruction.

near-infrared light irradiation and the subsequent aggregation of antibody and antigen proteins<sup>[7]</sup>. In this study, the aggregation of antibody-IR700 adducts by near-infrared light irradiation was revealed by nano-imaging, and the results were convincing enough to propose a new concept of cell death in NIR-PIT. Integrating these analytical chemical data, cell culture, and animal studies, the authors propose the following mechanism of action of NIR-PIT (Figure 6). The antibody-IR700 complex, which binds to cancer-specific antigens on cell membranes, suddenly aggregates upon NIR irradiation and acts as a physical cell membrane barrier by denaturing the cell membrane antigen to which it binds. This mechanism was considered to be a novel photocell death that is completely different from previous concepts. The clarification of this mechanism proved the superiority of NIR-PIT as a new anti-cancer modality, and led to its limited approval in Japan as a “fifth cancer treatment” ahead of other countries. In addition, as a further application, it is thought that it is possible to predict the effect of treatment in real time by measuring and evaluating near-infrared fluorescence before and after treatment as described in this paper, and devices are being developed that apply intraoperative fluorescence imaging devices. We look forward to the development of future applications.

\* Editorial note: This content is based on HORIBA’s investigation at the year of issue unless otherwise stated

## References

- [ 1 ] Mitsunaga M.; Ogawa M.; Kosaka N.; Rosenblum L. T.; Choyke P. L.; Kobayashi H. Cancer Cell-Selective in Vivo near Infrared Photoimmunotherapy Targeting Specific Membrane Molecules. *Nat. Med.* 2011, 17 (12), 1685–1691.
- [ 2 ] NCT03769506, <https://clinicaltrials.gov/ct2/show/NCT03769506>
- [ 3 ] <https://rakuten-med.com/jp/news/press-releases/2020/09/25/7637/>
- [ 4 ] Buss NA, Henderson SJ, McFarlane M, Shenton JM, de Haan L : Monoclonal antibody therapeutics : history and future. *Curr Opin Pharmacol* 12 : 615-622, 2012
- [ 5 ] Agostinis P, Berg K, Cengel KA, Foster TH, Girotti AW, Gollnick SO, Hahn SM, Hamblin MR, Juzeniene A, Kessel D, Korbelik M, Moan J, Mroz P, Nowis D, Piette J, Wilson BC, Golab J. Photodynamic therapy of cancer: an update. *CA Cancer J Clin.* 2011 Jul-Aug;61(4):250-81.
- [ 6 ] Sato K, Nakajima T, Choyke PL, Kobayashi H. Selective cell elimination in vitro and in vivo from tissues and tumors using antibodies conjugated with a near infrared phthalocyanine RSC *Adv.* 2015 Mar 3;5(32):25105-25114.
- [ 7 ] Sato K, Ando K, Okuyama S, Moriguchi S, Ogura T, Totoki S, Hanaoka H, Nagaya T, Kokawa R, Takakura H, Nishimura M, Hasegawa Y, Choyke PL, Ogawa M, Kobayashi H. Photoinduced Ligand Release from a Silicon Phthalocyanine Dye Conjugated with Monoclonal Antibodies: A Mechanism of Cancer Cell Cytotoxicity after Near-Infrared Photoimmunotherapy. *ACS Cent Sci.* 2018 Nov 28;4(11):1559-1569.

SATO Kazuhide

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Designated Assistant Professor  
Institute for Advanced Research, Respiratory Medicine / S-YLC,  
Nagoya University.

## Real-Time Monitoring and Control of Pharmaceutical Production Processes Using Spectroscopic Data

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Dr. Sanghong Kim

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To realize the continuous production in the pharmaceutical industry, it is essential to monitor the production process in real time, but it is often difficult. In addition, conventional methods for predicting pharmaceutical quality from near-infrared spectra have drawbacks, such as a decrease in prediction accuracy over time. In this study, a stable and accurate method for predicting the quality of pharmaceuticals from near-infrared spectra was developed using data science techniques. As a result, a new analysis method was established that can contribute to realize real-time quality management and control, and to improve the efficiency of the pharmaceutical production processes.

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### Introduction

In the pharmaceutical production process, strict quality control based on Good Manufacturing Practice (GMP) is required. In recent years, there has been a shift from batch processes to continuous processes and the introduction of Quality by Design (QbD) and Process Analytical Technology (PAT)<sup>[1]-[4]</sup> to improve the efficiency of pharmaceutical production processes. In order to achieve the above goals, it is necessary to monitor and control the quality of pharmaceuticals in real time. However, it is often difficult to measure the quality of a drug directly in real time. In order to solve this problem, this research developed and put into practical use a method for constructing a statistical model to estimate drug quality from variables that can be measured relatively easily in real time, such as near-infrared spectra.

### Main

#### Process Nonlinearity, Model Maintenance, Input Variable Selection

There are many studies on methods for predicting drug

quality from near-infrared spectra, but most of them use the classical method of partial least squares (PLS)<sup>[5]-[8]</sup>. However, PLS is known to have some problems such as inability to cope with nonlinearity and high maintenance load because the estimation accuracy decreases with time. In this research, a new method called locally-weighted PLS to improve PLS was developed. Locally-weighted PLS can cope with nonlinearity by selectively using historical data that are similar to those around the conditions for which quality prediction is required. In addition, since a new model is automatically constructed for each quality prediction, the estimation accuracy can be maintained for a long time. **Figure 1** shows the results of PLS and locally-weighted PLS prediction of the concentration of the active ingredients in the powder after mixing for a mixing process. Locally-weighted PLS reduced the root mean square error of prediction (RMSEP) from 1.84 to 1.13. This allows for more strict quality control and more efficient drug production processes.

In addition to the development of locally-weighted PLS, a method to select the variables that are useful for estimation was developed. It is said that various physicochemical


Figure 1 Prediction results by PLS (left) and locally-weighted PLS (right)

factors coexist in near-infrared spectra, and it is not always possible to make good predictions by selecting functional groups and their corresponding wavelength regions, or by selecting wavelength regions by referring to the spectra of active ingredients and other ingredients as pure substances. The wavelength regions selected by the developed method are shown in Figure 2 as the blue region. The target granule is composed of API (Active Pharmaceutical

Ingredient) and five other components. The peaks of API and other components do not exist in the selected region, which may seem useless for estimation. However, the RMSEP was reduced to 1.13 from 1.84 by using absorbance at those wavelengths compared to use the absorbance at wavelengths selected by conventional method: variable importance on the projection (VIP) [12].


Figure 2 Wavelength selection result of the proposed method and preprocessed absorbance of each component of a drug

The detailed information can be found in the paper [9].

### Input Variable Scaling

When using statistical methods to build a quality prediction model, it is necessary to pre-process the input variables. Although input variable scaling, a data preprocessing method in which the values of each input variable are multiplied by the scaling factor of the input variable, can have significant effect on the estimation performance of soft sensors, research on input variable scaling has not been actively conducted. Hence, this research focuses on input variable scaling. In past research, autoscaling was commonly used [11-13]. In addition, Pareto scaling, level scaling, Poisson scaling, range scaling, and VAST scaling [14] have been considered.

In this research, a new method of preprocessing near-infrared spectra was developed, and it was found that it can improve the prediction accuracy of the impurity concentration inside the device. In the case of spectral data, it is difficult to evaluate the importance of each variable because hundreds to thousands of variables need to be

handled simultaneously.

An example of the application of the proposed method to the prediction of magnesium stearate concentration is shown in Figure 3. Figure 3 shows the near-infrared spectra of magnesium stearate and methanol solutions of magnesium stearate with different concentrations. Unlike Section 2-1, it was known in advance that the blue region in the figure is important for concentration prediction because it is a solution system. The evaluation results are shown in Figure 4. From Figure 4, it can be seen that the proposed method is able to correctly determine the important regions for concentration prediction. Accordingly, the prediction error was reduced by 45%.

The detailed information can be found in the paper [15].

### Model and Parameter Selection

Research on statistical modeling methods has been actively conducted, however, the usefulness of the methods has usually been evaluated by using a single dataset in most of the research, and the robustness of the methods have not


Figure 3 Preprocessed absorbance of magnesium stearate and its methanol solution


Figure 4 Importance evaluation result of the proposed method

been evaluated. To propose a highly reliable method, it is necessary to evaluate whether the developed soft-sensor can stably estimate the difficult-to-measure variables for various processes and a wide range of operating conditions. In this research, partial least squares (PLS), locally weighted PLS, support vector regression (SVR), and artificial neural network (ANN), which are widely used for soft-sensor design, are compared using twelve real-world datasets. The result of comprehensive comparative study in 6 datasets is shown in Figure 5. In Figure 5,  $p$  denotes the ratio of the number of samples used to model construction to the number of all the samples in each dataset. The result showed that locally weighted PLS outperformed the other methods. At the same time, the drawback of locally weighted PLS has also been made clear. locally weighted PLS is a modified version of PLS which is a linear regression method, and locally weighted PLS can deal with nonlinearity. However, while the conventional parameter tuning method can adapt to the nonlinearity of the target process, it has also a risk of excessively adapting to variations in measured values due to the influence of noise. Thus, a new parameter tuning method is proposed to improve the accuracy while guaranteeing the minimum prediction accuracy, and showed that the prediction accuracy of active ingredient and impurity concentrations can be improved.

The detailed information can be found in the paper <sup>[16]</sup>.

## Conclusion

In this research, real-time monitoring methods of continuous pharmaceutical processes were investigated. Some of the developed technologies have already been put to practical use and have achieved social contribution. In addition to the content of the former section, the research

on measurement and automation using data science and process control technology for many processes in the pharmaceutical production process, including the granulation process, drying process, and tableting process have been conducted. This research is expected to make a significant contribution to society, as it is planned to be put to practical use by our joint research partners.


Figure 5 Result of comprehensive comparative study

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## References

- [ 1 ] FDA 2004. Pharmaceutical cGMPs for the 21st century—a risk-based approach final report.
- [ 2 ] ICH 2005. ICH harmonised tripartite guideline—pharmaceutical development Q8 (R2).
- [ 3 ] ICH 2005. ICH harmonised tripartite guideline—quality risk management Q9.
- [ 4 ] ICH 2008. ICH harmonised tripartite guideline—pharmaceutical quality system Q10.
- [ 5 ] Moes, J.J., Ruijken, M.M., Gout, E., Frijlink, H.W., Ugwoke, M.I., 2008. Application of process analytical technology in tablet process development using NIR spectroscopy: blend uniformity, content uniformity and coating thickness measurements. *Int. J. Pharm.* 357, 108–118.
- [ 6 ] Berthiaux, H., Mosorov, V., Tomczak, L., Gatamel, C., Demeyre, J., 2006. Principal component analysis for characterising homogeneity in powder mixing using image processing techniques. *Chem. Eng. Process.* 45, 397–403.
- [ 7 ] Virtanen, S., Antikainen, O., Yliruusi, J., 2007. Uniformity of poorly miscible powders determined by near infrared spectroscopy. *Int. J. Pharm.* 345, 108–115.
- [ 8 ] Wu, H., Tawakkul, M., White, M., Khan, M., 2009. Quality-by-design (QbD): an integrated multivariate approach for the component quantification in powder blends. *Int. J. Pharm.* 372, 39–48.
- [ 9 ] Eriksson, L., Johansson, E., Kettaneh-Wold, N., Wold, S., 2001. *Multi- and Megavariate Data Analysis. Principles and Applications.* Umetrics Academy.
- [10] Kim, S., Kano, M., Nakagawa, H., Hasebe S., 2011. Estimation of active pharmaceutical ingredients content using locally weighted partial least squares and statistical wavelength selection, *Int. J. Pharm.*, 421, 269–274
- [11] Engel J, Gerretzen J, Szymańska E, Jansen JJ, Downey G, Blanchet L, et al. 2013. Breaking with trends in pre-processing? *Trends Anal. Chem.* 50, 96–106.
- [12] van den Berg RA, Hoefsloot HCJ, Westerhuis JA, Smilde AK, van der Werf MJ., 2006. Centering, scaling, and transformations improving the biological information content of metabolomics data. *BMC Genomics*, 7, 142.
- [13] Todeschini R, Consonni V, Maiocchi A., 1999/ The k correlation index theory development and its application in chemometrics. *Chemom Intell Lab Syst* 46, 13–29.
- [14] Keun HC, Ebbels TMD, Antti H, Bollard ME, Beckonert O, Holmes E, et al., 2003. Improved analysis of multivariate data by variable stability scaling application to NMR-based metabolic profiling. *Anal Chim Acta*, 490 265–76.
- [15] Kim, S., Kano, M., Nakagawa, H., Hasebe S., 2015. Input variable scaling for statistical modeling, *Comput. & Chem. Eng.*, 74, 4, 59–65
- [16] Matsuyama, Y., Kim, S., Hasebe S., 2021 Robust parameter tuning method of LW-PLS and verification of its effectiveness by twelve industrial processes", *Comput. & Chem. Eng.* 146 107224

**Dr. Sanghong Kim**

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Associate Professor  
Department of Applied Physics and Chemical Engineering,  
Tokyo University of Agriculture and Technology.

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### The Role of Particle Design Studies in Developing Pharmaceutical Dosage Forms

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Emeritus Professor

Laboratory of Advanced Pharmaceutical Process Engineering  
Gifu Pharmaceutical University

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In order to realize patient centric therapy, studies on drug administration is progressed focusing on its route and type of dosage forms. In addition to guaranteeing efficacy and safety in drug administration, there are development of new dosage forms that take into consideration ease of administration, and development of DDS for the purpose of more effective drug delivery. In either study, particle design research in developing the dosage forms is important. Recent trends in dosage form design development and the particle design researchs to support them will be introduced.

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
Figure 1 Disintegration time of erythritol tablet containing composite particle prepared by spray-drying with silica and disintegrant  
Formulation of tablets: Erythritol + 10% composite particle + 2-10% disintegrant + 0.2% lubricant  
Disintegrant: C-PVP: cross povidone, CCS: croscarmellose sodium, L-HPC: low substituted hydroxy propyl cellulose

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Table 1 Premixed or co-processed excipients with mannitol for direct tableting of ODT

Trade name	Formulation
F-MELT	Mannitol, xylitol, microcrystalline cellulose, cross povidone, Metasilicate aluminic acid magnesium or Anhydrous calcium hydrogen phosphate
SmartEx	Mannitol, L-HPC, PVA
Ludiflash	Mannitol, Collidone CL-SF, Collicoat SR30D
Pearitol FLASH	Mannitol, corn starch
Granfiller D	Mannitol, carmelose, cross-povidone

L-HPC: Low-substituted hydroxy propyl cellulose  
PVA: poly vinyl alcohol

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
Figure 3 Confocal laser scanning micrographs of slice samples of upper part of the intestinal tract at 2hr after intragastrical administration of (A) Lip (MLVs), (B) CS-Lip, (C) ssCS-Lip in rats  
Lip: liposome, CS-Lip: chitosan coated liposome, ssCS-Lip: submicron sized CS-Lip  
The mean particle size of liposomes: (A) 7.56  $\mu\text{m}$  (B) 3.58  $\mu\text{m}$  (C) 0.28  $\mu\text{m}$ . Lipid formulation of liposomes: DSPC: DCP: Chol.=8: 2: 1. Red parts means presence of corresponding liposomes.

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
Figure 4 Delivery of coumarin-6 (fluorescence marker for liposomes) with liposome to posterior part of eye observed with fluorescence microscopy 30 min after instillation to the eye in mice  
 Lip: liposome, ssLip: submicron sized Lip, EPC: egg phosphatidyl choline, DSPC: di-stearoyl phosphatidyl choline  
 EPC liposomes are relatively softer compared with DSPC liposomes.

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- [1] 竹内洋文, ヒトに優しい固形製剤設計「医薬品製剤化方略と新技術」竹内洋文編, シーエムシー出版(2007)3-7.
- [2] Seager H., Drug-delivery products and the Zydis fast-dissolving dosage form., *J. Pharm. Pharmacol.*, 50(1998)375-382.
- [3] Takeuchi H., Fuji Sylsya Chemical Ltd., Jpn. Patent No. 5026675(2012)
- [4] Tanimura S, Tahara K, Takeuchi H. Spray-dried composite particles of erythritol and porous silica for orally disintegrating tablets prepared by direct tableting. *Powder Technol.* 286(2015)444-450.
- [5] 谷村信爾, 川添博義, 山本浩充, 竹内洋文. 噴霧乾燥HPC-Lを結合剤として用いた低成形性医薬品粉末の直接打錠. *粉体工学会誌.* 43(2006)648-652.
- [6] 林洗仁, 竹内淑子, 田原耕平, 内田実, 竹内洋文. 滑沢剤ショ糖脂肪酸エステルの乳糖錠剤における崩壊 挙動抑制とそのメカニズムに関する検討. *粉体工学会誌.* 49(2012)750-757.
- [7] Wening, K., J. Breitreutz J., Oral drug delivery in personalized medicine: Unmet needs and novel approaches. *Int. J. Pharm.* 404(2011)1-9.
- [8] Niese, S. Quodbach J., Development of a dosing device for individualized dosing of orodispersible warfarin films. *Int. J. Pharm.*, 561(2019)314-323.
- [9] Visser J. C., Woerdenbag H. J. Crediet S., Gerrits E., Lesschen M. A., Hinrichs W. L.J., Breitreutz J., Frijlink H. W., Orodispersible films inn individualized pharmacotherapy: The development of a formulation for pharmacy preparations. *Int. J. Pharm.*, 478(2015)155-163.
- [10] Matsumoto M., Arai K., Fujita A. Tsuji M., Hayakawa S., Sawakami K., Sawakami K, Shida S., Nohara Y., Sumiyaa K., Murata K., Takeuchi Y., Takeuchi H., Murata R., Characterization of Orally Disintegrating Films Containing Fluconazole, *Journal of Drug Interaction Research*, 43(2019) 55-59.
- [11] Takeuchi H., Yamakawa R., Nishimatsu T., Takeuchi Y., Hayakawa K., Maruyama N., *J. Drug Deliv. Sci. Technol.*, 23 (2013)471-475.
- [12] Takeuchi Y., Umemura K., Tahara K., Takeuchi H., Formulation design of hydroxypropyl cellulose films for use as orally disintegrating dosage forms. *J. Drug Deliv. Sci. Technol.*, 46(2018)93-100.
- [13] Takeuchi Y., Nishimatsu T., Tahara K., Takeuchi H., Novel use of insoluble particles as disintegration enhancers for orally disintegrating films. *J. Drug. Delivery Sci. Technol.*, 54,(2019)101310.
- [14] Takeuchi Y, Usui R, Ikezaki H, Tahara K, Takeuchi H. An advanced technique using an electronic taste-sensing system to evaluate the bitterness of orally disintegrating films and the evaluation of model films. *Int J Pharm.*, 531(2017)179-190.
- [15] Takeuchi Y., Kawamoto M, Tahara K, Takeuchi H. Design of a new disintegration test system for the evaluation of orally disintegrating films. *Int J Pharm.*, 553(2018)281-289.
- [16] Takeuchi H, Yamamoto H, Niwa T, Hino T, Kawashima Y. Mucoadhesion of Polymer-Coated Liposomes to Rat Intestine in Vitro. *Chem. Pharm. Bull.*, 42(1994)1954-1956.
- [17] Takeuchi H., Yamamoto H., Niwa T., Hino T., Kawashima Y., Enteral Absorption of Insulin in Rats from Mucoadhesive Chitosan-Coated Liposomes. *Pharm. Res.*, 13,(1996)896-901.
- [18] Takeuchi H, Yamamoto H, Kawashima Y. Mucoadhesive nanoparticulate systems for peptide drug delivery. *Adv. Drug Del. Rev.*, 47(2001)39-54.
- [19] Takeuchi H, Matsui Y, Sugihara H, Yamamoto H, Kawashima Y. Effectiveness of submicron-sized, chitosan-coated liposomes in oral administration of peptide drugs. *Int. J. Pharm.*, 303(2005)160-170.
- [20] Makhlof A, Werle M, Tozuka Y, Takeuchi H. A mucoadhesive nanoparticulate system for the simultaneous delivery of macromolecules and permeation enhancers to the intestinal mucosa. *J. Cont. Rel.*, 149(2011)81-88.
- [21] Murata M, Nakano K, Tahara K, Tozuka Y, Takeuchi H. Pulmonary delivery of elcatonin using surface-modified liposomes to improve systemic absorption: Polyvinyl alcohol with a hydrophobic anchor and chitosan oligosaccharide as effective surface modifiers. *Eur. J. Pharm. Biopharm.*, 80(2012) 340-346.
- [22] Murata M, Tahara K, Takeuchi H. Real-time in vivo imaging of surface-modified liposomes to evaluate their behavior after pulmonary administration. *Eur J Pharm Biopharm.* 86 (2014)115-119.
- [23] Hironaka K, Inokuchi Y, Tozuka Y, Shimazawa M, Hara H, Takeuchi H. Design and evaluation of a liposomal delivery system targeting the posterior segment of the eye. *J Contr. Rel.*, 136(2009)247-253.
- [24] Nakano K, Tozuka Y, Yamamoto H, Kawashima Y, Takeuchi H. A novel method for measuring rigidity of submicronsize liposomes with atomic force microscopy. *Int. J. Pharm.*, 355 (2008)203-209.
- [25] Hironaka K, Inokuchi Y, Fujisawa T, Shimazaki H, Akane M, Tozuka Y, Tsuruma K, Shimazawa M, Hara H, Takeuchi H. Edaravone-loaded liposomes for retinal protection against oxidative stress-induced retinal damage. *Eur. J. Pharm. Biopharm.*, 79(2011)119-125.
- [26] Shimazawa M, Inoue Y, Masuda T, Onodera R, Tahara K, Shimizu Y, Mibe Y, Tsuruma K, Takeuchi H, Hara H. Topical diclofenac-loaded liposomes ameliorate laser-induced choroidal neovascularization in mice and non-human primates. *Curr. Neurovasc Res.*, 14(2017)46-52.

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Current Status and Objectives on Technologies in Pharmaceutical Industry

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Professor

Research Organization of Science and Technology, Ritsumeikan University

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The Monozukuri, a manufacturing technology, that contributed to Japanese progress in the 20th century will also support Japanese prosperity in the 21st century. One of such Monozukuri is definitely a pharmaceutical technology that can offer health to the people all over the world. The technology of pharmaceutical Monozukuri that delivers valuable medicines to the world can be created with fusion among the technologies of drug discovery, drug delivery system and pharmaceutical manufacturing. Japan showed some contributions to development of those technologies until now. However, in each technology there are technical barriers to be overcome for the further progress and Japanese pharmaceutical industry can survive by providing some feasible solutions to those barriers. This paper describes the current status and objectives on the technologies of drug discovery, drug delivery system and pharmaceutical manufacturing

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
Figure 2 Pegylated liposomal doxorubicin. Left panel: schematic drawing of liposome; right panel: Cryo-TEM picture liposomes with characteristic rod precipitate (doxorubicin sulfate) and ovoid shape. (Quoted from Gabison, Patil and La-Beck)<sup>[1]</sup>

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
Figure 3 Scanning electron micrograph of LPGA microspheres containing Leuporelin acetate. Gifted by Dr. Hiroaki Okada.

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
Figure 4 Chemical structure of Tacrolimus. Mw:822.05


Figure 5 Time course of blood level of Tacrolimus after oral dose of crystalline powder or solid dispersion powder in dogs at a dose of 3.2 mg/kg. Gifted by Dr. Kazunari Yamashita.

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Table 1 Dosage forms developed for Tacrolimus and its first year approved all over the world

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- [1] Alberto Gabizon, Raphael Catane, Beatrice Uziely, Bela Kaufman, Tamar Safra, Rivka Cohen, Francis Martin, Anthony Huang, and Yechezkel Barenholz, "Prolonged Circulation Time and Enhanced Accumulation in Malignant Exudates of Doxorubicin Encapsulated in Polyethylene-glycol Coated Liposomes", *Cancer Research*, vol.154, pp.987-992, 1994.
- [2] H.Okada, "One- and three-month release injectable microspheres of the LH-RH superagonist leuporelin acetate", *Adv. Drug. Del. Rev.*, vol.28, pp.43-70, 1997
- [3] Kazunari Yamashita , Toshiomi Nakate, Kazuto Okimoto, Atsuo Ohike, Yuji Tokunaga, Rinta Ibuki, Kazutaka Higaki and Toshikiro Kimura, "Establishment of new preparation method for solid dispersion formulation of tacrolimus", *Int J Pharm.*, vol.28, pp.79-91, 2003

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### Analytical Chemistry for Biopharmaceutics

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Professor

Department of Bioengineering, School of Engineering,  
The University of Tokyo

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In the research and development of biopharmaceutics, including antibodies, importance of analytical chemistry has become more and more recognized. In particular, the proteins in biopharmaceutics exist as molecular colloid and have unique conformations regulated by their sequences, which means that precise analyses of the drug candidates from the viewpoints of both colloidal and conformational stabilities are required for the process. This article reviews recent progresses on the methods for evaluation of colloidal stability, and of conformational stability under higher concentration of the molecules.

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	Company A	News A	Center A	News B	Conference A	Paper A	Paper B
SEC	✓	✓	✓		✓	✓	✓
DLS	✓	✓	✓	✓	✓	✓	✓
SLS	✓		✓			✓	✓
AUC	✓	✓	✓	✓	✓	✓	✓
Nano Tracking Analysis (NTA)	✓	✓	✓		✓		✓
Resonant Mass Measurement (RMM)	✓	✓	✓		✓		
Coulter Counter	✓		✓			✓	✓
Field-Flow Fractionation (FFF)	✓				✓	✓	✓
Micro-Flow Imaging (MFI)	✓		✓	✓	✓	✓	✓
Laser Diffraction (LD)		✓		✓			
Light Obscuration (LO)	✓	✓				✓	✓
Optical Microscopy	✓	✓	✓			✓	✓
Infrared Microscopy/FT-IR	✓		✓				✓
Visual Inspection	✓	✓	✓			✓	✓

[Company A] <http://www.coriolis-pharma.com/service-and-key-expertise/aggregate-and-particle-analysis/protein-aggregation-analysis/>

[Center A] <https://mvsc.ku.edu/content/particle-characterization-submicron-subvisible-and-visible-particles>

[News A] <http://www.biopharminternational.com/analyzing-protein-aggregation-biopharmaceuticals-0>

[News B] <http://www.chemengonline.com/particle-sizing-technology-selection/?printmode=1>

[Conference A] <http://www.americanlaboratory.com/913-Technical-Articles/158757-Protein-Aggregation-at-PepTalk-2014/>

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- [1] A.L. Fink, Protein aggregation: folding aggregates, inclusion bodies and amyloid, *Fold. Des.* 3, R9-23(1998).
- [2] T. Ito, K. Tsumoto. Effects of subclass change on the structural stability of chimeric, humanized, and human antibodies under thermal stress, *Protein Sci.* 22, 1542-1551 (2013).
- [3] W. Fraunhofer, G. Winter, The use of asymmetrical flow field-flow fractionation in pharmaceuticals and biopharmaceuticals, *Eur. J. Pharm. Biopharm.* 58, 369-383(2004).
- [4] D.K. Sharma, P. Oma, M.J. Pollo, and M. Sukumar, *J. Pharm. Sci.* 99, 2628-2642(2010).
- [5] E. Krayukhina, K. Tsumoto, S. Uchiyama, K. Fukui. Effects of syringe material and silicone oil lubrication on the stability of pharmaceutical proteins, *J. Pharm. Sci.* 104, 527-535(2015).
- [6] S. Totoki, G. Yamamoto, K. Tsumoto, S. Uchiyama, K. Fukui, Quantitative laser diffraction method for the assessment of protein subvisible particles, *J. Pharm. Sci.* 104, 618-626(2015).
- [7] Kiyoshi M, Caaveiro JMM, Tada M, Tamura H, Tanaka T, Terao Y, Morante K, Harazono A, Hashii N, Shibata H, Kuroda D, Nagatoishi S, Oe S, Ide T, Tsumoto K, Ishii-Watabe A. *Sci. Rep.* 8, 3955(2018)
- [8] Kosuge H, Nagatoishi S, Kiyoshi M, Ishii-Watabe A, Tanaka T, Terao Y, Oe S, Ide T, Tsumoto K. *Biotechnol. Prog.* e3016 (2020).
- [9] Egashira Y, Nagatoishi S, Kiyoshi M, Ishii-Watabe A, Tsumoto K. *J. Biochem.* 163, 481-488(2018).
- [10] Z. Q. Wen, Raman Spectroscopy of Protein Pharmaceuticals, *J. Pharm. Sci.* 96, 2861-2878(2007).
- [11] A. Rygula, K. Majzner, K. M. Marzec, A. Kaczor, M. Pilarczyk, M. Baranska, Raman spectroscopy of proteins: a review, *J. Raman Spectrosc.* 44, 1061-1076(2013).
- [12] C. Ota, S. Noguchi, S. Nagatoishi, K. Tsumoto, Assessment of the Protein-Protein Interactions in a Highly Concentrated Antibody Solution by Using Raman Spectroscopy, *Pharm. Res.* 33, 956-969(2016).
- [13] C. Ota, S. Noguchi, K. Tsumoto, The molecular interaction of a protein in highly concentrated solution investigated by Raman spectroscopy, *Biopolymers* 103, 237-246(2015).
- [14] Rolinger L, Rüdert M, Hubbuch J, A critical review of recent trends, and a future perspective of optical spectroscopy as PAT in biopharmaceutical downstream processing, *Anal. Bioanal. Chem.* 412, 2047-2064(2020).
- [15] Murali K. Maruthamuthu, Scott R. Rudge, Arezoo M. Ardekani, Michael R. Ladisch, Mohit S. Verma, Process analytical technologies and data analytics for the manufacture of monoclonal antibodies, *Trends Biotechnol.* 38, 1169-1186(2020)
- [16] Yoann Le Basle, Philip Chennell, Nicolas Tokhadze, Alain Astier, Valérie Sautou, Physicochemical stability of monoclonal antibodies: A review, *J. Pharm. Sci.* 109, 169-190 (2020)
- [17] Anurag S Rathore, Saxena Nikita, Garima Thakur, Navnath Deore, Challenges in process control for continuous processing for production of monoclonal antibody products, *Curr. Opin. Chem. Engineer.* 31, 100671(2021).
- [18] Alaa A. Makki, Victor Massot, Hugh J. Byrne, Renaud Respaud, Dominique Bertrand, Elhadi Mohammed, Igor Chourpa, Franck Bonniera, Understanding the discrimination and quantification of monoclonal antibodies preparations using Raman spectroscopy, *J. Pharm. Biomed. Anal.* 194, 113734(2021)
- [19] Meng-Yao Li, Bruno Ebel, Cédric Paris, Fabien Chauchard, Emmanuel Guedon, Annie Marc, Real-time monitoring of antibody glycosylation site occupancy by in situ Raman spectroscopy during bioreactor CHO cell cultures, *Biotechnol. Prog.* 34, 486-493(2018).
- [20] Bethan S McAvan, Leo A Bowsher, Thomas Powell, John F O'Hara, Mariangela Spitali, Royston Goodacre, Andrew J Doig, Raman Spectroscopy to Monitor Post-Translational Modifications and Degradation in Monoclonal Antibody Therapeutics, *Anal. Chem.* 92, 10381-10389(2020).

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### Principles of Good Practice for Laboratory Medicine

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Professor

Department of Laboratory Medicine, Hamamatsu University School of Medicine.

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Laboratory medicine gives us indispensable for medicine and objective information, and the proportion of the influence of laboratory medicine upon clinical decisions is from 60 to 70%. Therefore, laboratory medicine should supply standardized results both at all times and all laboratories. For this purpose, total laboratory process consisted of pre-analytical, analytical and post-analytical process shall be precisely managed, and standardization or harmonization of laboratory medicine shall be promoted. Standardized laboratory medicine would be expected to lead to ultimate personalized medicine (precision medicine) by integration of artificial intelligence and big data, for example longitudinal personal laboratory results, genetic factors, lifestyle habit, biological information, imaging information, research outcomes including omics analysis and clinical findings.

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
Figure 1 Total Testing Process of Laboratory Diagnostics

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Table 1 Biological Variation Factors of Laboratory Data

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## Translational Raman Spectroscopic Approaches for Clinical Routine

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Prof. Dr. Jürgen POPP

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Leibniz Institute of Photonic Technology,  
Member of the Leibniz Research Alliance - Leibniz Health Technologies,  
Albert-Einstein-Str. 9, D-07745 Jena, Germany

Institute of Physical Chemistry & Abbe Center of Photonics  
Friedrich-Schiller University Jena,  
Helmholtzweg 4, D-07743 Jena, Germany

The article reviews our latest results on innovative technological and data analysis concepts for bringing Raman spectroscopy closer to clinical use in terms of diagnosing and targeted therapy of infectious diseases and cancer. The first part will report on clinical Raman solutions for rapid diagnosis of infectious diseases - in terms of rapidly identifying the infection causing pathogen its antibiotic resistance pattern and ideally also its immune response - being decisive parameters for a targeted antibiotic administration, which is crucial for the survival of patients. The second part will focus on nonlinear multi-modal Raman imaging for a fast and safe precise intraoperative tumor margin control, because reliable tumor margin recognition during an intervention is the key to effective tumor treatment.

#### Keywords

Raman spectroscopy, Artificial intelligence, Chemometrics, Microbial analysis, Infectious diseases, Cancer, Spectral histopathology

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### Introduction

Understanding the causes of diseases, recognizing them earlier and treating them more specifically - hopes that are associated with modern biomedicine - requires the determination of diagnostic, prognostic and predictive factors including their comprehensive evaluation in just a few steps or ideally in a single step. In this context, the sharp rise in cancer due to an ageing society and the rapid spread of life-threatening infectious diseases (due to in part unknown pathogens) and antibiotic-resistant germs, which is partly due to an increasing worldwide mobility

but also to the ill-considered administration of broad-spectrum antibiotics, should be mentioned in particular. An effective and early diagnosis and personalized therapy of cancer and infections requires new methods of differential diagnosis and represents an outstanding task of medicine. In principle, the following applies to all diseases: the earlier treatment begins, the better the chances of cure. There is therefore a great need for new diagnostic methods for a targeted early diagnosis of diseases in order to be able to use a targeted therapy as early as possible.

Raman spectroscopy plays a key role in the implementation

of these ambitious goals. The application of Raman methods to address biomedical research has grown rapidly over the past ten years and has advanced into a new era due to advances in instrumentation and most importantly due to an enhanced cross-disciplinary dialogue between spectroscopists and clinicians, which e.g. is fostered in Europe by the COST action Raman4clinics (<https://www.raman4clinics.eu/>).<sup>[1, 2]</sup>

This contribution will highlight our recent advances in translating Raman approaches towards routine clinical applications with focus on infectious diseases and cancer. Addressing this question requires new Raman instrumentation, which can be applied out of specialized labs in a clinical environment (e.g. operation theatre, bedside or in a doctor's practice).

### Point-of-care Raman spectroscopy for a rapid infectious disease diagnosis and treatment

Let's assume the following situation: a patient is at the doctor's office and after a short anamnesis the doctor decides that it is a bacterial infection and prescribes a certain antibiotic. Now one can ask the question, how does the doctor know that it is a bacterial infection, or in case is a bacterial infection, that the causing germs are not resistant to the prescribed antibiotic. While this approach was appropriate a few years ago, nowadays in the course of increasing antibiotic resistance, it would be very beneficial, that prior to prescribing an antibiotic or before taking an antibiotic, a precise diagnosis in terms of rapidly identifying the pathogen and characterizing its resistance profile takes place on which then a targeted therapy follows. In the following the great potential of optical technologies with focus on Raman spectroscopy to successfully address these unmet needs will be shown. At the same time, the hurdles that have to be overcome in order to move from an academic scientific process to a validated and verified product will be outlined.

First the question, what exactly is needed for an optimal and personalized treatment of infectious diseases will be addressed. Here, it is important to know whether the symptoms a patient - the so-called host response - is exhibiting are symptoms of an infectious disease or of another inflammatory disease. In case it is an infection, the next step is to find out what kind of pathogen is causing the infection, i.e. are the infection causing pathogens viruses, bacteria, or fungi? This knowledge is crucial for initiating an effective therapy. In the case of a bacterial infection, it is also necessary to know whether the bacterial pathogen is resistant to certain antibiotics or not.

### The host response:

In order to find out if the patient suffers from an infection or not it is important to consider the role of immune cells during an infection. When infectious agents enter the body, they interact with immune cells, i.e. among other with white blood cells including neutrophils, monocytes, eosinophils, T cells, B cells or the basophils. The interaction of pathogenic microorganisms, i.e. viruses, bacteria or fungi, with these immune cells leads to a certain cell response, causing molecular changes within the immune cells. Raman spectroscopy offers a very simple approach (requiring only a laser, a microscope, some filters and a spectrometer together with a sensitive camera) to monitor these molecular changes via the Raman spectral fingerprint. Raman spectroscopy requires very simple sample preparation steps. The isolation of white blood cells from whole blood can be achieved very easily via lysis of the red blood cells. After isolation, the white blood cells are placed under the microscope and irradiated with the Raman excitation laser. A notch filter separates the elastically scattered from the inelastically scattered Raman light, which is spectrally dispersed by a monochromator and detected by a camera (e.g. CCD camera). This simple analysis process, which has been established within the EU project Hemospec (<https://cordis.europa.eu/project/id/611682>), requires just 60 minutes from sampling to the final result. However, in order to achieve this white blood cell analysis two challenges had to be overcome: (1) since the Raman spectra of the different cell types or the activation of the different cells differ only minimally, a simple spectral analysis by eye is not possible and artificial intelligence methods are required to analyze the cellular Raman spectra (see below). (2) a large number of cells need to be examined and commercially available Raman setups are not able to measure several 1000 cells in a very short time. Therefore, a high-throughput Raman (HTR) setup has been developed within Hemospec and put to use in the clinic.<sup>[3, 4]</sup>

The investigation of the host response started with *in vitro* experiments, where isolated white blood cells here neutrophils from whole blood were spiked with various infectious agents such as bacteria and fungi. The individual cellular Raman signatures were analyzed by a PC-LDA analysis subsequent to a sophisticated data pre-treatment.<sup>[5]</sup> In detail, the Raman spectral data of infected neutrophils have been compared with those of uninfected neutrophils and with an accuracy of 90% the infected could be distinguished from the non-infected neutrophils.<sup>[5]</sup> The examination of the infected neutrophils revealed that it is possible to distinguish a fungal infection from a bacterial infection with an accuracy of 92%. Finally, neutrophils associated with Gram-positive bacteria such as *Staphylococcus aureus* could be separated from

Gram-negative *E. coli* with an accuracy of 84%, i.e. Raman spectroscopy also allows to distinguish Gram-positive bacteria from Gram-negative bacteria. This *in-vitro* microbial infection diagnostics study<sup>[5]</sup> (see Figure 1) demonstrates that pathogen-specific activation can be detected with Raman data and opens up new possibilities for clinical diagnostic applications in terms of timely personalized therapy. For this, however, it is necessary to transfer the concept of the *in-vitro* stimulated cells shown in Reference [5] to patient samples.

In a recently published clinical study<sup>[4]</sup> we were able to show that it is possible to distinguish between inflammation, infection and sepsis by analyzing the Raman spectroscopic fingerprints of leukocytes from the blood of hospitalised patients from an emergency department. Using the aforementioned HTR approach<sup>[3]</sup>, we were able to clearly distinguish patients with viral infection from patients with bacterial infection. This clinical study paves the way for translating this Raman spectroscopic host response approach as a first steps towards a targeted therapy, as it allows to quickly choose the appropriate therapeutic approach based on the host response.<sup>[5]</sup>

**Pathogen diagnostics and antibiotic resistance:**

The currently approved microbiological methods to identify the infection-causing pathogens require the cultivation of bacteria and are therefore time-consuming, which means that a targeted antibiotic therapy is usually performed too late. Therefore, broad-spectrum antibiotics are often used empirically before the appropriate antibiotics matched to the infection-causing pathogen are found. However, this approach leads to more and more pathogens becoming resistant to antibiotics. To treat an infection

successfully, the physician needs the information about the type of the pathogen and its resistance to antibiotics - as quickly as possible. As mentioned above, standard microbiological methods usually require more than 24 hours due to cultivation, which is often too late.

To overcome this unmet medical need we have developed the so-called RAMANBIOASSAY™ - an approach that allows the identification of the pathogens as well as the determination of the minimum inhibitory concentration (MIC) in a very short time. The unique feature of this RAMANBIOASSAY™ approach is the dramatically reduced diagnostic time: in only 2 to 3.5 hours and with a very small number of microbial pathogens (a few hundred bacteria are sufficient), the pathogens can be clearly identified and their resistance properties can be determined without a prior cultivation step. Based on this information a physician can adapt the therapeutic antibiotic treatment specifically to the pathogen in question.

The basis of the RAMANBIOASSAY™ approach is the combination of chip-based Raman microspectroscopy with classical imaging. The smart use of laser light in a Raman setup in combination with a light microscope enables a label-free, non-destructive and culture-independent optical and spectroscopic characterization and identification of bacteria down to the level of single cells. Each molecule generates an individual signature in the Raman spectrum, creating a specific molecular fingerprint for each bacterium. Statistical evaluation algorithms enable automated classification and identification of bacteria and resistance to antibiotics. Furthermore, the additional integration of chip-based enrichment methods into miniaturized structures covers the entire process chain


Figure 1 Summary of *in-vitro* Raman study to unravel the host response.[5] For details see text.


Figure 2 Schematic of the individual steps of the RAMANBIOASSAY™: The core is a chip onto which the bacteria from a patient sample are applied. Up to four different antibiotics are already present in four concentrations each. This means that the chip is ready for immediate use.

from sampling to final results (see Figure 2).

The Raman spectroscopic identification of the pathogen can be completed after 35 minutes. At the same time, the morphological changes of the pathogens in a short-term culture (on the RAMANBIOASSAY™ chip) together with the Raman spectra are used to determine to which antibiotic agent resistance exists and what the MIC is. This information is extracted from the morphological images and the Raman spectra using computer-assisted statistical image data analysis and chemometrics. By doing so the complete resistogram is thus already available after 2 to 3.5 hours and the physician can react much earlier than currently and does not have to resort to the application of a broad-spectrum antibiotics.

The potential of the RAMANBIOASSAY™ approach was demonstrated using the example of urinary tract infections with an identification of the pathogens, their resistance to antibiotics and their MIC from patient material.<sup>[6-11]</sup>

The uniqueness of the RAMANBIOASSAY™ approach lies in its independence from culture and the scalability of the approach. The independence from culture allows the analysis time to be reduced to 2 to approximately 3.5 hours after sample collection. The scalability enables universal use on the one hand in clinical microbiology in combination with pipetting robot compatible Raman chip platforms, which allows a high degree of automation in the clinical microbiology workflow, and on the other hand in the form of a closed cartridge-based approach with a miniaturized Raman reader for direct use in a physician's office or in a hospital without its own microbiology unit. Figure 3 shows the current status using the RAMANBIOASSAY™ chip together with the BioParticle Explorer (a manual Raman spectroscopy system

together with a microscope for morphological analysis) as well as future developments of the RAMANBIOASSAY™ for clinical microbiology or use in a physician's office.

Equally important as the development of clinically usable Raman devices (see Figure 3) is the development of tailored Raman spectroscopic analysis routines. As mentioned above the analysis of bacterial Raman spectra or changes within induced by e.g. antibiotics treatment cannot be analyzed by naked eye and require sophisticated artificial intelligence based spectral analysis routines.

In general, the success of Raman spectroscopy for medical diagnosis and therapy (and also other applications like e.g. in life sciences, process analytics, pharmacy or environmental analysis) is inherently connected with the development of customized Raman data evaluation algorithms for translating Raman measurement data (spectral data sets, image data, etc.) into qualitatively and quantitatively usable information for end users. In this context we have developed a universally applicable Raman data analysis software called RAMANMETRIX


Figure 3 Photograph of RAMANBIOASSAY™ chip (middle) together with the BioParticle Explorer(left) and future developments towards a cartridge-based system with miniaturized reader for doctor's office and clinics without clinical microbiology (right).

(see: <https://docs.ramanmetrix.eu/>). This software allows for a one-click data analysis of Raman spectroscopic data in a robust and reliable way.

### Non-linear spectroscopic multimodal instrumentation for intraoperative tumor identification

The diagnostic gold standard in most surgeries is the extraction of biopsies and their histopathological examination to confirm the tumor and the tumor borders. The histopathological examination is carried out by means of rapid sections of non-contrasted tissue sections during the operation and definitively also on tissue sections of fixed material. Depending on the number of removed rapid sections, the rapid section diagnostic takes about 20-30 minutes. Then the surgeon is typically informed by telephone about the result and can decide on his/her further approach. If, for example, a sample in the border region of the tumor still shows tumor tissue, the surgeon will resection here and, if necessary, provide a new border section for rapid section diagnostics until the rapid sections indicate complete tumor resection.

This procedure is time-consuming, investigator-dependent, and dependent on the size, number and quality of the removed tissue samples, which are, of course, limited. Furthermore, since the quality of intraoperatively prepared frozen sections is not comparable to that of embedded tissue, the results of frozen section diagnostics often differs from those of an examination on embedded sections, which is why subsequent diagnostic confirmation using embedded sections is still necessary. For example, in head and neck surgery, depending on tumor size and location, the incidence of residual tumor (R1), i.e., subsequent finding of tumor cells in the incision margin, is approximately 7.5 - 10%. For such cases, the patient has to undergo a new operation and often also postoperative radiotherapy, which is an enormous burden for the patient. Therefore, even within the primary surgery, where possible, more tissue is removed than would actually be necessary to reach the tumor limit. However, the more tissue that is removed, the greater the impact on the patient's healing process. Thus, new methods and approaches are urgently needed for a fast and reliable intraoperative diagnosis.

Over the last few years, we have been investigating a multimodal nonlinear imaging approach that has the potential to reliably assess tissue and the success of surgery or endoscopy directly in the operating theatre or endoscopy room. The approach combines three different nonlinear imaging techniques namely two-photon excited autofluorescence (TPEF), second harmonic generation

(SHG) and coherent anti-Stokes Raman scattering (CARS) displaying vibrations in the CH-stretch wavenumber region. The combination of these three modalities allows to determine the morphological and chemical composition (morphochemistry) of unfixed tissue sections in a label-free manner.<sup>[12]</sup> In order to translate the morphochemical information encoded in the multimodal spectroscopic images into medical relevant information photonic data science, i.e. machine and deep learning approaches are necessary.<sup>[13-16]</sup>

In the following we will briefly highlight the potential of multimodal nonlinear imaging in combination with innovative image analysis routines as a powerful tool for computer based spectral histopathology allowing for an automatic prediction of tissue types / disease and thus offering great potential to fulfill the aforementioned unmet medical needs in terms of reducing the time in an operative theatre due to instant feedback and smaller workload due to automatization.

In a recent study we investigated head and neck squamous cell carcinomas using the multimodal combination of CARS, TPEF and SHG. The analysis of the images by a machine learning classification model features a 90% accuracy compared to gold standard diagnosis of a blinded pathologist.<sup>[14, 15, 17, 18]</sup>

Furthermore, we could show that the utilization of deep learning approaches also allows for a pseudo-staining of multimodal images. The deep learning generated pseudo H&E images nicely represent the real H&E images and show that the combination of the three label-free non-linear imaging modalities CARS, TPEF and SHG yields information that can be translated into computational pseudo hematoxylin and eosin (HE) images.<sup>[13,16]</sup> For a clinical application compact and easy to use devices are needed. Thus, we have transferred the presented CARS/SHG/TPEF approach into a compact microscope suitable for clinical use in terms of a rapid *ex-vivo* tissue analysis.<sup>[12]</sup>

In order to further extend the applicability of this multimodal microscopy approach for *in vivo* tissue screening, various endoscopic probe concepts were also realized.<sup>[19-23]</sup> The core of all these setups are robust and alignment-free fiber laser concepts. In the following two concepts of endoscopic probes for multimodal nonlinear imaging, which have been researched or are still currently under investigation will be briefly summarized. The first is a rigid needle endoscope for neurosurgical applications.<sup>[21]</sup> This probe provides very good imaging in a compact design, there are no moving or electric parts in the probe head, but it is not flexible (which is fine for neurosurgery and some other fields of application, e.g., head and neck).

For applications requiring flexible endoscopes, a unique endoscopy concept using a double-core double-clad fiber and focus-combining micro-optical concept allowing for a background free, low-loss, high peak power laser delivery, and an efficient signal collection in backward direction has been successfully realized.<sup>[23]</sup>

## Conclusion

In conclusion this short review highlights the great potential of Raman spectroscopy in combination with innovative photonic data science concepts for clinical diagnosis and therapy.

In the first part of this contribution, we report on the culture-free isolation and identification of pathogens, their host-response and their antibiotic resistance by using a combination of Raman spectroscopy, chip-based sampling strategies as well as chemometric spectroscopic data analysis methods. The main advantages of this rather simple approach compared to conventional microbiological analysis methods is that the analysis result is available in just a few hours. This represents a major step forward, as antibiotic therapy can be started promptly and specifically tailored to the pathogen. It will be shown how this approach was transferred into an automated clinically applicable system the RAMANBIOASSAY™

The second part of this presentation reports on multimodal non-linear imaging solutions in the field of pathological cancer diagnostics. Here, the combination of CARS with SHG and TPEF imaging in combination with adapted image analysis methods represents a powerful *ex-vivo* and *in-vivo* approach for a label-free clinical intraoperative tissue diagnostic for tumor margin detection in terms of computer based spectral histopathology.

## Acknowledgement

*Financial support of the EU, the “Thüringer Ministerium für Wirtschaft, Wissenschaft und Digitale Gesellschaft”, the “Thüringer Aufbaubank”, the Federal Ministry of Education and Research, Germany (BMBF), the German Science Foundation, the Fonds der Chemischen Industrie and the Carl-Zeiss Foundation are greatly acknowledged.*

\* Editorial note: This content is based on HORIBA's investigation at the year of issue unless otherwise stated.

## References

- [1] C. Krafft, I. W. Schie, T. Meyer, M. Schmitt, and J. Popp, "Developments in spontaneous and coherent Raman scattering microscopic imaging for biomedical applications", *Chem. Soc. Rev.* 2016, 45, 1819-1849 and references cite therein.
- [2] C. Krafft, M. Schmitt, I. Schie, D. Cialla-May, C. Matthäus, T. Bocklitz, and J. Popp, "Label-free molecular imaging of biological cells and tissues by linear and non-linear Raman spectroscopic approaches", *Angew. Chem. Int. Ed.* 2017, 56, 4392-4430 and references cite therein.
- [3] A. S. Mondol, N. Töpfer, J. Rüger, U. Neugebauer, J. Popp and I. W. Schie, "New perspectives for viability studies with high-content analysis Raman spectroscopy (HCA-RS)", *Sci. Rep.* 2019, 9, 12653.
- [4] A. Ramoji, D. Thomas-Rüddel, O. Ryabchykov, M. Bauer, N. Arend, E. J. Giamarellos-Bourboulis, J. Eugen-Olsen, M. Kiehntopf, T. Bocklitz, J. Popp, F. Bloos, U. Neugebauer, "Leukocyte Activation Profile Assessed by Raman Spectroscopy Helps Diagnosing Infection and Sepsis", *Crit. Care Explor.* 2021, 3(5), e0394.
- [5] N. Arend, A. Pittner, A. Ramoji, A. S. Mondol, M. Dahms, J. Rüger, O. Kurzai, I. W. Schie, M. Bauer, J. Popp and U. Neugebauer, "Detection and Differentiation of Bacterial and Fungal Infection of Neutrophils from Peripheral Blood Using Raman Spectroscopy", *Anal. Chem.* 2020, 92, 10560-10568.
- [6] U.-C. Schröder, F. Bokeloh, U. Glaser, J. Popp, U. Neugebauer, M. O'Sullivan, J. Ducrée, K. Wolf and W. Pfister, "Rapid, culture-independent, optical diagnostics of centrifugally captured bacteria from urine samples", *Biomicrofluidics* 2015, 9, 044118.
- [7] U. C. Schröder, C. Beletes, C. Assmann, U. Glaser, U. Hubner, W. Pfister, W. Fritzsche, J. Popp and U. Neugebauer, "Detection of vancomycin resistances in enterococci within 3 (1/2) hours", *Sci. Rep.* 2015, 5, 8217.
- [8] U.-C. Schröder, J. Kirchhoff, U. Hübner, G. Mayer, U. Glaser, T. Henkel, W. Pfister, W. Fritzsche, J. Popp and U. Neugebauer, "On-Chip spectroscopic assessment of microbial susceptibility to antibiotics within 3.5 hours", *J. Biophotonics* 2017, 10, 1547-1557.
- [9] J. Kirchhoff, U. Glaser, J. A. Bohnert, M. W. Pletz, J. Popp and U. Neugebauer, "Simple ciprofloxacin resistance test and determination of minimal inhibitory concentration (MIC) within two hours using Raman spectroscopy", *Anal. Chem.* 2018, 90, 1811-1818.
- [10] S. Kloß, B. Kampe, S. Sachse, P. Rösch, E. Straube, W. Pfister, M. Kiehntopf and J. Popp, "Culture independent Raman spectroscopic identification of urinary tract infection pathogens - A proof of principle study", *Anal. Chem.* 2013, 85, 697-708.
- [11] S. Kloß, B. Lorenz, P. Rösch, S. Dees, I. Labugger and J. Popp, "A Raman spectroscopy-compatible isolation procedure of single bacteria cells from sputum samples", *Anal. Bioanal. Chem.*, 2015, 407, 8333-8341.
- [12] T. Meyer, M. Baumgartl, T. Gottschall, T. Pascher, A. Wuttig, C. Matthäus, B. F. M. Romeike, B. R. Brehm, J. Limpert, A. Tuennermann, O. Guntinas-Lichius, B. Dietzek, M. Schmitt and J. Popp, "A compact microscope setup for multimodal nonlinear imaging in clinics and its application to disease diagnostics", *Analyst* 2013, 138, 4048-4057.
- [13] T. Bocklitz, F.S. Salah, N. Vogler, S. Heuke, O. Chernavskaja, C. Schmidt, M. Waldner, F.R. Greten, R. Bräuer, M. Schmitt, A. Stallmach I. Petersen, J. Popp, "Pseudo-HE images derived from CARS/TPEF/SHG multimodal imaging in combination with Raman-spectroscopy as a pathological screening tool", *BMC Cancer*, 2016, 16, 534/1-1.
- [14] T. Meyer, D. Akimov, O. Guntinas-Lichius, M. Schmitt, B. Dietzek, F. v. Eggeling and J. Popp, "Multimodal nonlinear microscopic investigations on head and neck squamous cell carcinoma: Toward intraoperative imaging", *Head Neck* 2012, 35, E280-287.
- [15] E. Rodner, T. Bocklitz, F. von Eggeling, G. Ernst, O. Chernavskaja, J. Popp, J. Denzler and O. Guntinas-Lichius, "Fully convolutional networks in multimodal nonlinear microscopy images for automated detection of head and neck carcinoma: Pilot study", *Head Neck* 2019, 41, 116-121.
- [16] P. Pradhan, T. Meyer, M. Vieth, A. Stallmach, M. Waldner, M. Schmitt, J. Popp, T. Bocklitz, "Semantic Segmentation of Non-linear Multimodal Images for Disease Grading of Inflammatory Bowel Disease: A SegNet-based Application", *ICPRAM 2019*: 396-405; DOI:10.5220/0007314003960405.
- [17] T. Meyer, D. Akimov, O. Guntinas-Lichius, M. Schmitt, B. Dietzek, F. v. Eggeling and J. Popp, "Multimodal nonlinear microscopic investigations on head and neck squamous cell carcinoma: Toward intraoperative imaging", *Head Neck* 2012, 35, E280-287.
- [18] S. Heuke, O. Chernavskaja, T. Bocklitz, F. B. Legesse, T. Meyer, D. Akimov, O. Dirsch, G. Ernst, F. von Eggeling, I. Petersen, O. Guntinas-Lichius, M. Schmitt and J. Popp, "Multimodal nonlinear microscopy of head and neck carcinoma - toward surgery assisting frozen section analysis", *Head Neck* 2016, 38, 1545-1552.
- [19] A. Lukic, S. Dochow, O. Chernavskaja, I. Latka, C. Matthäus, A. Schwuchow, M. Schmitt and J. Popp, "Fiber probe for nonlinear imaging applications", *J. Biophotonics.*, 2016, 9, 138-143.
- [20] A. Lukic, S. Dochow, H. Bae, G. Matz, I. Latka, B. Messerschmidt, M. Schmitt and J. Popp, "Endoscopic fiber probe for nonlinear spectroscopic imaging", *Optica* 2017, 4, 496-501.
- [21] P. Zirak, G. Matz, B. Messerschmidt, T. Meyer, M. Schmitt, J. Popp, O. Uckermann, R. Galli, M. Kirsch, M. J. Winterhalder and A. Zumbusch, "Invited Article: A rigid coherent anti-Stokes Raman scattering endoscope with high resolution and a large field of view", *Apl. Photonics* 2018, 3, 092409.
- [22] J. Tragardh, T. Pikalek, M. Sery, T. Meyer, J. Popp and T. Cizmar, "Label-free CARS microscopy through a multimode fiber endoscope", *Opt. Express* 2019, 27, 30055-30066.
- [23] B. Messerschmidt, E. Pshenay-Severin, G. Matz, K. Reichwald, J. Bierlich, J. Kobelke, H. Bae, T. Meyer-Zedler, M. Schmitt, J. Popp, "CARS endomicroscopic imaging probe enabled by a double-core double-clad fiber and related focus-combining micro-optical concept", *Proceedings Volume 11636, Optical Biopsy XIX: Toward Real-Time Spectroscopic Imaging and Diagnosis*; 116360A (2021) <https://doi.org/10.1117/12.2578538>.

## Improving Speed and Safety of Biopharmaceutical Manufacturing with HORIBA Analytical Solutions

Innovative spectroscopic (A-TEEM and Raman) and optical techniques (ViewSizer)

offer faster, accurate and less costly alternatives to chromatographic quantitative analysis in biotechnology

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The transition of biotechnologies from research project to pharmaceutical products requires analytical solutions that ensure safety and efficacy. In this paper we describe several innovative HORIBA products and highlight how these novel approaches rapidly deliver key information to accelerate processes in biopharmaceutical manufacturing. The HORIBA product portfolio can transform a variety of biopharmaceutical QC processes including the rapid characterization of cell culture media and viral vectors, vaccine manufacturing and analysis of components needed for cell and gene therapies.

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### Introduction

The publication of the human genome in February of 2001 was a tremendous achievement, and is just one of the multiple pieces of fundamental knowledge that has fueled the rapid pace of discovery and innovation in biotechnology. The tangible benefits to human health of this and other advances, like CRISPR-Cas9 gene editing, are just now being felt in the mainstream, as there is natural delay between discovery and commercialization. In many instances this transition to commercialization of transformative

innovations is slowed by the lack of complementary analytical technologies to confirm safety and efficacy. The innovation of novel analytical tools (or repurposing of existing ones) to enable biopharmaceutical manufacturing may not receive the same recognition as the initial discoveries, but commercialization absolutely depends on it. In recognition of the need for innovations in analytical tools, governments around the world have created public-private partnerships to solve this problem.

NIIMBL (National Institute for Innovation in Manufacturing

Biopharmaceuticals)<sup>[1]</sup> in the US, for example, is focused on the development of “novel real-time analytical technologies for integrated continuous processing for process integration and intensification.” Their goals are to enable “flexible and adaptive manufacturing, de-risk innovation, lower manufacturing costs and accelerate development and approval.” NIIMBL want to speed commercialization of novel biopharmaceutical product categories, such as Cell and Gene Therapy, and also want to streamline processes for existing products that lack efficient analytical tools such as: monoclonal antibodies (mAbs), vaccines, antibody-drug conjugates (ADCs), bispecifics, and virus-like particles (VLPs). Another government funded consortium is the Cell and Gene Therapy Catapult (CGT Catapult) in the UK,<sup>[2]</sup> aimed at accelerating the development of process analytical technologies (PAT) for cell and gene therapy manufacturing. HORIBA has joined CGT Catapult with the purpose of assessing the A-TEEM fluorescence technology as a way of reducing batch failures and manufacturing costs.<sup>[3]</sup>

### Traditional Analytical Approaches and Their Limitations

Chromatography is the workhorse in biotechnology, configurable with a wide range of detectors, stationary and mobile phase options, it is used to solve many analytical challenges. It is used for quantitative analysis, purification, classification, in R&D and for QA/QC. As omnipresent and useful as it is, separations techniques have drawbacks - chromatography is almost exclusively a lab-based technique. Although there have been attempts to make process-ready chromatography platforms, the challenges presented by these physical separations approaches are difficult to overcome robustly. Separations also tend to be quite expensive on a per-measurement basis, requiring solvents or gases for the mobile phase, columns and standards, and waste disposal fees for expended solvents. It is also slow compared to optical techniques - a fast chromatographic method typically takes tens of minutes to run compared to seconds-to-minutes for spectroscopy. Additionally, there are fundamental limitations in what

chromatographic methods can detect, as they are limited to the chemistry exposed on the outside of molecules, with no visibility to what may be buried inside. Ultimately separations techniques confirm **expected** components, therefore, unexpected components often go unnoticed, being “invisible” to the developed chromatographic method.

Molecular spectroscopy approaches like Fluorescence, Raman, near-infrared (NIR) and Fourier-transform infrared spectroscopy (FT-IR), on the other hand, are green, rapid, robust, and non-destructive technologies that can sit near-, at-, or even in-line. These techniques are sensitive to molecular environment in total, not just the molecules exposed on the exterior of a molecule, and spectroscopic libraries can be used to identify “unexpected” components. These attributes have contributed to molecular spectroscopy approaches, largely Raman and NIR, being incorporated as standard Process Analytical tools (PAT).

There is however a continuum of capability across these optical techniques, and some struggle to detect low concentration components, or to easily distinguish very similar components. Vaccines, ADCs, bispecifics, Adeno-Associated Virus (AAV) characterization (serotype differentiation, full vs. empty capsid) are all samples that require the highest level of specificity and sensitivity, so that even Raman cannot always reliably provide the detailed characterizations needed. Multidimensional excitation-emission matrix (EEMs) fluorescence should be the technique of choice based on its sensitivity and specificity. However, standard EEM implementations have struggled with repeatability and reproducibility. HORIBA has solved this with a patented 2-in-1 approach to multidimensional fluorescence spectroscopy called Absorbance, Transmittance and Excitation Emission Matrix (A-TEEM).<sup>[4]</sup> Table 1 compares A-TEEM with other optical spectroscopy techniques.

Summarizing, we find that in support of biopharmaceutical manufacturing, standard separations approaches,

Table 1 Comparison of Molecular Spectroscopy Analytical Techniques for Biopharmaceutical Analysis

Technique	Sensitivity	Selectivity	Comments	LOD (PPM)
A-TEEM (2-in-1) Fluorescence	↑	↑	UV/Vis & Fluorescence, Quantitative across broad concentration range	<0.001
2D/EEM Fluorescence	↑	↑	Poor analytical quantification, molecular fingerprint is concentration dependent	<0.1
Raman	●	●	Struggles with low concentrations	25-150
FTIR	●	●	Water interferes with molecular fingerprint	100-2000
UV/VIS	●	↓	Low selectivity Low information content	0.3
NIR	↓	↓	Struggles with low concentrations Low selectivity	0.1-1

although effective, are slow and mostly confined to the laboratory, and unable to characterize unexpected sample components, or differentiate samples based on non-surface characteristics. Vibrational spectroscopic techniques such as Raman, NIR, and FT-IR, are rapid and effective for many samples, able to operate on the manufacturing line, and can characterize samples on overall molecular environments, not just surface variations; but have constraints for low concentrations or highly similar samples. Fluorescence A-TEEM shares the speed and efficiency of Raman, adds the high specificity and sensitivity of fluorescence EEMs, and finally implements reproducibility and repeatability that transform it into a true analytical tool. We'll go into a bit more detail about A-TEEM implementation and applications in the next section.

### Fluorescence spectroscopy progression - from EEM to A-TEEM™

Fluorescence spectroscopy is an exceptionally sensitive tool. In some sense this sensitivity presents challenges for its adoption as a routine analytical tool, as clearly assigning a spectral change to a corresponding sample change can be quite complicated. There have been several advances though, that have enabled even complex fluorescence profiles to be converted into actionable information. Multivariate analysis and fluorescence Excitation Emission Matrix (EEM) spectroscopy taken together have advanced the use of fluorescence spectroscopy for routine use. An EEM is a 3D scan, resulting in a contour plot of excitation wavelength vs. emission wavelength vs. fluorescence intensity. EEMs are used for a variety of applications and, coupled with chemometrics, have shown success in monitoring batch endpoint and quality of final product in a bio-fermentation process.<sup>[5]</sup>

Fluorescence EEM spectroscopy unites the sensitivity of chromatography with the attributes of spectroscopic implementation, such as speed, robustness, low cost of measurement, and is effective for low concentration samples, and in looking for very slight environmental or molecular structural changes.

### The Limitation of Standard EEM Spectroscopy

The molecular identification capabilities of a standard EEM spectroscopy has some fundamental limitations. The inner filter effect (IFE) distorts the measured fluorescence spectrum of a molecule due to re-absorption that occurs at higher concentrations. The standard correction protocol uses a secondary measurement from a different absorbance spectrometer to adjust the measured fluorescence signal. But as the measurement does not occur simultaneously and with the same exact volumes, this approach doesn't provide a perfect correction. The imperfect match between data and correction means that an

EEM for a single component won't necessarily be reproducible, and this severely limits the ability of a traditional EEM for robust component identification as well as quantification.

### A-TEEM™ Spectroscopy Overcomes These Limitations

To overcome this limitation, HORIBA Scientific developed a new technology, called A-TEEM™ spectroscopy,<sup>[4]</sup> giving a spectrofluorometer the ability to simultaneously acquire Absorbance, Transmittance and a fluorescence Excitation Emission Matrix (A-TEEM) of a particular sample. Because the absorption spectra used for IFE correction are collected at the same time on the same sample as the fluorescence measurement, the correction is reproducible, and the A-TEEM data can be used for accurate sample ID as well as highly sensitive and specific component quantification.

A-TEEM provides *true and accurate representations of fluorescent molecules over a broad concentration range* (typically up to ~2 absorbance units). Additionally, absorbance and color information for all molecules, including non-fluorescent ones, is acquired and can be used in multivariate methods for comprehensive multi-component molecular identification. Interestingly, color can be a very sensitive indicator of protein product quality, and this information is acquired as part of the measurement. The A-TEEM method is a robust and extremely sensitive analytical technique, comprising the best attributes of both molecular spectroscopy (speed, low per measurement cost, lab-to-line placement) and chromatography (sensitivity and selectivity). A-TEEM therefore has the potential to replace traditional instruments like HPLC, GC, LCMS, and GCMS as a simpler, and faster analytical tool, with significantly lower per-sample measurement costs. A-TEEM is also a green approach to quantification compared to chromatography, as the need for a mobile phase goes away. Based on simultaneous use of fluorescence and absorbance spectroscopy, A-TEEM is validatable using the standard approaches for these two techniques.<sup>[6]</sup> A-TEEM molecular fingerprints are extremely information-rich, and multivariate analysis techniques are used to simplify the data and answer specific questions, providing direct and unambiguous results.

The A-TEEM fluorescence method provides a complete and traceable optical fingerprint of liquid samples, overcoming limitations that previously hampered adoption of fluorescence EEMs. This technique is well suited to the characterization of vaccines and other biopharmaceutical samples, with not only the sensitivity and specificity comparable to chromatographic methods, but also the speed, cost savings, and reproducibility of vibrational spectro-

scopic approaches.

### Cell Culture Media - Where Bioproducts Begin

The fundamental process in creating biopharmaceutical products is the growth of cells in bioreactors, fed by cell culture media, a specially formulated broth which contains a balance of nutrients to optimize cell productivity. The cost of the media itself is a small fraction of the value of the final product, therefore assessing its quality prior to the start of a bio-fermentation is a cost-effective way to ensure end-product quality and quantity. Standard approaches to perform this critical QC step are chromatography and mass spectrometry, making these measurements time consuming, expensive, and designed only to confirm what is expected, not to detect the unexpected. There is a growing interest in the use of molecular spectroscopic approaches for media QC, and both Raman and EEMs have been tested.<sup>[7,8]</sup> As is often the case, the combination of these two spectroscopic approaches may be the best way to cover the breadth of components that need to be characterized for this application.

As expected, A-TEEM provides a robust approach to assess the quality of cell culture media prior to the addition of cells.<sup>[9]</sup> We investigated eight different commercially available cell culture media samples, with the goal to differentiate between them, and also to study the sensitivity of the technique to detect environmental degradation. A more in-depth analysis of specific cell media components is a topic for future studies, and not addressed in this overview. Five samples from each of the eight types of media were analyzed, each sample was measured in triplicate resulting in 120 measurements, and each measurement took ~60 seconds. All cell culture media was stored at 4°C before testing and allowed to equilibrate to ambient laboratory temperature prior to analysis. The PCA showed clear differentiation between the different categories of media. To explore the specificity of the technique, a separate analysis including data


Figure 1 3D principal component scores plot showing ability of A-TEEM to differentiate within a cell culture category. DMEM samples of varying composition (types 1-3) were analyzed.


Figure 2 A-TEEM characterization of fresh from degraded cell culture media samples

from only a single category of media was performed. Figure 1 shows three different DMEM samples with slight compositional differences were easily distinguishable.

The final test for A-TEEM characterization of cell culture media was to compare fresh samples with those degraded through exposure to light, and elevated (room) temperature. A-TEEM was able to track the degree of degradation over a period of 5 days, as shown in Figure 2. A-TEEM provides a rapid approach to differentiating eight commercial cell media formulations, classifying three similar formulations, and tracking degradation of a cell media sample against adverse storage conditions.

### Viral Vectors Characterization

Researchers have turned to nature to understand how to deliver genetic material into human cells, and viruses provide a perfect mechanism. Several specific viruses, Adeno-Associate Virus (AAV) and lentivirus, for example, have been found to have appropriate characteristics, allowing scientists to manipulate them for the purpose of delivering genetically engineered payloads into human cells. These so-called viral vectors are the work horses of cell and gene therapy, and are making inroads into other applications, such as vaccines. In developing and deploying viral vectors, there are multiple parameters that are important to characterize for a complete picture of product quality. Measurements need to discriminate: 1) whole, broken, empty, aggregates, and infectious or non-infectious viruses; 2) empty from full capsids; and 3) vector serotypes. We'll present results from HORIBA's multi-laser Nanoparticle Tracking Analysis (NTA) approach using the ViewSizer 3000™ technology showing how it can be used to determine infectious titer, and then briefly touch on as-yet unpublished results from the A-TEEM method on its ability discriminate empty from full capsids and discriminate between vector serotypes.

### Infectious Titer Determination with ViewSizer 3000

The ViewSizer 3000 is a nano-particle tracking (NTA) platform that incorporates three lasers to collect the most accurate distribution and concentration data across a wide range of sample sizes. The point of the three lasers is to compensate for a variety of known problems when looking at a polydisperse sample: large particles scatter too much light and will saturate the detector from one laser, whereas small particles will only weakly scatter, and will be hard to detect. With three different lasers with independent power control, particles size can be determined across a very broad range, with scattering signals optimized to provide the most accurate distribution and concentration results.

The Viral Plaque Assay is the most robust approach to determine Infectious Titer, and it is highly manual and very time consuming. Viral preparations are serially diluted and allowed to grow on the planar surface of cell culture monolayer. Over a period of time (hours to weeks) the viral cells proliferate, infecting the confluent cells and destroying them. Regions of dead cells appear as empty patches on the culture surface, and these regions or “plaques” are then manually counted to determine the number of Plaque Forming Units (PFU) in a given dilution. Multiplying the PFUs by the dilution factor determines the final infectious titer of the parent stock, expressed in PFU/mL.

To determine the utility of the ViewSizer 3000 approach, the following protocol was followed: samples were transferred to the ViewSizer 3000 cuvette, fitted with a magnetic stir bar, and twenty-five video segments were collected, with five seconds of stirring in between each to ensure completely independent sets of particles in each video. The measurements were recorded with the following parameters: frame rate: 30 frames/second; exposure: 15 ms; gain: 30; blue laser power: 210 mW; green laser power: 12 mW; and red laser power: 8 mW; temperature control: active, 22°C. Each run took between 15-20 minutes. Multiple users independently analyzed the samples using the same measurement and acquisition parameters described above. The data were then processed and plotted using the same settings. Particle sizes ranged from 50 nm to 1000 nm, with a significant enrichment of particles around 200 nm, as shown in Figure 3. The total viral particle concentration within the sample was  $1.6 \times 10^7$  particles/mL. The distribution tailing towards 1000 nm suggests the existence of host cell debris or other background particles contained within the cell culture growth media. The image extracted from the video recordings further validated the presence of larger viral aggregates.

Infectious titer data from multiple assay replicates performed by independent analysts were plotted as a function of measured particle concentration from the ViewSizer 3000 and fitted with a result shown in Figure 4. The result indicates an excellent  $R^2$  value of greater than 0.9. When the results are analyzed using only data from one of the three lasers available, correlation dropped to an  $R^2$  value of 0.6 because of its inability to properly count the aggregates. By using three lasers, all particles were counted and the improved correlation is presumably because infectious aggregates contribute to the titer and were quantified along with single virus particles. This result indicates that over the course of multiple experiments, the ViewSizer 3000 accurately and reproducibly determined the viral titer within a heterogeneous sample.

This demonstrates the ViewSizer 3000’s ability to assess not only intact viral particles and aggregates, but also potential cell debris and background noise generated during upstream processing. Since downstream processing depends strongly on the removal of these impurities to minimize further aggregation, quantifying infectious viruses reliably depends on understanding both the purity of the upstream sample and the downstream analytical precision. The ability to detect these aggregates and cell


Figure 3 Measurement result of a human viral vector sample. Note the distribution tailing towards 1000 nm suggests the presence of host cell debris and aggregates.


Figure 4  $R^2$  value of > 0.9 was achieved between titer and data collected from the ViewSizer 3000.

debris is also critical when examining the stability of the final product over time. As we saw with the roll out of Coronavirus vaccines, the stability of biologic preparations within a freezer is a crucial factor to understand. With ViewSizer we are able to compare fresh preparations of a sample to ones that have been in the freezer for a number of weeks. By monitoring the increasing number of aggregates/cell debris that are present in the stored samples versus time, it is possible to gain an understanding of how long a sample is stable in a freezer. Utilizing the multiple lasers of the system is a crucial factor in these studies, as the particles of interest are usually around 80-120 nm, while aggregates and debris can appear all the way up to a micron; suggesting that the ability to measure across that entire size range is an indispensable tool in understanding the stability of these biologic preparations.

In addition to the applications described above, there are several unpublished studies that point towards other very interesting applications of A-TEEM, including the differentiation of AAV serotypes (AAV2 from AAV9), and quantification of empty vs full capsid ratio.

#### Vaccine Characterization and Manufacturing QC

The NIIMBL Vaccines roadmap<sup>[10]</sup> noted that vaccine product release time is on the scale of weeks to months, rather than the one to two days that might be possible with the right analytical tools. Vaccines have historically been relatively simple formulations, with one or multiple vaccine component (live-attenuated vaccines, subunit, polysaccharide, conjugate, or toxoid) with a diluent, and perhaps an adjuvant (a component that enhances the immune response). With the events of the last 18 months, vaccine technology has made huge advances in the fight against the Coronavirus pandemic, and mRNA vaccines, and the use of viral vectors has proliferated. The ViewSizer has been used as an analytical tool in the analysis of adjuvants.<sup>[11]</sup> However, the focus of this section will be on describing an approach for QC for vaccine manufacture that might meet the NIIMBL goal of one to two days for product QC.

One of the struggles for standard vibrational spectroscopic approaches is that vaccines tend to be formulated at very low concentrations. Raman spectroscopy without resonant enhancement struggles to characterize protein samples under ~1 mg/ml. Although Surface Enhanced Raman (SERS) can be used to enhance signals, it is challenging to fully characterize the impact that SERS substrates might have on the vaccine components. A-TEEM has been shown to have very low limits of detection, recently calculated to be ~0.15 µg/mL for a well-characterized vaccine formulation.<sup>[12]</sup> For these same low

concentration vaccine samples, A-TEEM has shown the ability to readily detect post-translational modifications, as well as differentiate between samples with simple amino acid substitutions, even those of non-fluorescent amino acids. In addition to the low limits of detection, A-TEEM has also been shown to characterize complex mixtures, differentiating and quantifying components against strong backgrounds, most notably in wine.<sup>[13]</sup>

To demonstrate the suitability of A-TEEM to the analysis of vaccines, we tested four over-the-counter multicomponent canine vaccine products. We used the SOLO-JEC brand multi-component canine vaccines (Table 2) from Boehringer Ingelheim VetMedica that were available at a local retail store (Tractor Supply Store). We tested SOLO-JEC 5 (“protection against five common canine infectious diseases”); SOLO-JEC 6 (“delivers the same protection as SOLO-JEC 5 plus additional protection against coronavirus”); SOLO-JEC 9 (“delivers additional protection against four types of leptovirus that are known to infect dogs”); and SOLO-JEC 10 (“delivers the same protection as SOLO-JEC 9 plus additional protection against coronavirus”).

The vaccines were reconstituted according to the package directions, diluted with distilled water in a ratio of 1:61 and added to a cuvette for measurement. To ensure repeatability, six individual measurements were collected for each vaccine, and this was done on two separate days. The calibration data set therefore consists of 48 measurements. To test reproducibility, a third unique set of validation samples were collected on a different instrument by a different operator, with 5 repeat measurements for each vaccine formulation, resulting in a validation set with 20

Table 2 OTC multicomponent vaccine products tested for this study

Disease	Basic Program Solo-Jec 5&6		Lepto Included Solo-Jec 9&10	
	Solo-Jec 5	Solo-Jec 6	Solo-Jec 9	Solo-Jec 10
Distemper	•	•	•	•
Hepatitis		•	•	•
Adenovirus 1	•			
Adenovirus 2	•	•	•	•
Coronavirus		•		•
Parainfluenza	•	•	•	•
Parvovirus	•	•	•	•
Lepto (4 types)			•	•
Gentamicin	•	•	•	•
Amphotericin B	•			
Thimerosal		•	•	•
Adjuvant	•	•	•	•

Table 3 Samples used for Calibration and Validation

Vaccine	Calibration set (n=48)	Validation set (n=20)
Spectra 5	S5	U5
	S-5	
Spectra 6	S6	U6
	S-6	
Spectra 9	S9	U9
	S-9	
Spectra 10	S10	U10
	S-10	


Figure 5 PCA scores scatter plot showing clear differentiation between 4 vaccine products, and excellent reproducibility for multiple measurements of a single vaccine type.

measurements, as shown in Table 3. A simple principal component analysis (PCA) was able to easily differentiate between the 4 vaccine products, with tight clustering of the calibration data.

In order to classify “unknowns” from the validation samples, a discriminant analysis model was created using Extreme Gradient Boost - Discriminant analysis. A-TEEM was able to identify and validate “unknown” samples with 100 percent certainty, as shown in Figure 6.

A-TEEM is the only spectroscopic approach that can perform this level of analysis on low concentration and complex vaccine formulations. Standard chromatographic QC approaches for vaccine batch analysis take days to complete. The results presented here demonstrate the potential for an on-line, rapid, sensitive and reproducible spectroscopic approach for vaccine QC, offer a potential solution to meet the goals for vaccine release that NIIMBL has set.

### Conclusion - The Future of Biopharma Characterization

HORIBA Instruments Inc. has a unique portfolio of analytical solutions that offer rapid, cost-efficient alternatives to chromatography, manual titer, qPCR, Analytical Ultra Centrifugation and other traditional analytical approaches. From A-TEEM to ViewSizer and with Raman as a complementary approach, these tools have the potential to


Figure 6 100% success in correctly identifying the vaccine product for unknown samples

shorten manufacturing times, and improve product safety with the ability to identify unknown components. The applications presented here are not meant to be an exhaustive compilation of the value that HORIBA tools provide to life science, but to highlight how spectroscopic and optical approaches offer a robust solution, as biopharmaceutical production methods scale to meet growing demand, and more information is needed more quickly. Starting with cell culture media, the fundamental raw material for all bioprocessing, both Raman and A-TEEM have proven utility in characterizing the quality of this material. Viral vectors are used to deliver novel vaccines and the components that make cell and gene therapies successful; ViewSizer and A-TEEM can be used to speed up analytical characterization steps that are needed in their manufacture. Finally, vaccines production has taken on new urgency, and A-TEEM has the potential to reduce QC time from weeks to hours. Other contributions in this version of Readout will explore other aspects and applications, including characterization of up-and-coming biotechnologies, like exosomes, and how HORIBA technologies are critical for the characterization of bioprobes.

\* Editorial note: This content is based on HORIBA's investigation at the year of issue unless otherwise stated.

## References

- [ 1 ] NIIMBL. <https://niimbl.force.com/s/about-niimbl>.
- [ 2 ] About Cell and Gene Therapy Catapult. <https://ct.catapult.org.uk/about-us>.
- [ 3 ] Editors, B. P. I. (2021, June 18). Horiba Joins CGT Catapult Consortium to Accelerate Development of PAT. BioPharm International. <https://www.biopharminternational.com/view/horiba-joins-cgt-catapult-consortium-to-accelerate-development-of-pat>.
- [ 4 ] Gilmore, A. M. & Tong, X. (2014, December 2). System and Method for Fluorescence and Absorbance Analysis, US Patent 8,901,153.
- [ 5 ] Peter P. Mortensen; Rasmus, Bro. Real-time monitoring and chemical profiling of a cultivation process. Chemometrics and Intelligent Laboratory Systems, Volume 84, Issues 1-2, 1 December 2006, Pages 106-113.
- [ 6 ] Adam M. Gilmore. How to Collect National Institute of Standards and Technology (NIST) Traceable Fluorescence Excitation and Emission Spectra. National Center for Biotechnology Information. March 27, 2014.
- [ 7 ] Li et al., "Fluorescence Excitation-Emission Matrix (EEM) Spectroscopy for Rapid Identification and Quality Evaluation of Cell Culture Media Components," Applied Spectroscopy volume 65, Issue 11, 2011.
- [ 8 ] Li et al., "Rapid Characterization and Quality Control of Complex Cell Culture Media Solutions Using Raman Spectroscopy and Chemometrics." Biotech & Bioeng. Volume 2, Issue 107, 2010.
- [ 9 ] Alvin Togonon, Adam Gilmore, Ben Yang, Karoly Csatorday and Marinella Sandros. Monitoring Cell Culture Media Variability using A-TEEM™
- [10] The National Institute for Innovation in Manufacturing Biopharmaceuticals Vaccines Roadmap
- [11] P. Singh et al., "Particle size analyses of polydisperse liposome formulations with a novel multispectral advanced nanoparticle tracking technology", International Journal of Pharmaceutics, Vol. 566, 2019, pp 680-686.
- [12] Unpublished results
- [13] R. K. R. Ranaweera, "Authentication of the geographical origin of Australian Cabernet Sauvignon wines using spectrofluorometric and multi-element analyses with multivariate statistical modelling", Food Chemistry 335 (2021) 127592

**Linda Kidder YARLOTT, Ph. D.**

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HORIBA New Jersey Optical Spectroscopy Center  
HORIBA Instruments Incorporated

**Adam M. GILMORE, Ph. D.**

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HORIBA New Jersey Optical Spectroscopy Center  
HORIBA Instruments Incorporated

**Sean TRAVERS, Ph. D.**

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Technical Sales Representative and Consultant for  
HORIBA Scientific  
Anatom Technology

## Fluorescent Bioprobes for Life Science Applications

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Bioprobes, and especially fluorescent bioprobes, have become essential tools for the detection and monitoring of important biological species and their related biological processes in live samples. Key applications of bioprobes fall into two main categories. The first one is bioimaging, from the cellular level, up to the organ and even small-animal level, *in-vitro* but also under *in-vivo* conditions. The second is biosensing, where probes are employed to detect various entities like ions, reactive oxygen species (ROS) or macromolecules, to help detect cancer, to observe drug internalization, and much more. This featured article will detail the typical characterization needs for the development and the use of the most common families of bioprobes: 1) determination of optical emission profiles; 2) measurement of size; 3) assessment of elemental composition or molecular signature. The article will showcase how HORIBA Scientific instruments are used by chemists all over the world to obtain rapid feedback on the best direction to follow to improve their synthesis processes.

#### key words

Fluorescence, Bioprobes, Sensing, Bioimaging

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### Introduction

Bioprobes are fluorophores, fluorescent chemical compounds that re-emit light upon light excitation, intended

for Life Sciences applications. Bioprobes can be called many things, tags, labels, sensors, reporters, or sensors, but they always refer to same thing. They are sometimes used alone, as tracers in fluids, as dyes for staining, as

Figure 1 General families of fluorescent bioprobes.

indicators (when their fluorescence is affected by environmental factors such as polarity, temperature, viscosity, or pH), or bound to macromolecules, serving as a marker for affinity or bioactive reagents.

Bioprobes are employed in a variety of domains from biomedical research to clinical applications and diagnostics. The first application is biosensing, which is the detection of biomolecules using an analytical device (or biosensor) combining a biological component with a detector. Small point-of-care devices used for diagnostics often rely on fluorescent bioprobes for their optical readout. This is also the case for devices intended for the monitoring of anthropogenic pollutants in various water matrices and aquatic environments.

Bioprobes are widely used to stain tissues, cells, or other biomaterials in a variety of analytical methods, including fluorescence imaging and spectroscopy. For instance, for *in-vivo* bioimaging in cancer research, to help with guided surgery in combination with other multimodalities. For that purpose, many groups are currently working on the development of near-infrared (NIR) emitting tags which are better adapted to deep-tissue imaging, taking advantage of reduced interference from photon absorption and scattering in this optical window.

High-content screening (HCS) systems used in the pharmaceutical domain for drug discovery also rely on fluorescent probes. The most common analysis workflow involves labeling proteins with fluorescent tags, while changes in cell phenotype are measured using automated image analysis in multiwell plates, in order to optimize synthesis processes or to confirm if a drug has the expected effect on its target.

Flow cytometry, which is used for biomarker detection, cell counting and sorting, or microorganism identification, is based on labels for measuring the physical properties of cells or particles. In most cases, fluorescent bioprobes with characteristic excitation and emissions peaks are designed to attach to a specific biological structure, for example a surface protein, to allow sorting cells

in a fluidics setup, based on the fluorescence measured in different detection channels.

A new and rapidly growing field is the clinical application of bioprobes for photodynamic therapy (PDT), which shows great promises for cancer treatment by killing cell with the heat produced locally by the particles upon light absorption.

Here, research focuses on finding dyes with properties like good biocompatibility, specific targeting, or strong two-photon absorption.

Bioprobes can be grouped into 3 main categories, as illustrated in Figure 1. Organic dyes (e.g., fluorescein, rhodamine), biological fluorophores (e.g., green fluorescent protein, phycoerythrin, allophycocyanin) and fluorescent nanostructures or nanomaterials (e.g. quantum dots, AIEgens).

Organic dyes are small molecules mostly synthesized through chemistry, but which can also be found in nature, like in plants. Common tags employed for confocal fluorescence microscopy are found in this category.

Biological fluorophores on the other hand are large molecules, typically fluorescent proteins, the most famous of all being GFP (Green Fluorescent Protein). DNA, which falls into this category does not have intrinsic fluorescence but can be used as biotags once coupled to a fluorophore.

Finally, we have nanomaterials that are engineered with specific chemical functions or tailored structures to generate unique light emission. This class of material encompasses a wide range of objects, like metallic nanoclusters, core-shell fluorophores, carbon dots, supramolecule assemblies and many more.

## Characterization needs

### Simultaneous Absorbance and Fluorescence

When studying novel bioprobes with unknown characteristics

or with emission properties impacted by solvent effects, having a complete fluorescence profile of a given material over a large range of wavelengths is key to avoid reperforming additional measurements at a later stage. An Excitation Emission Matrix (EEM), which is basically a three-dimensional scan that results in a contour plot of excitation wavelength vs. emission wavelength vs. fluorescence intensity, can provide such information (Figure 2).

Of importance as well when characterizing concentrated solutions (typically above 0.1 to 0.2 absorbance units), which is often the case when synthesizing bioprobes, is to account for the Inner Filter Effect (IFE). The IFE is a well-known phenomenon that distorts the measured fluorescence spectrum of a molecule due to absorption within the optical path length of the liquid sample, traditionally measured in a cuvette. IFE can significantly affect fluorescence emission spectral profiles, distorting their general shape, shifting the spectral position of the peak maximums, and decreasing the emission intensities, preventing from performing true quantitative analysis.

IFE can be corrected using absorbance data, and the HORIBA Duetta two-in-one spectrometer is unique in that it combines absorbance and fluorescence measurements on the same platform, allowing real-time correction for the IFE, on the same sample at the same time, as illustrated in Figure 2. It also takes advantage of an extended wavelength detection up to 1100 nm to cover the useful


Figure 2 Fluorescence EEM profiles obtained with a HORIBA Duetta spectrometer, with (top) and without (bottom) inner filter effect correction.

NIR biological transparency window.<sup>[1]</sup>

### Quantum Yield

The fluorescence quantum yield (QY) is defined as the ratio of photons emitted to photons absorbed by a fluorophore, basically an indicator of its brightness. This parameter is crucial for fast bioimaging consideration, as samples are often light sensitive or present rapid time kinetics.

A reliable method for recording the QY is a comparative method which involves the use of well characterized reference standards with known properties.

A second approach compatible with liquid samples is based on lifetimes and different concentrations of a

Figure 3 Quanta-φ integrating sphere coupled to a FluoroMax-4 spectrometer by HORIBA for quantum field determination.


Figure 4 Interface window of the FelixFL software on the HORIBA Fluorolog-QM spectrometer with integrated quantum yield calculator.

quencher to calculate the QY of a molecule in accordance with the Stern-Volmer equation. HORIBA customers from the Swiss Federal Laboratories for Materials Science and Technology (EMPA) recently reported on this method for their work on polyimide-based photoemitters for oxygen biosensing.<sup>[2]</sup>

However, the most accepted method to determine the absolute QY is to employ an integrating sphere (Figure 3), a sphere coated with an all-reflective surface to capture all the light going in and out.

A recent example of the use of the HORIBA Quanta- $\phi$  integrating sphere by scientists from the University of Chile explains the synthesis of an anthracene derivative dye with solvatochromic properties and highlights the importance of designing bright bioprobes for microscopy experiments on living cells.<sup>[3]</sup> Accurate quantum yield determination can be made easier using a dedicated application embedded in most HORIBA software, as seen in Figure 4.

### Lifetime

The fluorescence lifetime is a photophysical property of bioprobes, providing complementary information to the simple steady-state monitoring of fluorescence intensity. The technique has gained popularity thanks to its high sensitivity to microenvironmental parameters and changes in molecular conformation such as temperature, viscosity, pH, and ion concentration. At the same time, lifetime measurements are largely independent of fluorophore concentration and do not require calibration steps like intensity-based experiments.

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
Figure 5 Lifetime decay (red) and fit (green) to a 5-exponential decay equation, expressed in photon counts. The lower graph shows the standard deviation error residuals from the fit to the fitted data.

which are widely used to monitor protein-protein interactions to elucidate cellular signaling, disease progression and drug efficacy, also benefit from a lifetime monitoring, especially for *in-vivo* applications. Figure 5 illustrates the lifetime decay for a solution of Qtracker 655 quantum dots, collected on a HORIBA DeltaPro spectrometer equipped with a photomultiplier and a DeltaDiode485L pulsed laser source.

### Polarization

Fluorescence anisotropy or fluorescence polarization is a measurement of the changing orientation of a molecule in space, with respect to the time between the absorption and emission events. Anisotropy is an excellent tool for understanding changes in macromolecule shape, as well as molecular binding.

Figure 6. exemplifies the confirmation of effective binding between an antibody and the bronchodilator drug Theophylline labeled with a fluorescein dye and mixed with an antiserum in pH 7.2 buffer, then incubated for 3 minutes. The free and the complexed drug were distinguished through polarization measurements.


Figure 6 Polarization emission scan for labeled theophylline with (red) and without (blue) antibody ( $\lambda_{exc} = 485$  nm). The emission maximum for fluorescein-labeled theophylline is marked by a gray line.


Figure 7 Fluorescein dissolved in sucrose solutions. With increased temperature, viscosity decreases, yielding faster rotation times and lower anisotropies.

Figure 7 shows additional polarization measurement of fluorescein dye dissolved in four aqueous solutions, performed with a HORIBA Fluorolog-QM spectrometer.

This indirect measure of the local viscosity gives information on sample aggregation, structural changes, molecular binding, and other mechanisms.

### Bioimaging

Fluorescent bioprobes, especially small organic molecules, have become indispensable tools in modern biology. For research on living cells, these probes need to fulfill specific requirements like high permeability and good intracellular solubility, good affinity and specificity for the target, photostability, as well as low toxicity.

Microscopes with epifluorescence illumination are routinely used to observe the penetration and distribution of labelled targets within cell organelles to provide information on their structure and functions. Figure 8 illustrates


Figure 8 Epifluorescence imaging of Bovine Pulmonary Artery Endothelial cells (BPAE Line) with Anti- $\alpha$ -Tubulin antibody, stained with BODIPY, Texas Red and DAPI fluorophores, using a HORIBA XploRA microscope.


Figure 9 Left: Decay associated spectra for the FUN-1 dye in yeast cells grown in different growth media, also showing the sum spectrum that relates to the steady state emission. Right: FLIM measurement of DASPMI stained *Saccharomyces cerevisiae* in a polysaccharide polymer grown.

cell staining with multiple probes, Texas Red as a marker of proteins, DAPI as an indicator of membrane viability, and BODIPY to reveal neutral lipids.

The influence of environmental parameters on intra- and intermolecular interactions at the cellular level can also be obtained with microscopy instruments by monitoring the fluorescence lifetime of fluorophores instead of their intensity.


HORIBA now offers a new concept of wide field imaging, called FLIMera, capable of video-rate lifetime determination, based on time-correlated single photon counting technique (TCSPC) realized independently in each pixel.<sup>[4]</sup>

Figure 9 highlights research done with the FLIMera camera by HORIBA application teams in collaboration with Glasgow Caledonian University (UK). *Saccharomyces cerevisiae*, a single-celled fungus microorganism, was stained with a viscosity sensitive mitochondrial fluorescent probe called DASPMI to assess its viability. The cells displaying a longer average lifetime for DASPMI are associated with cells having a higher microviscosity, which is indicative of cell death (Figure 9 Right).

In addition, the decay associated spectra were determined for yeast cells suspended in growth medium, to investigate the photophysics of the fluorescence emission (Figure 9 Left). The use of the two-color FUN-1 stain, which transports into cell vacuoles, revealed the condition of plasma membranes through the formation of cylindrical intravacuolar structures.<sup>[5]</sup>

### Research based on HORIBA instruments

HORIBA Scientific is proud to count world-leading research teams among its customers. This following section will highlight a selection, far from exhaustive, of their recent publications in the field of bioprobes.

### Photodynamic therapy

Photodynamic therapy refers to treatments which involve light-sensitive therapeutics which can be switched “on” with a light source to destroy abnormal cells and treat various types of diseases.

The development of multifunctional bioprobes for simultaneous specific near-infrared (NIR) imaging and phototherapy of tumors being of great significance. For this purpose, our Fluoromax-4 spectrometer users at the City University of Hong Kong (China) reported on the integration of protoporphyrin IX (PpIX) into nanoscale metal-organic frameworks (NMOFs).<sup>[6]</sup>

Another example can be found in a recent publication in the journal *Nanomedicine: Nanotechnology, Biology and Medicine*, where a fluorophore called indocyanine green was co-encapsulated with a CT (Computed Tomography) contrast agent within nanoliposomes for imaging guided treatment. For this work, a Fluorolog-3 instrument from HORIBA recorded the fluorescence spectra in the NIR range and helped with the formulation of the particles.<sup>[7]</sup>

### Bioluminescence

Recently, our European laboratory, based South of Paris (France) on the Paris-Saclay campus, was closely involved in research with the famous Institut Pasteur.<sup>[8]</sup> Time-dependent chemiluminescence properties of core-modified coelenterazine luciferin analogues were recorded in the range of 250 to 800 nm using a HORIBA Aqualog spectrometer. These small molecules generated in bioluminescent organisms show great promise for bioimaging and therapeutics, by eliminating the need for external optical illumination and avoiding autofluorescence background in existing fluorescence techniques.

A recent publication highlights the application of such probes for the detection and imaging of peroxynitrite, a reactive nitrogen species (RNS), in live cells and mice *in-vivo*.<sup>[9]</sup> A HORIBA Fluoromax-4 spectrofluorometer was used to measure the bioluminescence emission of the synthesized compounds.

### Ion sensing

As detailed in the introduction, a major application for bioprobes is biosensing, especially to monitor inorganic ions which are essential for cellular activity.

Researchers from the College of Pharmacy at Seoul National University (Korea) have developed novel fluorescence cellular iron ions sensors to better understand cell death processes.<sup>[10]</sup> Their work required advanced characterization tools. A HORIBA FluoroMax-4 spectrometer was used to study the fluorescence response of their che-


Figure 10 Fluorescence spectra of the widely used FITC dye (Fluorescein isothiocyanate) in different solvents and under varying pH conditions.

mosensors in the presence of various metal ions. The brightness of the probes was assessed by measuring the absolute quantum yield with an integrating sphere assembly. Finally, kinetics information was obtained through fluorescence decay curves and analyzed using the DAS6 decay analysis software.<sup>[11]</sup>

Another group from the College of Chemistry and Materials Science from Nanning (China) proposed an efficient fluorescent sensor to detect Ag<sup>+</sup> but also L-cysteine in complex biological fluids and living cells using nitrogen, boron, sulfur co-doped carbon dots. The fluorescence emission profile, quantum yield, as well as lifetime of those nanostructures were determined using a HORIBA QM-8075 QuantaMaster instrument.<sup>[12]</sup>

### Microenvironment

Bioprobe fluorescence spectra are very sensitive to local environment, including solvent polarity, pH and temperature. This sensitivity is often leveraged in cell physiology research as many processes (metabolism, cell proliferation, membrane potential, etc.) are affected by changes in these parameters. Figure 10 shows the evolution of FITC emission properties under different pH conditions.

Besides pH, other local environmental factors can be assessed with specific tailored probes. Scientists from the National Academy of Sciences of Belarus have for instance designed carboxyfluorescein bifluorophores as viscosity sensors in model lipid bilayers or internal membranes of live cells.<sup>[13]</sup> Indeed, viscosity is one of the key parameters that controls the diffusion rate of molecular species and affects the reaction speeds of diffusion-controlled processes at the microscopic scale.

For their research, they have used a HORIBA Fluorolog-3 spectrometer to measure the fluorescence spectra and degree of polarization of the molecules.

### A focus on quantum dots

An important class of fluorescent nanomaterials is nanoparticles of semiconductors called quantum dots or QDs. QDs are artificially designed nanostructures which display unique emission properties depending on their composition and shape. They offer several advantages over organic fluorophores, among them brighter emission, photostability, broader excitation spectra combined with sharper emission peaks to avoid signal overlap for multi-color labelling.

As quantum objects, QDs emission wavelengths directly relates to their dimension (Figure 11. a), therefore controlling their size at different steps of their production is key. Using a HORIBA SZ-100 analyzer relying on dynamic light scattering (DLS), CdTe QDs with a passivation coating layer were measured in aqueous solution in a microcuvette and revealed a monodisperse distribution centered around 9.9 nm with a standard deviation of 3.4 nm (Figure 11. d)

Because of their brightness, QDs emission profiles can readily be measured with any spectrofluorometer.

However, we have seen before that due to inner filter effects (IFE), the fluorescence of bioprobes may not coincide with the true fluorescence emission, even for low concentrations.<sup>[14]</sup> This is visible in Figure 11.b where InP/ZnS QDs were measured with the HORIBA Duetta instrument capable of simultaneous IFE correction. Besides the obvious difference in fluorescence intensity, we see that the position of the peak maximum is blue shifted by 15 nm.

Apart from controlling the expected optical emission properties, verifying the elemental composition of synthesized or purchased QDs is important to validate chemical synthesis routes, but also for safety and ecological consideration. This information can be easily obtained by using X-Ray Fluorescence (XRF) spectroscopy. Figure 11. c. shows a XRF spectrum of commercially available CdTe QDs functionalized with an organosulfur compound, characterized under solid form with a HORIBA XGT-9000 analyzer.

Typical QDs have the disadvantage of showing potential cytotoxicity for biomedical applications, as their core-


Figure 11 a) QDs of different sizes showing different emission colors. b) QDs fluorescence emission profile extracted from an Excitation Emission Matrix. c) QDs elemental composition determined by X-Ray Fluorescence spectroscopy. d) QDs size distribution obtained by dynamic light scattering.

shell structures often contain heavy metals in the form of binary chemical compounds such as lead sulfide (PbS), cadmium selenide (CdSe), gallium arsenide (GaAs), or zinc telluride (ZnTe), which can generate reactive oxygen species, which in turn damage cellular proteins, lipids, and DNA. Therefore, significant efforts have recently been made to find alternative biocompatible materials.

Such a structure, with low toxicity and chemical inertness can be found in the form of carbon dots. Researchers from the University of Madras (India) recently reported on carbon dot nanocomposites conjugated with chitosan, acting both as a nano-drug carrier (for dopamine encapsulation in this case) and as a delivery tracer through bioimaging.<sup>[15]</sup>

Besides photoluminescence characterization of the probes with a HORIBA Fluorolog-3 spectrometer, the authors also employed Raman spectroscopy with a HORIBA LabRAM HR microscope to confirm the intensity ratio of the D band (1378 cm<sup>-1</sup>) to the G band (1580 cm<sup>-1</sup>) to control the presence of defects or amorphous carbon.

Another approach to avoid QDs toxicity is to encapsulate them in order to prevent their degradation when internalized, and to control their release at the right target. Our FluoroMax user at the University of Leeds (UK) proposed “sweet” QDs capped with mannose, (a sugar monomer important in human metabolism), to probe multivalent interactions of HIV/Ebola receptors using a sensitive, ratiometric FRET assay. The detection scheme is based on the fluorescence intensity ratio measured at wavelengths of 626 nm and 554 nm.<sup>[16]</sup>

## Conclusion

We have seen that fluorescence bioprobes have become tools of choice for various applications in bioimaging and biosensing.

HORIBA Scientific offers a comprehensive suite of analytical tools capable of covering the variety of characterization needs required for the design and the optimization of chemical synthesis routes of such novel probes.

\* Editorial note: This content is based on HORIBA’s investigation at the year of issue unless otherwise stated.

## References

- [1] [https://www.horiba.com/en\\_en/products/detail/action/show/Product/duetta-1621/](https://www.horiba.com/en_en/products/detail/action/show/Product/duetta-1621/)
- [2] *Reversible and Broad-Range Oxygen Sensing Based on Purely Organic Long-Lived Photoemitters*. Applied Polymer Materials, 3 (2020).  
<https://pubs.acs.org/doi/10.1021/acsapm.1c00064>
- [3] *CAPRYDAA, an anthracene dye analog to LAURDAN: a comparative study using cuvette and microscopy*. Journal of Materials Chemistry B, 8 (2020).  
<https://doi.org/10.1039/C9TB01738K>
- [4] [https://www.horiba.com/en\\_en/products/detail/action/show/Product/flimera-1989/](https://www.horiba.com/en_en/products/detail/action/show/Product/flimera-1989/)
- [5] *Viability of Saccharomyces cerevisiae incorporated within silica and polysaccharide hosts monitored via time-resolved fluorescence*. Photochemical & Photobiological Sciences, 12 (2013).  
<https://doi.org/10.1039/C3PP50202C>
- [6] *Harnessing combinational phototherapy via post-synthetic PpIX conjugation on nanoscale metal-organic frameworks*. Journal of Materials Chemistry B, 31 (2019).  
<https://doi.org/10.1039/C9TB01154D>
- [7] *Co-encapsulating indocyanine green and CT contrast agent within nanoliposomes for trimodal imaging and near infrared phototherapy of cancer*. Nanomedicine, 29 (2020).  
<https://doi.org/10.1016/j.nano.2020.102269>
- [8] <https://research.pasteur.fr/fr/project/design-and-optimization-of-new-bioluminescent-systems-better-suited-to-in-cellulo-and-in-vivo-conditions/>
- [9] *A Bioluminescent Probe for Imaging Endogenous Peroxynitrite in Living Cells and Mice*. Analytical Chemistry, 90 (2018).  
<https://doi.org/10.1021/acs.analchem.8b00198>
- [10] *Novel turn-on fluorescent biosensors for selective detection of cellular Fe<sup>3+</sup> in lysosomes: Thiophene as a selectivity-tuning handle for Fe<sup>3+</sup> sensors*. Dyes and Pigments, 169 (2019).  
<https://doi.org/10.1016/j.dyepig.2019.05.008>
- [11] [https://www.horiba.com/en\\_en/products/detail/action/show/Product/das6-1377/](https://www.horiba.com/en_en/products/detail/action/show/Product/das6-1377/)
- [12] *Red emission nitrogen, boron, sulfur co-doped carbon dots for “on-off-on” fluorescent mode detection of Ag<sup>+</sup> ions and l-cysteine in complex biological fluids and living cells*. Analytica Chimica Acta, 1035 (2018).  
<https://doi.org/10.1016/j.aca.2018.06.051>
- [13] *Fluorescent Properties of Carboxyfluorescein Bifluorophores*. Journal of Fluorescence, 30 (2020).  
<https://doi.org/10.1007/s10895-020-02535-w>
- [14] [https://static.horiba.com/fileadmin/Horiba/Application/Materials/Material\\_Research/Quantum\\_Dots/How\\_IFE\\_Can\\_Affect\\_Fluorescence\\_Measurements.pdf](https://static.horiba.com/fileadmin/Horiba/Application/Materials/Material_Research/Quantum_Dots/How_IFE_Can_Affect_Fluorescence_Measurements.pdf)
- [15] *Luminescent chitosan/carbon dots as an effective nano-drug carrier for neurodegenerative diseases*. RCS Advances, 10 (2020).  
<http://xlink.rsc.org/?DOI=d0ra04599c>
- [16] [https://www.leeds.ac.uk/news/article/3841/quantum\\_dots\\_light\\_the\\_way\\_for\\_new\\_hiv\\_and\\_ebola\\_treatment](https://www.leeds.ac.uk/news/article/3841/quantum_dots_light_the_way_for_new_hiv_and_ebola_treatment)

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HORIBA Scientific, France  
Global Life Sciences Market Manager  
Ph. D.

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Providing Solutions for the Life Science Field Using Spectroscopic Analyzers

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We provide solutions using surface plasmon resonance imaging device, particle tracking analysis device, Raman spectroscopic device, etc. for research on drug discovery and formulation of antibody drugs, quality control, and research on medical treatment and drug discovery using exosomes, which have been attracting attention in recent years. I will introduce the outline of these products and measurement examples, I will also introduce leading research and development in the life science field, which we are working on in collaboration with external organizations.

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
b)

Figure 6 Exosomes detection results with ViewSizer3000  
a) Fluorescence imaging  
b) Particle size distribution comparison

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Introduction of “Rapica” for Sensitive and Rapid Detection of ATP at the Level of One Microbial Cell

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Since microorganisms are invisible to the naked eye, and they can be visually confirmed when forming colonies by incubating under the appropriate culture media and temperature. Cultivation of microorganisms takes time: 4 to 7 days for pharmaceutical water used in pharmaceuticals, and 14 days for sterility tests used in release testing for sterile pharmaceutical products. HORIBA has developed Rapica, a rapid microbiological test system that can detect extremely small amounts of ATP possessed by microorganisms in a short time without culturing. We introduce its features and applications in this article.

key words

RMM : Rapid microbiological method, ATP method

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Figure1 Appearance of Rapica, the rapid microbial detection system


Figure 2 Overview of the measurement principle

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History of HORIBA Medical Products to Contribute to the In Vitro Diagnostic Testing

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Medical and health care is one of the greatest topic on development on society development and sustaining. Most of the medical and health care treatment is based on diagnosis based on the result from various testing, analysis and interviews. In vitro diagnostic medical devices are critical testing source to achieve necessary care to the patient. This article describes how HORIBA medical have been supporting the medical field by providing IVD medical device.

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Figure 4 Microsemi CRP

Figure 5 HORIBA POCT Product: (Left) Yumizen M100 Banalyst, (Right) Antsense Duo

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## Disruptive Technologies in Haematology for POCT Market

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A blood count is one of the most common blood tests that provides valuable physiological information about the patient's state. The tests are generally performed in specialized laboratories using "Gold standard" high-throughput diagnostic instruments. These automatic haematology analysers are bulky, need qualified operators for sample and data manipulation, and have more maintenance-related requirements. They also must respond to a growing demand for new complex parameters, complexities that are not required for urgent screening at the Point-of-Care. In this article, we present disruptive technologies for sample preparation and cell detection allowing the development of new generic platforms in haematology for a wide range of Point-of-Care applications.

key words

blood count, sample preparation, cell detection

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### Introduction

Blood count is one of the first prescribed and most ordered tests by doctors as it generally constitutes the first step in the diagnosis and monitoring of many pathologies and treatments. It provides physicians with valuable information on the patient's physiological state, allowing them to make treatment decisions. Nowadays, most automated hematology laboratory analyzers are based on costly, bulky, reagent consuming and complex hardware systems for blood preparation and measurement. Moreover, highly qualified and specialized laboratory technicians are required to handle such instruments.

Today and tomorrow's world and society are shaping up new challenges for the healthcare systems such as the one generated by the Covid 19 pandemic. Therefore, for delocalized and personalized health care applications, there is an increasing need for automated, compact and transportable point-of-care blood cell counters, operable without requiring any training or maintenance. Such solutions would eliminate the need for patients to travel into lab testing facilities and shorten the time needed for doctors to make treatment decisions. Point-of-care devices can be suitably deployed at the bedside, private clinics, research laboratories, rural areas as well as in developing countries.


Figure 1 Scheme of the blood analyser principle based on innovative disruptive technologies

Using simple and smart technologies for blood preparation and cell measurement, low power and reagent consumption, these future low cost automated and miniaturized blood cell counters offer a genuine route towards a greater portability of hematology analyzers.

The work presented here is result of an “open innovation” collaborative program held since 2015 with a leading player in technologies development, the French Atomic Energy Commission (CEA), which is the biggest patent applicant in France. By combining their own expertise, both parties can accelerate the delivery of impactful technology solutions (Figure 1).

In this paper, we introduce first a new disruptive technology for cells detection, counting and differentiation, based on holographic lens-free imaging. Since all the optics and mechanical parts are removed, this low-cost technology, embedded in a device, should lead to ease-of-use with minimal operator training and remove the need for user maintenance. Secondly, we introduce an innovative microfluidic preparation cartridge for automated blood cell counters allowing low consumption and scaled-down integration. The whole blood preparation system embedded in these POCT devices answers different key challenges that need to be addressed regarding microfluidics technology: i) a large range of dilution, ii) dilution accuracy, iii) low cost per test and iv) maintenance free.

These two patented technologies<sup>[6, 7, 8]</sup> are well adapted to a POCT system, and, in the future, will enable us to propose a different approach to hematology, closer to the patient who needs it.

## Lens-Free Imaging

The Lens-Free Technology (LFT) is a simple imaging technique, developed in recent years and under continuing evolution and improvement. Based on in-line holography and numerical processing, LFT gives access to a reconstructed optical field revealing phase and absorption contrasts of the studied samples, which can be cells, particles, crystals, etc... It can generate highly contrasted images of biological objects on a large field-of-view with a micrometer-scale resolution using a very simplified optical setup without magnification. The components of a lens-free imaging system are uncomplicated, small and low-cost. In addition, the wide field of view, together with the high resolution of the sensor, leads to a count (> 10 000 cells) in a single-shot acquisition, consistent with the statistical performances required in haematology. The complexity of the system is transferred to digital processing requiring robust and powerful reconstruction algorithms. Hence, this technology allows the building of an automated, simple, cost-effective, robust, light weight and compact system meeting the requirements for screening point-of-care (POC) tests in human and veterinary markets (Figure 2).

In practice, the sample is simply positioned in a chamber between a light source and a high-resolution optical sensor, for example, a Complementary Metal Oxide Semiconductor (CMOS) (typically a resolution of around 10 MPix and 2  $\mu\text{m}$  pixel size). This sample is illuminated from above by the point-like source of light and the sensor records interferences patterns (holograms) created by a combination of the light diffracted by the objects and the transmitted light. The raw data are then computerized with a numerical reconstruction algorithm, and dedicated


Figure 2 The advantages of the LensFree Imaging system, compared to classical Flow Cytometry setup

image processing workflows are specifically designed for each new application. Briefly, the reconstruction step allows, starting from the raw diffraction patterns acquisition, the reconstruction of an in-focus absorption map. From the acquired defocused hologram, classical gradient autofocus algorithms were applied to compute the position of the object plane (i.e. the focus plane). In this selected plane, an iterative reconstruction algorithm from a single image based on specific norm minimization is used. On this reconstructed image, grey-level thresholding is performed to detect the more contrasted objects. A binary mask is created: a 0 value is associated to the background,

and a 1 value is assigned to the detected components. At this stage, there are between hundreds and thousands of components, depending on the sample. The detected components are sorted out based on morphometric characteristics or finer supplementary criteria using grey-level. This leads to the identification and characterization of most of the blood cells, depending on the considered application (Figure 3).

This new technique for measuring cells requires increased computation capacities. We implemented algorithms on an Nvidia Jetson Nano Graphical Processing Unit (GPU)


Figure 3 The steps of the algorithm for image reconstruction and cells detection and identification.


Figure 4 Summary of haematological analytical performances obtained with LFT in comparison with classical measurement systems embedded on current HORIBA Medical routine analysers

for massive parallelization. Increasing connectivity and the versatility of this sensor allow us to consider various and remote application fields.<sup>[2]</sup>

As a proof-of-concept, several applications have already been investigated with LFT, in the human and veterinary haematology field. We display, in Figure 4, performances obtained for red blood cells (RBC), white blood cells (WBC) and platelets (PLT) counting, as well as for lymphocyte (LYM) and neutrophil (NEU) identification, in comparison with results obtained with our reference HORIBA Medical analyzers (ABX Pentra 120 and Yumizen H2500).

Those preliminary results show good agreement between the LFT and conventional hematology analyzers on this set of parameters. No major discrepancies were observed among the different sample database evaluated (more than 100 normal and pathological samples). R correlation parameters are all located above 0.95.<sup>[1, 3, 4]</sup>

### Microfluidic Preparation

In this project, the main objective is to develop microfluidic cartridges allowing precise preparation of small samples with great accuracy and repeatability, and scaled-

down integration for the associated analyzer. However, developing such a system with a small footprint comes with several challenges regarding microfluidics technology such as: sample introduction, reagent management, efficient fluid sample mixing, low cost-per-test and more importantly, large yet precise dilutions (Figure 6a).

To reach the desired high dilution ratios of the small blood samples and still achieve high precisions, we opted for a technology based on a hyper-elastic membrane to create collapsible chambers with large volumes and high aspect ratios (Figure 6b). Optimal configuration of these chambers and channels was reached through extensive Computational Fluid Dynamic (CFD) simulation studies using the YALES2BIO software (Figure 5). This technology offers a high-volume precision as the amount of used liquids is directly controlled by the chamber size. Furthermore, the microfluidic cartridge can be easily replaced by a “Plug & Play” system, leading to a maintenance free system.

For a first application case as a proof-of concept, the cartridge was designed to handle two precise dilutions simultaneously, one at 1/1000 ratio and another one at 1/10 ratio. The first dilution of the whole blood is performed for counting the RBCs and the second dilution for


Figure 5 Comparison of sample injection on a dilution chamber. Experience confirms optimal results obtained by CFD computing.


Figure 6 (a) The advantages of the Microfluidic preparation, compared to classical technology, (b) Principle of FlowStretch Technology. Chambers and valves are actuated using pneumatic pressure (COC : Cyclic Olefin Copolymer)

counting the WBC and the PLTs, and for the Hemoglobin measurement. During this second dilution, a lysis of RBCs is performed. This lysis step is necessary for two reasons: first, RBCs are numerous and can hinder the proper visualization of WBCs. Second, the lysis of RBCs allows the release of hemoglobin in the solution, and in this way the measurement of the hemoglobin concentration can be performed.<sup>[5]</sup>

To achieve the dilution protocol, a specific architecture and an associated protocol is designed. The several zones for different preparations are visible on Figure 7(a) one sample preparation for RBC counting (Zone I), and another one for WBC/PLT counting and Hemoglobin measurement (Zone II). The whole blood is first sampled in two successive cali-brated meandershaped channels (Figure 7(b)). Each meander channel is connected to two stretchable chambers that can be filled with reagents. By successive actions of the two opposed chambers, the calibrated blood sample is efficiently and quickly diluted and

mixed with calibrated reagent (Figure 7(b)). The precise and repeatable dilution ratio is simply provided by the ratio volume of the meander channel and collapsible chamber. Thus, the resulting dilution ratio is generated solely by the geometry of the cartridge and does not depend on fluid properties. Repeatability performances obtained for blood parameters with microfluidic preparation are shown in Figure 8 (Left). In Figure 8 (Right), the correlation of Hemoglobin measurement obtained with mock-up coupling Microfluidic sample preparation and LensFree Imaging sensor and the one obtained with classical routine analyser is depicted.

Figure 8 (Left) shows a good agreement between the expected statistical performances and the actual Coefficient of Variation (CV) measured during repeatability standard operating protocols. For each parameter, the actual value is slightly higher than the target. This is a normal effect of the variability introduced by the hardware such as the pneumatic driving of the actuated


Figure 7 a) The architecture of the microfluidic cartridge with the two regions corresponding to the two performed dilutions. b) Representative images of the sampling and mixing steps in the meander shaped channel with the dilution chambers


Figure 8 (Left) Repeatability performances for blood parameters with microfluidic preparation - (Right) Performances on Hemoglobin Measurement with mock-up coupling Microfluidic sample preparation and LensFree Imaging sensor

membrane. Likewise, the correlation on the hemoglobin parameter (Figure 8 Right) shows a good agreement with our qualified internal reference HORIBA Medical Yumizen H2500 It suggests an accurate preparation of the blood sample, from its sampling to the lysis, the dilution and the measurement of concentration.

### Conclusion

These innovative technologies pave the way to the development of new generic platforms for a wide range of POCT applications. The technology leads that we are following are opening a new way of performing diagnostics. By reducing the hardware complexity and the footprint of the instruments, we may now envision smaller, more reliable and more versatile devices. Whether it is at the doctor's office, closer the patient's bed at the hospital, or somewhere remote in the field, we propose two key assets for building more efficient ways of delivering the patient the care he/she needs. We have demonstrated a proof-of-concept feasibility on essential parameters. We now need to work on the enhancement of the algorithms involved in the computation of the LFT results. Indeed, since our primary goal is to deliver impactful insights to the physicians, we need to ensure now that we can address pathological conditions as well as healthy ones.

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References

- [ 1 ] A. Ali-Cherif et al., International Conference on BioMedical Photonics | 16-17 March 2018, Towards a disruptive technology in Hematology: lens-free imaging for high precision red blood cell counting and sizing.
- [ 2 ] A. Ali-Cherif, et al., MIFOBIO 2018 (GDR ImaBio) | 5-11 Octobre 2018 : Imagerie holographique grand champ pour l'hématologie.
- [ 3 ] D. Isèbe, et al. "High accuracy platelet counting using LensFree imaging", Proc. SPIE 11075, Novel Biophotonics Techniques and Applications V, 110750O (22 July 2019);
- [ 4 ] T. Lavabre and al., Detection of circulating microfilariae in canine EDTA blood using lens-free technology: preliminary results, Journal of Veterinary Diagnostic Investigation, 2021, vol. 33(3) 572- 576.
- [ 5 ] D. Isebe, et al. Reusable microfluidic device for complete blood count applications, MicroTAS, October 2019
- [ 6 ] Allier, C. Method for identifying blood particles using a photodetector US10379027 filed on March 23, 2016 and issued on August 13, 2019.
- [ 7 ] Patent. Holographic method for characterizing a particle in a sample US10845286 filed on March 22, 2017 and issued on November 24, 2020
- [ 8 ] Patent. Microfluidic sample preparation device offering high repeatability US17/291139 filed on November 07, 2019.

**Damien ISEBE, Ph.D.**

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Innovation Project Manager  
Innovation Dept.  
HORIBA ABX SAS

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Development Engineer  
Innovation Dept.  
HORIBA ABX SAS

**Benoit THOUY, M.D.**

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Algorithm Engineer  
Innovation Dept.  
HORIBA ABX SAS

**Aurélien DAYNES, Ph.D.**

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Scientific Expert  
Innovation Dept.  
HORIBA ABX SAS

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Innovation & Technology Manager  
Innovation Dept.  
HORIBA ABX SAS

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