Readout

HORIBA Technical Reports

August 2022 No. **56**

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
HORIBA contributes with its analytical and measurement technology in three areas that lead megatrends: “energy and environment,” “biotechnology and healthcare,” and “advanced materials and semiconductors. In this issue, we introduce our analysis solutions for these three areas.

In search of a refreshing breeze blowing through the trees, I visited the highlands of Shinshu. As I strolled through the greenery, the sunlight filtering through the trees made me feel as if I had wandered into the world of an impressionist painting.

*-Photographer MATSUI Hideo-
(Member of Nikakai Association
of Photographers)*

Name of this Journal

This Journal is named “Readout” in the hope that “the products and technology we have created and developed will be read out and so become widely known”.

Analytical Solutions in Megatrends

Foreword

- 4** HORIBA's Global Analysis Services for Megatrends
KOMATANI Shintaro

Guest Forum

- 6** Frontiers and Standardization Trends of Particle Measurement and Analysis Technology
FUJIMOTO Toshiyuki

Feature Article

- 13** New Platform for Multimodal Spectroscopic Characterization of Semiconductors
Praveena MANIMUNDA, Francis NDI
- 18** In Operando Characterization of Li Composite Battery Electrode by Ellipsometry and Raman Spectroscopy. The Case of $\text{Li}_4\text{Ti}_5\text{O}_{12}$ Based Anode.
Valerie SILLER, Alex MORATA, Juan Carlos Gonzalez-ROSILLO, Marc Nuñez EROLES, Michel STCHAKOVSKY, Albert TARANCON
- 24** High-resolution Piezoresponse Force Microscopy Imaging of a Few Layered 2D Ferroelectric
Ana I. PEREZ-JIMENEZ
- 30** Particle Size Studies for Biotechnology and Life Science
- Application Proposal from HORIBA Analytical Solution Plaza -
SHINOZAKI Yoko, SAIHARA Koji, MITSUNARI Kyoko
- 36** Introduction of Non-destructive Elemental Analysis for In/Online Through Real Case Studies
AOYAMA Tomoki, UEDA Hideo

Facilities Introduction

- 40** HORIBA Techno Service Introduce “Analytical Solution Plaza”
HORIBA's Services that Connect “Analysis” and “Understanding (Finding out)”
SAKAGUCHI Mai
- 46** Introduction to HORIBA FRANCE Applications Laboratory: Missions and Resources
Alessia QUATELA
- 50** HORIBA Technical Centers at Premier Institutes in India
Dr. John Kiran ANTHONY, Dr. Rajeev GAUTAM, Dr. H.C. SUDEEKSHA
- 54** Introduction of Application Centers in HORIBA China
Jing SHEN

Topics

- 57** Kinki Region Invention Awards 2021
- 60** Award Winners of HORIBA Group IP World Cup 2021
- 64** Scientific Research Papers from HORIBA Group
- 72** HORIBA World Wide Network

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Frontiers and Standardization Trends of Particle Measurement and Analysis Technology

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This report summarizes the frontiers of measurement and analysis technologies related to particle size, which is a fundamental and extremely important Physico-chemical property of fine particles, with a focus on compliance with regulations in Europe, where leading technical standards are being introduced, as well as recent topics in standardization. In the case of particle size measurement in air, the calibration of detectors, which is directly related to the regulations, and in the case of particle size measurement in liquid, the number-based size distribution, which is the basis of the EC definition of nanomaterials, were introduced, including attempts to obtain international equivalence by using various measurement methods.

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Table 1 Requirements of particle detection efficiencies under European vehicle emissions controls

	導入年	粒子の検出効率
Euro5 ⁺	2011	50%±12%@23 nm 90%以上@41 nm (SPN23)
Euro6	2014	
Euro6a	2017	
Euro7	2025	65%±15%@10 nm 90%以上@15 nm (SPN10)


Figure 1 Requirements of particle detection efficiencies under European vehicle emissions control (indicated by bands) and typical detection efficiency curves of a CPC. Orange is for SPN23, blue is for SPN10 (planned).

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
Figure 2 Schematic representation of the CPC calibration setup.

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
Figure 4 Particle size distributions of bimodal reference sample evaluated by DLS and LD.
DLS (Dynamic Light Scattering), LD (Laser Diffraction)

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
Figure 5 Number based particle size distribution evaluated by PTA method for a bimodal reference sample.

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AFM(Atomic Force Microscopy),
 CLS(Centrifugal Liquid Sedimentation),
 CPC(Condensation Particle Counter),
 DMA (Differential Electrical Mobility Analyzer), DLS (Dynamic Light Scattering),
 ICP-MS (Inductively Coupled Plasma - Mass Spectrometry),
 LD(Laser Diffraction),
 PTA(Particle Tracking Analysis),
 SAXS(Small-Angle X-ray Scattering),


SEM(Scanning Electron Microscopy),
spICP-MS(Single Particle ICP-MS),
TEM(Transmission Electron Microscopy),
UV-Vis(Ultra Violet - Visible)

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
a) Carbon Nano Tube : Particle size distribution and electron microscope images.

b) Nanocellulose: Particle size distribution and correlation with AFM.

Centrifugal Nanoparticle Analyzer

New Platform for Multimodal Spectroscopic Characterization of Semiconductors

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Francis NDI

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Multimodal spectroscopy is the concept of combining several different spectroscopies onto one platform, thereby expanding the range of analytical capabilities available on that single platform. In advanced semiconductor fabrication, the need to streamline production and avoid unnecessary wafer handling makes it quite common to have several different metrologies at one wafer processing station. However, this concept is still not common at the material research or product & process development stages where it is more common to have several different instruments offering different analytical spectroscopies. Besides the obvious benefit of cost reduction, having multiple analytical spectroscopies offers the added benefit of sample co-location so that multiple complementary measurements can be made at the same location of the sample. The benefit of co-location is particularly important as feature sizes get smaller, from a few microns to nanometers in size. We show results from a multimodal microspectroscopy system from HORIBA for semiconductor characterization. The Standard Microscope Spectroscopy (SMS) system is a modular and flexible system capable of accommodating up to seven different spectroscopies, including Raman, Photoluminescence, Time-resolved Photoluminescence (Lifetime), Reflectance and Transmittance, etc. Furthermore, we show that even when there is a need to perform correlative measurements across different instruments, it is possible to use coordinate system transformation technologies, such as nanoGPS navYX™ from HORIBA to rapidly achieve sample co-location across different measurement platforms.

Key words

multimodality, wafer handling, speed, cost minimization, reproducibility

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Introduction

Semiconductor devices are some of the most complex engineered devices in today's world. Their manufacture often requires an even more complex sequence of steps, first to understand and tailor the material properties so that desired electrical, optical or mechanical properties are possible, followed by very intricate and controlled fabrication processes so that devices are reproducible and reliable, and at a reasonable cost. The above sequence of events demands that the practitioner is able to measure the material and device properties at every step of the design and production process to ensure that yields are optimal, and to remove defective material or components as early as possible in the production process, as every step adds cost and complexity.^[1] A vast number of techniques have been developed for these types of material, product and process characterization. Amongst these, optical or spectroscopic techniques are usually preferred. The preference for optical spectroscopy in semiconductor material and process characterization derives from several factors. Optical spectroscopic (OS) techniques are often non-contact, hence reducing the chances of contamination. Secondly, many OS techniques are comparatively fast and require little to no sample preparation.

Given the large number of material and process parameters that can affect the behavior of the final semiconductor product, it follows that control of the material design and fabrication process is necessarily a multimodal effort. This means that many measurement techniques are required to measure and control all the parameters that lead to desired device behavior. Instrument vendors have accordingly developed a wide variety of specialized optical spectroscopy instruments often specialized on each of the required techniques. In addition, as device features go from micro to nano, many of these instruments are based on a microscope or other submicron to nanometer measurement platform. For example, it is common in a semiconductor research or fabrication facility to have one micro-Raman spectroscopy instrument used to characterize crystallinity of epitaxial deposition or stress, and a separate photoluminescence instrument to measure wafer homogeneity, etc. Beyond the cost burden of having multiple instruments to perform these necessary measurements, the task itself, in going from instrument to instrument, has become quite challenging recently, as the features of interest become smaller and approach the micro to nanoscale.

In this article, we describe a novel approach to achieve multimodality on one microspectroscopy platform (Figure 1), enabling the practitioner to characterize various semiconductor samples using different spectroscopies

Figure 1 HORIBA multimodal, Standard Microscope Spectroscopy (SMS) system equipped with 12-inch wafer stage.

with the benefit of sample co-location. In this approach, various complementary measurements can be performed at the same micro location and in so doing, obtain deeper insights into the sample or process. In addition, when there is a necessity to perform micro or nano measurements across different instruments, we present a new coordinate transformation technology (nano-GPS) that enables fast and accurate location of nanostructures across different measurement instruments.

Example 1: LEDs

One of the most important considerations of LEDs is the emission wavelength. For example, most common Lidar sources have semiconductor materials designed to have light emission at 905 nm or 1550 nm. The primary material property that controls this parameter is the optical bandgap, which is readily measured using micro photoluminescence (Figure 2a). After the material properties are understood, a fabrication process engineer could be


Figure 2 (a) Photoluminescence (PL) map on a InP wafer. (b) PL spectrum (insert) and histogram showing intensity distribution- an indication of homogeneity. (c) High resolution PL map on a defect region (dotted square on (a)). (d) PL lifetime decay curves recorded at three different locations on a defect.


Figure 3 (a) Representative electroluminescence image of AllnGaP LED and (b) Optical micrograph and (c) PL intensity map recorded from a InGaN LED. Defect region appeared as cold spots on PL map.

interested in the uniformity of the epitaxial deposition on a large wafer, so that a Lidar laser die cut from one part of the wafer can be expected to perform and behave as one from any other part of the wafer. Once again, this property is readily characterized using a photoluminescence (PL) intensity distribution across the wafer (Figures 2a and b). If, for some reason, some of the dies from this wafer do not show optimal luminescence compared to others, the process engineer might be interested in understanding the nature of the defects causing this sub-optimal luminescence efficiency (cold spots). Figure 2c shows a high resolution PL map of a defect. Sometimes, the charge carriers are captured at the defect sites and don't emit light.^[2] Time-resolved PL sheds more light on defect properties.^[3] Figure 2d shows PL lifetime decay curves from three locations close to a defect. At the defect center, significantly shorter lifetime was observed.

As part of a QA process, the engineer might want to measure the performance of the device under conditions similar to actual use by measuring the electroluminescence of the device before it is packaged (Figure 3a). PL mapping can also be used to analyze the final device as well. Figure 3 (b-c) shows results from an InGaN LED sample, where the PL map detected local luminescence variations.

Example 2: Photovoltaics

Development and fabrication of Photovoltaics show a similar need for multimodal characterization. For example, at the material stage it is important that the material bandgap be engineered to optimally absorb the solar spectrum. Once again, PL is a good technique for determining that property (Figure 2a). The photovoltaic effect relies on the efficient movement of charge carriers either to the electrical load for use or to a battery for energy storage. Time-resolved PL is often used to characterize carrier dynamics (Figure 2d) or Raman spectra to


Figure 4 (a) Raman map and spectra (insert) of a possible contaminant flake on the surface of a semiconductor sample. (b) Photocurrent map of a sheet of Silicon PV material showing hot spots (possible defects). The red strip is a piece of the conducting electrode.

determine micro crystallinity (which in turn affects carrier dynamics (Figure 4a)). Finally, in the QA process of the solar cell device, one might be interested in measuring the overall device efficiency by measuring the spectral photocurrent response (Figure 4b).

Although contrived, the above measurements are typical occurrences in the design and manufacture of LED and

PV devices, and it is also common to see that all these measurements are typically made on different instruments. The novelty in this paper is showing that they can be made on one instrument (the SMS system from HORIBA), resulting in cost savings and adding convenience to the process.

Correlative spectroscopy on semiconductor materials:


Simple, fast and non-contact microspectroscopy techniques such as the ones described above, are usually preferred in semiconductor material characterization, but it is sometimes necessary to use other, more complex, techniques that are not easy to combine with the above. For example, defect characterization sometimes requires high spatial resolution (nanometers) that is only available on instruments such as a scanning electron microscope (SEM) or an atomic force microscope (AFM). In those instances, and due to the cost and complexity involved in using a specialized instrument such as an SEM, it is desirable to establish a correlative optical spectroscopy so that such a defect can be identified in the future using a simpler spectroscopic technique, rather than doing the measurement on an SEM repeatedly, which can be costly and slow. To achieve this correlation, it becomes necessary to identify a nanoscale feature in the SEM, and to also be able to identify the same feature under an optical microspectrometer which can be laborious and time-con-

suming without some type of coordinate system matching between the two instruments. To facilitate this process, automated nanoscale coordinate transformation technologies, such as nanoGPS navYX™ from HORIBA have been developed to enable fast, accurate and repeatable localization of nanoscale objects between different measurement systems (Figure 5a). Figure 5c shows the cathodoluminescence map recorded from a Ga₂O₃ nanostructure. Using nanoGPS, the sample was transferred to a HORIBA SMS system and a PL map was recorded successfully (Figure 5c).


Conclusion

In conclusion, speed, cost minimization and reproducibility are persistent drivers for decision making when it comes to instrumentation choices for semiconductor research and fabrication. In this paper, we introduce a novel modular and multimodal platform that enables the efficient combination of several complementary spectroscopic techniques relevant for semiconductor characterization on one platform. Furthermore, and for when it is necessary to measure across different platforms, we introduce a new coordinate transformation technology that enables the fast and accurate localization of nanoscale features across different measurement and metrology platforms. For further information about the different spectroscopies available on such a platform, please visit www.microspectroscopy.com and download the application handbook.

* Editorial note: This content is based on HORIBA's investigation at the year of issue unless otherwise stated.

(a)

(b)

(c)

Figure 5 (a) Schematic illustration of nanoGPS tag attached to the sample facilitates seamless coordinate transfer between two analytical systems. (b) Correlative cathodoluminescence (CL) and (c) PL measurements performed on a Ga₂O₃ nano structure.

References

- [1] Carlton Osburn, Henry Berger, Robert Donovan, Gary Jones, *The Effects of Contamination on Semiconductor Manufacturing Yield*, Journal of the IEST (1988) 31 (2): 45-57.
- [2] N. A. Modine, A. M. Armstrong, M. H. Crawford, and W. W. Chow, *Highly nonlinear defect-induced carrier recombination rates in semiconductors*. *Journal of Applied Physics* 114, 144502 (2013).
- [3] Michael A. Reshchikov, *Measurement and analysis of photoluminescence in GaN*. *Journal of Applied Physics* 129, 121101 (2021)

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Applications Scientist II
Optical Spectroscopy Division (OSD)
HORIBA Scientific

Francis NDI, Ph. D.

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Global Product Line Manager
Optical Spectroscopy Division (OSD)
HORIBA Scientific

In Operando Characterization of Li Composite Battery Electrode by Ellipsometry and Raman Spectroscopy. The Case of $\text{Li}_4\text{Ti}_5\text{O}_{12}$ Based Anode.

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Lithium titanium oxide thin films are increasingly popular anode materials for applications in lithium-ion microbatteries and hybrid supercapacitors, due to their improved safety, cost and cycle lifetime. Nanoengineering of their stoichiometry and crystallinity can modify their lithium-site occupancy and stretch the materials specific capacity. So far, research efforts have mainly focused on the pure spinel phase $\text{Li}_4\text{Ti}_5\text{O}_{12}$ (LTO) and only a small fraction is dedicated to a broader spectrum of titanium-based metal oxide thin films. In this work, Pulsed Laser Deposition (PLD) is used by alternating LTO and Li_2O ablations to create a wide landscape in the composition of titania-based micro-anodes. Their structures are studied by Raman Spectroscopy (RS) and their performances by In Operando Spectroscopic Ellipsometry (SE).

Key words

microbatteries, lithium titanium oxide, pulsed laser deposition, spectroscopic ellipsometry, Raman spectroscopy

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Introduction

The expected growth in the number of sensor nodes and other Internet-of-Things devices in the next decade imposes the development of sustainable micro-power

sources, whose manufacturing should preferably be compatible with Silicon processing methods for economic reasons. Titania-based electrodes are one among other alternative to metal Li matching these conditions, the most studied one being the spinel $\text{Li}_4\text{Ti}_5\text{O}_{12}$. The insertion


Figure 1 Multi-layer deposition in the PLD. A full cycle is completed after the ablation of each target with a fixed number of pulses. Two targets are alternated repeatedly in cycles until the desired film thickness are reached.


Figure 2 Raman spectra of LTO thin films with different amounts of additional Li_2O . Growth conditions are indicated with the deposition pulse ratio $\text{LTO} : \text{Li}_2\text{O}$. Raman references for TiO_2 Rutile and Anatase are from literature.^[1]

of extra lithium into the LTO structure would have a large impact on the accessible capacity of the battery assuming every Ti^{4+} is reduced to Ti^{3+} . Typically, LTO is used in the 1.5-3 V range because discharging to 0 V was reported to be unfavorable for the Li-ion diffusion beyond $\text{Li}_7\text{Ti}_5\text{O}_{12}$. Deeper insight about the structural consequences of cycling such alloys to very low potentials is one target of this work.

Thin film fabrication

The concept of the multi-layer PLD deposition with two different alternating targets is schematically shown in Figure 1.

In this work, the targets were $\text{Li}_4\text{Ti}_5\text{O}_{12}$ and Li_2O . Depositions have been carried out on 1cm^2 Si wafer previously covered with Pt layer to provide electrical accessibility of the anode thin films.


Figure 3 UVISEL (red) and Raman Superhead (green): 1: Polarizer, 2: Modulator, 3: Source, 4: Detector, 5: Laser, 6: Raman optics, 7 and insight: electrochemical cell for in operando studies.

Analyzing tools, HORIBA Raman SuperHead and UVISEL Spectroscopic Ellipsometer

Raman spectra were collected using a HORIBA Raman SuperHead, coupled to an iHR320 monochromator equipped with a Sincerity CCD detector and a laser at $\lambda = 532\text{ nm}$.

Figure 2 illustrates the Raman spectra of structural evolutions of the $\text{Li}_4\text{Ti}_5\text{O}_{12}$ spinel phase upon a staggered lithium addition. It reveals the main presence of anatase and rutile phases and no evidence of spinel LTO for lower lithium contents. With increasing the Li_2O addition, the anatase signal vanishes and some vibrational modes of LTO become evident. Rutile appears to be present at any lithium concentration, which was confirmed by X-ray diffraction.

Literature has reported that nano-scaled TiO_2 shapes have increased the capabilities of lithium uptake for titania-based anodes and promoted their electrochemical


Figure 4 Imaginary part of the permittivity ϵ_2 for LTO thin films for different LTO : Li_2O ablation ratios in (a). Vertical dashed lines indicate significant changes in absorption bands upon lithium addition, marked and color-coded according to X1 – X4. In (b) the evolution of the band gap for different LTO : Li_2O pulse ratio (as inset) are shown.

activity and performance in batteries^[2-5]. While this report is being written, the scanning probe microscope for characterization of such nano objects in operando (HORIBA Electrochemical Tip Enhanced RS) is running first experiments at IREC.

UVISEL phase modulation ellipsometer by HORIBA Scientific was used for optical characterization of the samples. The setup, combined with Raman Superhead is illustrated on Figure 3.

While Raman informs on the structural properties, UV-Visible ellipsometry reveals thickness and optical properties of thin films, the latter correlated to their electrical properties. Figure 4 displays optical properties as imaginary part of dielectric function in a) and bandgap deduced from Tauc plot in b) for LTO with different Li concentrations.

Absorption bands (X1 – X4) occurring in Figure 4 (a) can be attributed to interactions between the electronic states of the Li_{2s} orbital and the hybridized molecule orbitals of O_{2p} and Ti_{3d} forming the TiO_6 octahedra^[6]. The evolution of these bands is highly correlated to Li diffusion modes through the film and can be controlled by in operando ellipsometry.


Let us first state that a detailed ex situ SE study of the optical properties of LTO: Li_2O composites preceded the operando analysis. It allowed composition analysis by means of usual ellipsometric modeling of effective medium. Then, optical properties for layer evolution during cycling was compared to Figure 4, depending on the Li_2O content and optical absorption bands could be assigned to specific Li ions energetic levels (and corresponding occupancy inside the spinel structure) according to literature. These results set a solid ground for further

SE analysis in a dynamic environment.


Electrochemical performance of the films, discharge to very low potential

The electrochemical experiments have been achieved in 1 M LiPF_6 electrolyte and lithium acting as counter and reference electrode in a classic three-electrode setup.

From theory and literature^[7,8], the standard reaction potential for the intercalation of three Li^+ into spinel $\text{Li}_4\text{Ti}_5\text{O}_{12}$ is known to be 1.5 V vs. Li/Li^+ and follows a mechanism expressed in the following equation:

Under continued discharge below 0.6 V vs. Li/Li^+ , two additional lithium ions can be intercalated at tetrahedral 8a sites in $\text{Li}_7\text{Ti}_5\text{O}_{12}$ by reducing two remaining Ti^{4+} as described in the different following equation:

Cyclic voltammograms (CV) and chronopotentiometry measurements have been collected at different current densities in the potential range of 0.3 – 3.5 V vs. Li/Li^+ and are shown in Figure 5 for different LTO : Li_2O films. For the calculation of specific capacities, the thicknesses determined by SE have been implemented.

Constant current charge and discharge (chronopotentiometry) curves in Figure 5 (b, d, f) exhibit a steady increase in the specific capacities, strongly correlated to the addition of Li_2O during deposition. Excellent charge/discharge profiles are reached with well-defined plateaus for the


Figure 5 Cyclic voltammograms (left) and constant current charge/discharge curves (right) of LTO layers without additional Li_2O in (a,b) and with additional Li_2O in a pulse ratio of 2:1 (c,d) and 4:3 (e,f) for LTO : Li_2O , respectively. The electrochemical window ranges between 0.3 – 3.5 V vs. Li/Li^+ . The applied current densities are $13 \mu\text{A cm}^{-2}$, $25 \mu\text{A cm}^{-2}$ and $126 \mu\text{A cm}^{-2}$.

$\text{Ti}^{3+/4+}$ transition (Figure 5 (f)) under further Li_2O addition, resulting in the extraordinary high specific discharge capacity of 298 mAh g^{-1} at 0.5 V vs. Li/Li^+ for the 4:3 ablation ratio. An additional plateau appears to begin below 0.5 V vs. Li/Li^+ , which may be attributed to the extra lithium uptake described earlier by Yi et al.^[9] under the occupation of octahedral 16c and tetrahedral 8a sites.^[10]

To obtain further insights about phenomena taking place at these very low potentials, we performed operando SE during the first cycle on a pristine film with high Li-content (LTO: Li_2O ratio of 4:3).

The first output of our SE analysis is the thickness evolution of the lithiated layer with time for the three different regimes (open symbols in Figure 6a referred to the right y-axis). Remarkably, the thickness evolution of the entire film follows quite accurately the shape of the charge injection for the 3 regimes, which indicates a change in volume upon lithiation/de-lithiation. In the imaginary part of the dielectric function, two distinctive features evolving with time are present in all cases, namely, a lower band at $\sim 1.5 \text{ eV}$, and a higher band at $\sim 4 \text{ eV}$ (Figure 6b).

For regime I, our data are in excellent agreement with


Figure 6 (a) Constant current charge and discharge curves collected under a current density of $3.6 \mu\text{A cm}^{-2}$ (0.17 C-rate) monitored by operando SE. The inset shows the charge-discharge profiles of the operando experiment as a function of capacity. Right-axis (purple) shows the thickness evolution of the entire film with time. (b) Optical absorption as a function of the photon energy corresponding to the three different regimes described in the main text (arrows indicate the evolution with time along the experiment). Insets represent the lithium occupancy in the LTO lattice.

previous experimental results and first principle calculations predicting that during lithiation from $\text{Li}_4\text{Ti}_5\text{O}_{12}$ to $\text{Li}_7\text{Ti}_5\text{O}_{12}$, the intensity of the high energy band should decrease while the low energy band should increase^[11,12].

Beyond $\text{Li}_7\text{Ti}_5\text{O}_{12}$, i.e., regime II, the lithiation process is far less explored and, to the best of our knowledge, there is no experimental data to compare with. The increase in intensity of Dielectric function can be related to the role of oxygen in this regime upon further lithium insertions^[13-18]

Finally, during the charging process (regime III), both energy modes decrease their intensity. In terms of Li occupation, these observations indicate that only the Li inserted in regime II (at the 8a tetrahedral sites) was extracted upon charging, while the Li inserted in regime I (at 16c positions) could not leave the structure.

Overall, our results back up the relevant idea that the cycling range of LTO-based systems can be safely extended until the limit of $\sim 0.2 \text{ V}$, which is beyond the typical value of 1V reported in the literature.

Conclusion

In this work, Pulsed Laser Deposition (PLD) was used by alternating LTO and Li_2O ablations to create a widelandscape in the composition of titania-based

micro-anodes.

The alternating ablation of LTO and Li_2O targets in the PLD provides a good strategy to fabricate highly performing anodes. Stable cycling was achieved down to potentials of 0.2 V. The structural consequences of cycling in the low potential region during the first discharge were studied using in operando SE. These results provide guidelines to extend the useful cycling range of LTO-based anodes to increase their competitiveness versus Li-metal anodes.

Note

This HORIBA readout contribution is a summarized version of the original paper “Safe extended-range cycling of $\text{Li}_4\text{Ti}_5\text{O}_{12}$ -based anodes for ultra-high capacity thin film batteries” published in *Materials Today Energy* 25 (2022) 10 0979.

* Editorial note: This content is based on HORIBA’s investigation at the year of issue unless otherwise stated.

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Catalonia Institute for Energy Research (IREC)

High-resolution Piezoresponse Force Microscopy Imaging of a Few Layered 2D Ferroelectric

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Ana I. PEREZ-JIMENEZ

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Two-dimensional (2D) ferroelectrics are promising materials for future applications in advanced nanoelectronics. Piezoresponse force microscopy (PFM) is extensively used to investigate nanoscale-electromechanical properties. However, the underlying small mechanical deformations displayed by 2D ferroelectrics, ion conductivity, and side effects from tip-surface forces, complicate the high-resolution imaging of ferroelectric domains and domain walls in ambient conditions. Here, we demonstrate that the setup versatility of HORIBA AFM together with its scanning method “PFM-Top mode”, further improve signal sensitivity and spatial resolution when mapping low-dimensional ferroelectrics. This approach is demonstrated by the PFM imaging with up to 12 nm spatial resolution of ferroelectric domains, and domain walls, from a 16 nm thick CuInP_2S_6 flake.

Key words

PFM, thin-film ferroelectrics, high-resolution imaging, CIPS, domain walls

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Introduction

Ferroelectrics are crystals which structural lattice polarization can be reversed under an external electric field, and remain polarized even when the field is removed^[1]. This quality of ferroelectrics has motivated the scientific interest for their potential application in developing future electronic, optoelectronic, and photovoltaic nanodevices^[2,3]. However, the successful integration of classic ferroelectrics like thin-perovskites in nano-sized devices, remains

precluded by problems related to material degradation and ferroelectricity loss, when the thickness of the film shrinks to tens of nm^[4,5]. In this respect, 2D van der Waals (vdW) ferroelectrics offer a means to overcome these downsides. Mainly, 2D-vdW ferroelectrics are nanometer-to-atomic scale thick, structurally stable, and prone to be controlled or modulated^[5].

Up to now, several 2D vdW ferroelectric structures have been theoretically predicted but a few experimentally


Figure 1 PFM operation principle and experimental set up in (a) HORIBA AFM-Raman. (b) Basics of PFM: a conductive tip is brought into contact with the piezoelectric sample. An alternating voltage (Vac) is applied to the tip at a given amplitude and frequency by the function generator (generator 1), while the sample substrate (or back electrode) is kept grounded. The generated oscillating electric field leads to the localized mechanical deformation of the sample (i.e., either extension or contraction). In turns, the resulting strain induces the deflection of the tip cantilever. The reflected laser signal is read out by a position sensitive detector (PSD) that is connected to two lock-in amplifiers. The output signals from lock-in 2 are the PFM amplitude and phase, which contain the information about the strength of the vertical electromechanical response, and the direction of the polarization, respectively. In addition, lock-in 2 is the main amplifier of the feedback-controller loop. While, the output readings from lock-in 1 are the lateral PFM amplitude and phase signals. Schematic drawing of AFM heads (c) HE002 and corresponding optical registration system; and (d) HE001 and its optical registration system. (e) Simplified sketch of the point-contact mode operation of PFM-Top mode.

proven^[3,6]. To unequivocally verify ferroelectricity, the formation of well-defined ferroelectric domain patterns has to be observed^[7]. Ferroelectric domains are nanometer-scale structures with opposite polarization and separated by boundaries known as domain walls. Nevertheless, it is challenging to observe nanoscale domains on materials with intrinsic small-vertical deformations and so, weak ferroelectric responses. This is even more critical when characterizing domain walls. Therefore, to further the experimental discovery, and application of 2D ferroelectrics for advanced nanoelectronics, the electromechanical properties of domain patterns, and domain walls have to be characterized with nanoscale-spatial resolution. PFM is conventionally employed to explore nanoscale electromechanical properties. Since the mechanical deformation (i.e., crystal lattice extension or contraction) induced by the applied electric field is recorded by means of the lock-in technique, strains in the order of picometres can be detected^[1] (working principle in Figure 1b). Despite, when imaging piezoelectric responses at high-spatial resolutions, background noises considerably degrade the image signal-to-noise ratio (SNR). In addition, the tip-sample contact nature of PFM makes piezoresponse readings prone to artifacts generated by long-range electrostatic interactions. The use of stiff cantilevers (tens of N/m) helps to circumvent electrostatic-side

effects. Yet, at the cost of damaging either the tip or the sample. In this manner, the synergic effect of background noises and electrostatics interactions complicate the PFM imaging at high-spatial resolutions of 2D ferroelectrics.

Fortunately, HORIBA AFM offers practical solutions to these issues thanks to its versatile instrumental setup and scanning method PFM-Top mode (Figure 1). These advantages are explained and demonstrated below in the framework of working with HORIBA AFM-Raman (Figure 1a).

In the standard configuration of HORIBA AFM Raman, the Smart SPM (Figure 1a) is coupled to the AFM HE002 head (Figure 1c), which optical detection system is arranged such as the diode laser beam is reflected by the cantilever (at its rest position) to the PSD at an angle of 30°. The aim of this wide angle is to allow the top and side interface of the Smart SPM to the Omegascope (or Raman laser path, Figure 1a). Conversely, the reflected angle of the laser in the registration system of head HE001 is 60° (Figure 1d). With this narrower angle, the signal-collection efficiency (sensitivity) of the registration system is improved a tend fold compared to the sensitivity achieved with the wide angle: HE001's registration system noise <0.03nm vs HE002's registration system noise


Figure 2 Topography and height-line profile (a), amplitude (b) and phase (c) images of a CIPS flake with a thickness of 16 nm. Experimental settings: driving PFM voltage (V_{ac}) = 3 V; driving PFM frequency (f) = 114 KHz (contact resonance); Probe = Pt/Ir coated Si tip, spring constant (k) = 15 N/m; resolution = 256 x 256 pixels; humidity ~ 38%.

<0.1 nm. The superior sensitivity performance of HE001 makes this head suitable to conduct AFM and/or NanoRaman experiments in liquid media. Notably, users can take advantage of this detection scheme to achieve ambient experiments in which resolving very short cantilever deflections is a must. This is the case of high-resolution PFM studies of ultra-thin ferroelectric materials. Because unlike to most of the AFM methods in which the lateral resolution depends on the tip radius, in PFM, this merit is strongly influenced by the material's elasticity^[1].

On the other hand, when scanning via PFM-top mode, the piezoresponse signal is recorded when the tip is in contact with the sample, while the lateral motion of the tip to the next measuring point is performed in semicontact mode (Figure 1e). This operation is repeated on each pixel of the image. The key advantages of this point-contact operation compared to plane-contact mode are the reduction of the tip-sample contact time (in the order of milliseconds) and lateral forces. In turn, short contact times minimize tip and sample damage, loading force drifts and indirectly, electrostatics side effects because stiff cantilevers can safely be used. On the other hand, abrasive lateral forces are avoided because the lateral motion between measuring points is performed in semicontact mode. The absence of lateral forces allows to resolve more precisely the fine contours of the objects.

Explained the instrumental and scanning advantages, in the following sections we probe the capabilities of HORIBA AFM for the high-resolution PFM imaging on ambient conditions of the 2D ferroelectric CuInP_2S_6 (CIPS). CIPS is a vdW ferroelectric with unique properties such as negative piezoelectricity, negative capacitance behavior, and interchange between ferroelectricity and ionic conductivity^[8]. Yet, its reliable nanoscale imaging for fundamental research is scarce due to complications from low ferroelectric signals and ion conductivity effect (see below).

Piezoresponse imaging of a few layered CIPS

This section describes the PFM imaging of a 16 nm thick CIPS flake obtained by coupling the HE002 head to Smart SPM, and by the scanning method PFM-Top mode. The ambient topography, amplitude, and phase PFM images acquired simultaneously during scanning are shown in Figure 2a-c. The experimental conditions are described in the image's caption. The amplitude image (Figure 2b) reveals the existence of domain patterns with labyrinthine-like topology and uniform piezoresponse, but enclosed by domain walls exhibiting low signal strength (Figure 2b, dark-blue stripes). The phase image (Figure 2c) unveils a contrast between domains close to 180°, as expected for domains with opposite polarization. The formation of opposite polarized domain structures is direct evidence of material ferroelectricity and reliable PFM imaging. Notice that the PFM images depict clear signal contrast between domains and domain walls. However, a closer inspection of domain walls requires of higher spatially resolved images.

The PFM amplitude images obtained at different scanning areas are illustrated in Figure 3. At 3 μm^2 image size, the amplitude map depicts good SNR (Figure 3a) and consistent signal. At 2 μm^2 scan area (Figure 3b), the domain walls are well discerned, but the SNR is relatively lower compared to the amplitude map in Figure 3a. A further zoom of the scan area (Figure 3c) lead to the overwhelming of the piezoresponse signal by background noises. Such a poor PFM contrast precludes a robust evaluation of the lateral resolution, since background noises may broaden the width of the domain walls^[1]. A tentative solution to improve the SNR is to increase the driving voltage ($V_{ac} > 3\text{V}$). But this strategy is not recommended for CIPS since the polar structure and ferroelectricity can be severely damaged due to its high ionic conductivity^[9,10]. Besides, high voltages can damage the thin structure of the flake. Thus, to avoid sample damage and the alteration of the ferroelectric properties, we concluded that the sensitivity


Figure 3 PFM amplitude images taken at $3 \times 3 \mu\text{m}$ (a); $2 \times 2 \mu\text{m}$ (b); and $1 \times 0.4 \mu\text{m}$ (c) scan areas from a 16 nm-thick CIPS flake. Experimental settings: $V_{ac} = 3 \text{ V}$; $f = 114 \text{ KHz}$ (contact resonance); Probe= Pt/Ir coated Si tip, $k = 15 \text{ N/m}$; resolution = 256×256 pixels; humidity ~ 38%.

of the present approach was not well-suited for the high-resolution imaging of the sample.

Characterization of domain walls

To improve the image SNR and provide a most robust estimation of the achieved lateral resolution, the Smart SPM was coupled to the HE001 head. PFM Top mode was kept as the scanning method. Figure 4 shows the topography, amplitude, and phase high-resolved images ($0.5 \times 0.5 \mu\text{m}$) acquired from the same 16-nm thick CIPS flake. Notice that the PFM maps (Figure 4b, c) exhibit greater SNR compared to those recorded with the HE002 configuration, although a lower driving voltage ($V_{ac} = 2 \text{ V}$) was applied. Similarly, the domain patterns present clear 180° polarization contrast.

It is customary to estimate the PFM-lateral resolution by mathematically modelling the amplitude signal of a 180° domain boundary (Figure 4c, d). The resulting full width at half maximum (FWHM) is indicative of lateral spatial resolution^[1]. Here, domain wall amplitude profiles from different locations were fit to the Lorentz function (Figure 4e). According to the estimations, the highest lateral resolution is about 12 nm, which corresponds to the measure of the domain wall profile 2 (Figure 4e). It is important to clarify that the estimated lateral resolutions calculated from the domain wall profiles are not their actual width.


Figure 4 High-resolution AFM topography (a), amplitude (b), and phase (c) images of a 16 nm thick CIPS flake. Amplitude and phase profiles of a nearly 180° domain wall (d, black and blue dashed lines, respectively); domain walls amplitude profiles (e, black line) modeled by Lorentz function (e, blue dashed line). Experimental settings: $V_{ac} = 2$ V; $f = 114$ KHz (contact resonance); Probe = Pt/Ir coated Si tip, $k = 15$ N/m; resolution = 256 x 256 pixels; image size = 0.5 x 0.5 μm ; humidity $\sim 38\%$.

Conclusion

In this work we demonstrated the setup versatility of HORIBA AFM together with its scanning method “PFM-Top mode”.

The combined use of the HE001 head and PFM-Top scanning mode permitted the high-resolution piezoresponse imaging in ambient conditions and low driving voltages of a few layered CIPS flake at 12 nm spatial resolution. The reliability of the measurements is supported by phase contrast values between domains close to 180°, and the topography integrity of the sample, even though cantilevers with several N/m were used. Hence, confirming non-destructive operation of PFM-Top mode.

Perspectives

Besides the characterization of ultrathin films and domain walls, the present strategy could be implemented to the investigation, in either ambient or liquid environments, of delicate and weakly piezoelectric materials like polymers and biological systems. Certainly, PFM imaging and fundamental studies will be greatly enriched when colocalized with optical techniques like normal and polarized Raman Microscopy, Second-Harmonic Generation (SHG), and NanoRaman. Definitely, correlative analysis will be ideal approaches to fully exploit the capabilities of HORIBA AFM Raman.

* Editorial note: This content is based on HORIBA's investigation at the year of issue unless otherwise stated.

References

- [1] E. Soergel, *Journal of Physics D: Applied Physics*, 2011, 44, 464003.
- [2] M. Wu, *ACS Nano*, 2021, 15, 9229-9237.
- [3] Z. Guan, H. Hu, X. Shen, P. Xiang, N. Zhong, J. Chu and C. Duan, *Advanced Electronic Materials*, 2020, 6, 1900818.
- [4] R. Khosla and S. K. Sharma, *ACS Applied Electronic Materials*, 2021, 3, 2862-2897.
- [5] J. Chu, Y. Wang, X. Wang, K. Hu, G. Rao, C. Gong, C. Wu, H. Hong, X. Wang, K. Liu, C. Gao and J. Xiong, *Advanced Materials*, 2021, 33, 2004469.
- [6] M. Liu, T. Liao, Z. Sun, Y. Gu and L. Kou, *Physical Chemistry Chemical Physics*, 2021, 23, 21376-21384.
- [7] A. Gruverman, M. Alexe and D. Meier, *Nature Communications*, 2019, 10, 1661.
- [8] D.-D. Xu, R.-R. Ma, A.-P. Fu, Z. Guan, N. Zhong, H. Peng, P.-H. Xiang and C.-G. Duan, *Nature Communications*, 2021, 12, 655.
- [9] N. Balke, S. M. Neumayer, J. A. Brehm, M. A. Susner, B. J. Rodriguez, S. Jesse, S. V. Kalinin, S. T. Pantelides, M. A. McGuire and P. Maksymovych, *ACS Applied Materials & Interfaces*, 2018, 10, 27188-27194.
- [10] A. Belianinov, Q. He, A. Dziaugys, P. Maksymovych, E. Eliseev, A. Borisevich, A. Morozovska, J. Banys, Y. Vysochanskii and S. V. Kalinin, *Nano Letters*, 2015, 15, 3808-3814.

Ana I. PEREZ-JIMENEZ, Ph. D.

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Senior Application Scientist
HORIBA (China) Trading Co., Ltd

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Particle Size Studies for Biotechnology and Life Science
- Application Proposal from HORIBA Analytical Solution Plaza -

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Particle size analyzers are used for R&D and quality control purposes in a wide range of fields, including ceramics, battery materials, catalysts, pigments, cosmetics, foods, and pharmaceuticals. In various industrial processes, the size of particles, including powders, is an important factor in characterizing the functionality of products. The application of particle size analyzers has been remarkably expanding not only in the field of nanotechnology, but also in the fields of biotechnology and life sciences. Particularly in the process of drug development and manufacturing, drug modalities are diversifying, including antibody drugs, nucleic acid drugs, cellular drugs, and gene therapy, as well as small molecule drugs, which have been the mainstream in the past and there is a growing need to evaluate cell-derived products such as exosomes, liposomes and drug delivery systems (DDS). Biopharmaceuticals have more complex molecular structures and must be analyzed and evaluated while considering the degree of damage to the sample. Therefore, conventional analytical principles and methods may not be applicable and different analytical approaches need to be considered. For small molecule pharmaceuticals as well, there are increasing demands for evaluation using different scales, rapid testing methods, high throughput and it is essential to propose analyzers and applications that satisfy these demands.

The HORIBA Group possesses multiple particle size analysis technologies that are suited to different measurement targets and offers a wide range of solutions.

This paper focuses on dynamic light scattering technology, which has been gathering attention in recent years among particle size analysis technologies and introduces its applications in the fields of biotechnology and life science.

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
Figure 1 Analytical cases for small molecules and biopharmaceuticals.

Figure 2 nano Partica SZ-100V2.

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
Figure 3 Optical system of nano partica SZ-100V2.


Figure 4 Gel unit.

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

Figure 5 ViewSizer3000.

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Figure 6 Instrument configuration of ViewSizer3000.

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
Figure 7 Exosome whole particle image by ViewSizer3000.


Figure 8 Comparison of ViewSizer3000 results between exosome whole particles and fluorescently labeled particles.
 — : Exosome whole particles
 — : Fluorescence labeled particles

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
Figure 11 Particle size distribution depending on surface drug content.
— : 0.0% — : 0.5% — : 1.0% — : 5.0% — : 10.0%


Figure 12 Samples with various cellulose concentrations.

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
Figure 13 Mesh size distribution at various concentrations of CNF.

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Introduction of Non-destructive Elemental Analysis for In/Online Through Real Case Studies

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Automation and labor saving in the manufacturing process of industrial products have been adopted for the purpose of lower cost and higher quality. Inline analysis and online analysis are required actually to reduce labor costs and to perform full inspections, which used to be limited to sampling inspections. X-ray fluorescence analysis is a non-destructive and a non-contact method of elemental analysis, which was used mainly in laboratories. Now, it becomes to be applied to inline film thickness measurement in roll-to-roll processes, and online measurement of elements contained in oil, plating solution, wastewater, etc. In this report, examples of actual in-line and on-line measurements are presented.

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HORIBA Techno Service Introduce “Analytical Solution Plaza”

HORIBA's Services that Connect “Analysis” and “Understanding (Finding out)”

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HORIBA Techno Service Analytical Solutions Plaza has areas for each of HORIBA's focus markets such as “Energy & Environment”, “Advanced Materials & Semiconductors”, and “Biotechnology & Healthcare”. As a center connecting 18 application laboratories in Japan and overseas, HORIBA's core analytical technologies are concentrated here.

HORIBA promotes application proposals, contract analysis, and joint development tailored to customer needs. As a solution partner to our customers, we aim to provide high value-added services.

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Analytical & Testing Technology Department
Analytical & Testing Technology Department
Analytical Technology Division
HORIBA Techno Service Co., Ltd.

Facilities Introduction

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Introduction to HORIBA FRANCE Applications Laboratory: Missions and Resources

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Alessia QUATELA

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HORIBA FRANCE SAS Applications Center has its main analysis laboratory located in Palaiseau, Île de France. The principal role of this laboratory is to exhibit and demonstrate the HORIBA's scientific instruments to further increase the level of customer focus and application support. The laboratory offers application support and full access to our knowledge base and instruments. The laboratory is staffed by an experienced applications team of MSc. and Ph.D. scientists with a variety of backgrounds including Analytical Chemistry, Pharma, Semiconductors, etc.

Welcome, your applications laboratory visit starts here.

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Our Missions:

- ✓ Samples analysis for potential customers
- ✓ Paid analysis for all the systems
- ✓ Live or remote demonstrations of instruments' operations
- ✓ Support to existing customers by phone, e-mail or face-to-face meeting
- ✓ After-purchase training on customer site, scheduled training on HORIBA FRANCE site or by remote
- ✓ Development of new analytical methods and procedures for and with our customers
- ✓ Participation at scientific collaborations
- ✓ Participation at conferences and workshops with talks and posters
- ✓ Validation of new products (accessories, software...)
- ✓ Writing technical (Application Notes and Peer Reviewed Articles), commercial and marketing documents
- ✓ Customers' visits for application support with Sales Forces/Product Specialists/Product Managers

Material Resources:

The 250 m² laboratory hosts different techniques, from elemental analysis to molecular analysis, surface analysis and particle characterization: Atomic force microscopy (AFM), micro and macro-Raman, AFM-Raman, Fluorescence, C/S analyzer, X-Ray Fluorescence (XRF), Glow Discharge Optical Emission Spectroscopy (GD-OES), Inductive Coupled Plasma Optical Emission Spectrometry (ICP-OES), Laser diffraction (Particle Characterization Analysis-PCA), Spectroscopic ellipsometry, Surface Plasmon Resonance imaging (SPRi), etc. (Figure 1).

Here the potential customers have the possibility of testing the systems before the purchase, being trained by our experts (trainings are taught in English and/or in French) on the latest hardware and software releases, analyze their samples. The laboratory has also strong relationships with industry and academia and works in many collaborative projects:

Figure 1 HORIBA FRANCE Applications Laboratory.


Figure 2 Chemiluminescence spectrum of an investigated compound.

- local collaborations with research centers located on the Plateau de Saclay, a new scientific and technological cluster inspired by the Silicon Valley
- regional and national collaborations with laboratories from various communities (semiconductor, H₂, pharma, environment etc.)
- international collaborations as, just to mention one, the Harverstore project which received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 824072.

As an example of this scientific exchange, Dr. Y. Janin from Pasteur Institute contacted the Horiba application lab as he was looking for a system sensitive enough to detect the chemiluminescence/bioluminescence from luciferins molecules in the frame of the LuLISA project (*bioluminescence as a tool for human diagnostics, from allergy to Covid19*).

This research project aims developing diagnostic tools (i.e., serological tests) for multiple infectious diseases at local, regional or nation scales.

The use of the Aqualog was crucial to detect their weak signal and characterize this kind of molecules (Figure 2).

Two articles are the results of this fruitful collaboration:

- Gagnot G, Hervin V, Coutant EP, Goyard S, Jacob Y, Rose T, Hibti FE, Quatela A, Janin YL., *Core-Modified Coelenterazine Luciferin Analogues: Synthesis and Chemiluminescence Properties*, Chemistry 2021
- Gagnot G, Legrand P, Tadros A, Ezzhara-Hibti F, Quatela A, Janin YL., *On pyridopyrazinol chemistry, synthesis of chemiluminescent substances*, Synthesis 2021 Feb

Human Resources:

Technical support is provided by Ph. D.s, MSc.s, and Master students. Senior scientists and market/technical experts from other departments supervise and support the application team and the customers' projects (Table 1).

Table 1 Experts' background summary


Technique	Background
AFM/Raman	Biophysics, Advanced Materials, Material Science, Pharma, Electronics/Semiconductors, Nano-optics, Nano-Science
C/S analyzer-GD	Material Science
Ellipsometry	Material Science
Fluorescence	Material Science, Biophysics, Cosmetics
ICP	Analytical Chemistry
PCA	Physics, AgroFood, Material Science
PP-TOFMS	Analytical Chemistry and Material Science
SPRi	Biology, Agro-food
XRF	Material Science

Applications domains:

This wide background is daily employed to support our applications demand and the main applications trends (Figure 3).

The combination of our multidisciplinary experience and technical offer permits to propose to the customers a complete solution building trust, creating added value, and ensuring their satisfaction.

*Editorial note: This content is based on HORIBA's investigation at the year of issue unless otherwise stated.

Advanced Materials

- Energy
- Semiconductor
- Healthcare
- Quantum

Water

- Treatment & Recycling
- Quality control
- Drinking water
- μ -plastics

Life Sciences

- Agrofood & Botany
- Cosmetics
- Pharma
- Biology
- Medical

Energy

- PV & Li battery
- Fuel cell
- Biofuel
- Hydrogen
- Flow & μ -battery

Figure 3 Applications Mega Trends.

Alessia QUATELA, Ph. D.

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Head of the Applications Laboratory
HORIBA FRANCE SAS

Facilities Introduction

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HORIBA Technical Centers at Premier Institutes in India

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Dr. Rajeev GAUTAM

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Dr. H.C. SUDEEKSHA

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High-end instruments are indispensable tools in modern-day, cutting-edge research. While they are expected to add high value to the advanced research work, the understanding of potential of these sophisticated instruments for multiple applications plays a vital role in their wider usage. In a growing market like India, creating an effective awareness of such advanced equipment is the need of the hour to enhance HORIBA sales. In this direction, HORIBA India (HIN) Scientific had analysed the market situation to understand the existing distributor-driven sales picture and had realized that it was essential to create a suitable HORIBA-Academics-Industry platform to facilitate regular promotional activities to drive sales of all scientific products to both academic and industrial markets.

Key words

research lab, analytical facility, HORIBA-IISc, HORIBA-IIT Delhi

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研究所, 分析機器, HORIBA-IISc連携, HORIBA-IIT Delhi連携

Introduction

The Indian analytical market analysis and customer voice indicated us a different way of branding by partnering with a few premier institutes in India that had their brand significance in both academic and industrial domains. The Indian Institute of Science (IISc) located in Bangalore emerged as our preferred institute because of various aspects of branding it would bring to HORIBA: IISc is the top-most influential academy of India that was initially started by Jamsetji TATA before it was handed over to the central government of India for making it a highly funded, rank-1 institute of India that attained global importance.

Prof. C.V. Raman served as the first Indian director of this institute. With a timely support and guidance of Prof. Siva Umopathy, the well-known Raman spectroscopist from this institute, HORIBA India (HIN) could sign an MoU with IISc for setting up the HORIBA-IISc Technical Center (HITC) in 2018. The HITC was inaugurated by Dr. Jai Hakhu, Prof. S Umopathy, Prof. S Ramakrishnan, Dr. Rajeev Gautam, and Prof. G Mugesh (from R to L), followed by a grand customer connect (Figure 1).

Notably, it was the first-ever step taken by HORIBA in India to start a technical platform, with its own Application Scientists, for running a customer-outreach program to

strengthen the direct marketing in India and SAARC (HIN secured a global HORIBA BlackJack award for this idea in 2019).

Figure 1 Indian Institute of Science, India's topmost research institute and a glimpse of the inaugural ceremony of the HITC.

Figure 2 A view of the present HITC.

HITC strengthens the Industry-Academy relationship:

HITC is a state-of-the-art applications lab comprising LabRAM HR Evolution, Duetta, DeltaPro, and SZ-100 instruments, with the possibility of further expansion with other equipment for demonstration (Figure 2). In phase II, we are planning to showcase Aqualog, AFM-Raman, veterinary products, and Smart-Assay kits.

IISc scholars get free access to this analytical center and a chance to work with HORIBA scientists for their collaborative research. Through HITC, HORIBA gets the advantages of free access to all online literature, library, world-class faculty interaction and guidance. Our collaboration with Ph.D. and postdocs will create future potential customers. We are running various industrial projects to help sales of high-end techniques such as Raman, ViewSizer and Fluorescence. As IISc has several HORIBA customers on its campus, HIN can showcase

various other HORIBA techniques that are procured by IISc to prospective customers. Various industrial visitors at HITC get an opportunity to interact with both HIN Scientists as well as IISc faculty and HORIBA gets good chances for impressing these customers through faculty testimony and recommendation.

Welfare to the HORIBA group:

Since its inception, the HITC has been catering to the needs of potential academic and industrial customers from both India and SAARC countries. The center is used for HIN internal team training, COE interactions, distributor training, and marcom materials creation (application notes/publications). Excited by the successful operation of this much-needed Industry-Academy platform, the present director of IISc has happily allotted a bigger space for HITC in the Chemical Sciences Block that is

ideally located adjacent to the nanoscience center, the physical science block, the biological science block, and the aero science block. This gives us ample opportunities for bringing multiple researchers to the HORIBA community. We already started witnessing the results in terms of various collaborations with the faculty of IISc (who chair various purchase committees in India) to boost cross-segment activities as well.

World-class scientists and researchers from academia and industry frequently visit IISc, so it is much easier for us to invite them to HITC to discuss their active research problems and the new technology or solutions they are looking for. Customer feedback on our products and service received during these interactions will strengthen our technical capabilities and customer engagement programs. Various research scholars are being trained by our Applications Scientists Dr. John Kiran Anthony and Dr. H. C. Sudeeksha to regularly


Figure 3 Raman images of finished tablets showing distribution of API (green) and excipients (colors) (a) and polymorphs of API polymorphs (red and blue) (b).

Expansion and Way forward:

The presence of HITC in Bangalore has attracted customers from most parts of India, yet it is a bit distant for customers from North part of India. Encouraged by the success of HITC Bangalore, we have now signed an MOU with the Indian Institute of Technology (IIT) Delhi which is a prominent and influential technological institute located in Delhi (Figure 5). This center is currently hosting our LA-960V2 supervised by Application Scientist, Dr. Namrata Jain, who is offering hands-on trainings to potential customers and academic researchers. The center will soon be expanded as a technical hub for Scientific, Optical Smart Sensing, and Process and Environment divisions in the new campus of IIT D Sonipat.

HIN has witnessed several advantages of creating technical centers at these premier institutes in India. We plan to actively carry out our branding activities to enhance our

use all available instruments for their research work and to acknowledge HITC in their publications. Multiple projects in Raman, PCA and Fluorescence spectroscopy have been handled for pharmaceuticals, food, and battery industries. All potential customers are regularly invited for a demonstration at HITC, and this has served as a key factor in winning sales against our competitors! The customers see this joint venture between HORIBA and IISc as a confident technical platform, so that they consider HORIBA as their preferred partner. As a result, we see that many Contract Research Organizations (CROs) and 3rd Party Testing Labs have already started purchasing HORIBA equipment for their research and analysis work. This will open a potential market in pharma, biopharma, food, chemicals, metals, and cement in India for HORIBA products.


Figure 4 Polymorphs of API (blue, red, green) in skin cream (a) and distribution of fused graphene with silicon (red) graphene (blue) and silicon (green) on an electrode surface (b).

Figure 5 The MOU signing ceremony for HORIBA-IIT Delhi Technical Center.

market share and to create new domains of the market for the latest technologies of HORIBA.

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Segment Head & Applications Scientist
HORIBA Scientific India

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Corporate Officer
HORIBA, Ltd., Japan
President
HORIBA India

Dr. H.C. SUDEEKSHA

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Applications Scientist
HORIBA Scientific India

Facilities Introduction

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Introduction of Application Centers in HORIBA China

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Jing SHEN

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HORIBA Application centers in China have gone through different phases. All of them were dedicated to meet all requests from the Chinese market. At the same time, they are reference to show HORIBA capabilities on advanced techniques of optics and spectroscopy. In 2022, with the opening of Analytical Solution Plaza (ASP), the new era is coming.

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Introduction

In 2008, HORIBA Scientific settled its first demo system in Beijing. From that day, a new business mode started in China, with seeing is believing. In these early stages, customers were interested to see instruments specifications and reliable results from their own samples.

As, our business and customer base grew with time. In 2011, to better support our growth and our customer base, we expanded our capabilities in Beijing and opened a new application center in Shanghai. While, in Beijing, the application center was focused on demonstrating the performances of our Raman and Fluorescence systems, in Shanghai, our customers had access to most of our HORIBA Scientific products. At the same time, we grow the number of application engineer and their skills to provide regular training classes to our customers and education program to the next generation of customers. This

program based on optical technologies became popular in China under the name of HORIBA Optical School.

To increase our business, addressing industrial needs is required. Nowadays, we see an increasing demand from industry to provide solutions to different issues they are facing. Providing a solution means understanding the customer pains and having the capabilities to solve them in a short time.

Our company made a big strategic investment in China, by building a big facility in Shanghai, the HORIBA C-CUBE. In this new facility, a large area is dedicated to demonstrating all HORIBA technologies and to provide a series of solutions for current and future needs of the China Industry, named, Analytical Solution Plaza (ASP). The open ceremony is planned for this year. The ASP is an open space, designed to show in a glance for visitors, all HORIBA advanced techniques such as elemental

Figure 1 Application Center in Beijing.

Figure 2 Application Center in Shanghai.

analysis, molecular analysis, surface/coating analysis, particle size analysis etc. The ASP is divided into four markets area, Energy, Semiconductor, Environment, and Life Science, which are in line with both the China market trends and company megatrend directions. As example, clean energy is one of the key points to achieve the carbon peak and carbon neutrality. Lithium battery, hydrogen, solar cells are all playing a significant role. In ASP, the team will continue our mission to bring out solutions besides graphitization degree characterization for LiB materials, silicon crystallinity and thickness for solar cell, and so on.

These application solution examples will be shown briefly on digital panels in each relative area. Visitors can play

with the digital panels to understand HORIBA capabilities to solve their issues. The examples will grow with time by building collaborations with Key Opinion Leaders in the field of each application. To meet the local customization needs, the application team in ASP will work closely with Local Engineering team to fix the specifications and to evaluate the prototype. In such a way, customized solutions are provided to customers with China speed. With ASP, a close collaboration between HORIBA segments is expected to expand our offer and support our customers in semiconductors, automotive and environmental customers.

Let's light up a bright future with our ASP in HORIBA C-CUBE.

Figure 3 The entrance of ASP in HORIBA C-CUBE.

HCT application center solution case examples

Micro Fluorescence Solution for Perovskite based Solar Cell

Perovskite solar cells have become highly efficient nowadays. However, a number of challenges remain before they can become a competitive commercial technology, e.g. stability, scale, etc.

The film uniformity is one of the key specifications for perovskite solar cell produced by different processes. A good sample was acquired from one of the top labs for perovskite research in China. We made the investigation with Nanolog with CCD detector and microscope to see if there is any difference from the crystals and thus to check the film uniformity. Weak NIR emission was observed from the samples under microscope.

We made our first try to process the data with Labspec software. It is very interesting to get the distribution of peak position and peak width from the film on bare glass substrate. The peak position is different from different size and shape of the crystals. The FWHM distribution


Figure 4 Peak position distributions and Nanolog.

images further confirm that the perovskite film is composed of different defected or non-uniform crystal forms. Lifetime mapping will be our next goal for this kind of sample.

Solution with Multi Techniques for Archeology

Black-glazed porcelain, especially Jian (Tenmoku) wares, which are famous for their lustrous black glaze that exhibits unique colored patterns, is highly valued in China and other Asian countries.

In collaboration with SIOM, CAS, we had the chance to understand both the elemental and compositional information from black-glazed porcelain using XGT and Raman spectroscopy respectively.


Figure 5 Fe, Al, Si distribution and XGT-9000.

It is well known that the black-glazed porcelain is iron rich. From XGT images, we can see a clear Fe precipitation on the silvery appearance of the snow like patterns.

Besides, it is interesting to know the phase of iron oxide which can well reveal the firing conditions and firing procedure of ancient ceramics.

With Raman analysis, it is found that the snow like patterns are composed of hematite (α - Fe_2O_3) and the dark background are composed of ϵ - Fe_2O_3 wherever it appears in crystal or non-crystal shape under microscope.


Figure 6 α - Fe_2O_3 / ϵ - Fe_2O_3 spectra with LabRam Odyssey.

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Jing SHEN, Ph. D.

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Application Department Manager
HORIBA Scientific
HORIBA China Trading Co., Ltd.

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MAJOR BENEFITS OF THE WIDTH RATIO FILTER

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[Silver Award]

ANALYZER, ANALYSIS METHOD AND ANALYZER PROGRAM

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[Bronze Award]

Development of an innovative coating based on a mixture that significantly improve the grating resistance to High power short pulse laser (LIDT)

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Porosity and size analysis of porous microparticles by centrifugal sedimentation with and without density gradient

a: Yuichi Kato, Takushi Sugino, b: Takahiro Morimoto, Kazufumi Kobashi, Toshiya Okazaki, c: Tetsuji Yamaguchi, Tetsuya Mori

a: Nanomaterials Research Institute, National Institute of Advanced Industrial Science and Technology (AIST), Kansai Center

b: CNT-Application Research Center, AIST

c: HORIBA, Ltd.

Powder Technology Volume 407, July 2022, 117663

Determination of Cabernet Sauvignon wine quality parameters in Chile by Absorbance-Transmission and fluorescence Excitation Emission Matrix (A-TEEM) spectroscopy

a: Alvaro Gonzalez, Jorge Zincker, b: Adam Gilmore, Linxi Chen, ac: Doreen Schober

a: Center for Research and Innovation, Viña Concha y Toro, Ruta k-650 km 10, Penciahue, Chile

b: HORIBA Instruments Inc., 20 Knightsbridge Rd., Piscataway, NJ 08854, USA

c: Emiliana Organic Vineyards, Nueva Tajamar 481, Las Condes, Santiago, Chile

Food Chemistry Volume 392, 30 October 2022, 133101

Overall equipment effectiveness, efficiency and slide review analysis of high-end hematology analyzers

a: Sukesh Chandran Nair Dr, Pandiyan Murugan, Asady Sukanya Sukumar, Joy J. Mammen, Saravanan Mullai

b: Shubham Rastogi

a: Department of Transfusion Medicine & Immunohematology, Christian Medical College, Vellore, Tamil Nadu, India

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Practical Laboratory Medicine Volume 30, May 2022, e00275

Safe extended-range cycling of $\text{Li}_4\text{Ti}_5\text{O}_{12}$ -based anodes for ultra-high capacity thin-film batteries

a: Valerie Siller, Juan Carlos Gonzalez-Rosillo, Marc Nuñez Eroles, Alex Morata, b: Michel Stchakovsky, cde: Raul Arenal

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e: Fundacion ARAID

f: Catalan Institution for Research and Advanced Studies (ICREA)

Materials Today Energy Volume 25, April 2022, 10097

Perspectives for regulating 10 nm particle number emissions based on novel measurement methodologies

a: Z. Samaras, L. N.tziachristos, A.Kontses, Z.Toumasatos, A.Tsakis, b: M.Rieker, P.Kreutziger, c: E.Papaioannou, D.Zarvalis, P.Baltzopoulou, N.D.Vlachos, L.Chasapidis, D.Deloglou, E.Daskalos, d: W.F.van Dorp, e: M.Kousoulidou, f: J.Andersson, g: A. Bergmann, M.Bainschab, h: S.Hausberger, L.Landl, i: J.Keskinen, P.Karjalainen, S.Martikainen, j: A.Mamakos, k: Ch Haisch, l: O.Piacenza, G.Nicol, m: A.N.Bhave, K.F.Lee, no: J.Akroyd, nop: M.Kraft, q: M.Kazemimanesh, A.M.Boies, r: C.Focsa, D. Duca, Y. Carpentier, C.Pirim, s: J.A.Nobler, O.Lancry, S.Legendre, t: T.Tritscher, J.Spielvogel, H.G.Horn, u: A.Pérez, S.Paz, v: T.Lähde, A.Melas, cw: A.G.Konstandopoulos, x: S.Zinola, y: S.Di Iorio, F.Catapano, B.M.Vaglieco, z: H.Burtscher, aa: D.Zamora, ab: M.Maggiore

a: Aristotle University, Lab of Applied Thermodynamics

b: HORIBA Europe GmbH

c: Centre for Research and Technology-Hellas (CERTH)

d: Uniresearch B.V.

e: European Climate, Infrastructure and Environment Executive Agency

f: Ricardo Automotive & Industrial, Shoreham Technical Centre

g: University of Technology Graz, Institute of Electronic Sensor Systems

h: University of Technology Graz, Institute for Internal Combustion Engines and Thermodynamics

i: Tampere University, Aerosol Physics Laboratory

j: AVL List GmbH

k: Technical University of Munich

l: Centro Ricerche Fiat

m: CMCL Innovations

n: Department of Chemical Engineering and Biotechnology, University of Cambridge

o: CARES, Cambridge Centre for Advanced, Research and Education in Singapore

p: School of Chemical and Biomedical Engineering, Nanyang Technological University

q: University of Cambridge, Department of Engineering,

r: Univ. Lille

s: HORIBA France SAS

t: TSI GmbH

u: IDIADA Automotive Technology SA

v: European Commission, Joint Research Centre (JRC)

w: Department of Chemical Engineering, Aristotle University Thessaloniki

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z: Institute for Sensors and Electronics

aa: Mobility Ion Technologies SL (MION)

ab: Directorate General for Research and Innovation

Journal of Aerosol Science Volume 162, May 2022, 105957

Facile fabrication of 2D material multilayers and vdW heterostructures with multimodal microscopy and AFM characterization

ab: Siyan Dong, S. Shiva. P. Nathamgari, Xu Zhang, Jin Wook Hwang, Horacio D. Espinosa, c: Xiang Zhang, Pulickel M.Ajayan, d: Andrey Krayev

a: Department of Mechanical Engineering, Northwestern University

b: Theoretical and Applied Mechanics Program, Northwestern University

c: Department of Materials Science and NanoEngineering, Rice University

d: HORIBA Scientific, Novato

Materials Today Volume 52, January–February 2022, Pages 31-42

A taxonomy of validation strategies to ensure the safe operation of highly automated vehicles

ab: Felix Batsch, a: Stratis Kanarachos, Mike Blundella, b: Madeline Cheah, Roberto Ponticelli

a: Institute for Future Transport and Cities, Coventry University,

b: Horizon Scanning Department, HORIBA MIRA

Journal of Intelligent Transportation Systems Volume 26, Issue 1, 2022, Pages 14-33

Effects of aging time in hydrogen peroxide-glycine-based Cu CMP slurry

a: Juhwan Kim, Donggeon Kwak, c: Jungjae Park, Takayuki Kubota, ab: Taesung Kim

a: School of Mechanical Engineering, Sungkyunkwan University

b: SKKU Advanced Institute of Nanotechnology, Sungkyunkwan University

c: HORIBA STEC KOREA

Materials Science in Semiconductor Processing Volume 140, 15 March 2022, 106343

Multi-task manifold learning for small sample size datasets

a: Hideaki Ishibashi, Tetsuo Furukawa, ab: Kazushi Higa

a: Kyushu Institute of Technology

b: HORIBA Ltd.

Neurocomputing Volume 473, 7 February 2022, Pages 138-157

Evidence for 3-D network of P-centered pyramidal $P(\text{Se}_{1/2})_3$ and quasi-tetrahedral $\text{Se}=\text{P}(\text{Se}_{1/2})_3$ local structures and their 3-membered ring super structure counterparts decoupled from quasi 1D- ethylene-like $\text{P}_2\text{Se}_{2+x}$ ($x = 2, 1, 0$) chains in $\text{P}_x\text{Se}_{100-x}$ glasses

a: Matthew Burger, Aaron Welton, Maxwell McDonald, Ralph Chbeir, P. Boolchand, b: Soumendu Chakravart, c: Badriah S. Almutairi, d: S. Mamedov

a: Department of Electrical Engineering and Computer Science, University of Cincinnati

b: Department of Mechanical and Materials Engineering, University of Cincinnati

c: Department of Physics, College of Sciences, Princess Nourah bint Abdulrahman University

d: HORIBA Instruments Incorporated

Journal of Alloys and Compounds Volume 895, Part 2, 25 February 2022, 162645

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An Interlaboratory Comparison on the Characterization of a Sub-micrometer Polydisperse Particle Dispersion

a: Kurt D.Benkstein, Dean C. Ripple, b: Gurusamy Balakrishnan, c: Ashwinkumar Bhirde, d: Pascal Chalus, Atanas Koulov, Anja Matter, e: Tapan K.Das, Sambit R.Kar, f: Ngoc Do, g: David L.Duewer, h i: Nazar Filonov, j: Fook Chiong Cheong, Laura A. Philips, k: Patrick Garidel, Franziska Schleinzer, l: Nicole S.Gill, Heather Anne Wright, m n: Adam D.Grabarek, o: David G.Grier, Andrew D.Hollingsworth, p: Judith Hadley, q: Wesley W.Howard, Harshit Khansa, Yoen Joo Kim, r: Maciej Jarzębski, s: Wim Jiskoot, Stefan Romeijn, t: Vikram Kestens, Yannic Ramaye, u: Christine Probst, v: Theodore W.Randolph, Jared R.Snell, w: Miguel Saggi, x: Jan “Kuba” Tatarikiewicz, y: Dennis T.Yang

a: Biomolecular Measurement Division, National Institute of Standards and Technology, USA

b: Analytical Development and Attribute Science, Bristol Myers Squibb, USA

c: Office of Biotechnology Products, Office of Pharmaceutical Quality, Center for Drug Evaluation and Research, USA

d: Lonza AG, Drug Product Services, Switzerland

e: Biologics Development, Bristol Myers Squibb, USA

f: Spectradyne LLC, USA

g: Chemical Sciences Division, National Institute of Standards and Technology, USA

h: AlphaNanoTech, USA

i: Particle Metrix, USA

j: Spheryx, Inc., USA

k: Boehringer Ingelheim Pharma GmbH & Co. KG, Germany

l: Yokogawa Fluid Imaging Technologies, USA

m: Coriolis Pharma, Germany

n: Division of BioTherapeutics, Leiden Academic Centre for Drug Research, Leiden University, Netherlands

o: Department of Physics and Center for Soft Matter Research, New York University

p: Malvern Pananalytical, USA

q: Analytical Sciences, BioPharmaceuticals Development, R&D, AstraZeneca, USA

r: Department of Physics and Biophysics, Faculty of Food Science and Nutrition, Poznan University of Life Sciences, Poland

s: Division of BioTherapeutics, Leiden Academic Centre for Drug Research, Leiden University, Netherlands

t: European Commission, Joint Research Centre (JRC), Belgium

u: Luminex Corporation, USA

v: Department of Chemical and Biological Engineering, University of Colorado

w: Pharmaceutical Development, Genentech, USA

x: HORIBA Instruments Incorporated

y: Biopharmaceutical Research and Development, Lilly Research Laboratories, Eli Lilly and Company, USA

Journal of Pharmaceutical Sciences *Journal of Pharmaceutical Sciences* Volume 111, Issue 3, March 2022, Pages 699-709

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c: Société HORIBA Medical

Transfusion Clinique et Biologique Volume 28, Issue 4, Supplement, November 2021, Page S63

Development of Ag and Ag alloys-precipitated Ag₂O-TeO₂ glass and Ag₂O-TeO₂ glass/stainless steel reference electrodes for pH sensors

a: Tadanori Hashimoto, Keisuke Nakade, Atsushi Ishihara, b: Yuji Nishio

a: Division of Chemistry for Materials, Graduate School of Engineering, Mie University

b: HORIBA Advanced Techno, Co., Ltd.

Sensors and Actuators B: Chemical Volume 348, 1 December 2021, 130540

Bipolar self-doping in ultra-wide bandgap spinel ZnGa₂O₄

a: Z.Chi, Wu C.Sartel, I.Madaci, V.Sallet, Y.Dumont, E.Chikoidze, b: Fu-Gow Tarntair, R.H.Horng, c: M.Frégnaux, d: Wan-Yu, e: P.Chapon, f: A.Pérez-Tomás

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d: Department of Materials Science and Engineering, Da-Yeh University

e: HORIBA Jobin Yvon SAS

f: Catalan Institute of Nanoscience and Nanotechnology (ICN2)

Materials Today Physics Volume 20, September 2021, 100466

Yawn-like behavior in captive common bottlenose dolphins (*Tursiops truncatus*)

a: Akiko Enokizu, Tadamichi Morisaka, Motoi Yoshioka, b: Katsushi Murakami, Natsuko Sakurai, c: Nahoko Ueda

a: Graduate School of Bioresources, Mie University

b: Minamichita Beachland Aquarium

c: HORIBA Ltd.

Behavioural Processes Volume 189, August 2021, 104444

Spectrofluorometric analysis combined with machine learning for geographical and varietal authentication, and prediction of phenolic compound concentrations in red wine

a: Ranaweera K.R. Ranaweera, b: Adam M. Gilmore, ac: Dimitra L. Capone, Susan E.P. Bastian, David W. Jeffery

a: Department of Wine Science and Waite Research Institute, The University of Adelaide (UA)

b: HORIBA Instruments Inc.

c: Australian Research Council Training Centre for Innovative Wine Production

Food Chemistry Volume 361, 1 November 2021, 130149

Screening for potential interaction partners with surface plasmon resonance imaging coupled to MALDI mass spectrometry

a: Ulrike Anders, Renato Zenobi, ab: Maya Gulotti-Georgiev, b: Susann Zelger-Paulus, Roland K.O. Sigel, c: Fatima-Ezzahra Hibti, Chiraz Frydman, d: Detlev Suckau

a: Department of Chemistry and Applied Biosciences, ETH Zurich,

b: Department of Chemistry, University of Zurich

c: HORIBA France S.A.S

d: Bruker Daltonics

Analytical Biochemistry Volume 624, 1 July 2021, 114195

Magnetic field-enhanced agglutination as a readout for rapid serologic assays with human plasma

a: Nevzat Temurok, Martine Clot, Aurélien Daynès, b: Fanny Leon, Elena Pinchon, Vincent Foulongne, Jean-François Cantaloube, Philippe Vande Perre, Chantal Fournier-Wirth, Jean-Pierre Molès

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Talanta Volume 233, 1 October 2021, 122407

Time-resolved and excitation-emission matrix luminescence behaviour of boro-silicate glasses doped with Eu^{3+} ions for red luminescent application

a: S.H. Nandyala, A. Stamboulis, b: G. Hungerford, c: J.D. Santos, d: B.M. Walsh, e: L. Di Silvio

a: School of Metallurgy and Materials, University of Birmingham

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d: NASA Langley Research Center, Hampton

e: Faculty of Dentistry, Oral & Craniofacial Sciences, King's College London

Materials Research Bulletin Volume 140, August 2021, 111340

Fluorescence Anisotropy Eems: Adding a New Dimension to the Study of Protein Local Environment

Karen E.Gall, Alex Siemiarczuk

HORIBA Scientific, Piscataway

Biophysical Journal Volume 120, Issue 3, Supplement 1, 12 February 2021, Page 117a

Optimizing Model Calibrations for Natural Product Chemical Compositions with Absorbance-Transmittance Excitation-Emission (A-TEEM) Spectroscopy

Adam M.Gilmore

HORIBA Instruments Inc.

Biophysical Journal Volume 120, Issue 3, Supplement 1, 12 February 2021, Page 264a

Dynamic impact of cellulose and readily biodegradable substrate on oxygen transfer efficiency in sequencing batch reactors

a: Ahmed Shawki Ahmed, b: Ahmed Khalil, c: Yuichi Ito, d: Mark C.M.van Loosdrecht, e: Domenico Santoro, fg: Diego Rosso, ah: George Nakhla

a: Department of Civil and Environmental Engineering, Western University
b: Department of Mechanical and Materials Engineering, Western University
c: HORIBA Advanced Techno,Co. Ltd.
d: Dept. Biotechnology, Delft University of Technology
e: USP Technologies
f: Department of Civil and Environmental Engineering, University of California
g: Water-Energy Nexus Center, University of California
h: Department of Chemical and Biochemical Engineering, Western University
Water Research Volume 190, 15 February 2021, 116724

A label-free, direct solid-phase fluorimetric analysis of ochratoxin A in agricultural products with monoclonal antibody-immobilized monolith

ab: Koji Kikukawa, b: Miki Fukui, Ayano Oka, Ryohei Yamamoto, Atsushi Yamamoto, c: Mikiko Uchigashima, Takeshi Kono, d: Daisuke Kozaki, e: Shuji Kodama

a: Food Analysis Technology Center SUNATEC
b: College of Bioscience and Biotechnology, Chubu University
c: HORIBA, Ltd.
d: Department of Chemistry and Biotechnology, Faculty of Science and Technology, Kochi University
e: School of Science, Tokai University
Food Chemistry Volume 346, 1 June 2021, 128736

Sm³⁺-doped strontium barium borate phosphor for white light emission: Spectroscopic properties and Judd–Ofelt analysis

a: Sunil Thomas, Naser Qamhieh, Saleh T. Mahmoud, bc: Rani George, d: K.G. Gopchandran, e: Alessia Quatela

a: Department of Physics, UAE University
b: Department of Physics, St. Aloysius College
c: School of Pure and Applied Physics, Mahatma Gandhi University
d: Department of Optoelectronics, University of Kerala,
e: HORIBA France SAS
Spectrochimica Acta Part A: Molecular and Biomolecular Spectroscopy Volume 248, 5 March 2021, 119187

Authentication of the geographical origin of Australian Cabernet Sauvignon wines using spectrofluorometric and multi-element analyses with multivariate statistical modelling

a: Ranaweera K.R. Ranaweera, b: Adam M. Gilmore, ac: Dimitra L. Capone, Susan E.P. Bastian, David W. Jeffery

a: Department of Wine and Food Science, and Waite Research Institute, The University of Adelaide (UA)
b: HORIBA Instruments Inc, Piscataway
c: Australian Research Council Training Centre for Innovative Wine Production
Food Chemistry Volume 335, 15 January 2021, 127592

Green production of limonene diepoxide for potential biomedical applications

abc: Anderson J. Bonon, bc: Juliana O. Bahú, Rubens Maciel Filho, b: Bruno C. Klein, d: Dalmo Mandelli

a: HORIBA Scientific, Irvine

b: Laboratory of Optimization, Design and Advanced Control (LOPCA) – University of Campinas

c: INCT-Biofabris – University of Campinas, School of Chemical Engineering

d: Center of Natural and Human Sciences, Federal University of ABC

Catalysis Today Volumes 388–389, 1 April 2022, Pages 288-300

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HORIBA Advanced Technology Center

11-5 Hokotate-cho, Kamitoba, Minami-ku, Kyoto, 601-8116, Japan
Phone : 81-75-693-2300 Fax : 81-75-693-2350

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Phone : (39)-6-51-59-22-1 Fax : (39)-6-51-96-43-34

NETHERLANDS

HORIBA Europe GmbH, Netherlands Branch

Science Park Eindhoven 5080 (Industrial park
"Ekkersrijt") 5692 EA, Son, Netherlands

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RUSSIA

HORIBA OOO

Altufievskoe shosse, 13, building 5, 127106, Moscow,

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Phone : (34)-950-39-11-53

SWEDEN

HORIBA Europe GmbH, Sweden Branch (Gothenburg)

Grimboasen 10 A, S-417 05 Gothenburg, Sweden

Phone : (46)-10-161 1500 Fax : (46)-10-161 1503

HORIBA Europe GmbH, Sweden Branch (Sodertalje)

Sydhamnsvagen 55-57, SE- 15138, Sodertalje, Sweden

Phone : (46)-8-550-80701 Fax : (46)-8-550-80567

TURKEY

HORIBA Europe GmbH, Istanbul Office

Veysel Karani Mahallesi, Colakoglu Sokak No: 10,

Rings Rezidans D:23, PK: 34885, Sancaktepe /

Istanbul, Turkey

Phone : 90 216 572 1166 Fax : 90 216 572 1167

HORIBA World-Wide Network

UNITED KINGDOM

HORIBA UK Limited

Kyoto Close Moulton Park Northampton NN3 6FL
United Kingdom
Phone : 44 (0)1604 542500 Fax : 44 (0)1604 542699

HORIBA UK Finance Limited

Kyoto Close Moulton Park Northampton NN3 6FL
United Kingdom
Phone : 44 (0)1604 542500

HORIBA Jobin Yvon IBH Ltd.

133 Finnieston St. Glasgow G3 8HB, United Kingdom
Phone : (44)-141-229-67-89 Fax : (44)-141-229-67-90

HORIBA Test Automation Limited

Brook Court Whittington Hall Worcester WR5 2RX,
United Kingdom
Phone : (44)-1905-359359 Fax : (44)-1905 359332

Old HORIBA MIRA Certification Limited

Watling Street, Nuneaton, Warwickshire, CV10 0TU,
United Kingdom
Phone : (44)-24-7635-5000

HORIBA MIRA Limited

Watling Street, Nuneaton, Warwickshire, CV10 0TU,
United Kingdom
Phone : (44)-24-7635-5000

HORIBA MIRA Limited, Quatro Park

Unit 1, Quatro Park, Paycocke Road, Basildon,
Essex, SS14 3GH, United Kingdom
Phone : (44)-1268-290100

HORIBA MIRA Certification Limited

Watling Street, Nuneaton, Warwickshire, CV10 0TU,
United Kingdom
Phone : (44)-24-7635-5000

MIRA Land Limited

Watling Street, Nuneaton, Warwickshire, CV10 0TU,
United Kingdom
Phone : (44)-24-7635-5000

MIRA Service Limited

Watling Street, Nuneaton, Warwickshire, CV10 0TU,
United Kingdom
Phone : (44)-24-7635-5000

MIRA Technology Park Limited

Watling Street, Nuneaton, Warwickshire, CV10 0TU,
United Kingdom
Phone : (44)-24-7635-5000

CHINA

HORIBA INSTRUMENTS (SHANGHAI) CO., LTD

No.99, Chunxiu Rd, Anting Town, Jiading District,
Shanghai, China 201804
Phone : (86)-21-6952-2835 Fax : (86)-21-6952-2823

HORIBA Technology (Suzhou) Co., LTD.

No.1 building, Industry park, No.101 Chenmenjing
Rd, Taicang, Jiangsu, China (215400)
Phone : (86)-0512-3306-6388

HORIBA (China) Co., Ltd.

Room 1604, Building 1, No.185 Moyu Road, Anting
Town, Jiading District, Shanghai, China, 201805
Phone : (86)-21-6289-6060 Fax : (86)-21-6289-5553

HORIBA (China) Trading Co., Ltd.

Unit D, 1F, Building A, Synnex International Park,
1068 West Tianshan Road, 200335, Shanghai, China
Phone : (86)-21-6289-6060 Fax : (86)-21-6289-5553

HORIBA (China) Trading Co., Ltd., Beijing Branch

12F, Metropolis Tower, No.2, Haidian Dong 3 Street,
Beijing, 100080, China
Phone : (86)-10-8567-9966 Fax : (86)-10-8567-9066

HORIBA (China) Trading Co., Ltd., Guangzhou Branch

Room 1611/1612, Goldlion Digital Network Center,
138 Tiyy Road East, Guangzhou 510620, China
Phone : (86)-20-3878-1883 Fax : (86)-20-3878-1810

HORIBA Precision Instruments (Beijing) Co., Ltd.

Building 1, No.3 yuan, Xixing Road, Houshayu Town,
Shunyi District, Beijing, 101318 China
Phone : 86-10 8492 9402 Fax : 86-10 8492 7216

MIRA China Ltd.

Suite 501, Block B, Hongqiao Sunnyworld No. 1226
South Shenbin Road Shanghai, 201106, China
Phone : (86)-21-6220-6377 Fax : (86)-21-6220-6379

MIRA China Ltd. Xiangyang Workshop

A27-1, Jiahai industrial park, High-tech District,
Xiangyang, Hubei, 441004, China
Phone : (86)-710-2578-268

INDIA

HORIBA India Private Limited

246, OKHLA INDUSTRIAL ESTATE, PHASE 3 NEW
DELHI - 110020, India
Phone : 91-11-4646-5000 Fax : 91-11-4646-5020

HORIBA India Private Limited, Bangalore Office

No.55, 12th Main, Behind BDA Complex, 6th sector,
HSR Layout, Bangalore South, Bangalore - 560102,
INDIA

Phone : (91)-80-4127-3637

HORIBA India Private Limited, Chennai Office

No.9, 01 & 02 Floor, Ganapathy Colony, Thiru-Vi-Ka
Industrial Estate, Guindy, Chennai, 600032 India
Phone : (91)-44-42077899

HORIBA India Private Limited, Haridwar Factory

Plot No. 26, Sector-7, IIE, SIDCUL, Haridwar,
Uttarakhand - 249403, India
Phone : (91)-1334-239139

HORIBA India Private Limited, Kolkata Office

EK Tower/6th Floor/Office -4A, Action Area-II D,
Newtown, Pin Code 700161, India
Phone : (91)-90073-63356

HORIBA India Private Limited, Technical Center

D-225, Chakan MIDC Phase-II, Bhamboli Village,
Pune-410501, India
Phone : (91) 21-3567-6000

HORIBA India Private Limited, Nagpur Factory

Plot No B-3 & C-32 Butibori Phase 2, Industrial
Area, Dist Nagpur, Nagpur, 441122, India
Phone : (91) 71-0328-0200

INDONESIA

PT HORIBA Indonesia

Jl. Jalur Sutera Blok 20A, No. 16-17, Kel. Kunciran,
Kec. Pinang Tangerang - 15144, Indonesia
Phone : 62-21-3044 8525 Fax : 62-21-3044 8521

KOREA

HORIBA KOREA Ltd.

25, 94-Gil, Ijjik-Ro, Manan-Gu, Anyang-Si, Gyeonggi-
Do, 13901, Korea
Phone : (82)-31-296-7911 Fax : (82)-31-296-7913

HORIBA KOREA Ltd., Ulsan Office

613, Doosan We've the Zenith, 1877, Sinjeong-
Dong, Nam-Gu, Ulsan-Si, 44679, Korea
Phone : (82)-52-275-0122 Fax : (82)-52-276-0136

HORIBA MIRA Limited, Korea Branch Office

25, 94-Gil, Ijjik-Ro, Manan-Gu, Anyang-Si, Gyeonggi-
Do, 13901, Korea
Phone : (82)-70-4689-0680 Fax : (82)-31-296-7913

HORIBA STEC KOREA, Ltd.

98, Digital valley-ro Suji-gu, Yongin-si Gyeonggi-do
16878 Korea
Phone : (82)-31-8025-6500 Fax : (82)-31-8025-6599

PHILIPPINES

HORIBA INSTRUMENTS (SINGAPORE) PTE LTD., MANILA OFFICE

27/F Tower 2 Enterprise Centre, 6766 Ayala Avenue
cor Paseo de Roxas, Makati City, Philippines, 1226
Phone : 63 2 8885 8468

SINGAPORE

HORIBA Instruments (Singapore) Pte Ltd.

3 Changi Business Park Vista #01-01, Singapore
486051
Phone : (65)-6-745-8300 Fax : (65)-6-745-8155

HORIBA Instruments (Singapore) Pte Ltd., West Office

83 Science Park Drive, #02-02A, The Curie, Singapore
118258
Phone : (65)-6-908-9660

TAIWAN

HORIBA Taiwan, Inc.

8F.-8, No.38, Taiyuan St. Zhubei City, Hsinchu
County 30265, Taiwan (R.O.C.)
Phone : (886)-3-5600606 Fax : (886)-3-5600550

HORIBA Taiwan, Inc. Tainan Office

1F., No.117, Chenggong Rd., Shanhu Township,
Tainan County 741, Taiwan(R.O.C)
Phone : (886)6-583-4592 Fax : (886)6-583-2409

THAILAND

HORIBA Holding (Thailand) Limited

46/8 Rungrojthanakul Bld., 1st , 2 nd Floor,
Ratchadapisek Road., Huai Khwang Bangkok 10310,
Thailand
Phone : (66)-2-861-59-95 Fax : (66)-2-861-52-00

HORIBA (Thailand) Limited

46/8 Rungrojthanakul Bld., 1st , 2 nd Floor,
Ratchadapisek Road., Huai Khwang Bangkok 10310,
Thailand
Phone : 66 (0) 2 861 5995 ext.123 Fax : 66 (0) 2 861 5200

HORIBA (Thailand) Limited East Office

46/8 Rungrojthanakul Bld., 1st , 2 nd Floor,
Ratchadapisek Road., Huai Khwang Bangkok
10310, Thailand
Phone : 66 (0) 2 734 4434 Fax : 66 (0) 2 734 4438

VIETNAM

HORIBA Vietnam Company Limited

Lot 3 and 4, 16 Floor, Detech Tower II, No. 107
Nguyen Phong Sac Street, Dich Vong Hau Ward,
Cau Giay District, Hanoi, Vietnam
Phone : (84)-24-3795-8552 Fax : (84)-24-3795-8553

HORIBA Vietnam Company Limited, Ho Chi Minh Branch Office

The Office #01& 02, Level 23, AB Tower, 76A Le
Lai, Ben Thanh Ward, District 1, Ho Chi Minh City
Phone : 84 28 3829 9636 Fax : 84 28 3827 3068

Readout HORIBA Technical Reports August 2022 No.56

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