Product Introduction

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Development of Centrifugal Blood Analyzer “Yumizen Banalyst M120” for Point-of-Care Testing

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μ TAS, POCT(Point-of-Care Testing), HbA1c, CRP

Due to factors such as the increase in the number of patients with diabetes and revision of the medical care act, demand for point-of-care testing (POCT) and for greater accuracy and faster turnaround times in clinics is increasing. In response, we developed a blood analyzer that enables immediate testing at the point of care for analytes such as hemoglobin A1c, CRP, and high-sensitivity CRP, which are used for screening diabetes and infectious diseases, using only a small volume of blood. In this device, a dedicated reagent chip designed by applying μ TAS technology is loaded into the instrument, allowing processes equivalent to those performed by large-scale analyzers—from centrifugation and reagent reactions to optical measurement—to be reproduced within a small chip. This enables highly accurate measurements despite the analyzer’s small footprint. Furthermore, the analyzer reduces the measurement time by up to approximately 2 minutes and 30 seconds (about 33%) compared with the conventional analyzer. This paper describes the development process and key features of the centrifuge-based analyzer.

Keywords

μ TAS, POCT(Point-of-Care Testing), HbA1c, CRP

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
Figure 2 Functional sites of HbA1c reagent chip


Figure 3 Structure of YM-120

Table 1 Measurement time

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

Figure 4 Measurement procedure of YM-120.

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Result at the secondary reference laboratory [NGSP%]

The dotted line shows the $\pm 5\%$ of the result at the secondary reference laboratory.

Figure 5 NGSP Certification Test Results for HbA1c

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Figure 6 GATELINK connection configuration diagram

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Assistant Section Leader
Product Development Dept.
Bio & Healthcare Technology Division
HORIBA, Ltd.