Readout

HORIBA Technical Reports

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August 1996 ■ No.13

The Importance of Particle Size to the Performance of Abrasive Particle in the CMP Process

Michael C. Pohl, Duncan A. Griffiths

(Pages 41-47)

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The Importance of Particle Size to the Performance of Abrasive Particles in the CMP Process

Michael C. Pohl*, Duncan A. Griffiths*

1. Introduction

CMP is an area of technology which is growing at a very rapid pace. As Figure 1¹⁾ shows, it is projected to continue its phenomenal growth over the next five years. In spite of its growth, it is still a relatively new field of investigation²⁾. The mechanism of the process is understood to include oxidation and abrasion, but the details are a bit vague. The notion of applying an abrasive slurry to the wafer was quite radical and was expected to create some problems. One of the problems uncovered was the formation of micro-scratches in the wafer³⁾. From other applications for abrasives this was known to be caused by the presence of larger particles in the slurry⁴⁾. This problem suggested that particle size would be a critical area for study⁵⁾. This idea has been borne out by some recent experiments.


Fig.1 Consensus market forecast for CMP polishers

2. Abrasive materials for CMP slurries

A wide variety of particulate materials have been used in abrasive processes to planarize solid surfaces⁶⁾. These include Al_2O_3 , SiO_2 , SiC_2 , Si_3N_4 , diamonds and a host of other industrially important materials. In an abrasive procedure, the particle size is critical to controlling the rate of removal of the undesirable material. The problem is that at the point of total removal of the unwanted material, the etching must be halted very quickly. This need for etch rate control has led

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to extensive studies of the abrasive particles. Unpublished work has established a strong correlation between micro-scratching and the presence of stray larger particles. This knowledge has led to abrasives being specified by size according to U.S. Government Specification #C16195A. This same sort of classification has been applied to CMP slurries⁷⁾.

3. Particle size distribution analysis of the abrasive materials

A variety of techniques are available to size the abrasive particles used in the CMP slurries⁸⁾. Far and away the fastest growing and most widely used is laser light scattering⁹⁾. Most estimates have sales of these products growing at rates in excess of 15% per year. Its speed, ease of use, ready application to process control and other characteristics suggest its use to study CMP slurries. One can readily envision its use to study slurries immediately before and after contact with the wafer surface. Study of process operating parameters as well as checking incoming raw materials seem very useful applications. The information which could be obtained is very extensive.

The theory of operation of the laser light scattering instruments is very simple¹⁰⁾. Basically, different sized particles scatter light at different angles. If you measure the intensity of the light at a series of different angles, you can calculate the particle size distribution. At this point it becomes interesting because many different algorithms exist to calculate the particle size distribution. As a result, different results can be determined for the same sample, depending on the algorithm selected. Horiba has chosen to use a non-linear iterative regression analysis to solve this problem¹¹⁾. This selection has proven valuable in a number of well documented cases¹²⁾.

4. Measured data

One key area of validation is in the factory QA procedures performed on every instrument¹³⁾. Examples of this type of analysis are shown in Figures 2 and 3. These quartz and poly-styrene samples are considered to be valid standards by the particle sizing community. The close correspondence of the median size, 0.411 μ m for the 0.405 μ m PSL and the 1.080 μ m for the 1.150 μ m quartz suggest the accuracy of this instrument. Several overlays for the same two samples are presented in Figures 4 and 5. The close agreement in both cases document

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the reproducibility of this type of instrument. Both accuracy and reproducibility point to the proper selection of the algorithm.


Fig.2 Analysis of 0.405 μ m PSL standard


Fig.3 Analysis of BCR-66 quartz standard


Fig.4 Multiple runs of 0.405 μ m PSL standard

Data		Memory	Hist.	Cum.	SAMPLE NAME	ID#	Median(μ m)
Main	—	—	—	—	405 nm NIST-traceable standard	951013-145	0.405
M01	—	—	—	—	405 nm NIST-traceable standard	951130-950	0.416
M02	—	—	—	—	405 nm NIST-traceable standard	950830-003	0.411
M03	—	—	—	—	405 nm NIST-traceable standard	951201-004	0.414
M04	—	—	—	—	405 nm NIST-traceable standard	950929-006	0.416


Fig.5 Multiple runs of BCR-66 quartz standard

Data		Memory	Hist.	Cum.	SAMPLE NAME	ID#	Median(μ m)
Main	—	—	—	—	BCR 66 standard	951013-149	1.061
M01	—	—	—	—	BCR 66 standard	951130-958	1.061
M02	—	—	—	—	BCR 66 standard	951114-013	1.082
M03	—	—	—	—	BCR 66 standard	951121-033	1.073
M04	—	—	—	—	BCR 66 standard	951201-012	1.070

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In order to investigate typical applications of SiO₂ slurries, the chemical nature of the SiO₂ must be understood. In Figure 6¹⁴⁾ the procedure for preparing the raw SiO₂ is shown. The process for its production is shown below in Figure 7¹⁴⁾.


Fig.6 Aerosol flame(Schematic)


Fig.7 Production of aerosol(Flow chart)


Fig.8 TEM of aerosol 0 x 50


Fig.9 TEM of aerosol 380

Typical products of this process are shown in Figures 8 and 9. These two photographs reveal the complicated nature of the individual particles of SiO₂. The fundamental particles are truly small (~20nm), but they aggregate in the flame. These aggregates are the basic units present in the slurry which must be controlled. The shape of the aggregates make it a very difficult system to study and understand.

To make matters even worse, one must account for the behavior of the aggregates in solution. This added complexity is shown in Figure 10¹⁴⁾.


Fig.10 Schematic representation of the interaction between aerosol particles in liquids

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The figure clearly shows that the nature of the aggregates depends upon the energy state of the system containing the slurry. High shear energy retains the agglomerated particles as they were produced. Low energy allows the aggregates to agglomerate even further. Thus the performance would be characterized by this so called “network” of particles. This behavior makes performance prediction a very difficult procedure.

The investigation of CMP slurries must begin with the routine slurries which are being produced and used on an everyday basis. Figures 11 and 12 show some typical as received slurry samples.


Fig.11 Fumed silica A sample


Fig.12 Fumed silica B sample

While the median size does not show any significant difference, the volume percentage of particles greater than 1.0 μ m is drastically different. Since scratching by large particles is a major source of concern, these two slurries should behave quite differently. The experiments to validate this hypothesis remain to be performed.

Since the presence or absence of large particles is so clearly detected by this instrument, other applications were examined. A good example is presented in Figures 13,14 and 15.


Fig.13 Slurry prior to flushing


Fig.14 Slurry after one flushing procedure


Fig. 15 Slurry after multiple flushing procedures

This represents a procedure for monitoring the removal of a poorly performing slurry from a CMP delivery system. In Figure 13 the particle size of an unacceptable slurry is presented. Note the large quantity of particles greater than 1.0 μ m. This slurry was removed and the system was cleaned by a standard procedure. In Figures 14 good slurry was added, but analysis revealed that particles greater than 1.0 μ m remained. The entire system was thoroughly cleaned and additional slurry was added.

Figure 15 shows that all of the large particles were now removed. The system was now ready for further use.

5. Summary


These examples show the ability of particle size distribution analyzers to analyze slurry samples. Unfortunately, particle size as a physical test only provides half of the required information. As with any physical test, the most critical part is correlating the size to a performance criteria¹⁵⁾. In the case of CMP slurries it could be etch rate, residual slurry trapped, number of micro-scratches, etc. This is truly where the arduous work begins. One must prepare a large number of slurries and then measure the performance characteristics with them. The final result is a powerful predictive tool to specify performance. Presumably many such tests have been performed, but the proprietary nature of the work has prevented it from being published.

<Reference>

- 1) Michael A. Martinez, "A Year Later, CMP Market Has Grown Even Hotter," Solid State Technology, September 1995, p.44-48.
- 2) Beyer, et al., United States Patent 4,944,836, July 31,1990.
- 3) Geraldine C. Schwartz, "Planarization," paper presented at SemiCon West on July 14, 1995.
- 4) Jim Schlueter, "Identifying a Manufacturable Tungsten CMP Solution," SemiCon Southwest 95 on October 18, 1995.
- 5) Norton Diamond Compounds and Slurries, technical brochure from Norton Company, Abrasives Marketing Group, Worcester, Massachusetts 01606-2698.
- 6) International Energy Agency, "Development and Testing of Procedures for Characterization of Ceramic Powders," prepared by NIST for the DOE, September 1993.
- 7) Cab-O-Sil Fumed Silica Properties and Functions, technical brochure from Cabot Corporation Cab-O-Sil Division, Tuscola, Illinois 61953.
- 8) Terence Allen, Particle Size Measurement, Fourth Edition, Chapman and Hall, New York, 1990.
- 9) Michael C. Pohl, "Selecting a Particle Size Analyzer: Factors to Consider," Powder and Bulk Engineering, February 1990, p.58.
- 10)Tatsuo Igushi, "Laser Light Scattering Measurements," poster session presentation from Pittsburgh Conference in 1992 at Atlanta, Georgia.
- 11)S.Twomey, et al., "An Extension to the Chahine Method of Inverting the Radiative Transfer Equation," Journal of Atmospheric Sciences, July 1977, p.1085
- 12)Michael C. Pohl and Jeffrey Dickinson, "The Power of Particle Sizing Information When Combined With Visual Observation Techniques," technical presentation at Fine Powder Characterization course taught at Clemson University, March 1994.
- 13)Duncan Griffiths, "Quality Assurance Manual for Particle Characterization Analyzer," Horiba Instruments Inc., unpublished document.
- 14)"Technical Bulletin Pigment-Basic Characteristics of Aerosil," technical brochure from Degussa AG-Geschäftsbereich Anorganische Chemieprodukte, Frankfurt, Germany.
- 15)Michael C. Pohl, "Proper Use of Particle Size Distribution Information," technical paper presented at the Fine Particle Society Meeting in 1990 at Las Vegas, Nevada.

Michael C. Pohl, Ph.D.
HORIBA Instruments Inc.
Vice President

Duncan A. Griffiths
HORIBA Instruments Inc.
Product Specialist