# Readout

HORIBA Technical Reports

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### Optoelectronic Analysis of Blood Cells

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#### Abstract

In recent years, marked progresses have taken place in the automation of clinical laboratories and advances in measuring and information processing technologies have helped accelerate the move. However, since blood cell analysis deals with living cells which are prone to ready degeneration, much had to be done before the process could be performed without human intervention and the analytical results be made stable and accurate. This paper discusses recent trends and future courses of automated analysis of blood cells.

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

#### Methylene blue and its oxidation products


#### Thionine series

Methylene blue (CI 52015)	R = Me, R' = Me, R'' = Me, R''' = Me
Azure B (CI 52010)	R = Me, R' = Me, R'' = Me, R''' = H
Azure A (CI 52005)	R = Me, R' = Me, R'' = H, R''' = H
sym-dimethylthionine	R = Me, R' = H, R'' = Me, R''' = H
Azure C (CI 52002)	R = Me, R' = H, R'' = H, R''' = H
Thionine (CI 52000)	R = H, R' = H, R'' = H, R''' = H


#### Thionoline series

Methylene Violet Bernthsen	R = Me, R' = Me
Methyl thionoline	R = Me, R' = H
Trithionine	R = H, R' = H

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Constitutional formula of Methylene Blue and Eosin using for Romanowsky coloring

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WBC	RBC	Plt	Ret
Leukocytosis	Erythrocytosis	Thrombocytosis	Immatures
Leukopenia	Cold agglutinin	Thrombopenia	NRBC
Lymphocytosis	Anemia	Microcytosis	Reticulocytosis
Lymphopenia	Macrocytosis	Schizocyte	Reticulopenia
Neutrophilia	Microcytosis	Small cell	
Neutropenia	Anisocytosis	Plt aggregate	
Eosinophilia	Hypochromia		
Myelemia	Poikilocytosis		
Large immature cell	Pancytopenia		
Atypic. lymphocyto			
Left shift			
Nucleated RBC			
Monocytosis			
Basophilia			
Pancytopenia			

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- 1) 巽典之, 前田宏明, 津田泉, 臼井誠次, "自動血球計数の基礎知識", 厚生社(1991), p.19-44.
- 2) 巽典之. 三輪史朗編, "血液病学", 文光堂(1993), p.91-172.
- 3) National Committee for Clinical Laboratory Standards (NCCLS). "Method for reticulocyte counting NCCLS Document H16-P", NCCLS(1985), p.225-236.
- 4) J. A. Kopke, "Practical Laboratory Hematology", Churchill Livingstone (1991), p. 99-108
- 5) R. M. Schmit, "CRC Handbook Series in Clinical Laboratory Science", CRC Press (1980), p. 63-77.
- 6) National Committee for Clinical Laboratory Standards (NCCLS). "Reticulocyte counting by flow cytometry; proposed guideline NCCLS Document H44-P", NCCLS (1986), p.1-80.
- 7) J. V. Dacie, S. M. Lewis, "Practical Haematology, 7th Ed", Churchill Livingstone (1991), p.37-66.
- 8) B. H. Davis, "Immature reticulocyte fraction (IRF): by any name, a useful clinical parameter of erythropoietic activity", Laboratory Hematology, 2, 2-8 (1996).
- 9) I. Tsuda, N. Tatsumi, "Reticulocytes in human preserved blood as control material for automated reticulocyte counters", Am. J. Clin. Pathol. 94, 109-110 (1990).
- 10) National Committee for Clinical Laboratory Standards (NCCLS). "Romanowsky blood stains NCCLS Document H32-P", NCCLS (1986), p.1-80.
- 11) National Committee for Clinical Laboratory Standards (NCCLS), "Reference leukocyte differential count (proportional) and evaluation of instrumental methods NCCLS Document H20-A", NCCLS (1992), p.1-55.

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