Readout

HORIBA Technical Reports

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April 1998 ■ No.16

Material Analysis using the Emission Spectrometry

Y. LANG, D. ARNIAUD, P. CHAPON, Ch. DEREAD, O. ROGERIEUX

(Pages 45-56)

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Guest Article

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Material Analysis using the Emission Spectrometry

Y. LANG *, D. ARNIAUD *, P. CHAPON *, Ch. DERAED * and O. ROGERIEUX *

(* Instrument SA)

The first direct reading spectrometers were created by Dr HASLER for ALCAN ALUMINIUM. This happened approximatively in 1947. These « direct reading » spark spectrometers have replaced, in most of the applications, the emission photographic plates spectrometers. During nearly 3 decades, the world market has been dominated by 4 manufacturers : ARL (USA), BAIRD (USA) JARREL ASH (USA) and HILGER (U.K.). During this period, the emission spectrometry technique has mostly been dedicated to the spark (metal analysis) and marginally to D.C. arc (powders analysis).

In 1974, some manufacturers have adapted the ICP excitation technique to their simultaneous spectrometer. This powerful technique allows the analysis of liquids.

In 1977, JOBIN YVON created the first sequential computerized ICP instrument and in 1980, the first combined ICP : simultaneous and sequential spectrometer which gave an answer to the question : choice between the cadences of analysis and flexibility in the choice of elements to analyze.

The Emission Spectrometry Department of JOBIN YVON has been created in 1976 as part of JOBIN YVON (created in 1819 !). The Emission Spectrometry Department has known a very high growth (+ 24 % per year) during the 80's. During this period it has developed a complete range of Emission spectrometers dedicated to the elemental analysis of liquids, powders and solids.

At the present time, this range includes :

Spark spectrometers : • JY132 F	for the fast analysis of metals.
GDL spectrometers : • 5000 RF type • 10 000 RF type	dedicated to the thin films analysis in the semi-conductors field, analysis of wafers. dedicated to the thin films analysis for conductors and non-conductors samples. Main application : cars and steel industry.
ICP spectrometers : • JY 124 and JY 238 • PANORAMA and VHR PANORAMA	sequential instruments for analysis of liquid samples. combined instruments : simultaneous with sequential function for analysis of liquid samples.

We will cover the products of the Emission Spectrometry Department in 2 parts :

- analysis of liquid samples by ICP
- analysis of solid samples by spark and glow discharge.

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Analysis of Liquid Samples by ICP

1. Introduction

Inductively coupled plasma (ICP) sources combined with atomic emission spectrometers (AES) are based on the principle that by viewing the appropriate region in an Ar plasma tail flame, the atomic and ionic emission lines of analytes can be measured against very low background emission intensities. The very high temperatures (up to 10,000 K) in the ICP ensure more efficient atomisation than for lower temperature flames such as those used in atomic absorption spectrometry (AAS), and chemical interferences are consequently very small.

ICP-AES exhibit very wide linear response ranges of more than 5 orders of magnitude, making it possible to determine major, minor, trace and ultratrace elements in a single sample preparation. Potentially, up to 70 elements may be measured simultaneously in < 2 min. using < 2 ml of solution. Clearly, such powerful analytical technique has almost limitless applications. as i.e.:

Chemistry - Biochemistry, medicine, forestry and Agriculture - Geology, mining - Environmental protection - Metallurgy and metal production - Cement, glass and ceramics production - Machine and oil condition monitoring - Nuclear reactor engineering and power station - etc...

In this review, we summarise the characteristics of plasma spectrometry.

2. The ICP source

A plasma may be defined as any luminous volume of partially ionised gas. In ICP spectrometry the plasma is generated from radiofrequency (RF) magnetic fields induced by a copper coil which is

wound around the top of a glass torch, as shown schematically in Fig. 1. A sample is generally introduced as a solution, which is first nebulised to form a fine aerosol. The aerosol is transported into the centre of the ICP where it rapidly undergoes desolvation (removal of the solvent from the liquid sample), vaporisation to molecular level and dissociation into atoms, some of which are ionised. Both atoms and ions become excited in the plasma and as they revert to their ground states in the tail flame they emit light (photons). In ICP-AES, their characteristic emission in the tail flame is measured using an optical spectrometer. In ICP-MS, ions are extracted from the plasma into a mass spectrometer for analysis.


Fig. 1 Sample introduction system used in JY ICP-AES instruments

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2.1. Torches

The JY ICP torch consists of two concentric glass tubes fabricated from quartz, and a third one fabricated from alumina. The large injector tube diameter (3 mm) decreases the injection speed and increases the residence time of the sample in the plasma, thereby increasing the probability of emission, and, consequently signal amplitude.

2.2. Plasma formation

A water-cooled copper tube is coiled around the upper part of the plasma torch. The coil is connected to a 40.68 MHz solid state RF generator, which creates an oscillating RF magnetic field within the Ar flowing through the torch. Any small variations in plasma impedance are compensated by auto adjustment

Electron density and excitation temperature decrease with increasing frequency, but the argon continuum decreases with the square of the frequency. The overall result using a 40.68 MHz generator is an improvement in the signal-to-background ratio and, hence better limits of detection (LOD) by approximately a factor of four over 27 MHz.

Typical Detection Limits in $\mu\text{g/l}$ at 3 : As 1.2 / Ca 0.03 / Cd 0.09 / Cr 0.15 / Cu 0.18 / Fe 0.19 / Hg 1.2 / K 0.8 / Na 0.3 / P 1.5 / Pb 1.5 / Sb 1.5 / Sc 0.9 / Se 1.5 / Sn 1.3 / Tl 1.0
High frequency RF couples more effectively with the gas flow allowing lower argon flow rates to be used. The lower flow reduces turbulence in the plasma which decreases noise from the background continuum and improves precision. Plasma gas flow rates, generator power, start-up and shut-down are computer controlled and programmable. Safety interlocks are included for coolant water flow, generator power, argon pressure, exhaust flow and drain trap level.

2.3. Plasma gases

Ar is the preferred plasma gas because: (a) it is inert and therefore does not readily react chemically with samples; (b) it has a high first ionisation energy of 15.75 eV, causing effective ionisation of, and emission by, almost all other elements; (c) it is optically transparent; and (d) it has a

moderately low thermal conductivity, so that heat is retained in the plasma fireball enabling stable operation at moderate power inputs.

Inside the plasma, the "normal analytical zone" (NAZ) is characterised by an incredibly low intensity of background emission and is consequently nearly invisible until a sample is injected. Here, the emission spectrum exhibits a continuum with a relatively small number of discrete lines produced predominantly by Ar atoms. The large transfer optics of the JY spectrometers allow viewing of the entire 6 mm high NAZ at once. By viewing all of the emission zone at one time, optimization of conditions on one element gives excellent results on all elements.

Samples are injected into the centre of the plasma, and so pass through the centre of what is effectively a very high-temperature tube-furnace. There is little inter-mixing of the sample with gases in the plasma fireball. The sample is heated rapidly to ~ 8000 K by conduction, convection and radiation effects, and sample molecules undergo nearly instantaneous desolvation, vaporisation, dissociation, ionisation and excitation. The high temperature of an ICP ensures that dissociation is highly efficient and, unlike flame atomic absorption spectrometry (FAAS) and DCP-AES, chemical interferences (caused by the recombination of atoms and the formation or retention of stable molecular species in the flame) are negligible. Similarly, ionisation interferences are also generally insignificant.

3. Sample introduction

The purpose of the sample introduction system is to produce an aerosol from a sample and introduce it into the ICP where excitation occurs. The basic components consist of a nebulizer, spray chamber and torch all of which play a major role in the quality of the analytical results.

Sample introduction for plasma spectrometry is generally accomplished using solution nebulization. Sample solutions are aspirated by a nebulizer (Fig. 1) which shatters the liquid

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into fine droplets using an Ar gas stream of ~ 1 l/min. These droplets are directed into a spray chamber which removes the unsuitable larger (~ > 10 μm) material, and allows only the finest spray to pass into the plasma. The low injector Ar gas flows required in ICP spectrometry result in the efficiency of this nebulisation and sorting process being very low, typically only ~ 2% of the sample solution aspirated being transferred into the plasma.

No single introduction system is optimal for all types of samples. JOBIN-YVON has designed several systems in the form of individual cassettes permitting rapid interchange. Within a few minutes all hardware can be replaced without realignment. The following cassettes are available for the following applications: (1) Standard aqueous with cyclonic spray chamber and pneumatic nebulizer, (2) High salt content solutions and high solid suspensions, (3) HF containing solutions, (4) Ultrasonic nebulization, (5) Hydride generation and (6) Solid sample introduction by spark ablation.

<Patented Sheath Gas>


Mounted between the spray chamber and the injector, the sheath gas device enhances performance and stability by injecting a computer controlled, axial flow of argon around the aerosol prior to its insertion into the injector.

For samples with high salt content, e.g. 30 % NaCl, this minimizes contact of the sample with the injector wall and eliminates crystallization at the injector tip. The injector's uniquely large bore, 3 mm, also contributes to the ability to handle high solids such as slurries, and serves to improve detection limits by increasing the residence time of atoms in the plasma. An optional argon humidifier provides a moist nebulization gas which minimizes clogging of the nebulizer with high solids.

When alkali elements are to be determined, the sheath gas flow is automatically adjusted to optimize the con-

ditions for these easily ionized elements allowing a 7-10 fold improvement in their LODs.

- (a) A Paschen-Runge polychromator : The additional flat field mounting enables the determination of the alkali metals (K, Li, Na) which have their most sensitive Lines in the 500 -800-nm range.

- (b) A Czerny-Turner monochromator. contains only one (or two inter-changeable) photomultiplier tube(s) combined with a movable diffraction grating.


Fig. 2 Schematic representations of two major types of ICP-AES spectrometers.

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4. Inductively coupled plasma-atomic emission spectrometry (ICP-AES)

Quantification of elemental abundances in samples by ICP-AES is based on the ability to separate a complex emission spectrum into its component wavelengths, with sufficient sensitivity and resolution to precisely measure light intensities at the characteristic wavelength of each analyte. Schematics of typical ICP-AES systems are shown in Fig. 1 and Fig.2.

4.1 Spectrometer design

There are two major categories of ICP-AES spectrometer (Fig.2): polychromator (simultaneous), and monochromator (sequential) systems. Both are capable of measuring the full range of wavelengths used in AES, typically from the upper part of the vacuum ultraviolet (160 nm) to the limit of visible light (780 nm). Special MgF₂ optics and detector permits also the analysis at 134 nm.

Simultaneous spectrometers are direct reading designs (Fig. 2a). Light is focussed onto a fixed diffraction grating which splits the light into its component wavelengths. Exit slits located at pre-determined positions around a Roland circle are used to focus light of specific wavelengths onto individual photomultiplier tubes for measurement. The number of elements which can be determined simultaneously is limited by the physical space required to position photomultiplier tubes for adjacent wavelengths, up to 47 elements to be measured simultaneously. For the determination of elements such as P and S which have their optimum analytical wavelengths in the vacuum ultraviolet (< 200 nm), the spectrometer and its associated optics are purged with N₂. One major disadvantage of polychromators is their high cost, because they require a separate suite of electronics for every analytical wavelength installed. They are also relatively inflexible, since analytical wavelengths are best chosen before the instrument is assembled. These limitations are not as great as might be expected, because in reality there is often only limited choice of analytical wavelengths. A POLYSCAN function permits to scan the spectrum over 4.4 nm in front of each particular

photomultiplier. Each detector becomes then a monochromator and brings the flexibility to the polychromator. The main advantages of polychromator systems is their speed of analysis and low running costs, since they potentially enable the determination of 47 elements in 1-2 min., using only 1 - 2 ml of sample solution. In practise the number of elements which can be quantified is limited by the sample type and the sample preparation procedure employed, and a maximum of 25-30 elements is a more realistic expectation for simultaneous analysis. Finally, the simplicity of polychromators ensures reliable operation and optimum short and long-term reproducibility. Sequential systems are based on computer controlled scanning monochromators (Fig. 2b), which can be pre-programmed to drive rapidly from line-to-line, measuring sequentially at each analytical wavelength chosen. Such systems are inherently cheaper and more flexible than polychromators because they can measure any wavelength required. The main disadvantage of the monochromator is speed, since the analysis time is directly proportional to the number of elements being determined. A 20-30-element programme which requires 1-2 min. per sample on a polychromator will typically take 5-10 min. on a monochromator, depending on the integration times and number of the elements. The volume of sample solution consumed will also be increased accordingly. Combined spectrometers employing both a polychromator and a monochromator provide optimum performance and flexibility but are correspondingly expensive.

<Setting a new standard...>

Figures of merit for ICP systems were described in an article in Applied Spectroscopy, Vol. 49, No10, 1995 (12A-18A). When the JY 238 ULTRACE is compared to these figures, it has a perfect score. It is the only ICP today that can make this claim.

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Table 1 Figure of merit for ICP system

Figures of Merit	Measurements	Best Results	JY ULTRACE
Number of elements	Wavelength range	120-770 nm	120-800 nm
Selectivity	Ba(II) 230 nm line profile	5 pm	< 5 pm
Repeatability	RSD Mg(I) 285 nm line	0.2 % RSD	< 0.2 % RSD
Long-term stability	Warm-up time	15 min	< 10 min
	Long term	% RSD \leq 1 %	< 1 % RSD
Robustness	Mg(II)/Mg(I) ratio	\geq 10	> 10
Limits of detection	SBR Ni 1 ppm	\geq 30	> 30
	RSD background	< 0.3 % RSD	< 0.3 %
Accuracy	Calibration procedure	-Change in number of standards -Rejection of outliers -Possibility of weighting -No blank -Standard addition method	Yes Yes Yes Yes Yes

4.2. Emission spectra, interferences and line selection

ICP-AES is less susceptible to interferences than most other spectrometric techniques, but the coincidence or overlap of emission lines remains a serious analytical problem. Many elements have large numbers of individual lines, Fe alone, for example, exhibiting > 1000 separate atomic and ionic lines between 200 and 300 nm. Compilations of ICP-AES lines and spectra provide vital information on potential spectral overlaps.

ICP emission spectra are also characterised by background continua. The background spectrum is not only due to emission from the plasma itself but also originates from stray light within the spectrometer, recombination radiation and line broadening effects. Background enhancements caused by Al, Mg and Ca are particularly problematic.

Line selection and the identification and correction for spectral interferences is one of the most critical aspects of successful ICP-AES analysis. The use of holographic grating, long focal length and simple optical mounting reduces the level of background, of stray light and enhance the resolution. The actual lines chosen by the analyst will depend on a number of factors including the resolution,

sensitivity and other characteristics of the ICP-AES instrument employed, and the composition of the samples being analysed.

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Analysis of materials by Spark Emission

1. Introduction

Spark emission spectrometry is a rapid, precise and reliable method for elementary analysis of conductive solids.

Spark instruments are used on a systematic basis in the steel, metallurgy and foundry industries for developing, producing and providing quality control of metals and alloys.

Most elements are analysed directly and simultaneously including C, S and N at traces (from a few ppms) to major levels (several percents).

2. Source

A "spark" analysis consists in a succession of individual sparkings (each one with a duration of about 100 μ s) during a short period of time (usually less than 10s). Each sparking melts a small area of the sample to be analysed leading to the emission of a composite light.

For each element, a characteristic spectrum is emitted - formed by lines at specific wavelengths which correspond to the transitions between the various levels of energy of the cloud of electrons of this element. The intensity of the emitted lines allows, under well defined excitation conditions, the elements content to be measured.

3. Spectrometer

The optical system of a spark instrument consists of a polychromator which scatters the spectrum and isolates the analytical lines of the elements to be analyzed. In the JY 132F, a 50cm focal length Paschen Runge-type polychromator is used. It consists of an entrance slit, a con-

cave grating and exit slits. The slits, the center and the apex of the grating are located on a circle, known as Rowland's circle.


Fig.1 Principle of the polychromator

As a polychromator's performance basically depends on essential features such as brightness, stability, spectral purity and resolution, let's see how the various components of the system contribute to its performance, starting with the grating.

The diffraction gratings were initially engraved, then holographic. The last generation is that of the so called "blazed" holographic grating, whose mirroring angle is obtained by ionic machining.

Engraved gratings are cut line by line on special machines. Despite the very high precision of the engraving machines, these gratings may incorporate engraving errors which cause phantom spectra (ghosts) and stray light.

For holographic gratings, a polished support covered with photosensitive resin is exposed to the interference of two laser beams and processed chemically. This process allows strictly parallel and equally-spaced lines to be obtained.

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The grating is free of phantom spectra and stray light. However, the brightness achieved by the conventional holographic gratings is not as high as that of engraved gratings, as line profile is sinusoidal.

To achieve gains in terms of brightness, this profile is machined by ion bombardment, to achieve an echelle type profile according to the selected mirroring angle (blaze angle). Such blazed holographic gratings are much brighter than conventional holographic gratings.


Fig.2 Light efficiency holographic grating with ionic etching

For a given grating and in a given order, the resolving power is proportional to the focal length of the optical system and to the reciprocal of the entrance and exit slit width. Let us remind that it is exactly the opposite for the brightness.

To achieve the resolution necessary for analyses (considering that the spark lines are larger than the ICP lines), one can increase either the focal length of the system for a given grating, or the number of lines of the grating for a given focal length.


Fig.3 Line broadening in Spark compared to ICP

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When the system's focal length is increased, its mechanical stability is impaired and, for the same grating, brightness is reduced. Stability is a critical characteristic for equipment which must retain its calibration characteristics over a long period of time. Further, brightness is invaluable as it is synonymous with good detection limits.

Thus, a compact optical system with a 1/2 meter focal length offers excellent stability. Resolution is adequate for analyses if the system is equipped with a grating with a high resolution (3600 lines/mm) and high brightness (blazed holographic).

In an emission spectrometer, the exit slits of the polychromator isolate the lines corresponding to the elements to be analyzed. This is the spectrometer optical layout.

Among the many lines of an element, a selection of a few lines is adequate to meet most analytical requirements. This is the reason why the JY132F polychromator is equipped with a metallic mask of about 200 lines, allowing all possible elements in the usual matrices to be dosed at various concentrations.

When the mask is aligned, all its slits are aligned.

Only the slits of the optical programme of the instrument are kept opened but further on site extensions are easy.

It may be advantageous to explore the spectrum on either side of each line, either over a short distance, to check spectral alignment, or over a medium distance, to assess the influence of the emission background, or over a larger distance to use other analytical lines.

This is done by the Polyscan (JY patent) system which consists in the movement of the entrance slit along the Rowland system.

4. Instruments :


Fig.4 Principle of the Polyscan

The JY132F is a laboratory instrument for production control. Equipped with the JY Windows software it guarantees the quality and the traceability of the analysis, essential part of the ISO certification.

The JY Metaltest is a mobile instrument and can be brought on site to sort the scraps or, to identify final products and prevent mix ups before despatching.

Glow Discharge Optical Emission Spectrometry

1. Introduction

For many industrial applications, it is essential to know the chemical composition of the product, both at the surface and in the bulk. Whether the application is coated steels, e.g. for car bodies, or semiconductor wafers, e.g. for integrated circuit manufacture, the surface composition determines many important parameters such as appearance, corrosion resistance, adhesion, conductivity, etc, while the bulk composition is important for the stability and long life of the product.

Only one technique, Glow Discharge Optical Emission Spectrometry (GD-OES) can provide both the surface and bulk composition, quickly, cheaply, and with high sensitivity to all elements (including the gaseous elements), for almost all solid materials, including metals, metal alloy coatings, semiconductors, polymer coatings, powders, glass, ceramics, etc.


Fig.1 Depth profile analysis of a gold coated watch bracelet

As the name suggests, GD-OES combines a glow discharge source (GD) with an optical emission (OE) spectrometer. Glow discharges will be familiar to many readers from experiments at school or university. A glass tube containing two metal electrodes is filled with a low pressure of argon, then a potential is applied across the two electrodes. A plasma is created inside the tube and the tube is seen to glow.

2. Source

The source now used in GD-OES is quite different in detail from the glass tubes used in teaching laboratories. The cathode is replaced by the sample to be analysed while the positive electrode (anode) is a hollow metal tube.

When a potential is applied, argon ions are generated which will be attracted to the sample where they impact with sufficient energy to knock atoms off the surface.

These «sputtered» atoms then enter the plasma where they are excited by collisions. Then they de-excite by optical emission, hence the «glow». De-exciting atoms emit photons with characteristic wavelengths (e.g. Fe at 371.994 nm or H at 121.567 nm). By measuring the signals at these wavelengths, we can then measure the numbers of each type of atom coming from the sample.

The inside diameter of the tube, or hole, is typically 4 mm. The tube is positioned very close to the sample surface, typically 0.1-0.2 mm. This ensures that sputtering of the sample occurs only in the area of the sample opposite the hollow tube.

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When the operating conditions of the source are chosen carefully, for example, at specific values of argon pressure and applied rf power then uniform and rapid sputtering of the sample occurs.

(Nota : The potential applied to the sample is varying at high frequency for both conductive and non-conductive samples to be analysed with one source).

Sputtering rates are typically 1-10 μ m/min, equivalent to about one atomic layer every 6 ms. Since optical signals can be recorded at ms speeds for many minutes, it is possible, in GD-OES, from the one sample, to record signals from the first atomic layers of the surface down to depths of tens of micrometres.

Thus, in this rather complex process, the sample is sputtered


Fig.2 Principle of a glow discharge


Fig.3 Record of a crater shape and depth profile analysis

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and an optical signal is produced which measures the composition of the sputtered material.

3. Spectrometer

The optical signal emitted by the sputtered atoms leaves the source through a window and is focussed onto the entrance slits of one or more spectrometers (polychromator or monochromator).

If the instrument contains only one spectrometer, it is usually a polychromator with a number of fixed channels, one for each characteristic wavelength (element) of interest.

The principles of spectrometers have been described previously for Spark and ICP.


Fig.4 Quantitative depth profiles: Al on Zn on steel, for automotive

4. Instruments

2 instruments are proposed which differ by their optical resolution.

The JY 5000Rf features a 50cm poly and is an ideal instrument for production control. 2 dedicated versions of this instrument coupled with a robot have been delivered to one of the leading semi-conductor manufactures in Japan for the monitoring of wafers.


Fig.5 Example of quantified depth profile of a BPSG

The JY 10 000Rf (with a 1m polychromator) offers an unmatched resolution and is the instrument for complex applications and research.


Fig.6 Glow discharge optical emission spectrometer JY 10 000RF

Y. LANG
JY EMISSION - France
General Manager

D. ARNIAUD
JY Emission Division - USA
Vice president

P. CHAPON
JY Emission France
Spark & GDL Product Manager

Ch. DERAED
JY EMISSION - FRANCE
ICP Product Manager

O. ROGERIEUX
JY EMISSION FRANCE
R&D Manager