# Readout

HORIBA Technical Reports

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September 1998 ■ No.17

---

## Spectroscopic Phase-Modulated Ellipsometry ---Application to thin film metrology---

Ramdane BENFERHAT

(Pages 31-35)

---

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*Feature Article*  
特集論文

## *Spectroscopic Phase-Modulated Ellipsometry* *--- Application to thin film metrology ---*

Ramdane BENFERHAT

( Instruments S.A. JOBIN YVON, HORIBA Group )

### Abstract

Recent advances in thin film technology, especially in semiconductor one, are providing new opportunities and challenges for the development and application of spectroscopic ellipsometry in thin film metrology. Nowadays precise determination of film thicknesses, optical properties and surface morphology are extremely important to produce high quality devices .

The article shows why spectroscopic phase-modulated ellipsometry is extensively used for the characterization of materials (insulators, semiconductors and metals). Main application fields, particularly in semiconductor technology, of this popular optical technique are discussed. Through typical examples accuracy and sensitivity aspects are illustrated.

### 1. Introduction

Thin film materials pervade our everyday life. We are familiar with transparent conductors in LCD watches and computer displays, antireflection coatings, glass coatings for both color and energy efficiency, a large variety of microelectronic and optoelectronic devices, hard coating... Without thin film technology our way of life would not be the same.

Only about 25 years ago the variety of deposition and etching techniques for preparing thin films based devices was relatively limited. Now a certain level of process sophistication and system integration has been achieved and, consequently, more and more precise characterization of the thin films used is required in order to produce quality devices.

The sensitivity of ellipsometry to thin films and surfaces has been well known for over a century. However, it is only

in the last two decades that this non-destructive technique has become an automated, fast, sensitive and highly reliable tool for thin film structures characterization.

### 2. Spectroscopic ellipsometry : Theory and data analysis

Spectroscopic ellipsometry measures the change in the polarization state of probing light, introduced by its interaction (through reflection or, more rarely, transmission) with the sample under study, and as a function of wavelength. This change is usually described by two angles ( ) and ( ), which express the ratio of the complex reflection coefficients  $R_p$  and  $R_s$  of the sample (the indices  $\langle p \rangle$  and  $\langle s \rangle$ , stand for light polarized parallel (p) and perpendicular (s) to the plane of incidence - see **Fig. 1**):

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
Fig.1 Polarization of reflected light

$$= R_p/R_s = \tan \exp(i \quad )$$

The reflection coefficients are determined by the optical properties and composition of the substrate and the overlying layers, by their thicknesses and morphology, and by surface roughness. In order to deduce the unknown parameters of the sample under investigation, a model which optically describes the sample structure is built. This model contains initial estimates of the parameters sought. The parameters are then varied to generate a set of calculated  $E_p^{calc}$  and  $E_s^{calc}$  that fits best the measured data.

Data fitting consists in minimizing the mean-square deviation  $\chi^2$  between calculated and measured ellipsometric parameters ( $\Psi$  and  $\Delta$ ):

$$\chi^2 = \sum_{i=1..N} \sum_{j=1..M} [E_p^{calc}(i, j) - E_p^{exp}(i, j)]^2 + [E_s^{calc}(i, j) - E_s^{exp}(i, j)]^2 / (2N - P - 1) \dots (1)$$

where the sum is taken over all measured wavelengths  $\lambda_i$  and angles of incidence  $\theta_j$ ; P is the number of unknown material parameters and N, the number of experimental data points. The values of the varied parameters obtained after a  $\chi^2$  minimum has been achieved are then taken to be the best statistical estimates of the material parameters describing the sample properties.

### 3. Spectroscopic Phase-Modulated Ellipsometry (SPME): Principle of operation and instrumentation

The optical set-up of SPME in the UV-visible range is displayed in Fig.2. The excitation head consists of a light source, a polarizer and a photoelastic modulator. After reflection on the sample the emerging light beam goes through an analyzer and a monochromator before being detected by the detector. The last three optical elements form the detection head.

The photoelastic modulator consists of a fused silica block glued to a piezoelectric crystal quartz bar oscillating at the frequency of 50 kHz. This generates a periodic phase shift  $\delta(t) = A \sin(\omega t)$  between two orthogonal amplitude components of the transmitted beam.

The detected intensity takes the general form :

$$I(t) = I_0 + I_s \sin(\omega t) + I_c \cos(2\omega t)$$

Where  $I_s$  and  $I_c$  are function of ( $\Psi, \Delta$ ).

In the UVISEL system, the Fourier analysis of the signal is carried out by a Digital Signal Processor. This fast microprocessor system (cycle time < 100 ns) is especially dedicated to fast fourier transform computation. Finally the DC and the first two harmonics of the signal are continuously transmitted to a PC. The intensity ratios then

directly yield the ellipsometric parameters  $\Psi$  and  $\Delta$  characterizing the sample under study.

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
Fig.2 Schematic diagram of the Spectroscopic phase modulated ellipsometer

## 4. Applications

### 4.1 Flat panel display (FPD)


The application of spectroscopic ellipsometry in the field of the flat panel display often involves thin film thickness uniformity and optical properties measurements of a large-area multilayer structures.

#### (1) SiO<sub>2</sub> / Glass

Silica films are currently deposited as initial layers on glass substrates before further film deposits. Therefore, the thickness uniformity of the silica layer is of crucial importance for the quality of the whole multilayer structure. Phase-modulated ellipsometry is particularly well adapted to the characterization of transparent materials (e.g. silica) since it gives a precise absolute measurement of the ellipsometric phase. Thus, a typical thickness mapping of a large-area silica-deposited glass substrate by SPME is characterized by an Å-scale precision (see Fig. 3, 4).

#### (2) a-Si and poly-Si / Glass

Amorphous and polycrystalline silicon (Fig. 5) are the building materials of the actual thin film transistor (TFT) technology. The former is easy to deposit and process while the latter exhibits superior electrical properties (higher carrier mobility). The degree of crystallinity, thickness

Fig.3 n,k of SiO<sub>2</sub>/Glass

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
Fig.5 Pseudo <n,k> POLY-Si/Glass


Fig.6 Pseudo <n> of SnO<sub>2</sub>/Glass

uniformity, surface roughness and oxidation are the clue parameters to be monitored for a high quality TFT fabrication.

Owing to its high sensitivity and the use of advanced optical constant models SPME is capable of determination of the optical and morphological properties of both a-Si and poly-

Si layers . With its help the deposition processes can be better understood and controlled and, consequently, optimized .

(3)Transparent Conducting Oxides (ITO and SnO<sub>2</sub>)  
Transparent conducting oxides (especially indium tin oxide, ITO) are widely used as transparent electrodes since they combine high electrical conductivity with good optical transmissivity. Because of the particular deposition conditions these films are known to exhibit an optical constants gradient and non-negligible surface roughness thus making their optical characterization delicate( Fig.6). SPME inherently possesses the required sensitivity in order to successfully characterize such films. Thus, the use of a graded model of the optical constants of the layers enables the determination of essential material parameters such as index profile, surface thickness and roughness . Moreover, through modeling the of the free-carrier absorption in the films SPME provides valuable information on their electrical resistivity.

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#### (4) Multilayers

An essential step in the FTP fabrication process is the deposition of multilayered thin film coatings. The determination of the thickness and the optical properties of each individual layer represents a classical problem in ellipsometry and is successfully solved by SPME. Multilayer structures comprising a great number of layers, gradients, periodicity and substrate backside reflection are currently used in modeling SPME measurement results. Once the parameters of the structure having been determined at a given (x, y) point of the FPD, a fully automated mapping procedure can be performed over the entire surface of the panel.

#### 4.2 Lithography

Photolithography is one of the basic processes in the production of microelectronic devices. The actual integrated circuits (particularly high-density memories) require smaller and smaller patterns thus pushing the wavelength of the lithography light to the deep UV (DUV). Therefore, the probing light of the characterization optical technique (e.g. ellipsometry) must follow this trend and, as a consequence, large spectral range DUV (down to 190 nm) SPME has become available. The list of materials studied in this « extreme », spectral region includes various photoresists, antireflective coatings, semiconductors (poly-Si), insulators, and metals. SPME is capable of their optical characterization down from DUV up to near IR. The precise determination of the optical constants of these materials using SPME is of great importance for the optimization of the photolithographic process and the microdevice fabrication in general.

### 5. Conclusion

Spectroscopic phase-modulated ellipsometry, being a non-destructive, fast, precise and easy-to-use optical technique is on the way of becoming an invaluable tool for the thin film materials characterization. The application fields cited above are only a part of a whole list of various actual and potential applications of SPME to microelectronics or elsewhere. With the increasing complexity and lowering the size of the samples SPME has an ensured future owing to its high instrumental sensitivity and accuracy, and advanced modeling software capabilities.

#### References :

- D.E Aspnes and A.A Studna, Applied Opt. 14,220 (1975)
- S.N Jasperson and S.E Schatterely, Rev.Sci.Instrument 40,761 (1969)

**Dr. Ramdane BENFERHAT**

Ellipsometry Product Manager  
Instruments S.A. JOBIN YVON, HORIBA Group