Feature Article

Pittcon解説

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Portable X-ray Fluorescence Analyzer MESA-50 with New Developed Optics and New Spectrum Processing

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We have developed a portable XRF which has almost same performance as conventional desktop XRF. RoHS, soil and particulate matter are typical application examples for rapid and accurate analysis with portable XRF. Detection limits of several ppms have been achieved with Pb, As, Cd, Br, Se and Hg in plastics and soil. New developed optics has achieved highly accurate and rapid measurement. In addition, new spectrum processing method has reduced the effect of overlaps between peaks. Detail of new developed optics and spectrum processing method are presented.

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Measurement condition	
Sample	Cu alloy (ZBY921-925)
Voltage	50 kV
Current	Auto
Filter	5 element Filter
Process time	P2 (2 μ s)
Measurement time	100 s

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Table 3 MESA-50の仕様

Model name	MESA-50
X-ray tube	50 kV
Detector	SDD
Measurable elements	Al-U
Signal processor	Digital Pulse Processor
X-ray filter changer	4 position (including open)
X-ray bench (optional)	Atmosphere
Mass and Size	12 kg/202 (W) × 280 (L) × 178 (H)
Radiation area	7 mm, 3 mm, 1.2 mm (Automatic changer)
Sample image	CCD camera
Software	Reporting, Judge OK/NG
Correction algorithm	FPM Thickness, Material, Pb/As peak separation, Residual component estimation
Power supply	Ni-MH battery or AC adapter (100-240 V 50/60 Hz)
Battery life	6 hours (Ni-MH battery)

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