# Feature Article

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Introduction of Oil Measurement Applications using Oil Content Analyzer

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Oil is a commonly used material which is essential for mankind's activity. Thus, it is necessary to measure and control oil content for its effective use. HORIBA's OCMA-500 series and OCMA-550 series have achieved easy and accurate measurements by extracting oil to a solvent. In this article, we introduce some applications of oil extraction and measurement using OCMA-550 series, discussing the significance of each case. Specifically four difficult samples that require special technique for the extraction, which are oil in soil, residual oil on metal parts after being cut, residual oil on surgical instruments and oil in rice (rice balls from a store), are shown with detailed explanation regarding the setting of extraction time and removal of detergent-like substances.

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Table 1

Model	OCMA-500	OCMA-505	OCMA-550	OCMA-555
Solvent	S-316	H-997	S-316	H-997
Extraction	Automatic	Automatic	Manual operation	Manual operation
Country of sale	Overseas	Mainly Japan*	Overseas	Mainly Japan*

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Table 2 The example of the oil measurement technique in soil.

	Measurement time	Advantage	Fault
Hydrogen flame ionization gas chromatograph (GC-FID)	30 min.	Distinction of the grade of crude by a carbon number is possible.	Knowledge is required to read data. A hydrogen cylinder is required.
Extractive substance in normal-hexane	2 hours~	Regulating method (JIS)	Detection of volatile oil is impossible
Infrared absorption spectrophotometry	3 min.~	Operation is simple.	Distinction of a grade of crude is impossible.


Figure 2 Operation flow of measurement of the oil in soil

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
Figure 3 Example of extraction of the oil in soil.

Table 3 The measurement result of the oil extracted from soil.

Extraction object	S-316	H-997	
Soil extract : Measured value (mg/L)	17.7	16.4	A
The amount of oil in the soil 0.5 g : Calculated value (mg)	0.44	0.41	A' = A / (1000/25)
Soil extract which 1.7 mg of B heavy oil is added : Measured value (mg/L)	77.6	69.4	B
The amount of oil in the soil 0.5 g which 1.7 mg of B heavy oil is added : Calculated value (mg)	1.94	1.74	B' = B / (1000/25)
The extraction rate ( $\rho$ ) over B heavy oil addition (%)	88.09	77.94	$\rho = (B' - A') / 1.7$

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
Figure 4 The example of oil extraction of parts (a) and the cut parts (b)

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Table 4 The remains oil measurement result of the parts by OCMA-550 and OCMA-555.

	S-316	H-997
Washing parts (mg/PC)	0.0	0.1
Not washing parts (mg/PC)	2.0	3.0

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
Figure 5 The example of immersion of the edge of a cutter (a) and tweezers (b)

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Table 5 The cleanliness measurement result when OCMA-550 and OCMA-555 are used

	With Washing	With no washing
Tweezers : Extract by S-316(mg/L)	0.0	93.0
The amount of oil adhering to tweezers. (mg)	0.0	9.3
Cutter edge : Extract by H-997(mg/L)	0.0	6.0
The amount of oil adhering to a cutter edge. (mg)	0.0	0.6

Table 6 The ingredient of polished rice.(Per 100g)

Moisture	60.0 g
Protein	2.5 g
Lipid	0.3 g
Carbohydrate	37.1 g
Ash	0.1 g

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
Figure 6 Operation flow of measurement of the oil in rice.

Table 7 The measurement result of the oil extracted from the grain of rice.

The candidate for extraction	S-316	H-997	
Polished-rice extract : Measured value (mg/L)	20.4	22.1	C
The amount of oil in 1 g of the polished rice : Calculated value. (mg)	0.51	0.55	$C' = C / (1000/25)$
rice-ball extract : Measured value (mg/L)	119	144	D
The amount of oil in 1 g of the rice ball : Calculated value. (mg/g)	2.98	3.60	$D' = D / (1000/25)$

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
Figure 7 Example of extraction of the oil from rice.

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