# Selected Article

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Evaluation methods of carbon-nanotube dispersed aqueous solution

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Carbon nanotube is well known nano-material with many superior features as thermo stable, toughness, low gravity, unique electrical properties, and biocompatibility and applied to practical use in many fields. Especially, for electronics application such as transparent conductive film, antistatic coating, touch panel display, and FET (Field Effect Transistor), dispersion technique of CNT in solution is now becoming more important. Aggregation degree, length, semiconductor and metallic features and conductivity of CNT are important properties on its practical application. For its bulk material, elemental analysis and thermogravimetric analysis are used, and for nano-structure observation, TEM and AFM are often used. Otherwise, for dispersed solution of CNT, analysis method used with light probe or spectroscopic methods are more useful. In this report, the analysis method will be introduced through reporting the analysis results of CNT dispersed solution which was prepared by HORIBA.

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Figure 1 Structure of SWCNT Carbon hexagonal direction are expressed by chiral angle  $\theta$  against the tube axis of longitudinal direction expressed as long dashed short dashed lines which correspond to  $\theta = 0$  degree and  $\theta = 30$  degree.

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Figure 3 NIR-PL of SWCNT by Alcohol CVD method (a) EEM map (b) Chirality distribution. A pair of numbers is chiral indices of SWCNT. Data courtesy of Prof. S. Maruyama, the University of Tokyo.

Figure 4 Comparison of Raman spectra on the process of purification SWCNT, which was dispersed in 20% cycloamylose aqueous solution, was purified by centrifugation. F-1 and F-2 were fraction of before and after purification respectively. F-3 was residue fraction after purification. Data courtesy of Prof. S. Kitamura, Osaka Prefecture University.

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Figure 5 Nanoparticle Analyzer nano Partica SZ-100

Figure 6 Time variation of particle size distribution of SWCNT in 0.2% carboxymethyl cellulose aqueous solution by ultrasonication Time values at left-side of figures indicate process time of sonication.

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Figure 7 Particle size distribution of SWCNT/SDBS measured by DLS (Dynamic Light Scattering)method

Figure 8 D and G Raman bands of SWCNT/SDBS He-Ne Laser(633 nm) for excitation.

Figure 9 RBM Raman bands of SWCNT/SDBS He-Ne Laser (633 nm) for excitation. A pair of numbers in parenthesis is chiral indices of SWCNT.

Figure 10 Absorption Spectra of SWCNT/SDBS A pair of numbers in parenthesis is chiral indices of SWCNT

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Figure 11 NIR-PL EEM map of SWCNT/SDBS A pair of number in parenthesis is chiral indices of SWCNT

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Figure 12 Peak shift of NIR PL bands of SWCNT/SDBS NIR-PL bands were assigned as (10,2), (9,4) and (8,6) for each. Excitation wavelength was 720nm. The left figure is the initial spectrum. Every bands shifted about 8nm longer wavelength after 214 days.

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