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Expansion of Determination Scheme for SI-traceable Organic Standard Materials
Using Standard Materials Calibration System

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The various measurement values such as length, mass, amount of substance are bases of economic activity. In order to ensure international consistency of measurement values, it is necessary to assure metrological traceability to international system of units (SI), whereas most of organic standard materials are not assured SI-traceability. One of the reasons is that the supply of the SI-traceable standard materials is not enough to satisfy the requirement of the quantitative analysis. New system, Standard Materials Calibration System, was developed in collaboration with HORIBA STEC Co., Ltd. and National Institute of Advanced Industrial Science and Technology (AIST). In this report, we studied new determination method which was easy and quick to assure SI-traceability to various analyte using this system.

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
Figure 1 Schematic diagram of Standard Materials Calibration System(reactor)

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
Figure 2 The chromatogram of standard solution A


Figure 3 The calibration curve of standard solution A on Standard Materials Calibration System


Figure 6 The calibration curve of standard solution B on Standard Materials Calibration System


Figure 4 The chromatogram of the internal standard materials solution for crude oil analysis (ASTM[®] D5307)

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Table 1 The specific values of the internal standard materials for crude oil analysis (ASTM[®] D5307) by Standard Materials Calibration System

Substance	Mean value mmol/kg	Expanded uncertainty ($k=2$) mmol/kg
Tetradecane	4.10	0.05
Pentadecane	3.82	0.04
Hexadecane	3.60	0.04
Heptadecane	3.36	0.04


Figure 5 The chromatogram of standard solution B

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
Figure 7 Comparison between two measurement results. The unit of the values (mean \pm expanded uncertainty, $k=2$) is mmol/kg

Table 2 The specific values of aromatic hydrocarbons solution determined by Standard Materials Calibration System

Substance	Mean value mmol/kg	Expanded uncertainty (k=2) mmol/kg
Benzene- <i>d</i> ₆	7.61	0.16
Ethylbenzene- <i>d</i> ₁₀	5.68	0.09
Naphthalene- <i>d</i> ₈	2.53	0.05

Table 3 The specific values of fatty acid methyl esters (FAMES) solution determined by Standard Materials Calibration System

Substance	Mean value mmol/kg	Expanded uncertainty (k=2) mmol/kg
Methyl hexanoate	7.59	0.17
Methyl heptanoate	6.72	0.15
Methyl octanoate	6.38	0.14
Methyl nonanoate	5.78	0.13
Methyl decanoate	5.36	0.12

Table 4 Applicable fields and compounds of Standard Materials Calibration System

Applicable fields	Emission gas, Drug, Environment Food, Flavor, Fuel
compounds	Alkanes, Alkenes, Alcohols Esters, Ethers, Ketones, Terpenes Aldehydes (derivatives), Aromatics Fatty acid methyl esters (FAMES)

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