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Finding the Time for Fluorescence
Its Measurement and Applications in Life Science

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We summarise how developments in technology have brought fluorescence lifetime spectroscopy from being the preserve of the specialist to becoming a major tool for research across many science and engineering disciplines. We highlight the advantages which fluorescence lifetime measurements can bring, not only to underpin research, but also through application in helping to solve real-world problems. We illustrate this with recent examples in cancer and Alzheimer's research, which are aimed at improving disease understanding, diagnosis and therapeutics.

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
Figure 1 Model 5000F coaxial nanosecond flashlamp (a) and DeltaDiode picosecond LED/LD (b).

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
Figure 2 Model 5000W SAFE and its 16 channel TCSPC multiplexing module.

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
Figure 3 DeltaFlex dedicated lifetime system (a) and a hybrid FluoroMax-TCSPC lifetime system (b).

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
Figure 4 Fluorescence decay of melanin. In the presence of 4 μM Cu²⁺ ions the decay becomes a single exponential of τ ~ 6 ns at 480 nm emission, suggesting a simplified structure.


Figure 5 Relative amplitude of two decay parameters (p_1/p_2) when a 3-exponential decay model $i(t) = \sum p_n \exp(-t/\tau_n)$, for $n=3$, is used to describes the tyrosine fluorescence decay of Aβ during the early stages of aggregation. F_1/F_0 describes the fluorescence intensity of the sheet-sensing probe ThT during this time and it can be seen to provide little information on the cytotoxic oligomer formation which occurs as a precursor to Alzheimer's disease.


Figure 6 GNR hairpin sensor sprung by binding to c-DNA and the consequential increase in the average fluorescence lifetime of oligonucleotide labelled with 6-carboxyfluorescein as hybridization occurs. Excitation was at 474 nm using the NanoLED laser diode.

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- [1] *A single-photon counting fluorescence decay-time spectrometer*. D J S Birch and R E Imhof. J. Phys. E, **10**, 1044-9, 1977.
- [2] *Coaxial nanosecond flashlamp*. D J S Birch and RE Imhof, Rev. Sci. Instrum. **52**, 1206-12, 1981.
- [3] *A new sub-nanosecond LED at 280 nm: application to protein fluorescence*. C D McGuinness, K Sagoo, D McLoskey and D J S Birch. Meas. Sci. Technol. **15**, L19-22, 2004.
- [4] *Fast time-correlated single-photon counting fluorescence lifetime acquisition using a 100 MHz semiconductor excitation source*. D. McLoskey, D. Campbell, A. Allison and G. Hungerford. Meas. Sci. Technol. **22**, 067001, 2011.
- [5] *Pulse fluorometry using simultaneous acquisition of fluorescence and excitation*. D J S Birch, R E Imhof and A Dutch, Rev. Sci. Instrum. **55**, 1255-64, 1984.
- [6] *A Multiplexed single photon instrument for routine measurement of time-resolved fluorescence anisotropy*. D J S Birch, A S Holmes, J R Gilchrist, R E Imhof, S M Al-Alawi and B Nadolski, J. Phys. E. Sci. Instrum. **20**, 471-3, 1987.
- [7] *Multiplexed array fluorometry*. D J S Birch, A S Holmes, R E Imhof, B Z Nadolski and K Suhling, J. Phys. E.: Sci. Instrum. **21** 415-7, 1988.
- [8] *Multiplexed single-photon counting 1: A time-correlated fluorescence lifetime camera*. D McLoskey, D J S Birch, A Sanderson, K Suhling, E Welch and P J Hicks, Rev. Sci. Instrum. **67**, 2228-37, 1996.
- [9] *Eumelanin Kinetics and Sheet Structure*. J Sutter, T Bidláková, J Karolin and D J S Birch. App. Phys. Letts., **100**, 113701 (4 pages), 2012.
- [10] *Metal ion influence on eumelanin fluorescence and structure*. J-U Sutter and D J S Birch. Methods Appl. Fluoresc. **2**, 024005, 2014.
- [11] *Early detection of amyloid aggregation using intrinsic fluorescence*. O J Rolinski, M Amaro and D J S Birch. Biosensors and Bioelectronics. **25** (10), 2249-52, 2010.
- [12] *Radiative decay engineering: biophysical and biomedical applications*. J. R. Lakowicz, Anal. Biochem. **298**(1), 1-24, 2001.
- [13] *Hairpin DNA-functionalized gold nanorods for mRNA detection in homogenous solution*. W Guoke, Y Jun, W Jinliang, G Peng, D J S Birch and Yu Chen. J. Biomed. Opt. **21** (9), 097001, 2016.

**David J S Birch**

President
HORIBA JOBIN YVON IBH Ltd.
Professor of the University of Strathclyde
Ph.D.

**David McLoskey**

Managing Director
HORIBA JOBIN YVON IBH Ltd.
Ph.D.