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Life Science Analytical Tools for Real Life Problems

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When light interacts with matter (i.e. molecule/sample) many interesting phenomena occur. Light may get absorbed by your sample and depending on the sample properties you will observe light emission, reflection and scattering. HORIBA Scientific builds tools that measure these phenomena allowing our users to gain meaningful information about their samples. This information is used to solve outstanding problems in research or guide scientists to monitor product development in various fields. Life Science, in particular, is an attractive market because, over the years, it has experienced exponential growth, and continues to do so. All the divisions that engage in the scientific study of living organisms including plants, animals and human beings fall under the life science umbrella. This article will highlight the key contributions our products offer in the life science sub-divisions - pharmaceutical and health sciences - and will compare their advantages over existing technologies.

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
Figure 1 A schematic comparison between small molecule (TAXOL) and protein (antibody IMMUNOGLOBULIN) drugs.

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
Figure 2 Raman Analysis of lysozyme (200 mg/mL) in 20 mM Citrate-PBS buffer at pH 4.0 before and after the addition of 20 % ethanol. Sample was excited with a 532 nm laser and a grating of 1800 lines/mm was used on the XploRA, a 200 mm focal length instrument.^[1]


Figure 3 Comparison of the Trp band at 1555 cm^{-1} full width half maximum (FWHM) width (top left) and I_{856}/I_{830} tyrosyl ratios (top right) at concentrations ranging from 10 to 200 mg/mL. A schematic representation of the microenvironment around IgG molecules changes as the concentration increases (bottom (i)-(iii)).^[2]


Figure 4 Monitoring cell proliferation over 4 days using ATEEM. The EEM data were analyzed using PARAFAC analysis. As a result, four spectral components were successfully extracted and assigned to tyrosine, tryptophan, NAD(P)H and Riboflavin.^[3]

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
Figure 5 A schematic representation of the three different configuration of peptides (ligands) used on the microarray chip: linear (L), branched (B) and Tandem (T).^[5] (Reproduced with permission from ref.5, Copyright 2017, Elsevier)

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

Figure 6 SPRi study of alpha lactalbumin (LAC) antibody binding to the microarray peptide pre-functionalized with various orientation of alpha lactalbumin IgE epitopes LAC1_L, LAC1_B, LAC1_T. The SPRi association (A-C; association rates were deduced from 3 min of association) and dissociation curves (D-F, dissociation rates were deduced from 4 h of dissociation) after the injection of 8 different concentrations (0.98-500 nM) of the antibody (colored plots).^[5] (Reproduced with permission from ref.5, Copyright 2017, Elsevier)

Raman

SPRi

A-TEEM

Figure 7

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- [1] Marinella G. Sandros and Fran Adar. Assessing biotherapeutics stability using Raman Spectroscopy. (Application note) <http://www.horiba.com/scientific/products/horiba-life-science-solutions/application-notes/>
- [2] Chikashi Ota, Shintaro Noguchi, Satoru Nagatoishi and Kouhei Tsumoto. Assessment of the Protein-Protein Interactions in a Highly Concentrated Antibody Solution by Using Raman Spectroscopy, *Pharm. Res.*, 33, 956-969(2016).
- [3] Yuichi Kitagawa, Takumi Moriyama, Daisuke Irikura, and Yasushi Nakata. Holistic Analysis of Mammalian Cell Proliferation using Fluorescence Spectroscopy. (Application note) <http://www.horiba.com/scientific/products/horiba-life-science-solutions/application-notes/>
- [4] D.W.Zabriskie and A.E.Humphrey, Estimation of Fermentation Biomass Concentration by Measuring Culture Fluorescence. *Appl. Environ. Microbiol.*, 35, 337-343(1978).
- [5] Alessandro Gori, Marina Cretich, Renzo Vanna, Laura Sola, Paola Gagni, Giulia Bruni, Marta Liprino, Furio Gramatica, Samuele Burastero, Marcella Chiari *Analytica Chimica Acta* 983, 189-197(2017).

Marinella SANDROS, Ph.D.

Business Development Manager for Life Sciences
SPRi Product Manager
HORIBA Instruments Inc.

Fran ADAR, Ph.D.

Worldwide Raman Applications Manager
HORIBA Instruments Inc