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
Figure 1 Daily measurements of UVA(A) and linear relationships between UVA and TOC(B) for corresponding raw(RAW) and finished(FIN) water samples. For panel B the linear equation for the RAW samples was $UVA = 0.196 \times (TOC)$, adjusted $r^2 = 0.975$ and for the FIN samples $UVA = 0.099 \times (TOC)$, adjusted $r^2 = 0.950$. The linear fit predictions (red lines) are shown compared the 95% confidence intervals (green) and prediction intervals (blue).

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
Figure 2 Comparison of typical Aqualog EEMs(top) and UV-VIS absorbance spectra(bottom) for corresponding raw (left) and finished (right) water samples measured on the same day.

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
Figure 3 Excitation-emission contours for PARAFAC model components c1 (A), c2(B) and c3(C). The data set included $n=1484$ samples with duplicate daily measurements for each raw and finished sample. The model accounted for 97.3% of the variance, the split-half validation match was 98.6% and the core consistency was 94%.

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
Figure 4 Normalized concentration PARAFAC loadings for c1-c3 for the daily raw (RAW) and finished (FIN) samples.

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
Figure 5 Comparison of the PARAFAC component concentration ratio c2 : c1 (A) to SUVA (B) for the daily raw (RAW) and finished (FIN) water samples. SUVA was calculated using the data in Figure 1 using the formula in EPA Method 415.3; $SUVA = 100 \times UVA (cm^{-1}) / TOC (mg/L)$.


Figure 6 Box plot analysis of the means, maxima, minima and ranges for SUVA (left) and c2 : c1 (right) for the raw and finished water samples corresponding to the data in Figure 5.

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

Figure 7 Model linear correlation of the TOC calculated using the SUVA and c2 : c1 concentration ratio (A) and the comparison of the model TOC calculations to actual daily TOC measurements for the raw (RAW) and finished (FIN) water samples (B). The linear relationship in panel A equates to $TOC = [UVA / (c2 : c1)] / \text{Slope}$, where the slope=0.160.

Table 1 DBPR2 removal rules for TOC as a function of alkalinity in drinking water treatment systems defined by Subpart H that use conventional filtration.^[1]

Source Water TOC (mg/L)	Source Water Alaklinity (mg/L) as CaCO ₃		
	0-60	>60 to 120	>120
>2 to 4	0.35	0.25	0.15
>4 to 8	0.45	0.35	0.25
>8	0.50	0.40	0.30

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
Figure 8 The daily alkalinity values for the raw water and calculated required removal fraction (A) and the Aqualog model TOC values for the raw water (RAW) and finished (FIN) and the calculated target removal (RULE) (B).

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
Figure 9 Multiple linear regression model of the SDS-THMFP for the raw and finished water samples (A) and comparison of the model and actual periodic SDS-THMFP measurements (B). The red line indicates the DBPR2 MCL for TTHM and the green symbols represent the actual distribution site values for TTHM.

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- [1] OFFICE OF WATER (2010) Comprehensive Disinfectants and Disinfection Byproducts Rules (Stage 1 and Stage 2): Quick Reference Guide. USEPA; <http://water.epa.gov/drink> [Accessed 20 April 2009]
- [2] POTTER, B.B. and WIMSATT J.C. (2003) EPA method 415.3- Determination of total organic carbon and specific UV absorbance at 254 nm in source water and drinking water. Rev. 1. USEPA
- [3] MURPHY, K. R. et. al. (2013) Fluorescence spectroscopy and multi-way techniques. PARAFAC. Anal. Methods. 5: pp. 6557-6566.
- [4] CLESCERL, L. S. et. al. (Eds) (1999) Standard Methods for the Examination of Water and Wastewater (20th Ed), American Public Health Assn.

Adam M. GILMORE, Ph.D.

HORIBA Instruments Inc.