

## August Slides

### Slide 1

Anaemia film shows  
Anisocytosis(++)  
Microcytes(++)  
Macrocyte(+)  
Dacryocytes(++)  
Thrombocytopenia(+)

### Slide 2

Perioperative medical  
resuscitation unit.  
Anisocytosis.  
Thrombocytopenia.  
Lymphopenia.

### Slide 3

Anaemia.  
Leukocytosis.  
Monocytosis.  
Discrete associated  
myeloma.  
PMS? CMML?

### Slide 4

Monomorphic  
hyperlymphocytosis  
associated with numerous  
smudge.  
Suspected CLL

### Slide 5

Nothing to report.

### Slide 6

See on the right.



## This issue


June Slides [P.1](#)  
Monthly Digital Case Study [P.1](#)  
Cell Quiz [P.2](#)  
Cold Agglutinin Disease [P.2-3](#)

## Monthly Digital Case Study Presentation August 2023, Slide 6


### FBC Results

#### Following incubation at 37°C

WBC 6.93 ( $10^3/\text{mm}^3$ )  
RBC 3.89 ( $10^6/\text{mm}^3$ )  
HGB 13.8 (g/dL)  
HCT 36.3 (%)  
MCV 93 (fL)  
MCH 35.6 (pg)  
MCHC 38.1 \* (g/dL)  
PLT 65 ( $10^3/\text{mm}^3$ )  
Neutrophils 78.4%  
Lymphocytes 12.3%  
Monocytes 7.3%  
Eosinophils 0.7%  
Basophils 1.4%



RBC auto agglutination



Platelet clumping

### Clinical Details


Male 74 years old

### Slide Information

Presence of agglutinins (visible at room temperature).  
Analyses performed after placing the EDTA tube in an incubator at 37°C (for 30 minutes). Note that the MCHC result remains high. Platelet result flag: presence of aggregates? False thrombocytopenia?



# Cell Quiz

Can you name the cell below?



## Last month's quiz:

The slides below are from the same patient. What feature can be seen in both slides?



### Right answer:

Both slides show platelet clumping. Platelet clumping is most commonly caused by pre-analytic errors such as improper or delayed sample mixing, clotted sample, over-filled or underfilled tubes, and time delay between sample collection and testing. It can also be an ex vivo phenomenon due to EDTA-dependent antibodies. Occasionally this may occur in due to infections, drugs, cold agglutinins or immunoglobulins.

## Expert's Comment

Residual presence of few red blood cells clumps. It's visible on the "RBCs wall", despite the passage of the sample in an oven at 37°C for 30 minutes.

The result of red blood cells measured in impedance method seems plausible because it is not associated with a flag and is comparable to the result of red blood cells measured in optical methods.

Platelet aggregates(+++): DO NOT REPORT the result of platelets.

Add a comment of the type "underestimation of the number of platelets: cluster of PLT".

## Cold Agglutinin Disease


Cold Agglutinin Disease (CAD) is a form of autoimmune haemolytic anaemia characterised by red cell autoagglutination (clumping of red cells ranging from a few red cells to numerous red cells) at cold temperatures (< 30°C). The incidence is fairly rare and has an estimated incidence of about 1 per million people per year and comprises 15% of all autoimmune haemolytic anaemia cases.

As most haematology analysers use the aperture impedance principles which measure the number of cells and volume as they pass through the aperture, if there is red cell autoagglutination then several red cells will be counted as one cell. This causes the red cell count to be inappropriately low with a low Hct and high MCV. The other calculated parameters are also affected with the MCHC potentially being above the physiological limit (> 36 g/dL depending upon type of analyser). Any MCHC above this value should be investigated as the cause could be presence of spherocytes, hyper lipidaemia (causing an artificial increase in the Hb due to increase in turbidity) or RBC auto agglutination. Just looking at the sample may be sufficient to detect hyper lipidaemia as the blood will have a pinkish colouration and if extreme RBC auto agglutination is present then the red cells may show separation from the plasma and appear as if clotted or if only mild then small discrete agglutinates may be seen if the sample is gently mixed end over end.


If RBC auto agglutination is suspected, then re-testing after warming the sample to 37°C for a period of time will show an improvement in the MCHC to more acceptable levels, but in extreme cases the MCHC may never return to acceptable limits regardless of the time left at 37°C. A blood film should be made from the sample at room temperature and at 37°C (use slides and spreaders pre-warmed to 37°C).

RBC agglutination in blood films is seen as irregular clumps of red cells and can be differentiated from the stack like appearance of red cells in the case of Rouleaux (see below).

A) RBC Agglutination



B) Rouleaux formation



As can be seen above, picture A shows a discrete clump of red cells sticking together, whereas in B) the cells appear to be in a line or stacked.

RBC CAD was first described in 1903 by Landsteiner (Nobel Prize winner for his discovery of the major blood group system and the development of the ABO system of blood typing which heralded the safe procedure of routine blood transfusion). Cold agglutination is now known to be caused by the production of IgM autoantibodies directed against either the I or I antigen. Since IgM is a large molecule, once bound onto red cells it can bridge the gap between the RBC overcoming the natural repulsive forces between cells therefore allowing spontaneous in-vitro agglutination. Bound IgM activates the complement pathway causing C3b to coat the cell, these cells lose surface membrane by receptor specific macrophages in the liver which results in extra vascular haemolysis. The degree of haemolysis is not due to the concentration of the IgM but the thermal amplitude of the IgM. The degree of anaemia is varied in CAD and can range from severe and life threatening, to mild/moderate.

CAD may either be primary or secondary to another condition with a median age of onset of 72. CAD has been seen across a wide age group secondary to infective causes, such as M pneumonia or Infectious Mononucleosis (Epstein – Barr virus). The CAD presenting as an acute onset and is usually transient. In adult cases, CAD is a well-recognised entity that is often caused by a bone marrow clonal B-Cell lymphoproliferative disorder such as Lymphoma, CLL, IgM monoclonal gammopathy or Waldenstrom's macroglobulinaemia.

Additional testing includes direct antiglobulin test (Coombs test), IgM cold agglutinin titre test, reticulocytes, LDH, and Haptoglobin.

### Cold Agglutinin Titre Test Explained

Serial dilutions of the patient's plasma are mixed with a 5% solution of group O cells, incubated at 4°C for a defined time period after which each dilution is checked for RBC agglutination. The titre is the highest dilution of the plasma at which RBC auto agglutination occurs. The tubes are then incubated at 37°C and are examined for RBC agglutination if the agglutination was caused by a cold acting antibody, then the agglutination will disappear. A titre greater than 1:64 is considered positive.

## Yumizen H500 haematology analyser wins the "Best New Clinical Instrumentation of 2022"

We are proud to announce that our [new Yumizen H500 benchtop haematology analyser](#) has won the **Best New Clinical Instrumentation Award for 2022** at the **Scientists' Choice Awards®**. The instrument is designed for small laboratories and is ideal for use in a variety of clinical laboratory environments and point-of-care (POC) settings or for anyone looking to upgrade from a 3-part to a 5-part differential technology solution. Requiring a sample of **20 µL** whole blood and an analysis time of just **60 seconds**, this makes the Yumizen H500 a highly flexible solution for **rapid blood analysis** in a multitude of situations, including **point-of-care** settings for **paediatric care, oncology units** and **drug monitoring services**. [Read more](#).



## QSP 2.0

Available as a single use license and a site license which allows up to 50 concurrent users.

To find out more, [contact us](#).

### Bibliography

Review Article – Cold Agglutinin Disease: Blood, 15 August 2013, Volume 23, Number 7

Transfusion Medicine and Hemostasi (Third Edition), Chapter 51c- Autoimmune Haemolytic Anaemia, Nancy L. Van Buren

Cold Agglutinin Titer Cold Agglutinin Disease News

### Editorial Team

Kelly Duffy  
Andrew Fisher

### About us

HORIBA UK Limited  
Kyoto Close  
Moulton Park  
Northampton, UK  
NN3 6FL

HORIBA Medical  
Parc Euromédecine, 390  
Rue du Caducée, 34790,  
France

[www.horiba.com/medical](http://www.horiba.com/medical)

