



**Patient Demographics:** Male, 6 years old



**Symptoms:** Fever and generalized weakness

**Diagnosis:** Acute Lymphoblastic Leukemia (ALL) – B-cell type, confirmed by immunophenotyping and molecular cytogenetics

**Microscopic Review:** Blasts: 88%, Neutrophils: 3%, Lymphocytes: 9%. Peripheral blood smear revealed small to large heterogeneous blast cells with irregular nuclear shapes, frequent clefting, and inconspicuous to prominent nucleoli.

**How Yumizen H500 Assisted:** The HORIBA Yumizen H500 raised an alarm for “Blast” suspicion and flagged suspected pathologies, prompting a review for large immature cells and atypical lymphocytes. On the DIFF scattergram, a lymphoblast population was observed in the atypical lymphocyte area, corroborated by extended immature WBC differential values (ALY & IML). This advanced feature, available on a mid-range hematology analyzer, enables mid-segment laboratories to support early diagnosis and intervention—crucial for managing challenging cases like childhood leukemia at the primary care level.

CBC			
			Range
WBC	34.16 *H	10 <sup>3</sup> /mm <sup>3</sup>	4.00 - 11.00
RBC	3.57 l	10 <sup>6</sup> /mm <sup>3</sup>	3.90 - 5.80
HGB	9.6 l	g/dL	11.5 - 16.7
HCT	28.7 L	%	35.0 - 49.0
MCV	80.3 *	fL	75.0 - 97.0
MCH	26.8	pg	26.5 - 33.0
MCHC	33.4	g/dL	32.0 - 36.0
RDW-CV	14.9 *	%	12.0 - 18.0
RDW-SD	37.0 *L	fL	37.0 - 56.0
MIC	9.5	%	0.0 - 20.0
MAC	2.9	%	2.0 - 10.0
PLT	15 *L	10 <sup>3</sup> /mm <sup>3</sup>	150 - 450
MPV	11.0 *h	fL	7.4 - 11.0
PCT	0.016 *l	%	0.150 - 0.400
PDW	7.9 *l	fL	11.0 - 20.0
P-LCC	6 *l	10 <sup>3</sup> /mm <sup>3</sup>	44 - 140
P-LCR	38.4 *	%	18.0 - 50.0



### Recommended actions

Slide Review

#### Alarms

##### LMNE

Abnormal differentiation  
WBC Interference

##### PLT

RBC PLT interference  
Schistocytes/Macro PLT?

##### Synthesis WBC


LIC?  
ALY?  
PLT aggregates or NRBC?  
Blasts?

##### Susp. Pathologies

Thrombopenia  
Leukocytosis  
Lymphocytosis  
Large Immature Cells  
Atypical Lymphocytes

NLR : 0.07

DIF				
	%	Range	10 <sup>3</sup> /mm <sup>3</sup>	Range
NEU	6.8 *l	40.0 - 75.0	2.32 *	1.50 - 7.00
LYM	90.9 *h	15.0 - 45.0	31.08 *H	1.25 - 4.00
MON	2.0 *l	4.0 - 13.0	0.67 *	0.20 - 0.80
EOS	0.3 *l	0.5 - 7.0	0.09 *	0.00 - 0.40
BAS	0.0 *	0.0 - 2.0	0.00 *	0.00 - 0.10
IMG	0.9 *	0.0 - 2.0	0.31 *	0.00 - 0.50
IMM	0.3 *	0.0 - 0.5	0.11 *h	0.00 - 0.10
IML	1.2 *h	0.0 - 0.2	0.41 *h	0.00 - 0.05
ALY	51.3 *H	0.0 - 2.5	17.51 *H	0.00 - 0.20
LIC	3.1 *H	0.0 - 3.0	1.06 *H	0.00 - 0.20



Lymphoblastic population

