

## INTRODUCTION

The routine D-dimer quantification to exclude venous thromboembolism has led to the development of many assays, the usefulness of which depends on their reliability and performance.

## AIM


To evaluate the analytical performances of the Yumizen G DDi 2 assay (HORIBA Medical) performed with the Yumizen G800 analyzer and to compare it with other available D-Dimer assays: Vidas D-dimer Exclusion II for VIDAS<sup>®</sup> 3 (BioMérieux), STA<sup>®</sup>-Liatest<sup>®</sup> D-Di Plus for STA-R Max<sup>®</sup> (Diagnostica Stago), Innovance<sup>®</sup> D-dimer for Sysmex<sup>®</sup> CS-2100i (Siemens Healthineers) HemosIL<sup>®</sup> D-dimer HS500 for ACL TOP 700 (Werfen).

## METHOD

Within-run and between-run imprecision were evaluated using low- and high quality control plasma samples. The limit of detection, limit of quantification and linearity were determined according to the Clinical and Laboratory Standards Institute guidelines<sup>1</sup>. Interference due to hemolysis, icterus, lipemia, rheumatoid factor (RF) or heterophilic antibodies (HAMA) was evaluated by spiking plasma samples with hemolysate, bilirubin, Intralipid<sup>®</sup>, RF, or HAMA<sup>2</sup>. The measurements obtained with the different D-dimer assays were compared with the Bland-Altman plot method.

## RESULTS

The coefficients of variation values of within- and between-run were in accordance with the specifications by Ricos and colleagues: <3% for D-dimer values >1000 ng/mL fibrinogen-equivalent units (FEU) and <6% for D-dimer values close to the threshold of 500 ng/mL FEU. The assay linearity was very good for a broad range of concentrations, up to 32700 ng/mL FEU. Hemolysis and icterus did not have any effect up to 10 g/L hemoglobin and 300 mg/L bilirubin. Lipemia seemed to generate an underestimation of D-dimer concentration when Intralipid<sup>®</sup> concentration was > 5 g/L (Figure 1). RF and HAMA did not have any effect. The Bland-Altman analysis showed strong agreements between the Yumizen G DDi 2 assay and other assays performed with 66 fresh plasma samples with a wide range of D-dimer concentrations. Discordant values were found in D-dimer values > 2000 ng/mL FEU, without clinical impact.




**Figure 1.** Bias (%) in D-dimer measurement in function of hemoglobin (1A), bilirubin (1B) and lipid (1C) levels. Dotted lines represent the limit of acceptable bias. Each color corresponds to a plasma sample, and the different points correspond to the different concentrations of hemoglobin, bilirubin and intralipid spiked in the samples.

## CONCLUSIONS

Its analytical performances and main technical features indicate that the new Yumizen G DDi 2 is suitable for the rapid quantification of D-dimers in clinical laboratories.

## REFERENCES

- Clinical and Laboratory Institute (CLSI).** Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; approved guideline. 2nd Edition. CLSI Document EP17-A2 2012. <https://clsi.org/standards/products/method-evaluation/documents/ep17/> (accessed December 9, 2020).
- Nougier C, Jousset E, Sobas F, Pousseur V, Négrier C.** Effects of hemolysis, bilirubin, and lipemia interference on coagulation tests detected by two analytical systems. *Int J Lab Hematol* 2020;42:88–94.



**Figure 2.** Bland Altman plots of the D-dimer values obtained with the Yumizen G DDi 2 assay and the Yumizen G800 analyzer and other available assays. The X axis represents the mean of the measurements and the Y axis represents the difference between the measurements obtained with the two systems. Continuous and dotted lines represent the bias and the lower and upper limits of agreement with the 95% confidence intervals, respectively.