ISSUE 26 June

2022

CSP Monthly Slide PROGRAM

Last Month's Slides

Slide 1 See Case Study

Slide 2 Nothing to report

Slide 3

Leukocytosis. Neutrophilia (hyposegmented, degranulated and highly vacuolated neutrophils). Myelemia (Immature Granular cells degranulated). Blastosis. Expert's comments: AML post MDS probable. Difficult +++ for the cellular differential. Do not hesitate to reclassify "atypical" cells as blasts. Significant dysgranulopoiesis

Slide 4 Nothing to report

Slide 5 Nothing to report

Slide 6 Nothing to report

Monthly Digital Case Study May 2022 Slide 1

Presentation


Haematology clinic patient. Clinical haematology unit. Deterioration of general condition. Myelemia. Reticulocytosis.

FBC Results

WBC 27.0 * (10^3/mm3)
RBC 2.4*(10^6/mm3)
HGB 7.5*(g/dL)
HCT 23.9* (%)
MCV 100(fL)
MCH 31.2(pg)
MCHC 31.3(g/dL)
PLT 227(10^3/mm3)
RDW-SD 80fL
Reticulocytes 225 x10^9/L


Slide review

Anisocytosis(RDWsda70fL). Microcytes(++). Macrocyte(++). Schistocytes(++). Echinocytes(++). Prescence of Schistocytes. Check for Thrombotic Microangiopathy (TMA).



Schistocytes





This issue

Cell Quiz P.2

Last Month's Slides P.1

Monthly Digital Case Study P.1-2

Previous 12 Months' Issues P.3

Echinocyte

Macrocytic Polychromatic Cell



HORIBA

Cell Quiz:




What term best describes the above neutrophil

- a) Hypersegmented
- b) Hypogranular
- c) Pelger Huet form

Last Month's Cell Quiz:

An elderly patient presents with extensive red rash (erythroderma), enlarged spleen and lymph nodes. Blood film shows numbers of abnormal lymphocytes like the one above. What condition may be present?



a) Glandular Feverb) Sezary Syndromec) Allergic reaction

Answer:


b) Sezary Syndrome. Sezary Syndrome is a rare aggressive form of cutaneous T cell Lymphoma with patients often presenting with a red severely itchy rash (erythroderma). Sezary Cells have a convoluted or cerebriform nucleus with tightly intertwined nuclear lobes.

Slide 1 Continued

Schistocytes or schizocytes (from the Greek word schisto, broken or cleft, corresponding verb schizo) are circulating fragments of red blood cells (RBC) or RBCs from which cytoplasmic fragments have been lost. Schistocytes are absent in the blood film of normal healthy individuals. Schistocytes are formed by extrinsic mechanical destruction caused by the passage of the TBC through fibrin strand. Fibrin stands are formed at the site of damage to a blood vessel and a clot begins to form.

Schistocytes are a very important prognostic marker for Thrombotic Macroangiopathic Anaemia (TMA) and if present in a blood film, further investigation must be performed, particularly if they are the dominant RBC morphological feature. In 2011, ICSH published guidelines on the identification, diagnostic value, and quantitation of schistocytes.

A suggested flow chart from the SCH guidelines is detailed below:



TMA includes 2 major syndromes: Thrombotic Thrombocytopenia Purpura* (TTP) and Haemolytic Uraemic Syndrome (HUS) Schistocytes may also be found in patients with malfunctioning prosthetic valves, HELLP syndrome, malignant hypertension, and metastatic cancer. RBC fragments similar to Schistocytes can also be found in genetic or acquired RBC disorders (RBC membrane defects, thalassaemia, megaloblastic anaemia, primary myelofibrosis, and burns.

Schistocytes can be classified as Keratocyte, helmet, or triangle. Microspherocytes are also part of the schistocyte family in TMA.

A schistocyte count should be performed on a well stained blood film with the number of schistocytes counted per 1000 RBC.

HORIBA

* For more information, please see QSP issue 9, January 2021.

The first issue of the QSP Newsletter was published in May 2020. To celebrate the 2nd anniversary we looked back to the last 12 issues:

Issue 25 |May 2022: Case Study - Slide 3, 27% Erythroblasts, Target Cells, Howell Jolly Bodies. Article - Malaria testing

Issue 24 | April 2022: Case Study - Slide 2, High WBC, smudge cells - lymphoproliferative syndrome Article - Malaria overview

Issue 23 | March 2022: Case Study - Slide 2, Blast Cells indicative of AML Article - Rare Lymphoid Clonal Abnormalities (hairy cell leukaemia, SLVL)

Issue 22 | February 2022: Case Study - Slide 3, Atypical Lymphocytes - Infectious Mononucleosis Article - An interesting case of Thrombocytopenia (May-Hegglin Anomaly)

Issue 21 | January 2022: Case Study - Slide 5, Rouleaux formation Article - Red Cell Morphology terminology part 2 - Variation in shape

Issue 20 | December 2021: Case Study - Slide 6, Vacuolated and hypergranular Neutrophils Septic Shock Article - Red Cell Morphology Part 1 – variation in size and colour

Issue 19 | November 2021: Case Study - Slide 3, High WBC, RBC target cells, Howell Jolly Body, Cabot's Ring Article - Monocytes

Issue 18 October 2021: Case Study – Slide 2, Low Platelet count due to Platelet aggregation Article - Platelet aggregation detection on HORIBA Yumizen H1500/H2500, Eosinophils

Issue 17| September 2021: Case Study - Slide 2, Neutrophilia, Anisocytosis Article – Basophils

Issue 16 | August 2021: Case Study - Slide 1, Burkitt's Lymphoma Article – Granulocyte Development

Issue 15| July 2021: Case Study - Slide 3, Lymphoproliferative Sezary Syndrome Article - Paroxysmal Nocturnal Haemoglobinuria

Issue 14 | June 2021: Case Study - Slide 1, target cells, spherocytes Article - Platelets and counting methodology

To find a full list of QSP Newsletters visit our website.

<u>QSP 2.0</u>

Available as a single use license and a site license which allows up to 50 concurrent users. To find out more, <u>contact us.</u>

Bibliography

ICSH recommendations for identification, diagnostic value, and quantitation of schistocytes

Editorial Team Kelly Duffy Andrew Fisher

About us HORIBA UK Limited Kyoto Close Moulton Park Northampton, UK NN3 6FL HORIBA Medical Parc Euromédecine, 390 Rue du Caducée, 34790, France

www.horiba.com/medical



HORIBA