### Last Month's Slides



### This issue

Last Month's Slides P.1

Monthly Case study P.1

Platelet morphology P.2-3

Monthly Quiz P.2

Top Tips P.3

### July 2020 Slide Summaries

### Slide 1

Mostly normal slide with some platelet aggregates and giant platelet

### Slide 2

Generally normal film

### Slide 3

Lymphoproliferative disorder, probable CLL. Smear cells and abnormal lymphocytes

### Slide 4

Myelodysplastic blood film with circulating megakaryocytes, Blasts Multiple cell abnormalities. Erythroblasts, Cabot's rings in RBCs Basophilia

### Slide 5

Lymphoblasts, abnormal lymphoid cells. Probable ALL

### Slide 6

Leukoerythroblastic blood film. Myelodysplastic changes, approx. 50% blasts, auer rods, monocytoid blasts. Probable M4 Acute

### Monthly Digital Case study July 2020 Slide 5

### **Presentation**

Child (5 years old)

Presenting with weakness and jaundice

### **FBC Results**

WBC	2.63 (10^3/mm3)
RBC	2.37 (10 <sup>6</sup> /mm3)
HGB	65 (g/L)
HCT	20.3 (%)
MCV	85 (fL)
MCH	27.4 (pg)
MCHC	32.1 (g/dL)
PLT	104 (10 <sup>3</sup> /mm3)

### Neutrophils 7.6% Lymphocytes 87.3% Monocytes 1.7% Blasts 3.4%


### Slide review

Leukopaenia/neutropaenia confirmed

Small blasts present with minimal cytoplasm – lymphoid in appearance Abnormal lymphoid cells (counted with lymphocytes) with cytoplasmic projections. Almost certainly some of these abnormal cells are lymphoblasts making the differentiation between these abnormal cells and obvious blasts cells subjective

### **Diagnosis**


Probable Acute lymphocytic Leukaemia (immunophenotyping required to complete the diagnosis)






Monthly Morphology Quiz

Look closely at this red cell:




What is unusual about it and what could this Indicate?

### Last month's cells:



The blood film shows neutrophils with Dohle bodies in the cytoplasm

These are remnants of rough endothelial reticulum and can accompany toxic left-shift



# Platelet morphology in peripheral blood

An overview of laboratory findings


### Introduction

Platelets are cell fragments circulating in the peripheral that are involved in secondary haemostasis. Aggregated platelets, platelet polymorphism, platelet satellitism and the presence of platelet precursors, megakaryocytes, in the peripheral blood are useful indicators of a variety of conditions.

### **Platelet production**


Platelets are produced in the bone marrow by fragmentation of the tips of cytoplasmic extensions of cells called megakaryocytes. Each cell produces approximately 1000 to 5000 platelets. They are released into the blood stream through the endothelium of the vascular niches of the bone marrow.

There is a 10-day cycle for the production and release of platelets:



### **Platelet aggregates**

Platelet aggregates are very common and can give rise to falsely low platelet counts and occasional interference with other parameters. They can be cause by slow or difficult venepuncture but, in some individuals, platelets may be sensitive to EDTA anticoagulant and invariably clump when blood samples are taken. In these instances, a sample taken into tri-sodium citrate can give a correct result after correction for anticoagulant dilution factor.





# Platelet morphology in peripheral blood

### **Continued from page 2**

### **Platelet Satellitism**

Platelet satellitism (illustration above) is a rare phenomenon in EDTA blood where the platelets congregate around a white cell (usually a neutrophil). It is generally considered to be an artifact, but it has been observed in some conditions including lymphoproliferative disorders, lupus, vasculitis and liver disease but a cause has not been established.

### **Giant platelets**

The presence of giant platelets may indicate that there is increased platelet production as younger platelets tend to be larger. It is also associated with Idiopathic thrombocytopaenic Purpura (ITP). In this case, the platelets may not function correctly as they are unable to stick to the injured blood-vessel wall. Giant platelets also occur in hereditary platelets disorders such as Bernard-Soulier Syndrome.

### Megakaryocytes

Megakaryocytes are normally not seen in the peripheral blood and their presence is indicative of bone marrow such as Myelodysplastic Symdrome, Myeloid Leukaemia, Polycythaemia Rubra Vera and myelofibrosis.


Megakaryocytes in the peripheral blood are often 'bare' – ie. with minimal cytoplasm (illustrations of giant platelets and megakaryocyte from July 04 slide).

### This Month's Top Morphology Tip

### The Hunt for Malaria

Common errors in looking for malaria in thin films are caused by platelets overlying a red blood cell, if in doubt, rack in and out of focus on the particular cell, if the platelet is on top of the red cell, it should come into focus before the red cell, which will still be slightly blurred. If it is a malaria parasite within the red cell, clear focus will be on the same plane.





### **Other News**

QSP 2.0 Available now!

Options for a single operator or site license which allows up to 10 concurrent users

### **Bibliography**

QSP July2020

Hoffbrand's Essential Haematology 7<sup>th</sup> edition Wiley Blackwell

### **Editorial Team**

Kelly Duffy Mandy Campbell Shubham Rastogi

### **About us**

HORIBA UK Limited Kyoto Close Moulton Park Northampton, UK NN3 6FL

HORIBA Medical Parc Euromédecine, 390 Rue du Caducée, 34790, France

www.horiba.com/medical

