

Our bodies are constantly exposed to pathogens present in the environment, including viruses, bacteria, parasites, and fungi. The immune response is responsible for defending us against them through a wide range of mechanisms. It is mediated by white blood cells (WBCs), also called leukocytes, which include neutrophils, eosinophils (acidophiles), basophils, lymphocytes, and monocytes. Innate and adaptive immunities work together to protect us (1).

Phagocytic WBCs (in particular, neutrophils) are part of the innate branch of the immune system, along with tissue macrophages, that have a number of different functions such as triggering the immune response and cell to cell communication. Lymphocytes mediate the adaptive immune response via antigen-specific activities.

Plasma cells, not normally seen in the peripheral blood, are differentiated B-lymphocytes. They are non-proliferating and much larger than B cells and can secrete large amounts of antibodies (2).

Blood Smear

Seen microscopically on a blood smear, lymphocytes have a regular, spherical nucleus and a high nucleo-cytoplasmic ratio. The cytoplasm, which varies in volume, is stained a pale blue with Romanowski staining and is either clear or with scanty, small pink granules. The presence of granules in a larger cell with more cytoplasm may indicate that the cell is a T-lymphocyte but immunophenotyping is required to further differentiate the different lymphocyte subtypes (3).

Leukocyte differential count

A complete blood count (CBC) is a highly automated laboratory test which generates data for all blood cells. Part of the CBC is the WBC differential which identifies and counts the different types of normal blood cells and also provides additional measurements and flags regarding atypical, abnormal and immature WBCs.

In the Yumizen hematology analyzers, white blood cells can be visualized on the Diff scattergram or LMNE matrix (Figure 1). In this matrix optical extinction (Y-axis) is plotted against volumes (X-axis) to form an image with several clusters of cellular events. The matrix has thresholds, some fixed and some mobile, which allows the different populations to be determined and to detect abnormalities.




Figure 1 : Lymphocytes (LYM) and ALYs on LMNE matrix

This precise technology allows Yumizen hematology analyzers to define two additional leukocyte populations of relevance to lymphocytic cells: Large Immature Cells (including immature lymphoid cells – IML) and Atypical lymphocytes (ALY), reported in percentage as well as absolute counts.

	Female	Male	Bibliographic data (male and female)
Lower reference limit (ALY%)	0.20	0.20	0.00
Upper reference limit (ALY%)	1.90	1.70	2.50
Lower reference limit (ALY#)	0.01	0.01	0.00
Upper reference limit (ALY#)	0.14	0.09	0.20

Table 1: Lower and upper reference limits determined by the study for ALY parameter compared to values from bibliographic data (5, 6).

What is their clinical significance?

The presence of few atypical lymphocytes in most blood samples is probably of little clinical significance. However, the ALY parameter can help the clinician to consolidate a diagnosis or to use it for a prognostic purpose. Atypical lymphocytes are most commonly found in Infectious Mononucleosis, Toxoplasmosis, Hepatitis and other viral infections.

Atypical lymphocytes have also been observed in the peripheral blood of patients in a large number of other clinical situations, including immune reactions to transplantation and immunization, collagen diseases and other autoimmune disorders, malignant disease, drug reactions, certain bacterial infections, medications and stress (7). In Acute Lymphoid Leukemia (ALL) with small blasts, these lymphoblasts are frequently seen in the ALY area. Figure 4 shows an example of T-cell ALL where the lymphoblast population is especially present in the atypical lymphocyte area of the LMNE matrix.




Figure 4: Lymphoblast population in acute lymphoid leukemia T (ALL T)

In chronic lymphoid leukemia, atypical lymphocytes and immature lymphocytic cells also can be observed on the ALY area of the DIFF scattergram (8).

Small compact blast cell population from acute myeloid leukemia can feature in this region owing to their smaller size. In this case, the compact blast population can be observed on the LMNE matrix between the lymphocyte, atypical lymphocyte and monocyte areas and are flagged accordingly. For monitoring the progress of a given infection, ALYs can be used as a prognostic tool. For example, atypical lymphocytes circulating in blood have been reported in COVID-19 patients. COVID-19 reactive lymphocytes were found to be related to a better evolution and prognosis of the disease. Their presence suggests an abundant production of virus-specific T cells, thus indicating the better prognosis of patients with these cells in the circulation (9, 10).

Comments

Today's high-throughput analyzers are designed to provide the most comprehensive information about a sample. This data is based on threshold values and complex algorithms, and it must be subsequently carefully interpreted by the clinician. In this context, the ALY parameter provides information that may complement the routine blood parameters and add to the patient's clinical picture. It may also be useful in disease management or for prognostic purposes.

Yumizen Chase the Case #2





Patient demography: 32 years, Male

Suspicion: Lymphoma suspected

Clinical Sign: Splenomegaly and enlarged lymph nodes

Other results: Immunophenotyping: The analyzed population corresponds to mature B lymphocytes with clonality for light chain immunoglobulin Lambda Immunophenotyping compatible with Chronic Lymphoproliferative Disease B.

Diagnosis: Burkitt's lymphoma



HORIBA Team:
Fábio Oliveira, Daniele Pacheco, Shubham Rastogi
To find out more, [contact us](#)
HORIBA Medical Brazil & HORIBA Medical France

Dr. Jaqueline Valente
Application Specialist
Rio de Janeiro - Brazil

References

1. H.C. Allen, P. Sharma. *Histology, Plasma Cells*, January 28, 2021
2. R. Warrington, W. Watson, H. Kim & F. Antonetti, *An introduction to immunology and immunopathology, Allergy, Asthma & Clinical Immunology volume 7, Article number: S1 (2011)*
3. A. Adewoyin and B. Nwogoh, *Peripheral blood film - a review, Ann Ib Postgrad Med. 2014 Dec; 12(2): 71–79.*
4. M.W. Simon, *The Atypical Lymphocyte, International Pediatrics, January 200318(1):20-22*
5. C. Sultan / M. Gouault - Helman / M. Imbert, *AIDE MEMOIRE D'HEMATOLOGIE. Service Central d'Hématologie de l'Hôpital Henri Mondor, Faculté de médecine de Créteil (Paris XII). 1998 Troussard X, Vol S,*
6. E. Cornet, V. Bardet, JP. Couaillac, C. Fossat, JC. Luce, E Maldonado, V. Siguret, J. Tichet, O. Lantieri, J. Corberand. *Full blood count normal reference values for adults in France. Journal of Clinical Pathology (2014) 67 (4): 341-4.*
7. T. Shiftan, J. Mendelsohn, *The circulating "atypical" lymphocyte, Human Pathology, 1978 Jan;9(1):51-61*
8. *Horiba Medical: Clinical Case Studies. Yumizen H1500/H2500*
9. A. Merino, J. Rodellar et al., *Atypical lymphoid cells circulating in blood in COVID-19 infection: morphology, immunophenotype and prognosis value. Journal of Clinical Pathology, 2020 Dec : jclinpath-2020-207087.*
10. D. Gérard, S. Henry, B. Thomas, *SARS-CoV-2: a new aetiology for atypical lymphocytes, Images In Haematology, 20 April 2020*

Thanks, note:

“The authors thank Dr Parag Dharap, MD, Dr. Dharap's Diagnostic Centre, Mumbai, India for manuscript review and for general advice”