

# A-TEEM Multidimensional Fluorescence Spectroscopy as a Powerful Tool for Unveiling Drug Release Kinetics in Liposomal Delivery Systems

Hakam Alaqabani<sup>1,\*</sup>, Jeffrey A. Julien<sup>2</sup>, Yvonne Perrie<sup>1</sup>

<sup>1</sup> Strathclyde Institute of Pharmacy and Biomedical Sciences, University of Strathclyde, 161 Cathedral St, Glasgow G4 0R, UK

<sup>2</sup> HORIBA Instruments Inc., 20 Knightsbridge Rd, Piscataway, NJ 08854, USA




## Introduction

Recent advancements in microfluidic technology have revolutionized drug delivery by enabling precise fabrication and manipulation of liposomes for controlled drug release. Liposomes, versatile carriers for therapeutic agents, have shown substantial efficacy in enhancing drug delivery efficiency. This study pioneers using Absorbance-Transmission Excitation Emission Matrix (A-TEEM) spectroscopy to investigate doxorubicin (DOX) release from liposomes. A-TEEM combines absorbance and fluorescence spectroscopy, providing a contour plot that utilizes the UV/Vis absorbance spectrum and full fluorescence excitation and emission spectra to perform an automated correction of the Inner Filter Effect (IFE). Unlike traditional fluorescence EEMs, A-TEEM provides a true molecular fingerprint, enabling quantitative analysis of responsive components.

This study highlights A-TEEM spectroscopy as a novel, powerful method for analysing drug release dynamics from liposomal formulations. A-TEEM offers a sensitive, rapid, non-destructive, and cost-effective way to quantify drug release from liposomal nanoparticles, avoiding shear stress and potential damage associated with established centrifugation methods. Initially, DOX release was detected and quantified in various buffered media at different temperatures in DSPC/cholesterol/DSPE-PEG2000 liposomes but has now expanded to include cellular analyses in MDA-MB-231 breast cancer cell lines.

## Material and method

Liposomes were synthesized using the NanoAssemblr™ Benchtop and a staggered herringbone micromixer. The formulations consisted of Distearoylphosphatidylcholine (DSPC), cholesterol, and N-(carbonyl-methoxy (polyethylene glycol)-2000)-1,2-distearoyl-sn-glycero-3-phosphoethanolamine (DSPE-PEG2000-SH), in a w/w ratio of 3:1:1 labeled with 1% w/w DiR'. DiIC18(7) (1,1'-Diocadecyl-3,3',3'-Tetramethylindotricarbocyanine Iodide) Additionally. Purification was achieved via Tangential Flow Filtration, and loading with 12.5% w/w doxorubicin using the pH gradient method. Following this, the liposomes underwent characterization to assess their size, polydispersity index (PDI), and zeta potential. Additionally, UV-Vis from, FTIR analyses were carried out to further evaluate their properties, and the drug release was monitored over 144 hours at varying temperatures and different media using Absorbance-Transmission Excitation Emission Matrix (A-TEEM) spectroscopy from HORIBA.



## Conclusion

- A-TEEM has effectively measured drug release from liposomes and from within cells, proving its capability to detect drugs in buffer solutions and within cells.
- This versatility and precision make A-TEEM an invaluable tool for studying drug release mechanisms, ensuring accurate monitoring and understanding of therapeutic agents' behaviour within biological systems

## Acknowledgments

We would like to express our gratitude to the Centre for Process Analytics and Control Technology (CPACT) at the University of Strathclyde for their invaluable support and collaboration.

| Feature                 | Principle  | Strengths   | Weaknesses  | Data Output                                  | Suitability for  |
|-------------------------|--|---|---|--|--|
| Dialysis Centrifugation | Separates drug release media from dosage form using a semi-permeable membrane              | <ul style="list-style-type: none"> <li>Simple setup.</li> <li>Suitable for small sample volumes</li> </ul>  | <ul style="list-style-type: none"> <li>Time-consuming process (40 min per sample)</li> <li>Costly</li> <li>Yields varying outcomes</li> </ul> | Cumulative amount of drug released over time | In vitro studies for initial screening   |
| Dissolution Apparatus   | Flow-through system continuously pumps dissolution media past the dosage form              | <ul style="list-style-type: none"> <li>Mimics in vivo conditions with continuous flow.</li> <li>Suitable for a wide range of dosage forms</li> </ul>                                      | <ul style="list-style-type: none"> <li>Potential for shear stress affecting release kinetics,</li> <li>limited to sink conditions</li> </ul>  | Dissolution profile (concentration vs. time) | In vitro studies for formulation development, USP compliance testing                                       |
| A-TEEM (Spectroscopy)   | Measures drug release by monitoring changes in absorbance or fluorescence of released drug | <ul style="list-style-type: none"> <li>Non-destructive analysis</li> <li>provides information on both drug release and degradation</li> <li>Time effective with rapid analysis</li> </ul> | <ul style="list-style-type: none"> <li>May not be suitable for all drugs</li> </ul>   | Dissolution profile (concentration vs. time) | In vitro studies for detailed analysis of release mechanisms, compatibility studies, and stability studies |

| Sample  | Scatter  | Fluorescence |
|---|----------|--------------|
| Mean [nm]   | 90       | 88           |
| Mode [nm]   | 72.5     | 72.5         |
| D10 [nm]  | 64       | 51           |
| D50 [nm]  | 80       | 78           |
| D90 [nm]  | 133      | 131          |
| Concentration [particles/ml]                          | 8.75E+12 | 8.93E+12     |
| Bin Size [nm]   | 5        | 5            |
| Fluorescence Labelling Efficiency - Concentration (%) | -        | 102          |

