

Aqualog[®] For Water Treatment Organics Monitoring Featuring Datastream Dashboard



www.aqualog.com

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Aqualog[®]

For Drinking Water Treatment Plants

Validated, Automated Organic Analysis

The Aqualog monitors Dissolved Organic Carbon (DOC) concentration and composition, Disinfection By-Products (DBPs), DBP precursors, aromatic petroleum hydrocarbons, algae, and algal pigments associated with toxins, taste and odor issues. It comprises a patented two-in-one spectroscopic instrument providing absolute molecular fingerprints and a software package called "Aqualog Datastream" that automatically interprets the water sample measurements. Reported are high, low and Maximum Contaminant Level (MCL) threshold values of water quality parameters, providing timely actionable data to the water treatment plant operator. This is accomplished with simple push button operation and minimal sample preparation, and consumables.

Early Warning Sentinel

The Aqualog is an early warning sentinel for water treatment plants since it provides nearly instantaneous determination of DOC and composition, as well as Trihalomethanes (THM) formation potential, among many other key parameters. The early alerts facilitate prompt follow-up measurements according to established protocols, SOPs and regulations.

Saves Money

The DOC and enhanced THM formation predictive features of Aqualog have been documented to potentially save, on average, 5-10% of the annual chemical budget of a typical drinking water treatment plant. Chemical dosing needs to be applied only when the monitored levels are predicted to rise above pre-determined thresholds to ensure spending remains within the chemical budget for a given water treatment plant.

Rapid Return on Investment

The Aqualog provides rapid return on investment. With the typical annual chemical expense savings, the purchase of an Aqualog pays for itself within the first three to six months, and thereafter, those savings accrue every year. *For large municipalities, and large water companies, the Aqualog can save millions of dollars.*



It's All About Organics

The Aqualog is a novel and valuable optical tool in the organics laboratory for water treatment plants. It is superior to conventional instruments in speed, sensitivity, and selectivity for organics.


A primary concern for Drinking Water Treatment Plants (DWTPs), is humic and fulvic acids, due to their reactivity with halogenated disinfectants and tendency to form toxic DBPs, including THMs and haloacetic acids (HAAs).

Together, humic and fulvic acids comprise the majority of the DOC of natural organic matter in most surface water sources.

Humic acids are high molecular weight, aromatic compounds with multiple phenolic and carboxylic moieties linked together. They have a significant negative surface charge which gives them high affinity to positively charged coagulant compounds.


Fulvic acids, which are soluble in water at any pH, are relatively lower in molecular weight than humic acids, and have lower relative affinity to coagulants.


Protein-like compounds are another significant component of natural organic matter present in most surface water sources. They are also known to be associated with the presence of municipal wastewater and microbially available substrates. Compared to humic and fulvic acids, the protein-like compounds have a lower affinity to coagulants.



Humic Acids

- High MW
- Hydrophobic
- Highly Aromatic
- Absorbed by Coagulants
- Main SUVA (Specific UV Absorbance) Component)





Lower MWHydro/Transphilic

Fulvic Acids

Less Aromatic

Tryptophan and Tyrosine

- Protein-like
- Associated with Wastewater Effluent
- Biopolymers

The Gold Standard for Water Treatment

A sampling of Aqualog customers around the world

Water Companies

American Water Chelsea Technologies Group Doosan Heavy Industries and Construction Eskom (South Africa) Hazen and Sawyer Kurita Water (Japan) Public Utilities Board (PUB) of Singapore Sabesp (Brazil) Suez (Worldwide) WET Labs

Municipal Water Facilities

City of Åkron City of Philadelphia Water Department City of Sandusky City of Wheeling Water Department Denver Water Hampton Roads Sanitation District Las Vegas Valley Water District Louisville Water Company Metropolitan Water Company Metropolitan Water Company Orange County Water District Umgeni Water - Amanzi (South Africa) West Basin Municipal Water District Water Recycling Facility


Environmental Research

US Environmental Protection Agency (US EPA) US Geological Survey (USGS) Water Science Center US Naval Research Laboratory National Institute of Standards and Technology (NIST) National Aeronautics and Space Administration (NASA) National Oceanic and Atmospheric Administration (NOAA) Woods Hole Oceanographic Institution Stroud Water Research Center **Trussell Technologies** Vietnam Environment Administration National Laboratory for Civil Engineering (Portugal) Korea Institute of Civil Engineering and Building Technology (South Korea) Arizona State University Chinese Academy of Sciences (China) Colorado School of Mines Columbia University Florida International University Florida State University Georgia Institute of Technology Harbin Institute of Technology (China) Indiana University Kangwon National University


Kangwon National University (South Korea) King Abdullah University of

King Abdullah University of Science and Technology (Saudi Arabia) Kobe University (Japan) Louisiana State University Michigan Technological University New Mexico State University Northeastern University Oregon State University Rutgers University San Diego State University Seattle University Sejong University (South Korea) Sichuan University (China) Southwest University (China) Swedish University of Agricultural Sciences (Sweden) The Ohio State University The University of Vermont Tongji University (China) Umeå University (Sweden) University of Alaska University of Alberta (Canada) University of East Anglia (UK) University of Extremadura (Spain) University of Maryland, Center for Environmental Science University of Massachusetts at Amherst University of Michigan University of Minnesota University of Montana University of New Orleans University of Science of Technology of China (China) University of South Africa (South Africa) University of Western Ontario (Canada)


Aqualog in Use at a Drinking Water Treatment Plant




This illustration shows the Aqualog in use at a DWTP. It is found in the analytical lab where water is collected from various points along the water treatment process. Typically, Aqualog collects and measures water from three points in a WTP: The untreated water (Raw); after the sedimentation tank (Settled); and in the final treated effluent (Finished). The charts on the next page demonstrate the many parameters that are automatically reported and tracked over time with the Datastream Dashboard HTML interface.



Humic/Fulvic ratio decreases due to humic acid coagulation in transitioning from Raw to Settled/Finished state.




A254 is an indicator of aromatic organic (humic compounds) concentration as it decreases from Raw to Settled/Finished state.




DOC decreasing from Raw to Settled/Finished state.


Aqualog offers unique capabilities for monitoring Natural Organic Matter (NOM). Many DOC components are DBP precursors. Halogenated disinfectants, such as chlorine, can react with the DBP precursors to form unwanted DBPs,



SUVA, an indicator of aromatic organic concentration, decreases upon coagulation treatment. This is primarily due to the relative decrease in humic acid content and consequent reduction in aromaticity.



Trihalomethane Formation Potential (THMFP) decreases from a Raw to Settled/Finished state due to removal of the disinfection by-product precursors. The MCL indicates the EPA regulated maximum contamination limit of 80 μg/L allowed for THMs.



Aqualog model prediction of THM species in finished water. A good correlation $(R^2 = 0.972)$ was observed with the independently measured data.

which include toxic substances, such as THMs and HAAs. Because these substances are potentially carcinogenic and are regulated by the US EPA, their formation should be controlled by properly managing and optimizing the water treatment process.



Automation Accessories


Aqualog can be used with grab samples for manual measurements in an optical cell, as shown in the picture above, however, it is best configured with Sipper accessories described below.

AQ-Sipper:

HORIBA Sipper for automated sample extraction and sample measurement with Aqualog software. Includes sipper accessory and sample tray with thermostated cuvette holder and leak sensor. Extracts from a single sipper tube. External water bath not included.


AQ-Sipper-4S:

Optional four-position sampler for AQ-Sipper. Allows the AQ-Sipper to automatically extract from up to four different sipper tubes (raw, settled, and finished). Multiple units may be used in parallel.




A stable baseline absorbance for over 90 days of operation shows excellent sipper performance in a real plant condition. The flow cell prevents air bubbles and biofouling in the sample cell. Only quarterly cleaning is required for the flow cell, thanks to the automatic cleaning protocol.

Fingerprints of Wastewater, Algae and Aromatic Hydrocarbons




Blue-green Algae Fingerprint



Green, bluegreen and brown algae cells and pigments, associated with toxins, taste and odor issues, can be easily measured with the Aqualog because they each have characteristic A-TEEM fingerprints.

The Aqualog exhibits highsensitivity to aromatic hydrocarbons, and oils. These can also be identified and quantified through their distinct A-TEEM fingerprints.





*A-TEEM

Absorbance-Transmission and Fluorescence Excitation Emission Matrix

Instant Water Quality Reports

The Datastream Dashboard is the industrial enhancement to the HORIBA Aqualog. With push button operation, the Aqualog Datastream Dashboard facilitates completely automated analysis and reporting of a wide range of organic matter parameters that are critical for managing and optimizing the drinking water treatment process. The Datastream parameters have been selected to specifically target disinfection by-product issues, organic membrane fouling, algal issues, and other contamination components.

The Dashboard provides the latest readings, time series and tables or trend analyses, % removals, thresholds and MCLs for all of these parameters. It reports on fit statistics and residuals evaluation for system performance monitoring, contamination detection and early warning alerts.

Water treatment laboratories can also upload their independently measured data into the HTML Dashboard to simultaneously analyze pH, alkalinity, turbidity, Cl_2 and other key parameters.

Expert installation and model calibration included.

Vater Treatment PI	iant Latest San	nple					
Parameter	Units	Update	Raw	Settled	Finished	% Removal	
Aqualog Data							
Aqualog_TOC		2017-11-03	2.7	1.2	1.1	59 %	
Aqualog_A254	Au/cm	2017-11-03	0.14	0.037	0.033	77 %	
Aqualog_SUVA	L/mg-m	2017-11-03	5.2	2.9 Treatment Plant Time Server	2.9	44 %	_
Aqualog_THMFP	ug/l	2017-11-03	1.5e+2	Contraction of the second s			
ulvic	RSU	2017-11-03	0.015	• 10			O Raw O Setter O Finith O Balter
numic	RSU	2017-11-03	0.016 200	00 F	•		ADDVE CIUME
гур	RSU	2017-11-03	0.0053			•	
121	Ratio	2017-11-03	1.0	Reveal Part			
gofit	%	2017-11-03	99	Read That		100 m	• •

Aqualog Datastream Dashboard is powered by Solo_Predictor software from Eigenvecto Research, Incorporated



Aqualog Method vs. Conventional Methods

Parameters for Monitoring DBPs Patented Aqualog A-TEEM Method **Conventional Methods Dissolved Organic Carbon (DOC)** Provides DOC composition information **TOC Analyzer** not available from TOC Single parameter, no information on organic composition or contaminations UVA@ 254nm (UV254) Aqualog measures both DOC and **UV-Vis Spectrometer** UV254, thus also providing SUVA DOC UV254 only (Patented method) Disinfection By-products (DBPs) No sample prep, less expensive, real-GC/GC-MS time analysis Lots of sample prep, much more expensive, not suitable for real-time Natural Organic Matter (NOM) Aqualog measures DOC, DBPs, and LC-OCD Precursors to DBPs organic foulants For research and chemical optimization, Less expensive, fast, and real-time more expensive, slow (~2 hrs/sample) analysis and not suitable for real-time

Date 00

We have been using the Aqualog instrument to monitor organic matters in seawater to better evaluate the membrane fouling of a reverse osmosis membrane at seawater desalination plants. Compared to TOC measurements, it provides a sensitive quantification, as well as properties of organic matters, to membrane fouling.

> Kwanghee Shin, Ph.D. Senior Research Engineer Corporate R&D Institute Doosan Heavy Industries & Construction

Our HORIBA Aqualog has provided Hazen and Sawyer and our clients with a state-of-the-art tool for managing multiple challenges in the production of drinking water. We have used the instrument for rapid response to the spill events and environmental events that may negatively impact community water supplies, providing a way to monitor for contamination in a timely fashion.

> Ben Stanford, Ph.D. Senior Director/Water R&D American Water (Formerly with Hazen and Sawyer)

Aqualog with Datastream Dashboard Software Specifications

Scan this QR code for more information



Dashboard Software Speci	d.me/aqualog		
Reported Parameters (up to 60 user selectable parameters)	Specifications	Notes	
Dissolved Organic Carbon Concentration (DOC)	30µg/l to 20mg/l	Requires filtration (0.45µm)	
A254	1cm path length		
SUVA	L DOC mg ⁻¹ A254m ⁻¹		
Simulated Distribution System, Trihalomethane Formation Potential (SDS THMFP)	10 to 500µg/l	US EPA MCL = 80µg/l	
Parallel Factor Analysis Component Scores	Up to 7	Can include algal, oil/PAH, tracer dyes, and other components	
Residuals (Q)		Detect contaminants and measurement issues	
% Variance Accounted for			
Absorbance Spectrum	200 - 800nm; 1cm path length	Any wavelength coordinate or ratio can be analyzed	
Excitation Emission Matrix (EEM) Regions	I-V plus Algal (Blue-green/Brown/Green)	Custom EEM regions also available; Based on interpolated EEM processing	
Total Fluorescence	Sum of EEM Regions I-V		
Humic Index			
Fluorescence Index			
Independent Treatment Plant Data (Additional Parameters Available)			
DOC	mg/l		
Alkalinity	mg/l		
Chlorine Residual	mg/l		
THM/SDSTHMFP	µg/l		
A254	Adjustable Path Length Specifications		
рН			
	32.72kg (72 lbs)	Aqualog	
Weight	9kg (19.8 lbs)	Sipper WS-10	
	4.2kg (9.26 lbs)	4 Sample Changer WS-10-S	
	LWH (618 x 435 x 336mm); (24 x 17 x 13")	Aqualog	
Dimensions	LWH (250 x 250 x 308mm); (9.84 x 9.84 x 12.13")	Sipper WS-10	
	LWH (250 x 125 x 308mm); (9.84 x 4.92 x 12.13")	4 Sample Changer WS-10-S	

Beyond Water Treatment Plants

Aqualog was designed for quantitative and predictive water analysis, and it is ideal for the task. However Aqualog, with its unique A-TEEM molecular fingerprinting benefits, has proven itself to be an invaluable tool in a wide variety of other industrial QC/QA, as well as environmental and academic research applications.

Please contact the Fluorescence Division of HORIBA Scientific to learn more about the ever growing list of exciting applications where Aqualog provides unique benefits.

Molecular Weight Natural Organic Matter Powdered Activated Carbon Simulated Distribution System

Trihalomethanes

Total Organic Carbon

Standard Operating Procedure(s) Specific Ultraviolet Absorbance

Trihalomethanes Formation Potential

Acronyms		
A-TEEM	Absorbance - Transmission and Fluorescence	MW
	Excitation and Emission Matrix	NOM
DBPs	Disinfection By-product(s)	PAC
DOC	Dissolved Organic Carbon	SDS
DWTP	Drinking Water Treatment Plant	SOP(s)
EEM	Excitation and Emission Matrix	SUVA
EPA	Environmental Protection Agency	THM
HAA	Haloacetic Acids	THMFP
MCL	Maximum Contaminant Level	TOC

www.aqualog.com info.sci@horiba.com

HORIBA Scientific has a policy of continuous product development, and reserves the right to amend part numbers, descriptions and specifications without prior notice.

The Future of Fluorescence

A-TEEM™

MOLECULAR

FINGERPRINTING

HORIBA

USA:	HORIBA Instruments Inc., 20 Knightsbridge Rd., Piscataway, NJ 08854 - Toll-free: +1-866-562-4698 - Tel: +1 732 494 8660 - Fax: +1 732 549 5125 - Email: info.sci@horiba.com
France:	HORIBA Jobin Yvon S.A.S., 16-18 rue du Canal, 91165 Longjumeau cedex - Tel: +33 (0)1 69 74 72 00 - Fax: +33 (0)1 69 09 07 21 - Email: info-sci.fr@horiba.com
Japan:	HORIBA Ltd., Tokyo Branch Office, 2-6, KandaAwaji-cho, Chiyoda-ku, Tokyo 101-0063, Japan - Tel: +81-(0)3 6206 4721 - Fax: +81 (0)3 6206 4730 - Email: info-sci.jp@horiba.com
Germany	: HORIBA Jobin Yvon GmbH, Neuhofstraße 9, 64625 Bensheim - Tel: +49 (0)6251 8475 0 - Fax: +49 (0)6251 8475 20 - Email: info-sci.de@horiba.com
Italy:	HORIBA Jobin Yvon Srl., Via Cesare Pavese 21, 20090 Opera (Milano) - Tel: +39 06 51 59 22 1 - Fax: +39 2 5760 0876 - Email: info-sci.it@horiba.com
UK:	HORIBA UK Ltd., Kyoto Close, Moulton Park, Northampton NN3 6FL - Tel: +44 (0)1604 542 500 - Fax: +44 (0)1604 542 699 - Email: info-sci.uk@horiba.com
China:	HORIBA (China) Trading Co. Ltd., Unit D 1F, Bldg A, Srynnex International Park, No. 1068 West Tianshan Road, Shanghai 200335 - Tel: +86 (0)21 6289 6060 - Fax: +86 (0)21 6289 5553
	Email: info-sci.cn@horiba.com
Brazil:	HORIBA Instruments Brasil Ltda., Rua Presbítero Plínio Alves de Souza, 645, Loteamento Polo Multivias, Bairro Medeiros, Jundiaí / SP, CEP 13.212-181 - Tel: +55 (0)11 2923 5400
	Fax: +55 (0)11 2923 5490 - Email: infocientifica.br@horiba.com
Other:	Tel: +1 732 494 8660 - Email: info.sci@horiba.com
Printed June	
Printed June	2018