

# Duetta

Fluorescence and Absorbance Spectrometer





# **A Superior Spectrometer Concept**

- 2-in-1 Simultaneous Fluorescence and Absorbance Spectrometer
- UV-Vis-NIR Fluorescence Detection Wavelength Range from 250 to 1,100 nm
- Full 3-D Fluorescence EEM Acquisition in Less Than One Second
- Best in Class Fluorescence Sensitivity Specification of 6,000:1 RMS for Water Raman
- Automatic Correction for Primary and Secondary Inner Filter Effects (IFE)
- High Fidelity Molecular Fingerprinting with Unique A-TEEM™
   (Absorbance-Transmittance Excitation Emission Matrix) Technology
- Millisecond CCD Detection of Entire Fluorescence Spectrum

### Fluorescence and Absorbance in one Spectrometer!

Duetta can be used as a fluorometer, as a UV-Vis-NIR spectrometer to measure absorbance, or as an instrument that measures true molecular fingerprints, which require the acquisition of fluorescence and absorbance, correcting for IFE in real time.

### Fluorescence in the Blink of an Eye!

Duetta's ultra-fast, state-of-the-art CCD technology allows it to outperform, by orders of magnitude, any competitor using PMTs. It's the fastest bench-top fluorometer in the world, and the only fluorometer that can acquire a full spectrum in less than 100 milliseconds! This state-of-the-art CCD technology also extends NIR wavelength detection to 1,100 nm, well beyond the limit of a standard PMT fluorometer.

### **Highest Fluorescence Sensitivity in its Class**

Having the highest level of fluorescence sensitivity means that Duetta's optimized design lets you measure lower concentrations of a sample, and provides better quality data.

### **EzSpec™ Touch-Screen Software**

HORIBA debuts Duetta's new generation software interface with touch-screen utility, and a selection of Apps for routine analysis.








The Future of Fluorescence



# **A New Spectrometer Design**



### Simultaneous Fluorescence and Absorbance

Duetta's two specialized sample detectors collect fluorescence at right angles to the excitation light, and absorbance straight through the sample.



Ergonomic, Hands Free, Sample Access

# 

### Absolute Fluorescence Spectra with Real Time Inner Filter Effect (IFE) Correction


IFE limits the linearity of the fluorescence signal at higher sample concentrations due to primary and secondary re-absorption of fluorescence. Duetta's unique design applies real time IFE correction to obtain absolute fluorescence spectra at higher sample concentrations.

### **CCDs Collect the Entire Emission Spectrum**

Duetta's CCD detection takes a snapshot of the entire fluorescence emission spectrum from 250 to 1,100 nm, so, not only is it much faster, but it also gives much more information than the competition.

### Extended NIR Wavelength Range up to 1,100 nm

Near Infrared applications are growing rapidly, and Duetta is ready to meet the challenge with fluorescence detection from 250 to 1,100 nm, well beyond the reach of standard PMTs.



NIR detection comparing Duetta with PMT.

## Easy Lamp Replacement

Lamp cartridges are pre-aligned at the factory, guaranteeing perfect alignment upon snap-in replacement. No tools or service are required.



### **Smart Sample Accessory Trays**


Switching from one type of sample holder to another could not be easier.

Without any tools or cables, Duetta's spill-proof sample accessory trays are identified by EzSpec software, enable electrical power for the appropriate module (e.g. 4-position Peltier holder), and read an enabled sensor (e.g. temperature), if available.

### **Fully Automated Order Sorting Filters**


Duetta includes excitation and emission cut-off filters for automated elimination of second order grating effects, which otherwise interfere with the fluorescence signal of interest.


### **EEMs Acquired in Seconds**



Excitation Emission Matrix (EEM) acquired in 1 second, shows no second order scattering.










Pyrene in a solution of A-B-A triblock copolymer micelles is solvated in the hydrophobic core solvent environment. 3D EEMs can be shown in contour, waterfall plot, or 2D overlay mode.

BSA protein measured from 5 °C to 60 °C and the peak intensity and peak wavelength are plotted vs. temperature of the protein solution. Measured using the SampleSnap-4Pelt 4-position Peltier temperature-controlled cuvette holder.









Polymer Aggregation and Structure. Temperature dependent spectra showing polymer micellization. A. Kinetic spectral scans of native BSA protein with rapidly added ANS (3x10<sup>-6</sup> M). Spectra taken every 100 ms after ANS addition. As ANS binds to native BSA, the BSA Trp emission decreases and the ANS emission increases as a result of FRET, due to proximity of excited Trp and ANS.

B. Kinetic spectral scans of SDS-denatured BSA after addition of ANS. Much higher ANS concentration (4x10<sup>-5</sup> M) is required to affect binding to denatured BSA; no FRET is observed due to increased distances between ANS and Trp.

Emission spectra of a series of glass materials containing different lanthanides. Measured using the SampleSnap-Uni solid sample holder.






Comparison of bovine serum albumin protein and the same protein denatured by sodium dodecylsulfate (SDS). The emission peak narrows and shifts 15 nm to the blue.

Quenching: Fluorescence monitoring of the titration of 1 M potassium iodide (KI) into a solution of 2-aminopyridine, a common small molecule drug pre-cursor.

# **EzSpec Software**

### **Apps at your Fingertips**

From the EzSpec<sup>™</sup> launch screen, each App icon intuitively helps you identify and select the analysis you want to perform with a simple point and touch.



### **Standard Apps:**

- Build/use a Concentration Curve
- Protein A280
- DNA/RNA Purity (A260/A280)
- Pass/Fail

### **Upcoming Apps:**

- Autonomous Acquisition
- Concentration from Absorbance
- FRFT
- Ion Kinetics
- CIE Colorimetry
- Enzyme Kinetics

EzSpec is now translated into 7 languages: Chinese, English, French, Japanese, Portuguese, Russian, and Spanish.

EzSpec makes your most recent methods and data easily accessible.



### **Widest Range of Spectroscopic Acquisition Modes:**

- Fluorescence emission spectra plus Abs & IFE\* \*\*
- Fluorescence excitation spectra plus Abs & IFE\* \*\*
- Capture fluorescence value plus Abs & IFE\* \*\*
- Fluorescence EEM plus Abs & IFE\* \*\*
- Fluorescence kinetics (emission spectra\* \*\* or single point intensities)
- Absorbance and % Transmittance (spectra or kinetics)
- \* Unique to Duetta \*\*Anisotropy available for any fluorescence mode on Duetta-Bio



### **Acquire Menu:**



### **PDF Report Generator and Data Export:**

- Includes data and results of fitting/processing
- Instrument and user identification information
- Optional and customizable signature fields
- Full data and method traceability
- Data export to ASCII or SPC

- Full method details
- 2D and 3D graphs
- Customizable logo
- Graphical image export

# **Accessories**

# Optional Sample Holder Trays

### SampleSnap-WJ: Thermostated and Micro-Stirred Cuvette Sample Tray

Self-identifying thermostatable, water jacketed cuvette holder (1 cm) with magnetic micro-stir motor. The sample tray can be installed in seconds without the need of a tool, and provides power for micro-stirrer. Requires, but does not include, external water bath.

### SampleSnap-1Pelt: Single Cuvette Peltier Holder Sample Tray

Self-identifying single cuvette Peltier holder (1 cm) for automated temperature control and ramping from -15 °C to +105 °C, with magnetic micro-stirrer. Sample tray can be installed in seconds without the need of a tool, and provides power and control of the Peltier device, and reads temperature. Includes external water circulator for heat sink.

### SampleSnap-4Pelt: Four Position Peltier Cuvette Sample Tray

Self-identifying four position Peltier holder (1 cm) for automated temperature control and ramping from -15 °C to +110 °C, with magnetic micro-stirrers under all four positions. Sample tray can be installed in seconds without the need of a tool, and provides power and control of the Peltier device, and reads temperature. Includes external water circulator for heat sink.

### SampleSnap-UNI: Universal Variable Angle Sample Tray

Self-identifying sample holder with 360 degree variable angle goniometer for empirical determination of optimum sample angle. Sample tray can be installed in seconds without need of a tool. Requires, but does not include, solid sample holder accessory (SampleSnap-SS), powder sample holder accessory (SampleSnap-PWD) or cuvette holder accessory (SampleSnap-VACVT).

### Sample Accessories for SampleSnap-UNI: Universal Variable Angle Sample Tray

### SampleSnap-SS: Solid Sample Accessory

Solid sample accessory for universal sample holder to be used with flat or solid samples such as thin films, cover slips, paper, gemstones, and more. Requires, but does not include, SampleSnap-UNI.

### SampleSnap-PWD: Powder Sample Accessory

Powder sample accessory for universal sample holder. Requires, but does not include, SampleSnap-UNI.

### • SampleSnap-VACVT: Variable Angle Cuvette Sample Accessory

Cuvette sample accessory for universal sample holder. Requires, but does not include, SampleSnap-UNI.

### Additional Accessories

K-210: Six pack of quartz windows for powdered sample holder

K-151: 10 mm quartz cuvette

LampSnap: Replacement Duetta light source assembly with pre-aligned Xenon arc lamp

K-AFILTER: Replacement Duetta air intake filter

\* Duetta includes a standard room temperature single cuvette sample holder tray.



Duetta **Duetta-Bio Acquisition Modes** Absorbance and Fluorescence Spectrometer Fluorescence Sensitivity Water Raman SNR >6,000:1 RMS, 350 nm excitation, 5 nm slits **Spectral Acquisition Rate** 510.000 nm/min **EEM Acquisition Rate** As fast as 1 second (sample and wavelength dependent) **A-TEEM Acquisition Rate** As fast as 30 seconds (sample and wavelength dependent) Fluorescence Detector CCD / Spectrograph Fluorescence Detector Range 250 to 1.100 nm Fluorescence Bandwidth 1, 2, 3, 5, 10, 20 nm (excitation and emission) **Light Source** 75 W Xenon arc lamp. Dedicated cartridge for snap-in replacement **Excitation/Absorbance Wavelength Range** 250 to 1.000 nm **Absorbance Detector** Silicon Photodiode 250 to 1.000 nm **Absorbance Detector Range Absorbance Bandwidth** 1, 2, 3, 5, 10, 20 nm **Absorbance Range** 0 to 2 A **Absorbance Accuracy** +/- 0.02 A +/- 1 nm **Wavelength Accuracy Computer Requirements** Windows® 7/8/10 compatible (Windows10 recommended), PC with 1 USB minimum Software EzSpec Dimensions (W x D x H) 17 in × 20 in × 15.4 inches; 43 cm x 51 cm x 39 cm

45 pounds; 20.4 kg

Yes (280 to 750 nm)

Not Available

Weight

Polarizers (Time of Purchase)

# FRET PLQY FRET PLQY FRET PLQY FURA-2 PV LRET FLIM Our Roots Grow Deep Same Normal Socientific Scientific FLIM SAM Normal Socientific 1973 John Von Open USA Office Socientific SPEX FLUORESCENCE 1973 John Von Spex FLUORESCENCE 1974 John Von Spex FLUORESCENCE 1975 John Von Spex FLUORESCENCE 1975 John Von Spex FLUORESCENCE 1975 John Von Spex FLUORESCENCE 1976 John Von Spex FLUORESCENCE 1977 John Von John

HORIBA has a long history of global leadership in high performance fluorescence research instrumentation.

Now, with the introduction of Duetta, HORIBA has perfected an affordable spectrofluorometer that takes advantage of decades of experience in fluorescence. (Think SPEX, Jobin Yvon, SLM, IBH, PTI and HORIBA).



### http://www.horiba.com e-mail: info.sci@horiba.com

HORIBA, Ltd. Japan Head Office

2 Miyanohigashi, Kisshoin, Minami-ku, Kyoto, Japan Phone: 81 (75) 313-8123 Fax: 81 (75) 321-5725

HORIBA (China) Trading Co., Ltd.

Unit D, 1F, Building A, Synnex International Park, 1068 West Tianshan Road, Shanghai, 200335, China Phone: 86 (21) 6289-6060 Fax: 86 (21) 6289-5553

Beiiing Branch

12F, Metropolis Tower, No.2, Haidian Dong 3 Street, Beijing, 100080, China

Phone: 86 (10) 8567-9966 Fax: 86 (10) 8567-9066

Guangzhou Branch

Room 1611 / 1612, Goldlion Digital Network Center, 138 Tiyu Road East, Guangzhou, 510620, China Phone: 86 (20) 3878-1883 Fax: 86 (20) 3878-1810

HORIBA Instruments (Singapore) Pte Ltd. Singapore

3 Changi Business Park Vista #01-01, Akzonobel House, Singapore 486051

Phone: 65 (6) 745-8300 Fax: 65 (6) 745-8155

HORIBA KOREA Ltd.

Seoul Branch

China

10, Dogok-Ro, 6-Gil, Gangnam-Gu, Seoul-Si, 06259, Korea Phone: 82 (2) 753-7911 Fax: 82 (2) 756-4972

**HORIBA India Private Limited** 

246, Okhla Industrial Estate, Phase 3 New Delhi-110020, India Phone: 91 (11) 4646-5000 Fax: 91 (11) 4646-5020

**Bangalore Office** 

No.55, 12th Main, Behind BDA Complex, 6th sector, HSR Layout, Bangalore South, Bangalore-560102, India Phone: 91 (80) 4127-3637

HORIBA Vietnam Co., Ltd.

Unit 6, 10 Floor, CMC Tower, Duy Tan Street, Dich Vong Hau Ward, Cau Giay District, Hanoi, Vietnam Phone: 84 (24) 3795-8552 Fax: 84 (24) 3795-8553

PT HORIBA Indonesia

JI. Jalur Sutera Blok 20A, No.16-17, Kel. Kunciran, Kec. Pinang Tangerang-15144, Indonesia Phone: 62 (21) 3044-8525 Fax: 62 (21) 3044-8521

nc Printed in USA, Duetta, 4/19

HORIBA Instruments Incorporated USA HO

Brazil

9755 Research Drive, Irvine, CA 92618, U.S.A. Phone: 1 (949) 250-4811 Fax: 1 (949) 250-0924

**Piscataway Office** 

20 Knightsbridge Rd., Piscataway, NJ 08854, U.S.A. Phone: 1 (732) 494-8660 Fax: 1 (732) 549-5125

HORIBA Instruments Brasil, Ltda.

Rua:Presbitero Plinio Alves de Souza, 645, Loteamento Polo Multivias Barirro Medeiros-Jundiai Sao Paulo-CEP 13.212-181 Brazil

Phone: 55 (11) 2923-5400 Fax: 55 (11) 2923-5490

HORIBA (Thailand) Limited Thailand

393, 395, 397, 399, 401, 403 Latya Road, Somdetchaopraya, Klongsan, Bangkok, 10600, Thailand Phone: 66 (0) 2-861-5995 ext.123 Fax: 66 (0) 2-861-5200

East Office 850 / 7 Soi Lat Krabang 30 / 5, Lat Krabang Road, Lat Krabang,

890 / 7 Sol Lat Mabang 30 / 5, Lat Mabang Hoad, Lat Mabang Bangkok 10520, Thailand Phone: 66 (0) 2-734-4434 Fax: 66 (0) 2-734-4438 HORIBA FRANCE SAS

16-18, rue du Canal, 91165, Longjumeau Cedex, France Phone: 33 (1) 69-74-72-00 Fax: 33 (1) 69-09-07-21

HORIBA Jobin Yvon GmbH Germany

Neuhofstrasse 9, D\_64625, Bensheim, Germany Phone: 49 (89) 6251 84 75-0 Fax: 49 (89) 6251 9475-20

HORIBA ITALIA SRL Italy

France

Via Luca Gaurico 209-00143, Roma, Italy Phone: 39 (6) 51-59-22-1 Fax: 39 (6) 51-96-43-34

HORIBA UK Limited UK

Northampton Office

Kyoto Close Moulton Park, Northampton NN3 6FL, UK Phone: 44 (1604) 542-500 Fax: 44 (1604) 542-699

Duetta, EzSpec and A-TEEM are trademarks of HORIBA Instruments Inc. Windows is a registered trademark of Microsoft, Inc.

HORIBA Scientific has a policy of continuous product development, and reserves the right to amend part numbers, descriptions and specifications without prior notice.



Korea

India

Vietnam

Indonesia