


Raman spectroscopy produces characteristic spectral bands that correspond to the vibrational modes of molecules within a sample. These bands provide valuable information about molecular composition and structure.

While the exact positions of these bands can vary depending on factors like sample composition and environment, below are some of the more common Raman band assignments.

### Common Raman Bands



Functional Group/Vibration	Region	Raman
Lattice vibrations in crystals, LA modes	10 - 200 cm <sup>-1</sup>	strong
δ(CC) aliphatic chains	250 - 400 cm <sup>-1</sup>	strong
ν(Se-Se)	290 - 330 cm <sup>-1</sup>	strong
ν(S-S)	430 - 550 cm <sup>-1</sup>	strong
ν(Si-O-Si)	450 - 550 cm <sup>-1</sup>	strong
ν(Xmetal-O)	150 - 450 cm <sup>-1</sup>	strong
ν(C-I)	480 - 660 cm <sup>-1</sup>	strong
ν(C-Br)	500 - 700 cm <sup>-1</sup>	strong
ν(C-Cl)	550 - 800 cm <sup>-1</sup>	strong
ν(C-S) aliphatic	630 - 790 cm <sup>-1</sup>	strong
ν(C-S) aromatic	1080 - 1100 cm <sup>-1</sup>	strong
ν(O-O)	845 - 900 cm <sup>-1</sup>	strong
ν(C-O-C)	800 - 970 cm <sup>-1</sup>	strong
ν(C-O-C) asym	1060 - 1150 cm <sup>-1</sup>	strong

$\nu(\text{CC})$ alicyclic, aliphatic chain vibrations	600 - 1300 $\text{cm}^{-1}$	medium	medium
$\nu(\text{C}=\text{S})$	1000 - 1250 $\text{cm}^{-1}$	strong	weak
$\nu(\text{CC})$ aromatic ring chain vibrations	*1580, 1600 $\text{cm}^{-1}$	strong	medium
	*1450, 1500 $\text{cm}^{-1}$	medium	medium
	1000 $\text{cm}^{-1}$	strong/medium	weak
$\delta(\text{CH}_3)$	1380 $\text{cm}^{-1}$	medium	strong
$\delta(\text{CH}_2)$ $\delta(\text{CH}_3)$ asym	1400 - 1470 $\text{cm}^{-1}$	medium	medium
$\delta(\text{CH}_2)$ $\delta(\text{CH}_3)$ asym	1400 - 1470 $\text{cm}^{-1}$	medium	medium
$\nu(\text{C}-\text{NO}_2)$	1340 - 1380 $\text{cm}^{-1}$	strong	medium
$\nu(\text{C}-\text{NO}_2)$ asym	1530 - 1590 $\text{cm}^{-1}$	medium	strong
$\nu(\text{N}=\text{N})$ aromatic	1410 - 1440 $\text{cm}^{-1}$	medium	-
$\nu(\text{N}=\text{N})$ aliphatic	1550 - 1580 $\text{cm}^{-1}$	medium	-
$\delta(\text{H}_2\text{O})$	~1640 $\text{cm}^{-1}$	weak broad	strong
$\nu(\text{C}=\text{N})$	1610 - 1680 $\text{cm}^{-1}$	strong	medium
$\nu(\text{C}=\text{C})$	1500 - 1900 $\text{cm}^{-1}$	strong	weak
$\nu(\text{C}=\text{O})$	1680 - 1820 $\text{cm}^{-1}$	medium	strong
$\nu(\text{C}\equiv\text{C})$	2100 - 2250 $\text{cm}^{-1}$	strong	weak
$\nu(\text{C}\equiv\text{N})$	2220 - 2255 $\text{cm}^{-1}$	medium	strong
$\nu(-\text{S}-\text{H})$	2550 - 2600 $\text{cm}^{-1}$	strong	weak
$\nu(\text{C}-\text{H})$	2800 - 3000 $\text{cm}^{-1}$	strong	strong
$\nu(=\text{C}-\text{H})$	3000 - 3100 $\text{cm}^{-1}$	strong	medium
$\nu(\equiv\text{C}-\text{H})$	3300 $\text{cm}^{-1}$	weak	strong
$\nu(\text{N}-\text{H})$	3300 - 3500 $\text{cm}^{-1}$	medium	medium
$\nu(\text{O}-\text{H})$	3100 - 3650 $\text{cm}^{-1}$	weak	strong

It is important to note that Raman bands can vary depending on factors such as molecular symmetry, conjugation, and the presence of functional groups. Interpretation of Raman spectra often involves comparing observed bands to reference data and considering the specific molecular context.

#### Reference

Socrates, G. (2001). Infrared and Raman Characteristic Group Frequencies, Third Edition. John Wiley and Sons.

[info.sci@horiba.com](mailto:info.sci@horiba.com)

**USA:** +1 732 494 8660  
**UK:** +44 (0)1604 542 500  
**China:** +86 (0)21 6289 6060  
**Taiwan:** +886 3 5600606

[ramanacademy.com](http://ramanacademy.com)

**France:** +33 (0)1 69 74 72 00  
**Italy:** +39 06 51 59 22 1  
**India:** +91 80 41273637  
**Brazil:** +55 (0)11 2923 5400

[horiba.com/raman](http://horiba.com/raman)

**Germany:** +49 (0) 6251 8475 0  
**Japan:** +81(75)313-8121  
**Singapore:** +65 (0)6 745 8300  
**Other:** +33 (0)1 69 74 72 00