



## Optical Chopper

ACH-C, ACH-C-OPEN

ELEMENTAL ANALYSIS
FLUORESCENCE
GRATINGS & OEM SPECTROMETERS
OPTICAL COMPONENTS
FORENSICS
PARTICLE CHARACTERIZATION
RAMAN
SPECTROSCOPIC ELLIPSOMETRY
SPR IMAGING

### Direct-digital-synthesis optical chopper

The optical chopper from HORIBA Scientific is a microprocessor-based control system that utilizes direct-digital-synthesis to deliver precise optical chopping rates from 4 Hz to 5000 Hz. An external controller with large five-digit LED readout is included to enable setting of the chopping rate from the front panel. Utilizing a phase-locked-loop system, the chopping rate may also be synchronized to a user-supplied external clock.

The precision etched blade may be mounted directly to the spectrometer entrance and protected from inadvertent damage in the enclosed housing (ACH-C). A 3-slot blade and 30-slot blade are available for low and high repetition rates. An open variety (ACH-C-OPEN) is also available for mounting to an optical table in a free space configuration. This variety uses a single blade with 30 slots on the outer ring and 3 slots on the inner ring, eliminating the need to swap blades.


### Features and Benefits

- External controller- large five-digit readout enables digital entry of the desired chopping rate directly from the front panel
- Superior speed and accuracy- 4 Hz to 5000 Hz speed with accuracy of 0.001 Hz
- Simplicity of operation- set the repetition rate with the internal clock via the front panel on the controller or synchronize to an external clock

### Control

The optical chopper can be controlled locally via the settings on the front panel of the controller or remotely via a bidirectional RS-232 port.


The front panel includes controls for switching between high frequency (40 Hz to 5000 Hz) and low frequency (5 Hz to 500 Hz) with a manual dial to adjust the frequency in variable size increments. Once the frequency is set, a front panel BNC connector (Sync Output) can be connected to a lock-in amplifier, oscilloscope, or other instrument for synchronization.



## Specifications

Part number	ACH-C	ACH-C-OPEN
Frequency range	<b>3-slot blade, low speed:</b> 4 – 500 Hz, <b>30-slot blade, high speed:</b> 40 – 5000 Hz	
Frequency accuracy	±0.0025% of setting	
Frequency stability (internal clock)	< 10 ppm/°C	
Frequency-setting resolution	<b>High-speed range:</b> 0.1 Hz, <b>Low-speed range:</b> 1 Hz	
Phase jitter	<b>3-slot blade, low speed</b> 0.1% peak-to-peak, <b>30-slot blade, high speed</b> 1% peak-to-peak	
Settling time to lock	<b>Full-scale change:</b> < 3 s, <b>10% change:</b> < 1 s	
External frequency input	TTL/CMOS compatible, 4 – 5000 Hz	
Display	5-digit LED, green (565 nm), 14 mm (0.56")	
RS-232 port	9600 baud, N-8-1, three-wire	
Aperture size	<b>30-blade:</b> 15 mm high x 4.5 mm wide at center, <b>3-blade:</b> 15 mm high x 15 mm wide at center	
Operating temperature	0 – 40°C	
Power requirement	95 – 260 V AC, 50 – 60 Hz, 15 VA max.	
Controller dimensions (H x W x D)	69 mm (2.7") x 178 mm (7") x 231 mm (9.1")	
Head dimensions (H x W x D)	114.3 mm (4.5") x 114.3 mm (4.5") x 51 mm (2")	
Weight	1.36 kg (3 lbs)	

## Mechanical Diagram



**OPTICAL BUILDING BLOCKS**



[info.sci@horiba.com](mailto:info.sci@horiba.com) [www.horiba.com/opticalbuildingblocks](http://www.horiba.com/opticalbuildingblocks)

**HORIBA**  
Scientific

**USA:** +1 732 494 8660  
**UK:** +44 (0)20 8204 8142  
**China:** +86 (0)21 6289 6060

**France:** +33 (0)1 69 74 72 00  
**Italy:** +39 2 5760 3050  
**Brazil:** +55 (0)11 5545 1500

**Germany:** +49 (0)89 4623 17-0  
**Japan:** +81 (0)3 6206 4721  
**Other:** +1 732 494 8660