




## Tunable KiloArc™ Output Curves

You can customize a Tunable KiloArc™ Illuminator to suit your specific needs. The lamp and grating selected primarily determine the performance of the illuminator. However for any given lamp and grating the slit adjustment (or bandwidth) selected will also affect the intensity output of the unit.

Below are a series of intensity output curves for different configurations of the Tunable KiloArc™ Illuminators. There are far too many variables to provide output curves for all possible combinations for these illuminators. If you are interested in an illuminator with a different grating than those shown below, then you can compare the grating curves of that grating with one listed below to get a pretty good idea of how your particular illuminator will perform.



## Optical Performance Specifications

<b>Optical Power</b>	> 300 m W (grating, bandpass & wavelength dependent)
<b>Spot Size at Slit Exit</b>	10 mm (slit dependent)
<b>Diverging Beam angle (full)</b>	14.4 degrees
<b>Numerical Aperture (N.A.)</b>	0.12
<b>Short Term Optical Noise*</b>	from 0.15 to 0.2% RMS

\* 1,000 points/s, 1 s duration, 1.5 KHz detector bandwidth

## KiloArc™ Other Specifications

<b>Input</b>	210–240 V AC 50/60 Hz
<b>Starting</b>	45 kV starting pulse
<b>Power Rating</b>	800–1200 watts (adjustable) — recommended 800–1000 watts
<b>Lamp Module Type</b>	1000 W Xenon, 1000 W Mercury/Xenon (proprietary to OBB)
<b>Lamp Life</b>	Typically 1,500 hrs
<b>Focusing Optics</b>	High efficiency f/4 ellipsoid reflector
<b>Power Precision</b>	0.04% (0.4 watts)
<b>Output Volts Compliance</b>	17–23 VDC
<b>Output Current Limit</b>	70 A rms
<b>Height</b>	329 mm (12.9 inches)
<b>Width</b>	375 mm (14.8 inches)
<b>Length</b>	489 mm (19.3 inches)
<b>Weight</b>	31 kg (68 pounds)
<b>Window Diameter (D)</b>	127 mm (5.0 inches)
<b>Center Beam Line Height (without feet)</b>	128 mm (5.0 inches)

## Monochromator Dimensions

<b>Length</b>	241 mm long (9.5 inches)
<b>Width</b>	272 mm wide (10.7 inches)
<b>Height</b>	115 mm high (4.5 inches)
<b>Weight</b>	7.3 kg (16 lbs)

## Monochromator Digital Slit Assembly Dimensions

<b>Length</b>	10 mm thick
<b>Width</b>	85 mm wide
<b>Height</b>	180 mm high of which 113 mm (4.4 inches) is above the lid of the monochromator

## Optional Motor Controller Specifications

<b>Power Supply</b>	Universal power supply included
<b>TTL Output</b>	Synchronization TTL output each time motor stops
<b>Stepper Motor</b>	Two phase motor 1 A per phase, 200 steps per revolution, 1.8 degree per step
<b>Maximum Motor Speed</b>	1200 RPM with zero torque
<b>Maximum Speed with Mono</b>	15,000 nm per minute
<b>Stepping Motor Voltage</b>	5–12 V
<b>Stepping Modes</b>	Full, Half and Micro steps: 1/8, 1/16, 1/32, 1/64 computer selectable
<b>Slew Rate</b>	1 to 62,500 micro steps per second
<b>Calibration</b>	Auto calibration of wavelength
<b>Ramping</b>	Linear ramping rate for heavy duty, fast, precision operation


## Higher Resolution Gratings Can Provide Higher Intensities

The data curves above are all from illuminators that had a 1,200 g/mm grating. You will note from these curves that doubling the bandpass, or slit size, typically results in a factor of four increase in intensity. This is true as long as the slit size is the same or equal to the optical spot size. Therefore if you used a 2,400 g/mm grating you would be doubling the slit size from a 1,200 g/mm grating to maintain an equivalent bandpass. Thus a Tunable KiloArc™ Illuminator equipped with a xenon lamp and a 2,400 g/mm grating blazed at 500 nm will give you up to 1 watt of optical power in a 20 nm bandpass. However the 2,400 g/mm grating will only mechanically scan up to 600 nm, see monochromator for more details.

## Mounting Information


The KiloArc™ can be mounted rigidly by removing the adjustable feet that come with the illuminator and exposing the four tap holes (6-32) shown on the bottom view below. There are also four mounting holes (6-32) on the front face of the KiloArc™, as shown on the front view below. These front holes are for attaching OBB's adapter tube but they can be used for any other mechanical assembly.

### Front View




# Mounting Information Continued


## Left View



## Right View



## Bottom View



## Top View

