



Minimize sample and dispersant volume requirements with all the convenience of a fully-automated circulation system.

Unique Features

- Minimizes sample dispersant volume requirements.
- Permits use of any organic solvent.
- Built-in ultrasonic probe available for dispersing agglomerated materials.
- Auto-fill pump and auto-drain are standard features.
- Normal operation from the PC software.
- LA-950 Method Expert software compatible.

MiniFlow

Small Volume Circulation System

Quick, Accurate Operation:

The navigation function documents the automatic measurement procedure.

High Performance without Compromise:

The software controls all sequences such as fill, rinse and drain just like the normal flow system. No manual operation is required! Exchange the MiniFlow system with the regular flow cell in seconds, without any tools.

Easy Maintenance:

The measurement cell is easily removed for cleaning.



partica LA-950 Accessory “MiniFlow”

Potential Applications

- Limited or precious sample applications.
- Any sample requiring organic solvent dispersion.
- Oil dispersions of pigments, polymers and resin products.
- Sample handling systems to minimize operator exposure.
- Minimize hazardous solvent use to limit waste disposal.

Specifications

Measurement range: 0.01~1000 μm

Measurement time: Approximately 60 seconds from filling & measurement to rinsing.

Measurement method: Circulation system with the flow cell.

Sample amount: Approximately 0.2mg~1.0g (Depending on the sample's particle size and density).

Dispersing liquid volume: 35~55 ml

Power supply: 120V AC+10% or 230V AC+10%, AC 50/60Hz+1%, 120VA

Circulation System Details:


- **Ultrasonic probe:** 40kHz, 7 steps of power control, Titanium alloy.
- **Centrifugal pump:** 15 steps of speed control
- **Auto-fill pump:** Teflon diaphragm pump
- **Other materials contact in fluids:** Tempax glass, SUS316, Titanium alloy (Ultrasonic probe), Teflon, Perfluoro Elastomers

Operational environment condition

Temperature: 15~35°C

Relative humidity: 85% or less, no condensation.

External dimensions



Email: labinfo@horiba.com **Website:** www.horiba.com/particle

USA: 1-800-446-7422 **France:** +33 (0)1 64 54 13 00 **Japan:** +81 (0)3 38618231

© 2017 HORIBA Instruments, Inc.
For further information on this document
or our products, please contact us.