

New Nanoparticle Analyzer nanoPartica SZ-100V2 series



Unravel the nano-universe with HORIBA's highest level of nanoparticle analysis*

This new model from nanoPartica series with even high sensitivity consolidates measurement of three major elements that characterize nanoparticles into a single unit: particle size, zeta potential and molecular weight!


*Compared to conventional HORIBA products

High-power Laser Lineup (100mW)


Measurement of Dilute Samples

A new high-power laser in addition to dual optics improved the measurement sensitivity to low-concentration samples, which is about 15 times more powerful than the previous model (SZ-100).

This enables highly accurate and reproducible measurement, even for dilute samples or samples with weak scattered light intensity.




Standard polystyrene polymer particles



Measurement of 20 nm standard polystyrene polymer particles diluted to 0.2 ppm

Measurement of single nanoparticles

Equipped with HORIBA's unique high-precision and high-speed correlator and low stray-light 90° optics to enable highly accurate measurement of single nanoparticles.



Diameter measurement of a 2 nm gold colloidal particle* with a 100 mW high-output laser


*Sample: Provided by Dr. Tai, National Institute of Advanced Industrial Science and Technology. TEM diameter: 1.8±0.3 nm

Wide Range of Applications

Measurement of NIST SRM 1980 α -FeO O H Zeta Potential Measurement Results

The SZ-100 measures zeta potential of particles using the iontophoresis laser doppler method so that both average value and the zeta potential distribution can be obtained.


Additionally, using the pH controller enables easy detection of the isoelectric point.



Sample concentration: 50ppm, pH = 2.5, mobility (standard): $2.53 \pm 0.12 \mu\text{m} \cdot \text{cm}/\text{V} \cdot \text{s}$
Measured result: Mobility = $2.53 \mu\text{m} \cdot \text{cm}/\text{V} \cdot \text{s}$, zeta potential = 32.9mV

Measurement of Antibody Pharmaceutical Temperature Control

Chronologically measured the rate of oligomerization (aggregation) in immunoglobulins G (IgG) at 60°C.



*Measured in cooperation with Tsumoto Laboratory, (University of Tokyo)

Expanding Application in a Wide Range of Fields

Functional Nanomaterials	Polymers	Bio/Life Sciences	Semiconductor	Ceramics	Environment and Agriculture	Gel Materials
Metal Nano Colloidal Particles Catalysts Carbon Nanotubes	Cellulose Nanofibers Electrolytes Adhesives	Antibody Pharmaceuticals Nano-capsules Viruses Proteins	CMP Slurries	Titanium Oxide (Titania) Silica Aluminum Oxide (Alumina)	Fine bubbles	3D Printer Materials Medical Materials (Artificial Cartilage, etc.) Automobile Materials (Sealants, Vibration Proofing, etc.)

nanoPartica SZ-100V2 series: lineup and main specifications

	Type	Particle size	Zeta potential	Molecular weight	High-power laser
SZ-100	S2	○		○	
	Z2	○	○	○	
	HS2	○		○	○
	HZ2	○	○	○	○

Measurement principle	Particle size: Photon correlation Zeta potential: Laser doppler electrophoresis Molecular weight: Static light scattering Debye plotting
Measurement range	Particle size: 0.3 nm – 10 μm (diameter) Zeta potential: -500 mV – +500 mV (Particle : 2.0 nm - 100 μm*1) Molecular weight: 1000 - 2 x 10 ⁷ (Debye plotting), 540 - 2 x 10 ⁷ (MHS *2)
Particle size measurement angle	90° and 173° (depending on sample concentration)
Sample cells	Cuvette cells (particle size, molecular weight), cells with electrodes (zeta potential)
Measurement optics	Light source: 532 nm 10 mW or 100 mW semiconductor excitation solid-state laser Detector: Photomultiplier tube (PMT)
Set/adjustable holder temperature	0 - 90°C (up to 70°C for plastic cells and cells with electrodes)
Laser class	1

*1: Depends on sample

*2: Calculated with the Mark-Howink-Sakurada equation (depends on sample)

Options and Accessories

- pH Controller
- Various cells
- 21CFR part11 software
- IQ/OQ/PQ compatibility

Customized Accessories

- Flat surface zeta-potential measurement cells
- Fluorescence Filter

(It is necessary to confirm the specification to HORIBA group)



Two kind of zeta potential measurement cells
(For zeta potential and particle size measurement, 100 μL,
Aqueous sample and Organic solvent)

Various sample cells

HORIBA
Scientific

Bulletin: HRE-3948 (E)

labinfo@horiba.com

Americas: +1 (800) 446-7422
China: +86 (21) 6289 6060
France: +33 (1) 69 74 72 00
Germany: +49 (89) 62 51 84 750

India: +91 (80) 4127 3637
Japan: +81 (75) 313-8121
Korea: +82 (31) 296 7911
Singapore: +65 (6) 745-8300

www.horiba.com/particle

Taiwan: +886 (3) 560-0606
Thailand: +66 (0) 2 861 5995
Vietnam: +84 (24) 3795 8552
Other: +81 (75) 313-8121