# HORIBA Explore the future



HORIBA Instruments Particle Analysis

Jeffrey Bodycomb, Ph.D.

# Getting a Charge Out of Nanoparticles: Zeta Potential

May 3, 2018



# **Controlling Suspensions**

Suspension behavior is significantly affected by surface interactions. And that is controlled by size and chemistry of the interface.





# **Small size -> Large surface area**





# **Controlling Suspensions**

Suspension behavior is significantly affected by surface interactions. And that is controlled by size and chemistry of the interface.





### **Fine Suspended Particles**

# Fine particles will tend to flocculate to reduce surface energy.





# Why don't particles flocculate?

## Most particles in aqueous suspension have a surface charge and therefore repel each other; they never touch.



# Charge often determines suspension behavior!



### What is Zeta Potential?

 Zeta potential is the charge on a particle at the shear plane.





# **Effect of Liquid**

Charge at shear plane depends on liquid environment. The red and blue lines here corresponds to different salt concentration.





#### **Quiz Time**

Your colleague hands you a bucket of dry particles and asks "What is the zeta potential of these particles?"

Before determining zeta potential of a particle, what do you need to know?





#### **Quiz Time**

Your colleague hands you a bucket of dry particles and asks "What is the zeta potential of these particles?"

Before determining zeta potential of a particle, what do you need to know?



#### A. Liquid environment (ionic strength, pH, nature of ions)





## **How do Surfaces Acquire Charge?**

### **Ionization of surface groups**





### **How do Surfaces Acquire Charge?**





#### **How do Surfaces Acquire Charge?**

# Specific adsorption of ions, e.g. ionic surfactants

 $-SO_3^- Na^+$  Anionic surfactant (SDS)

Na<sup>+</sup>





# **Why Zeta Potential?**

- Good way of evaluating <u>electrostatic</u> stabilization of suspensions
- Can use to predict interactions





#### **How to Measure Zeta Potential:**

- Acoustic techniques (use sound to probe particle response)
- It is much more popular to use <u>light</u> <u>scattering</u> to probe motion of particles due to an applied electric field. This technique is known as electrophoretic light scattering.
- Used for determining electrophoretic mobility, zeta potential.



# **Other Light Scattering Techniques**

 Static Light Scattering: over a duration of ~1 second. Used for determining particle size (diameters greater than 10 nm), polymer molecular weight, 2<sup>nd</sup> virial coefficient, R<sub>g</sub>.

 Dynamic Light Scattering: use scattered light to probe random motion of particles. If this motion is due to Brownian motion the technique can be used to determine particle size.



### How to Measure? With the SZ-100

 Single compact unit that performs size, zeta potential, and molecular weight measurements: the SZ-100





#### How to determine zeta potential

- Apply an electric field and probe response of particles to applied field.
- You need to see Doppler shift in scattered light due to particle motion with respect to fixed electrodes.





# **Optical System**

# Use optical mixing to extract motion of particles relative to electrodes.





# **Optical Mixing**



Sum of waves has easy to see oscillation



## **Data Analysis**

Analyze observed spectrum.



Frequency



# **Doppler Shift Calculations**

$$\omega_D = \vec{q} \bullet \vec{V}$$

$$\mu_e = \frac{V}{E}$$

- ω<sub>D</sub> = frequency (Doppler) shift, measured
- q = scattering vector
  - (4 $\pi$ n/ $\lambda$ )sin( $\theta$ /2), known
- **V** = particle velocity
- E = electric field strength , known
- μ<sub>e</sub>= electrophoretic mobility (desired result)



# **Zeta Potential Calculation**

- Need to use a model to obtain zeta potential (desired quantity) from mobility (measured quantity.
- Most common is Smoluchovski (shown here)

 $\mu_e = \frac{\varepsilon_r \varepsilon_0 \zeta}{\eta(T)}$ 

- $\varepsilon_r$  = relative permittivity (dielectric constant)
- $\varepsilon_0$  = permittivity of vacuum
- $\zeta$  = zeta potential
- $\eta$  = viscosity (function of temperature)



# **Application (Zeta Potential)**

#### Ceramics; Ludox<sup>R</sup> Silica

Silica	L	Ludox Silica <sup>R</sup> TM-50 with 0.01M_KCI										
Sample Preparat	tion 1	100 ppm										
<b>Conditions</b>	·											
<b>Temperature: 25</b>	C degr	'ee		100								
Solvent: Water Refractive Index: 1.333 Distribution base:			90									
			80									
				00								
ocattering in	y			70-			1					
	Resu	lts	., u.)	60								
obility	-3.0	2	ity (a.	50-								
ata Potential	_31	8	tens	40								
nV)	-51.		Ľ	30								
				20								
				10-								
				0								
					-150	-100	-50	0	50	100	150	200
Zeta potential (mV)												

From GRACE's catalogue



#### **Isoelectric point**



# X-axis can also be Ca++ or other ion concentration. But, pH is most common.



# Why is IEP important?

- When particle charge is important for stability (which is often), the IEP is a quick way of describing the region at which particle charge is low. And those are conditions under which flocculation will occur.
- Will depend on <u>surface (and liquid)</u>, not bulk material composition.



# Why Measure IEP?

- It serves as a method of characterizing sample.
- It tells you amount of reagent to add (or regimes to avoid) when processing the sample.



# **Grades of TiO**<sub>2</sub>



#### Care needed when dispersing!



# **Effect of Surface on IEP of TiO<sub>2</sub>**

Bulk %	coating	IEP (pH units)			
SiO <sub>2</sub>	Al <sub>2</sub> O <sub>3</sub>				
		6.8			
	4.5	<b>8.4</b> (R900)			
6.5	3.5	<b>5.8</b> (R960)			
8.0	8.0	<b>4.6</b> (R931)			

Bulk percentages (elemental analysis) of each chemical coating not reliable indicator of how the surface will behave in solution

#### IEP of $Al_2O_3$ at pH 9, IEP of $SiO_2$ at pH 2~3

Imperative to check ZP vs pH profile for any material prior to use



#### Using Zeta Potential to Predict/Control Particle Interactions

- Note size maximum at IEP due to flocculation.
- If you mix Al<sub>2</sub>O<sub>3</sub> with Fe<sub>2</sub>O<sub>3</sub>, what happens?
  - At pH 2, both are positive: no interaction
  - At pH 6,  $Al_2O_3$  his positive and  $Fe_2O_3$  is negative: particles stick together.
  - At pH 9, both are negatively charged, no interaction.



Data from Berg et al., *Nanotoxicology*, Dec. 2009; 3(4): 276-283

Zeta



# **Application (Zeta Potential)**

#### Iso Electric Point of Coffee Mate

	Results
lso-electric point	pH 4.0

If you want to bind Ludox <sup>a</sup> HS (negatively charged at <sup>a</sup> all "allowed\*" pH values) to<sup>o</sup> Coffee Mate, which pH (between 2 and 12) should you choose?



pH 2, pH 4, pH 6

\*allowed for this example



# **Application (Zeta Potential)**

#### Iso Electric Point of Coffee Mate

	Results
lso-electric point	pH 4.0

If you want to bind Ludox do HS (negatively charged at do HS (negative) to hS (negative) to hS (negative) to hS (negative) to hS (ne



Coffee mate has a positive charge at pH 2. Since it is positively charged, it will be attracted to the negatively charged Ludox.

\*allowed for this example



# **Can compress all of this information**

- Coffee Mate IEP: pH 4
- Al<sub>2</sub>O<sub>3</sub> IEP: pH 7
- Fe<sub>2</sub>O<sub>3</sub> IEP: pH 4



# **Practical Charge: Béarnase Sauce**

#### Monique's\* Recipe for Béarnase Sauce

Mix <u>vinegar</u>, shallots, tarragon and ground pepper. Add ½ glass of white wine.

Boil until practically all the liquid has evaporated. Cool until nearly <u>cold</u>. Add <u>egg-yolks</u>, one-byone while stirring vigorously. Add ½ glass of white wine and mix well.

Heat the mixture on a "bain-marie" while <u>stirring all the time</u>. When the sauce becomes creamy, cool slowly again continuously stirring. When the pot temperature is such that you can touch it with your hands, <u>add clarified butter</u>, a <u>small amount at a time</u>. The <u>temperature of the clarified</u> <u>butter should be about the same as the temperature of the pot mixture. Continue stirring</u>. When all the butter is incorporated, sieve the mixture and add one spoonful of cut chervill. Keep the sauce tepid on the bain-marie. <u>Never heat it up again</u>. The consistency of the sauce should be similar to that of mayonnaise .

This tastes great, but sometimes it fails

Packet Sauce Instructions

Add butter and milk to pan. Heat to boiling, add mix with stirring.

Simmer and stir until sauce has creamy consistency.

Q: Why does a homemade sauce taste better?

Q: Why does a packet sauce mixture always work?



### **Particle Science has answers!**



Emulsion Particle Size determines sauce mouth feel/texture – large particles taste better.



Egg Yolk Zeta Potential determines sauce stability

Vinegar will adjust pH of sauce. If you get the wrong egg, the IEP is at pH 4.8 and the emulsion droplets coalesce. You need to adjust the vinegar on a case by case basis.



#### **Environmental**



# **Clay IEP**



# To flocculate clay so it settles, pH must be quite low. You will need a lot of acid.



# **Clay IEP with other ions**



# To flocculate clay so it settles, choose alum at 0.01 g alum/g clay. Too much or too little and floculation is not ideal.



### **Zeta Practical Tips**



### Impurities

- Zeta is a <u>surface property</u>.
- Result is sensitive to surface active impurities.
  - Soaps/detergents
  - Specific ions (e.g., Cl<sup>-</sup>, TSPP)
  - Grease/oil/fingerprints (hydrophobic materials will go to surfaces of aqueous suspension)
- Keep everything extremely clean
- Keep surfactant additives in mind when interpreting data.



#### **Electrolyte**

- Recall that potential is a function of ionic strength.
- Pure water has an ionic strength of  $\sim 10^{-7}$  M.
- A little bit of CO<sub>2</sub> from the air can raise ionic strength by a factor of 10 or 100 to 10<sup>-5</sup> M.
- Use 1 mM electrolyte instead of no electrolyte to keep electrolyte levels (and therefore results) consistent from sample to sample.
- This doesn't apply for samples that already have substantial electrolyte.



# **Electrolyte Continued**

- Titration is popular.
- Remember that acid and base will add to system ionic strength. pH 3 corresponds to 10<sup>-3</sup> M electrolyte.
- Adding acid or base will increase ion concentration.
- Start with a 10 mM (10<sup>-2</sup> M) salt (KNO<sub>3</sub>) concentration to keep acid/base concentration from affecting results.



# **Concluding Comments**

 Use measured surface charge (zeta potential) to predict colloidal stability

 Use measured surface charge to <u>control</u> particle-particle interactions





#### Ask a question at <a href="mailto:labinfo@horiba.com">labinfo@horiba.com</a>

#### Keep reading the monthly <u>HORIBA Particle</u> e-mail newsletter!

Visit the <u>Download Center</u> to find the video and slides from this webinar.

Jeff Bodycomb, Ph.D.

- P: 866-562-4698
- E: jeff.bodycomb@horiba.com





#### Thank you



# **Thank you**

감사합니다

Cảm ơn

HORIBA

ありがとうございました

Dziękuję

Grazie

Merci

谢谢

நன்ற

ขอบคุณครับ

धन्यवाद

**Obrigado** 

Σας ευχαριστούμε



Tack ska ni ha

Большое спасибо Gracias