


An Overview of the Concept, Measurement, Use and Application of Zeta Potential



**David Fairhurst, Ph.D.
Colloid Consultants, Ltd**


Fundamental Parameters that control the Nature and Behavior of all Particulate Suspensions



How Particle Surfaces Acquire a Charge in Water

(a) Differential ion solubility


Net positive surface charge



Net negative surface charge

How Particle Surfaces Acquire A Charge in Water



(b) Direct ionization of surface groups



(c) Isomorphous ion substitution


How Particle Surfaces Acquire a Charge in Water

(d) Specific ion adsorption



(e) Anisotropic crystals

Origin of Charge in Clays

Isomorphous Substitution		FACE
Lattice Imperfections		-ve
Broken Bonds		EDGE
Exposed Structural OH		+ve

In neutral water, *net* charge will usually be *negative*

Particle Association: **F – F E – F E – E**

Many structures are possible

Particle Charges of Various Surfaces in Neutral Water

Positive


Ferric Hydroxide
Aluminium Hydroxide
Chromium Hydroxide
Thorium Oxide
Zirconium Oxide
Basic Dyes
(Methylene Blue)
Base Proteins
(Protamines, Histones)


Negative

Silicon Dioxide
Au, Ag, Pt, S, Se
As₂S₃, PbS, CuS
Acidic Dyes
(Congo Red)
Acid Protein
(Casein, BSA)
Viruses, Microbes
Air bubbles

Charge in non-aqueous media often opposite in sign!
(Electron Donor - Acceptor Theory)

The Electric Double Layer





Ψ_0 - cannot be measured
 Ψ_d - mathematical concept

ζ - experimental parameter
 $\zeta \approx \Psi_d$

$$\Psi = \zeta \exp [- \kappa x]$$

The Debye-Hückel parameter, κ , defined as:

$$\kappa = [2e^2 N_A c z^2 / \epsilon \epsilon_0 k_b T]^{1/2}$$

The Debye length, κ^{-1} is a measure of the “electric double layer thickness”


For single symmetrical electrolyte:

$$\kappa^{-1} = 0.3041 / Z C^{1/2}$$


c is the concentration of electrolyte of valence, z

The electric potential *depends* (through κ) *on the ionic composition of the medium*. If κ is increased (i.e. the electric double layer is “compressed”) then the potential must decrease.

Effect of addition of electrolyte on the zeta potential



Effect of specific adsorption of an anion on the zeta potential of a cationic surface



Zeta Potential is the “effectiveness” of the surface charge in solution

Depends upon:

- Fundamental “surface” sites – how many, what type
- Solution conditions – temperature, pH, electrolyte concentration

**Useless to quote a zeta potential value
without specifying suspension conditions**

Calculation of the zeta potential

ζ is not determined directly

Most common technique:
microelectrophoresis (ELS/PALS)

Electrophoretic mobility,

$$U_E = V_p / E_x$$

V_p is the particle velocity ($\mu\text{m/s}$) and E_x is the applied electric field (Volt/cm)

Relation between ζ and U_E is non-linear: $U_E = 2\varepsilon\varepsilon_0 \zeta F(\kappa a) / 3\eta$


The Henry coefficient $F(\kappa a)$ is a complex function of ζ

Simplest solution: use electrophoretic mobility, U_E
as the measurement metric




SZ-100
nano particle analyzer


Effect of Electrolyte Concentration on Particle Charge



Zeta Potential of Corundum (Al_2O_3) in Solution of Various Electrolytes




Effect of pH on Particle Charge



Sign of Zeta potential: $\text{pH(iep)} - \text{pH(solution)}$

Aqueous Isoelectric Points



Isoelectric Points of some Oxides

Oxide

pH value of I.E.P.

Silicon Dioxide	2
Manganese Dioxide	3
Zirconium Dioxide	4
Titanium Dioxide (Rutile)	6
Chromium Oxide	7
Iron Oxide	8
Aluminium Oxide	9
Lead Oxide	10
Cadmium Oxide	11
Magnesium Oxide	12

Force of Repulsion

$$V_R \approx D a \zeta^2 \exp(-KH)$$

D is a constant related to the permittivity (dielectric constant) of the material.

a is the particle or droplet radius.

ζ is a measure of the surface potential (charge).


K is proportional to the ionic strength (“conductivity”).

H is the distance between particle surfaces.

For a fixed medium, particle size and zeta potential: repulsive force decreases as the ionic strength increases

For a fixed medium, particle size and ionic strength: repulsive force becomes larger with increase in zeta potential

Zeta potential and stability



Electrostatic Stabilization


Material	ZP (mV)
O/W emulsions	>15
Polymer Latices	>20
Metal oxides	>40
Metal Sols	>70

Critical Coagulation Concentration

$$ccc \approx \zeta^4 / z^2$$

z is electrolyte *counterion* valence

Effect of Zeta Potential on Suspension Properties



Stabilization


Any material added into solution can affect suspension stability

Water soluble polymers – “thickeners”, viscosity modifiers.

Presence in solution affects Repulsive Potential via the **DIELECTRIC** term:

$$V_R \propto \mathbf{D} \, a \, \zeta^2 \quad [\text{Geometric term}]$$

Surface Modification



Care needed when dispersing!

Effect of Surface Modification on the IEP of TiO_2

Bulk%coating		IEP (pH units)
SiO_2	Al_2O_3	
--	--	6.8
--	4.5	8.4 (R900)
6.5	3.5	5.8 (R960)
8.0	8.0	4.6 (R931)

Bulk percentages (elemental analysis) of each chemical coating not reliable indicator of how the surface will behave in solution

Imperative to check ZP vs pH profile for any material prior to use

Surface Modification

Typical “coatings” on TiO_2

Inorganic

Metal oxides

Silicones

Organic


Fatty acids

Organosilanes

Check the material MSDS!

Care needed in choice of dispersing aids!

Zeta Potential of Non-oxides



Surface impurities and contamination

In Conclusion

Zeta potential (ζ) measurement very useful technique

- ❑ provides information about the material surface-solution interface
- ❑ knowledge of ζ used to predict and control stability of suspensions/emulsions
- ❑ Measurement of ζ often key to understanding dispersion and aggregation processes
- ❑ The presence/or absence of surface charged moieties on materials (revealed by their ζ) directly affect their performance and processing characteristics in suspension
- ❑ The sign and magnitude of ζ affects process control, quality control and product specification
- ❑ at simplest level: help maintain a more consistent product
- ❑ at complex level: can improve product quality and performance



Q&A

Ask a question at labinfo@horiba.com

Keep reading the monthly **HORIBA Particle** e-mail newsletter!

Visit the [Download Center](#) to find the video and slides from this webinar.

Jeff Bodycomb, Ph.D.
P: 866-562-4698
E: jeff.bodycomb@horiba.com

Thank you

Thank-you