




Particle Characterization of Pigments

Mark Bumiller
Mark.bumiller@horiba.com



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



© 2012 HORIBA, Ltd. All rights reserved.

Outline

- Why is particle size important?
- Organic pigments using laser diffraction
- Inorganic pigments using laser diffraction
- Ink pigments
- TiO_2 study using laser diffraction
- TiO_2 particle size and zeta potential using dynamic light scattering
- Conclusions


Pigments

A pigment is a material that changes the color of reflected or transmitted light as the result of wavelength-selective absorption



Properties Dependent on Particle Size

- Hue/Tint Strength
- Hiding/Transparency
- Gloss/Flatting and Film Appearance
- Flocculation
- Viscosity
- Stability
- Weather resistance



Organic Pigments


- Commonly powder form
- Dry milling procedures monitored by dry laser diffraction method
- Particle size effects dispersibility of pigment in application media
- Wet laser diffraction method to monitor dispersion of organic pigments in low viscosity liquid media

Organic Pigments Study



Samples*:

- 2583 Diarylide Yellow 83
- 1274 Hansa Yellow 74
- 2574 Hansa Yellow 74
- 5576 Phthalo Blue 15:3




1. Select the small, high-dispersion nozzle for the fine particle sizes expected.
2. Set the automatic measurement conditions to proper testing range (96-94%T for the laser).
3. Set dispersion air pressure to Low.
4. To prevent nozzle clogging, remove large agglomerates from sample using 1mm sieve
5. Load enough filtered sample onto feeder chute to perform multiple measurements.
6. Take 3 consecutive measurements using the Auto Measurement function.


*Thank-you to Lansco Colors, Warwick, RI for supplying the samples for this study. www.pigments.com

Dry Organic Pigment Results

2583 Diarylide	Yellow 83
Sample ID	Results
Run 1	18.386
Run 2	18.082
Run 3	18.032
Mean	18.17
COV %	1.05




1274 Hansa	Yellow 74
Sample ID	Results
Run 1	50.125
Run 2	50.096
Run 3	50.442
Mean	50.220
COV %	0.38




Dry Organic Pigment Results


2574 Hansa	Yellow 74
Sample ID	Results
Run 1	36.490
Run 2	36.618
Run 3	36.603
Mean	36.570
COV%	0.19




5576 Phthalo	Blue 15:3
Sample ID	Results
Run 1	25.114
Run 2	25.546
Run 3	25.531
Mean	25.397
COV %	0.97



Inorganic: Iron Oxide, LA-950 Wet



Inorganic: Chromium Oxide Green




Example data

Median: 1.51 μm
Mean: 1.73 μm
D(10%): 0.89 μm
D(90%): 2.87 μm


- Micro 90 surfactant to wet powder
- 0.1% sodium pyrophosphate

Inorganic: Alumina Al_2O_3




Example data

File Name: 6472B
 Median: 0.300 μm
 Mean: 1.481 μm
 S. D.: 4.950 μm
 D(10%): 0.194 μm
 D(90%): 0.643 μm




Ink Pigment



Data Name	Graph Type	Sample Name	Median Size	Std.Dev.
200511011226001		COLOR WONDER SPRAY INK	0.10583(μm)	0.0120(μm)
200511011227002		COLOR WONDER SPRAY INK	0.10583(μm)	0.0120(μm)
200511011227003		COLOR WONDER SPRAY INK	0.10574(μm)	0.0120(μm)

TiO₂ Commercial Data Sheet

Figure 1. Particle Size Distribution (Horiba Laser)



Product Description


DuPont™ RPD Vantage® rutile paper dry titanium dioxide delivers improved opacifying value for coating and wet end paper applications.

Key Features


- Reduced TiO₂ consumption
- Excellent brightness
- Superior pigment dispersion
- Highest opacity/hiding power

Particle Size

Particle size and distribution are important contributors to light scattering by diffraction. **Figure 1** shows the RPD Vantage® particle size distribution compared against anatase. The lower median particle size (D50) and the narrower geometric standard deviation (GSD) of RPD Vantage® combined with the higher refractive index of rutile result in the highest scattering efficiency.




Crystal structure of Anatase titanium dioxide




Crystal structure of Rutile titanium dioxide

TiO₂ : Effect of Size

Relative scattering power
rutile TiO₂ vs. size



Opacity of TiO₂
mean 0.25 and 10 μm




TiO₂ Measurement Procedure

- Step 1: Buy a HORIBA LA-950
- Decide wet or dry analysis
- Wet: disperse sample using surfactant and ultrasound
- Method development using Method Expert
 - RI, concentration, ultrasound
- Dry: Use smallest nozzle and highest air pressure

TiO₂: Effect of Refractive Index


Measure sample once
Fix real RI, vary imaginary
Minimize R parameter (error calculation)

partica LA-950
Method Expert



Use 1.0 for imaginary RI

TiO₂ for Cosmetics



% under 100 nm
“nanoparticles”


Mean Size	: 0.21941(μm)
Median Size	: 0.15503(μm)
R Parameter	: 0.2560
D(v,0.1)	: 0.09093(μm)
Q(v,0.5)	: 0.15503(μm)
D(v,0.9)	: 0.36061(μm)
Cumulative % on Diameter	: (1)0.100 (μm)- 14.473(%)

TiO₂: Wet vs. Dry

Mean Size : 0.42307(μm)
 Mode Size : 0.3635(μm)
 Geo.Std.Dev. : 1.7940(μm)
 R Parameter : 6.8993E-2
 D(v,0.1) : 0.16831(μm)
 D(v,0.5) : 0.35496(μm)
 D(v,0.9) : 0.71036(μm)

Only 8 nm difference!
 Dry dispersion below 350 nm median

Wet




Mean Size : 0.44614(μm)
 Mode Size : 0.3615(μm)
 Geo.Std.Dev. : 1.8683(μm)
 R Parameter : 1.3147E-1
 D(v,0.1) : 0.16820(μm)
 D(v,0.5) : 0.34653(μm)
 D(v,0.9) : 0.74838(μm)

Dry

TiO₂: Dry Analysis of 2 Samples


Mean Size : 0.44614(μm)
 Mode Size : 0.3615(μm)
 Geo.Std.Dev. : 1.8683(μm)
 R Parameter : 1.3147E-1
 D(v,0.1) : 0.16820(μm)
 D(v,0.5) : 0.34653(μm)
 D(v,0.9) : 0.74838(μm)

Able to reproducibly distinguish
 2 samples with 12 nm difference
 @ D(v,0.5)




Mean Size : 0.45286(μm)
 Mode Size : 0.3183(μm)
 Geo.Std.Dev. : 1.9285(μm)
 R Parameter : 8.0177E-2
 D(v,0.1) : 0.16169(μm)
 D(v,0.5) : 0.33394(μm)
 D(v,0.9) : 0.77890(μm)

Nano TiO₂*



Data from HORIBA LB-500




* PREPARATION AND UV CHARACTERIZATION OF TiO₂ NANOPARTICLES SYNTHESIZED BY SANSS


Ching-Song Jwo¹, Der-Chi Tien², Tun-Ping Teng³, Ho Chang⁴,
Tsing-Tshih Tsung⁴, Chih-Yu Liao⁵ and Chi-Hsiang Lin¹

Rev. Adv. Mater. Sci. 10 (2005) 283-288

TiO₂ Size & Zeta Potential on SZ-100



Size, note: Intensity distribution



Zeta potential @ pH = 7

Conclusions

- Particle size important physical parameter requiring analysis for control
- Laser diffraction excellent tool down to 30 nm
- Dynamic light scattering for sub micron particle size + zeta potential
- HORIBA has products and support for your application requirements

To Learn More: www.horiba.com/particle

HORIBA

HORIBA Scientific

Chisissez un pays ou une zone

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific | All Segment Product Browser

Products | Markets & Industries | Service & Support | News & Events | Newsletter | About Us | Employment | Feedback

Particle Characterization

Home → Scientific → Products → Particle Characterization

Particle Characterization

HORIBA designs, manufactures, and supplies state of the art particle characterization instruments.

Every instrument across the five business segments must meet stringent requirements before the HORIBA name is attached. The Particle Characterization group of analyzers has incorporated this principle into each new design since entering the business in 1979. Relentless innovation united with high performance to attain the ultimate goal: a new standard in usability.

Particle Characterization Products

HORIBA offers instruments for particle size, particle shape, zeta potential, and surface area analysis. Measurable particle size range is from 1 nanometer to 30 millimeters, at concentrations ranging from 1 ppm to 50 vol% with shape determination available starting at 1 micrometer. A range of analytical techniques are employed including laser diffraction (Mie Theory), dynamic light scattering, acoustic and electroacoustic spectroscopy, and dynamic and static image analysis. (measuring both particle size and shape information).

HORIBA's advanced designs and powerful software, combined with flexible sample handling systems are available to meet every analysis need. These instruments can incorporate small volume pumping systems for precious materials, high throughput automation, dry powder dispersers and temperature controlled flow systems in order to provide the user with the best possible solution with none of the trade-offs that might otherwise be necessary.

Particle Size

- Laser diffraction
 - LA-950V2
 - LA-300
- Dynamic light scattering
 - SZ-100

Request Information

Quick Request

E-mail us

Particle Size Essentials

eBook

Newsletter

Sign Up

Download Center

More Info

Member login

Please enter your e-mail address and password in order to login on to www.horiba.com. We may ask you later to answer additional questions if you are the first time to this area.

E-Mail

Address:

Password:

LA-950 Wet and Dry Sizing of Titanium Dioxide AN204 [Download](#)

Talk to us, ask questions

Receive news of updates

View application/technical notes (>130), webinars (>40), white papers.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

© 2012 HORIBA, Ltd. All rights reserved.

Thank you

ありがとうございました

ขอบคุณครับ

谢谢

اشكر

Gracias

Grazie

Σας ευχαριστούμε

धन्यवाद

நன்றி

Tacka dig

Danke

Merci

감사합니다

Большое спасибо

Obrigado

おかしく おかしく

Omoshiro Okashiku