


Carbon Nanotube Analysis





Particle Analysis: Jeffrey Bodycomb, Ph.D. Fluorescence: Adam Gilmore, Ph.D.



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

© 2012 HORIBA, Ltd. All rights reserved

Complementary Size and Aspect Ratio Analysis of SWCNTs: Photoluminescence Excitation-Emission Mapping



Adam M. Gilmore, Fluorescence Product Manager

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

© 2012 HORIBA, Ltd. All rights reserve

Me and My Wife





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific




© 2012 HORIBA, Ltd. All rights reserved.

クカーボンナノチューブ

1981年カーボンナノチュープに211の発見ます。 CSTのナノテクノロジーは着くほど通いスピードで研究用品が進んでかり、 今年単時代の成果を提供ロナノナクノロジーと数件し、 次応用に見図すら刻帯にあていると言えます。

カーボンテノチューブ・3つの構造



ナノマアリアルの条件 ・ボルイエアの時代になった時期 ・ボル・ロークスイデアの後 ・モンの時に、たちになりたちになったのかでたらなてアークメイデアの後 ・モンの時代、たま時になっていたけができたかがらなっ、

04

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

08

© 2012 HORIBA, Ltd. All rights reserved.


Carbon

Nanotube

SWCNTs are Graphene Roll-ups







Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Why are SWCNTs Significant?

- 2 Major SWCNT families:
- 2/3 =Semiconducting:
 - photo- and electro-luminescence (PL/EL),
 - field-effect transistors (FETs)
 - Precise size and bandgap selection for device engineering
 - Bright, non-blinking PL/EL for chemical and biological sensing
 - NIR emission- fits ideally in biological window
 - Faciliate dense transistor networks
- 1/3 =Metallic:
 - high electric and thermal conductivity, efficient connectors
 - Faciliate transparent conductive films (TCFs)
 - Enhanced efficiency photovoltaic materials

© 2012 HORIBA, Ltd. All rights reserved





SWCNTs Hot News Flashes





- Effects of Gamma radiation characterized, show promise for medical apps and sterilization
- CNTs can divert heat from current flow in devices
- CNTs enhance photoacoustic imaging of tumours
- SWCNTs may replace ITO in solar cells
- Strain paint-stress changes shape and absorbance-emission spectra-aircraft etc..
- Sub-10 nm SWCNT transistors more efficient than predicted by models

SWCNTs: More Big News in 2012





- Discovery of how sonication shortens and damages SWCNTs
- Not all tubes are cylindrical they can be flattened, like graphene, as long as their diameter is wide enough (> 2 nm)
- Chiral selection possible by growth on stainless steel wires
- Selective dispersion methods improved for chiral selection.







- Def: Subsequent emission of light from a material caused by light it had previously absorbed.
- Semiconductor PL: Light emission from around the semiconductor material's bandgap energy level excited by absorption of light energy above the material's bandgap energy level.

© 2012 HORIBA, Ltd. All rights reserved




Quantum Confinement Effect






- The major principle defining the relationship between a nanoparticles' physical diameter and it's bandgap energy
- The 'Bohr-Exciton Radius' for a bulk semiconductor material can be defined is the physical distance (in nm) between the electron and it's hole across the bandgap.
- When the diameter of a nanoparticle of a material is smaller than it's 'Bohr-Exciton Radius' the nanoparticle's bandgap is inversely related to the diameter of the nanoparticle due to 'quantum confinement'.

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



© 2012 HORIBA, Ltd. All rights reserved

Suspension of SWNTs for PL





- 1. Nonfluorescent SWNT-SDS aggregates are sonicated
- 2. Single SWNTs are suspended
- 3. Centrifugation separates fluorescent SWCNTs

M. J. O'Connell et al., Science 297 (2002) 593 S. M. Bachilo *et al.,* Science 298 (2002) 2361.


Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific




finture.O

PL Characterization of SWNTs



HORIBA

Explore the future



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA


© 2012 HORIBA, Ltd. All rights reserved.

The iHR-320

Any 2 Detectors

- CCDs
- InGaAs Arrays
- PIN Diodes
- Photomultipliers
 - •Steady State- 10 ms
 - •TCSPC-50 ps
 - •MCS-20 ns
- Microchannel Plate PMTs
 •5 ps TCSPC

Kinematic Grating Turret 200 nm – 3000 nm



Explore the future



Single Walled Carbon Nanotubes (SWCNTs)

- SWCNTs can be characterized using photoluminesence spectroscopy for:
 - Diameter
 - Helical twisting (chirality)
 - Length
 - Bundling-(SWCNT to SWCNT interactions)




PL Yields Multidimensional Data



- absorbance and emission peaks correlate according to quantum confinement rules, smaller SWCNTs higher energy
- Helical twisting (chirality):
 - absorbance and emission peak energies are also influenced secondarily, and nonlinearly according to the SWCNT families
- Length:
 - intensity of absorbance extinciton and emisison intensity correlate.
- Bundling:
 - energy transfer from smaller (donor) to larger (acceptor) diameter SWCNTs influences relationship between absorbance and emission peak intensities.


Energy Transfer in Bundles





Length Dependent SWCNT PL





Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

SWCNT Absorbance Spectra

A.4 Results

Length-Dependent Op




Figure 2. Optical ab

Key

X Wavelength (nm)

Y Absorbance (arb. unit)

Figure A.2 — <u>The</u> absorption spectra of SWCNT samples dispersed in 1% (mass fraction) SC-D₂O solvent.

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA





ICLES

The Nanolog®-EXT



- Deeper NIR Excitation and Emission ranges are needed for larger diameter SWCNTs
- Larger diameter SWCNTs with small bandgaps are important for device manufacturing.


Essential Components:

- 1. Tunable Mainframe CW TiS Laser
- 2. Extended InGaAs Array (1100-2200 nm)



© 2012 HORIBA, Ltd. All rights reserved

Newport 3900S Power Ti-S Laser



- Tuning from 700-1000 nm
- Useful for SWCNTs from 1.3 – 2.2 nm
- The only light source compatible with extended InGaAs array for large diameter SWCNTs



Single-axis DC motor controller/driver

Model: SMC100CC Availability: 1 Week

The SMC100CC is a single axis motion controller/driver for DC servo motors up to 48 VDC at 1.5 A rms. It provides a very compact and low-cost solution for driving most of Newport's stages, including the popular... Read more

HORIBA Tuning Software with SMCC100 Stepper Motor Controller

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific


HORIBA

© 2012 HORIBA, Ltd. All rights reserved

3900S 5W Pump Power









© 2012 HORIBA, Ltd. All rights reserved.





Large Dia. SWCNT Absorbance




ORIBA

Scientific


Large Diameter SWCNTs EEM







Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



-®



HORIBA Scientific

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

© 2012 HORIBA, Ltd. All rights reserved.


HORIBA Scientific



Nanosizer®







Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

© 2012 HORIBA, Ltd. All rights reserved.



finture O

SWCNTs: TOP TEN STORIES

NEWSREEL 2011-2012

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA


© 2012 HORIBA, Ltd. All rights reserved.





Engineers build first sub-10-nm carbon nanotube transistor

February 1, 2012 by Lisa Zyga feature




9-nm CNT transistor with electron microscope images. Image credit: Franklin, et al. ©2012 American Chemical Society

(PhysOrg.com) -- Engineers have built the first carbon nanotube (CNT) transistor with a channel length below 10 nm, a size that is considered a requirement for computing technology in the next decade. Not only can the tiny transistor sufficiently control current, it does so significantly better than predicted by theory. It even outperforms the best competing silicon transistors at this scale, demonstrating a superior current density at a very low operating voltage. By Cameron Chai

Using carbon nanotubes, Rice University scientists, Bruce Weisman and Satish Nagarajaiah, have developed a new class of paint named 'strain paint,' which is capable of detecting strain in airplanes, bridges and buildings.

The researchers believe that their strain paint will be helpful in detecting deformations in structures such as airplane wings. This composite coating can be read with the help of a handheld infrared spectrometer. The study findings have been reported online in Nano Letters, a journal of the American Chemical Society.

Using this novel paint, it is possible to detect the signs of deformation in a material much earlier than the impact becomes detectable to the naked eye, and most importantly without contacting the structure. Moreover, the nanotube-based system is capable of measuring strain along any direction and at any spot.



An illustration shows how polarized light from a laser and a near-infrared spectrometer could read levels of strain in a material coated with nanotubeinfused paint invented at Rice University. (credit: Bruce Weisman/Rice University)





Explore the future





Rice University Researchers Reveal Results of Study on Flattened Nanotubes


Published on June 22, 2012 at 4:20 AM

By Will Soutter

A study carried out at the Richard E. Smalley Institute for Nanoscale Science and Technology, Rice University, revealed the potential of flattened carbon nanotubes.

Labeled as closed-edged graphene nanoribbons, they are the result of carbon nanotubes collapsing during growth. The nanoribbons demonstrate the properties of both graphene ribbons and nanotubes and hence could have a host of applications in the fields of materials and electronics. Collapsed nanotubes have the chemistry of both graphene in the middle and carbon-60 molecules (buckyballs) on the sides.

Researchers led by Robert Hauge found that the two portions can be separated by addition of functional groups on the sides. With the sides acting as insulators, the top and bottom layers are isolated and do not interact with the exception of excited-state or van der Waals-type interaction. Hauge believes that it is this process that generates new electronic and physical properties.



Flattened nanotubes (Credit: Ksenia Bets/Rice University)

HORIBA

Kequest

Quote

9

More Info



Science News

... from universities, journals, and other research organizations

Tiny Bubbles Snap Carbon Nanotubes Like Twigs

ScienceDaily (July 9, 2012) — What's 100 times stronger than steel, weighs one-sixth as much and can be snapped like a twig by a tiny air bubble? The answer is a carbon nanotube -- and a new study by Rice University scientists details exactly how the much-studied nanomaterials snap when subjected to ultrasonic vibrations in a liquid.

See Also:

Matter & Energy

- Nanotechnology
- Graphene
- Engineering

Computers & Math

- Computer Science
- Distributed
- Computing
- Computer Modeling


Reference

- Nanowire
- Carbon nanotube
- Fullerene
- Nanoparticle

"We find that the old saying 'I will break but not bend' does not hold at the micro- and nanoscale," said Rice engineering researcher Matteo Pasquali, the lead scientist on the study, which appears this month in the Proceedings of the National Academy of Sciences.

Carbon nanotubes -- hollow tubes of pure carbon about as wide as a strand of DNA -- are one of the moststudied materials in nanotechnology. For well over a decade, scientists have used ultrasonic vibrations to separate and prepare nanotubes in the lab. In the new study, Pasquali and colleagues show how this process works -- and why it's a detriment to long nanotubes. That's important for researchers who want





The mechanism by which carbon nanotubes break or bend under the influence of bubbles during sonication is the topic of a new paper led by researchers at Rice University. The team found that short nanotubes are drawn end-first into collapsing bubbles, stretching them, while longer ones are more prone to breakage. (Credit: Pasquali Lab/Rice University)



Letter to the Editor

Chiral-selective growth of single-walled carbon nanotubes on stainless steel wires

Maoshuai He^{a,} A^{a,} Marita Niemelä^a, Juha Lehtonen^a, Ville Viltanen^c, Jani Sainio^c, Hua Jiang^d, Esko I. Kauppinen^d, Marita Niemelä^a, Juha Lehtonen^a

^a Department of Biotechnology and Chemical Technology, School of Chemical Technology, Aalto University, P.O. Box 16100, FI-00076 Aalto, Finland

^b A.M. Prokhorov General Physics Institute RAS, 38 Vavilov Street, 119991 Moscow, Russia

° Department of Applied Physics, School of Science, Aalto University, P.O. Box 11100, FI-00076 Aalto, Finland

^d Department of Applied Physics and Center for New Materials, School of Science, Aalto University, P.O. Box 15100, FI-00076 Aalto, Finland

Received 9 March 2012. Accepted 6 May 2012. Available online 14 May 2012.

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





NANO RESEARCH

Volume 4, Number 4 (2011), 393-404, DOI: 10.1007/s12274-011-0094-0



Abstract

There is currently great interest in the potential use of carbon nanotubes as delivery vessels for nanotherapeutics and other medical applications. However, no data are available on the effects of sterilization methods on the properties of nanotube dispersions, the form in which most medical applications will be processed. Here we show the effects of gamma irradiation from a ⁶⁰Co source on the dispersion and optical properties of single-wall carbon nanotubes in aqueous dispersion. Samples of different length-refined populations were sealed in ampoules and exposed to a dose of approximately 28 kGy, a level sufficient to ensure sterility of the dispersions. In contrast to literature results for solid-phase nanotube samples, the effects of gamma irradiation on the dispersion and optical properties of the nanotube samples were found to be minimal. Based on these results, gamma irradiation appears sufficiently non-destructive to be industrially useful for the sterilization of nanotube dispersions.

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





Carbon nanotube circuits could outsource their heat to a separate device

When most materials carry an electric current, the motion of the electrons ...

by Matthew Francis - Apr 9 2012, 1:40pm EDT

PHYSICAL SCIENCES SCIENCE 6

Electron micrograph of two experimental setups demonstrating remote Joule heating. A carbon nanotube (thin gray line) is connected between electrodes (dark gray).

Photograph by courtesy of John Cumings, University of Maryland

The phenomenon is familiar: if you run an electric current through a wire, the wire heats up. Known as "Joule heating" (for James Joule, the physicist-brewer who quantified it in the 19th century), the cause is usually very simple: the electrons carrying the current transfer some of their energy to the atoms in the wire, and the increased vibration of the atoms is measured as a rise in temperature. While it's very useful in some applications, Joule heating can often be a problem, especially in electronic devices like computer processors, where excess thermal energy can cook the chips.

A new experiment performed at the University of Maryland has produced Joule heating where the current is separated from the heat. Kamal H. Baloch, Norvik Voskanian, Merijntje Bronsgeest, and John Cumings determined that a current flowing in a carbon nanotube can transfer thermal energy into the material the tube is sitting on, a process they dub "remote Joule heating." In other words, there is a separation between the electric flow—confined to the nanotube—and increased heat, which ends up in the substrate, even though it carries no current. Using electron thermal microscopy, the researchers determined that as much as 84 percent of the power in the nanotube was transferred to the substrate, pointing to a possible new way to manage excess heat in electronics.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific


Latest News Web Date: June 1, 2012

Carbon Nanotubes Highlight Tumors

Medical Imaging: The materials could aid cancer detection using a safe, cheap imaging method

By Prachi Patel

Department: Science & Technology | Collection: Life Sciences News Channels: Materials SCENE, Biological SCENE, Analytical SCENE, Nano SCENE Keywords: photoacoustic imaging, contrast agent, carbon nanotubes, cancer



Doctors rely on magnetic resonance imaging and X-ray computed tomography to spot tumors and monitor cancer treatment. But some researchers think a cheaper and safer alternative could come from photoacoustic imaging, a noninvasive technique that produces images based on sound. Scientists have now developed **photoacoustic contrast agents based on carbon nanotubes** that home in on tumors in mice, highlighting the tissue in scans (*ACS Nano*, DOI: 10.1021/nn204352r).

To create an image of tissue using photoacoustics, researchers first shine visible or near-infrared light on a stretch of skin. The light heats up tissue and blood beneath the skin. As the tissue warms, it expands and contracts, generating sound waves. Based on features of these waves, computer software then recreates an image of the tissue. These images' resolution matches those of MRI or CT scans, but unlike those standard methods, photoacoustic imaging doesn't require expensive equipment or harmful X-rays.

© 2012 HORIBA, Ltd. All rights reserved.



-IORIBA

he finture.O

Scientific

3

В

.

Email 🔤 Print



Science News

... from universities, journals, and other research organizations

Researchers Use Carbon Nanotubes to Make Solar Cells Affordable, Flexible

ScienceDaily (Sep. 27, 2011) — Researchers from Northwestern University have developed a carbonbased material that could revolutionize the way solar power is harvested. The new solar cell material -- a transparent conductor made of carbon nanotubes -- provides an alternative to current technology, which is mechanically brittle and reliant on a relatively rare mineral.

See Also:

Matter & Energy

- Solar Energy
- Graphene
- Nanotechnology

Earth & Climate

- Renewable Energy
- Energy and the Environment
- Geomagnetic Storms

Reference

- Indium
- Metallurgy
- Solar cell
- Riomase

Due to Earth's abundance of carbon, carbon nanotubes have the potential to boost the long-term viability of solar power by providing a costefficient option as demand for the technology increases. In addition, the material's mechanical flexibility could allow solar cells to be integrated into fabrics and clothing, enabling portable energy supplies that could impact everything from personal electronics to military operations.

The research, headed by Mark C. Hersam, professor of materials science and engineering and professor of chemistry, and Tobin J. Marks, Vladimir N. Ipatieff Professor of Catalytic Chemistry

Ads by Google

Merrill Edge® Trading — Get 30 Free Online Trades Per Month w/Merrill Edge. Start Trading Now. www.MerrillEdge.com

Solar Defect Analysis — Solar Wafer Defect Analysis using semiconductor metrology www.silybwafers.com

Lowest cost wind power — New technology recently introduced OTCBB stock ticker MMMW www.massmegawatts.com

Affordable Accountants — Get Competitive Accounting Bids. Compare, Read Reviews. It's Free! MiNeeds.com/Accountants

Best in Residential Solar — Is Your Utility Bill Over \$50? Sunrun Solar Can Help You Save! www.SunrunHome.com



Flinders News - Latest news from Flinders University

Posts Tagged 'carbon nanotubes'

Solar cell turns windows into generators

Posted on: March 19th, 2012 by Marketing and Communications



Imagine a world where the windows of high-rise office buildings are powerful energy producers, offering its inhabitants much more than some fresh air, light and a view.

For the past four years a team of researchers from Flinders University has been working to make this dream a reality – and now the notion of solar-powered windows could be coming to a not too distant future near you.

As part of his just-completed PhD, Dr Mark Bissett (pictured) from the School of <u>Chemical</u> and <u>Physical Sciences</u> has developed a revolutionary solar cell using <u>carbon nanotubes</u>.

A promising alternative to traditional silicon-based solar cells, carbon nanotubes are cheaper to make and more efficient to use than their energy-sapping, silicon counterparts.

"Solar power is actually the most expensive type of renewable energy – in fact the silicon solar cells we see on peoples' roofs are very expensive to produce and they also use a lot of electricity to purify," Dr Bissett said.


"The overall efficiency of silicon solar cells are about 10 per cent and even when they're

operating at optimal efficiency it could take eight to 15 years to make back the energy that it took to produce them in the first place because they're produced using fossil fuels," he said.

Dr Bissett said the new, low-cost carbon nanotubes are transparent, meaning they can be "sprayed" onto windows without blocking light, and they are also flexible so they can be weaved into a range of materials including fabric – a concept that is already being explored by advertising companies.

Explore the future


Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific







Carbon Nanotube Analysis





Particle Analysis: Jeffrey Bodycomb, Ph.D. Fluorescence: Adam Gilmore, Ph.D.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific


HORIBA

What is Dynamic Light Scattering?

- Dynamic light scattering refers to measurement and interpretation of light scattering data on a <u>microsecond</u> time scale.
- Dynamic light scattering can be used to determine
 - Particle/molecular size
 - Size distribution
 - Relaxations in complex fluids


Other Light Scattering Techniques


- Static Light Scattering: over a duration of ~1 second. Used for determining particle size (diameters greater than 10 nm), polymer molecular weight, 2nd virial coefficient, R_g.
- Electrophoretic Light Scattering: use Doppler shift in scattered light to probe motion of particles due to an applied electric field. Used for determining electrophoretic mobility, zeta potential.



Brownian Motion

Particles in suspension undergo Brownian motion due to solvent molecule bombardment in random thermal motion.






- Brownian Motion
 - Random
 - Related to Size
 - Related to viscosity
 - Related to temperature

© 2012 HORIBA, Ltd. All rights reserved

DLS signal



Random motion of particles leads to random fluctuations in signal (due to changing constructive/destructive interference of scattered light.




time (microseconds)

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Correlation Function




Random fluctuations are interpreted in terms of the autocorrelation function (ACF).



Gamma to Size





 $D_h = \frac{k_B T}{3\pi\eta(T)D_t}$

Note effect of temperature!


- Γ decay constant D_t diffusion coefficient
- q scattering vector
- n refractive index
- λ wavelength
- θ scattering angle
- D_h hydrodynamic diameter
- η viscosity
- k_B Boltzman's constant

Explore the future



What is Hydrodynamic Size?

DLS gives the diameter of a sphere that moves (diffuses) the same way as your sample.





BA

Hydrodynamic Size




The instrument reports the size of sphere that moves (diffuses) like your particle.


This size will include any stabilizers bound to the molecule (even if they are not seen by TEM).

Gold Colloids

Technique	Size nm
Atomic Force Microscopy	8.5 ± 0.3
Scanning Electron Microscopy	$9.9~\pm~0.1$
Transmission Electron Microscopy	8.9 ± 0.1
Dynamic Light Scattering	13.5 ± 0.1

SEM (above) and TEM (below) images for RM 8011





© 2012 HORIBA, Ltd. All rights reserved

Why DLS?





- Non-invasive measurement
- Requires only <u>small quantities</u> of sample
- Good for <u>detecting trace amounts</u> of aggregate
- Good technique for <u>macro-molecular</u> <u>sizing</u>



New Nanoparticle Analyzer

Single compact unit that performs size, zeta potential, and molecular weight measurements.



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



RIBA


Carbon Nanotubes



Nanotubes are not spheres....hence the "tube" in the name.

Nanotube

Model for DLS




Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Motion of Carbon Nanotubes

Two types of motion: Translation and rotation



Here I ignore rotation

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

BA

What to do?



- Recognize that big tubes diffuse slowly and recognize that DLS is a fast, easy size indicator.
- Try to match DLS data to model using some other information

$$D_{t} = \frac{k_{B}T}{3\pi\eta L} \left[\ln(L/d) + 0.32 \right]$$
 L= length d=diameter

Nair, N., Kim, W., Braatz, R, and Strano, M. Langmuir 2008, 24, 1790-1795.



Well Maasurement of Nanotubes

knew aspect ratio (1000) and diameter (0.7 to 0.9 nm) in advance

- Well dispersed: few if any aggregates
- Expect value between 97 and 125 nm using equation on previous slide

Obtain a nice match.

	Z-avg. Diameter, nm
Repeat 1	115.0
Repeat 2	104.5
Repeat 3	105.3
Repeat 4	109.5
Repeat 5	115.2
Repeat 6	106.2
Average	109.3
Standard Deviation	4.8
Coefficient of Variation	4.4%

How to Measure with DLS?

- Ensure that the sample is well dispersed
- Ensure that the tubes are not colliding (dilute sample)
- Use results as an indicator of tube size (or aggregation if you are thinking about dispersion quality)

Don't blindly turn D_h into L





Other NanoSamples

Nanometals

- Nanogold
- Nano Iron Oxide
- Nanorods

Colloidal particles



Questions? labinfo@horiba.com www.horiba.com/us/particle

Adam.Gilmore@horiba.com Jeff.Bodycomb@horiba.com





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific