

Particle Characterization of Trona



Ian Treviranus

ian.treviranus@horiba.com

www.horiba.com/us/particle

More Information



HORIBA
Scientific

■ Best Practices/Training

- Setting Attainable Size Specifications TR007
- Understanding Laser Diffraction PSA Results TR008
- Troubleshooting Laser Diffraction Data TR010
- Help! How Can I Trust My Size Results? TR015
- Refractive index selection, sampling, dispersion, system verification, method development and more

■ Technology

- BET Flowing Gas Surface Area TE005
- Find the Best Analyzer for Your Application TE006
- Intro to Laser Diffraction TE010

What is trona?

Why are power plants using trona?

Impact of size and surface area

Featured technologies

Interpreting particle size results


Q&A

What is trona?



HORIBA
Scientific

■ Google says...



Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

© 2013 HORIBA, Ltd. All rights reserved.

What is trona?

- Naturally formed sodium sesquicarbonate
- $\text{Na}_2\text{CO}_3 * \text{NaHCO}_3 * 2\text{H}_2\text{O}$



What is trona?

Why are power plants using trona?

Impact of size and surface area

Featured technologies

Interpreting particle size results

Q&A

Why is trona useful?

- EPA regulates SO_x emissions
- Multiple strategies to remove SO_x
- Dry Sorbent Injection attractive option
- Trona
 - Capital cost low relative trad. scrubbers
 - Operational cost proportional to sorbent cost
 - SO_x removal proportional to reactivity
 - Reactivity directly proportional to surface area
 - Surface area directly prop. to particle size

Balance SO_x removal vs. cost

■ Unmilled trona

- Coarser PSD → Smaller SSA
- Need to inject more sorbent to hit removal target

■ Milled trona

- Finer PSD → Larger SSA
- Need to inject less sorbent to hit removal target
 - Added cost → either material cost or on-site milling
 - Finer powders typically bring flow and handling challenges

What is trona?

Why are power plants using trona?

Impact of size and surface area

Featured technologies

Interpreting particle size results



Q&A

Calcination aids removal



- Trona has low natural SSA
- Trona injected into flue gas undergoes calcination and produces sodium carbonate with higher SSA
- $2(\text{Na}_2\text{CO}_3 \cdot \text{NaHCO}_3 \cdot 2\text{H}_2\text{O}) + \text{heat} \rightarrow 3\text{Na}_2\text{CO}_3 + 5\text{H}_2\text{O} + \text{CO}_2$
- Evolution of water and CO_2 creates micropores \rightarrow higher SSA

Effect of Calcination




Measure surface area to screen incoming material

Unmilled trona



HORIBA
Scientific




Measure particle size to screen incoming material

Milled trona



HORIBA
Scientific




Measure particle size to screen incoming material and monitor milling end point, efficiency, etc.

Unmilled vs. milled trona



HORIBA
Scientific



What is trona?

Why are power plants using trona?

Impact of size and surface area

Featured technologies

Interpreting particle size results

Q&A

Featured technologies



HORIBA
Scientific

- **LA-950 & LA-300**

Laser Diffraction

- **SZ-100**

Dynamic Light Scattering & Zeta Potential

- **CAMSIZER & CAMSIZER XT**

Dynamic Image Analysis

- **PSA300**

Static Image Analysis


- **SA-9600**

Flowing Gas BET Surface Area

Size Range by Technique





HORIBA
Scientific



Laser Diffraction



HORIBA
Scientific



- Converts scattered light to particle size distribution
- Quick, repeatable
- Powders, suspensions
- Most common technique

LA-950: Laser Diffraction



HORIBA
Scientific

- Lowest total cost of ownership
- Measures in less than 60 seconds
- One-button operation
- Wet and dry measurement
- Ultra durable
- Ninth generation
- 10 nanometer – 3 mm



LA-300: Laser Diffraction



HORIBA
Scientific

- Unique portable design
 - Shippable in Pelican case
- Ultra durable
- One-button operation
- Wet measurement
- 100 nm – 600 μm



SA-9600 Surface Area Analyzer

- 0.1 - >2000 m²/g
- Single or multi-point analysis
- One or three station systems
- Single point: up to 30 analyses per hour
- High value route to quick and easy surface area analysis



What is trona?

Why are power plants using trona?


Impact of size and surface area

Featured technologies

Interpreting particle size results

Q&A

Size Terminology




The most common designation is **micrometers** or **microns**. When very small, in colloid region, measured in **nanometers**, with electron microscopes or by dynamic light scattering.

The Basics



HORIBA
Scientific

Which is the most meaningful size?



different
size definitions




different
results

The Basics



HORIBA
Scientific

What sizes can be measured?



The Basics




HORIBA
Scientific

Particle



Particle Distribution






The Basics



HORIBA
Scientific

Particle Size Particle Size Distribution



The Basics





HORIBA
Scientific

Laser Diffraction

Assumes hard, spherical shape model


q% = amount
of each size
by volume



Central Values



HORIBA
Scientific



Mean

Weighted Average
Center of Gravity

Median

50% Point

Mode


Peak of the distribution
Most common value

Remember: D[4,3] is sensitive to large particles

Mean vs. Median



HORIBA
Scientific




15%
difference

Classic example:
CMP slurries

Most Common Statistics

half are smaller than this diameter half are larger than this diameter




Bimodal Distribution



HORIBA
Scientific


Which numbers to use for specifications?
D50 still an option, but some prefer finer details



Multimodal Report



HORIBA
Scientific



Automate COV Calculation



- Coefficient of Variation indicates precision of multiple measurements
- ISO 13320 makes recommendations using COV

Summary Report				
Export Summary	Print Summary	Edit Layout	Best Fit Columns	Hide Selected
File Name	Sample Name	D(v,0.1)	D(v,0.5)	D(v,0.9)
200811061138068.NGB	Zircoa Slurry	0.065	0.107	0.185
200811061140069.NGB	Zircoa Slurry	0.071	0.145	11.896
200811061144070.NGB	Zircoa Slurry	0.069	0.129	3.838
Average		0.068	0.127	5.306
Std. Dev.		0.003	0.019	5.992
CV (%)		4.471	15.023	112.921
ISO 13320-1 (20.0, 15.0, 20.0)		PASSED	FAILED	FAILED

Result Verification

Verification Setting

Parameter: Median Size

Specification: ISO 13320-1

Standard Value: 50 (μm)

Tolerance: ± 5 (μm)

Certified range of values

D(v,0.5) >= 10μm ± 0 %

D(v,0.5) < 10μm ± 0 %

Result Display Setting

Pass:

Color: [Blue]

Text: OK

Fail:

Color: [Red]

Text: NG

OK Cancel


Distribution Graph		Data Table	Result Data
Mean Size	:	0.18408(μm)	
Variance	:	1.8988E-3(μm ²)	
Median Size	:	0.17730(μm)	
Mode Size	:	0.1649(μm)	
Std.Dev.	:	0.0436(μm)	
Chi Square	:	4.162519	
R Parameter	:	3.7379E-1	
Diameter on Cumulative %	:	(2)10.00 (%) - 0.1345(μm)	
	:	(9)90.00 (%) - 0.2450(μm)	
Cumulative % on Diameter	:	(1)850.0 (μm) - 100.000(%)	
	:	(2)600.0 (μm) - 100.000(%)	
	:	(3)425.0 (μm) - 100.000(%)	
	:	(4)300.0 (μm) - 100.000(%)	
	:	(5)212.0 (μm) - 100.000(%)	
	:	(6)150.0 (μm) - 100.000(%)	
	:	(7)106.0 (μm) - 100.000(%)	
	:	(8)75.00 (μm) - 100.000(%)	
	:	(9)53.00 (μm) - 100.000(%)	
	:	(10)38.00 (μm) - 100.000(%)	
Verification	:	1.OK 4.3% [D(v,0.5) 0.170 (μm)(± 6.000%)]	
	:	2.OK 3.5% [D(v,0.1) 0.130 (μm)(± 10.00%)]	
	:	3.OK 6.5% [D(v,0.9) 0.230 (μm)(± 10.00%)]	

Data Name	Graph Type	Transmittance(R)	Median Size	R Parameter
andy1	[Red Box]	88.3(%)	0.17730(μm)	0.373795
200801181026014	[Green Box]	81.1(%)	9.35329(μm)	0.069234
andy1	[Blue Box]	88.3(%)	0.17730(μm)	0.373795

Specification on X or Y Axis



HORIBA
Scientific




Specification with Error



HORIBA
Scientific

Must tighten internal spec by lab error %
Then product always within performance specification



http://www.spcpress.com/pdf/Manufacturing_Specification.pdf, By David Wheeler

Specification with Error



- Must tighten internal spec by lab error
- Therefore minimize lab error makes life easier
- How to minimize error?
 - Get sampling right
 - Structured method development
 - Eye on the goal: reproducibility

What is trona?

Why are power plants using trona?

Impact of size and surface area

Featured technologies

Interpreting particle size results

Q&A

Thank you

ありがとうございました

ขอบคุณครับ

谢谢

اشكر

Gracias

Grazie

Σας ευχαριστούμε

धन्यवाद

Tacka dig

Danke

Merci

நன்றி

감사합니다

Большое спасибо

Obrigado

おかしく

Omoshiro Okashiku

www.horiba.com/particle



HORIBA
Scientific

Ian Treviranus
Product Line Manager
ian.treviranus@horiba.com

← Talk to us, ask questions
labinfo@horiba.com

← Receive news of updates

← View application & technical notes (170+), webinars (60+), white papers.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

© 2013 HORIBA, Ltd. All rights reserved.