



Julie Chen Nguyen Julie.Nguyen@horiba.com

Emerging Food Trend: Plant Based Proteins, healthy fat, sugar

July 21, 2016

© 2016 HORIBA, Ltd. All rights reserved. Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA



Previously...

"Particle Technologies for Food & Beverage" February, 2014 (AP006) http://www.horiba.com/scientific/products/particle-characterization/downloadcenter/webinars

Covered applications:

- 1. Traditional mayonnaise
- 2. Milk homogenization
- 3. Flavor emulsions
- 4. Flavor powder
- 5. Wheat flour
- 6. Coffee beans
- 7. Instance coffee
- 8. Sugar crystals
- 9. Chocolate
- 10. Pinto bean paste





What we'll talk about today...

- The Trend
- FDA Guidance for Industry
- Applications
 - 1. Mayonnaise
 - 2. Cow's milk
 - 3. Plant based milk
 - 4. Sugar





The Trend...that drives sales



http://www.npr.org/sections/thesalt/2016/05/10/477514200/why-the-fdais-reevaluating-the-nutty-definition-of-healthy-food

© 2016 HORIBA, Ltd. All rights reserved.

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





The Trend: Clean Label

- Definition: "...natural ingredients with no artificial ingredients and chemicals*..."
- Not regulated
- Ingredients people can recognize and pronounce
 e.g. Monosodium Glutamate (MSG), Calcium phosphate, Potassium Bromate
- Food scientists' tool box is shrinking
- Now what?

* Source: http://www.clean-label.de/index.php?page=that-is-clean-label





FDA Guidance for Industry: Nutrition Facts Label

- Removed calories from fat
- Vitamin D and potassium are required
- Changed serving size
- Declare added sugar
- Applies to vending machines too
- Compliance date: July 28, 2018
- Annual sales < \$10M: July 28, 2019
- Adjust!

Nutrit Serving Size 2/3		Fac	cts	Nutrition Fa	cts
Servings Per Co		out 8		8 servings per container	
				Serving size 2/3 cup	o (55g)
Amount Per Servi	.		5.1.70		
Calories 230	Ca	lories fron	n Fat 72	Amount per serving	200
		% Dail	y Value*	Calories 2	230
Total Fat 8g			12%	9/ D-1	ly Veluet
Saturated Fat	: 1g		5%		ly Value* 10%
Trans Fat 0g				Total Fat 8g	
Cholesterol 0	0		0%	Saturated Fat 1g	5%
Sodium 160mg	,	7.0	7%	Trans Fat 0g	
Total Carbohydrate 37g 12%			12%	Cholesterol Omg	0%
Dietary Fiber	49		10%	Sodium 160mg	7%
Sugars 1g Protein 3g				Total Carbohydrate 37g	13%
Protein 5g				Dietary Fiber 4g	14%
Vitamin A			10%	Total Sugars 12g	
Vitamin C			8%	Includes 10g Added Sugars	20%
Calcium			20%	Protein 3g	
Iron			45%		
* Percent Daily Values are based on a 2,000 calorie diet. Your daily value may be higher or lower depending on your calorie needs.				Vitamin D 2mcg	10%
			ding on	Calcium 260mg	20%
·	Calories:	2,000	2,500	Iron 8mg	45%
Total Fat Sat Fat Cholesterol	Less than Less than Less than	65g 20g 300ma	80g 25g 200mg	Potassium 235mg	6%
Sodium Total Carbohydrate Dietary Fiber	Less than Less than	300mg 2,400mg 300g 25g	300mg 2,400mg 375g 30g	* The % Daily Value (DV) tells you how much a serving of food contributes to a daily diet. 2 a day is used for general nutrition advice.	

http://www.fda.gov/Food/GuidanceRegulation/GuidanceDocumentsRegulat oryInformation/LabelingNutrition/ucm385663.htm

© 2016 HORIBA, Ltd. All rights reserved.

Explore the future





"Whether you want to call it disruptive or exponential innovation, there's going to be a transformation in the food industry."

-Barb Renner Vice chairman and US consumer products practice leader

© 2016 HORIBA, Ltd. All rights reserved.

Explore the future


Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



Mayonnaise: Oil in Water (O/W) Emulsion

HORIBA Scientific

- Oil (dispersed phase) + vinegar (continuous phase) + egg yolk (emulsifier) + salt (taste)
- 4 ingredients or 15+ ingredients depending on target cost.



- Avoid canola oil and stick with "healthy fat" trend such as extra virgin olive oil, avocado, and almonds.*
- Samples produced with EVOO showed the lowest consistency and firmness when compared to other oils.**
- Physiochemical properties

References:



^{*}Hyman, Mark. *Eat Fat, Get Thin.* New York: Little Brown and Company, 2016. Print pg.77 *Key TJ, Appleby PN, Davey GK, Allen NE, Spencer EA, Travis RC. Mortality in British Vegetarians: review and prelimiary result from EPIC-Oxford. Am J Clin Nutr 2003 Sep; 78 **Carla Di Mattiaa. Physical properties, microstructure and stability of extravirgin olive oil based mayonnaise. InsideFood Symposium, 9-12 April 2013, Leuven, Belgium



Mayonnaise: Physiochemical Properties

- Physiochemical Properties:
 - 1. Emulsion stability

(choosing an appropriate plant-based emulsifier to mimic egg)


- 2. Rheological properties (spreadibility)
- 3. Sensory characteristics (taste, color, odor, consistency, texture, appearance, and overall acceptability)
- 4. Particle size distribution
- 5. pH
- 6. Cholesterol content
- 7. Microstructure




Mayonnaise: FDA Guideline

- The FDA requires that "mayonnaise" contain 65% vegetable oil and at least one egg yolk-containing ingredient (21 CFR 169.140(c)).
- "Mayo" is ok. "Mayonnaise" is not ok.
- Be careful if the "cholesterol free" claim meets 21 CFR 101.62(d).

Mayonnaise: An overview






© 2016 HORIBA, Ltd. All rights reserved.

Explore the future

HORIBA



Mayonnaise: Plant-based Protein vs. Traditional



*source: www.amazon.com

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific




Undersize (%)



Mayonnaise: An overview





© 2016 HORIBA, Ltd. All rights reserved.

Explore the future

HORIBA



Cow's Milk: Oil in Water (O/W) Emulsion

HORIBA Scientific

- 3-4% fat (dispersed phase) + 87% water (continuous phase) + 3.5% protein [casein micelles] (emulsifier) + 5% lactose + other essential nutrients
- 2016 "In the absence of any evidence for the superior effects of low fat dairy, and some evidence that there may be better benefits of whole fat dairy products for diabetes, why are we recommending only low fat diary?" Dariush Mozaffarian of Dean of Tufts University of Nutrition
- Particle size of fat determines stability, shelf life, taste, and mouth feel
- Homogenization to reduce fat droplet size
 - 1. Conventional 2 stage Gualin homogenizer
 - 2. Commonly used high-shear fluid processor/high pressure processor
 - 3. Typically aim for 0.2-2um particle size range




Cow's Milk: Fat size measurement methods

- 1980 Emulsion Quality Analyzer (EQA) measures the amount of light passing through <u>diluted</u> milk sample
 - Dilute milk with EQA solution (5% tetracetic acid, 2% sodium hydroxide, 2% sodium hexametaphosphate in water)
 - Sample prep by adding 10ml diluent to 250ml DI water, heat to 80-100F
 - Transfer 1ml to EAQ sample cell
 - Get the absorbance index reading, go to the chart, locate the fat % content, and finally read the fat globule size diameter.
- Laser diffraction particle size analyzer measures the angular intensity of the particles scattered from sample
 - Click fill to fill the analyzer with DI water
 - Add milk to approximately 5-10% concentration (aka obscuration)
 - Click measure to read size diameter



Cow's Milk: An Overview

Data name	Graph type	D(v,0.1)	D(v,0.5)	D(v,0.9)	Mean size
Popular Whole Milk avg	\Box —	0.22319 (µm)	0.55434 (µm)	1.20575 (µm)	0.65445 (µm)
Popular 2 Percent Milk avg		0.11486 (µm)	0.24057 (µm)	1.02227 (µm)	0.40825 (µm)
Popular 1 Percent Milk avg	\Box —	0.10564 (µm)	0.19535 (µm)	0.55576 (µm)	0.27659 (µm)
Popular Fat Free Milk avg		0.08744 (µm)	0.13889 (µm)	0.21124 (µm)	0.14495 (µm)
Non Homogenized Milk avg		2.50955 (µm)	3.95841 (µm)	6.37796 (µm)	4.24636 (µm)





© 2016 HORIBA, Ltd. All rights reserved.

Explore the future



Cow's Milk: Store Brand – the Success Story






Explore the future



Cow's Milk: A1 and A2 beta-casein in cow's milk



- Each cow carries two copies of the gene encoding beta-casein, A1 or A2g
- A2 mimics human breast milk "better" for digestive systems
- "A2 beta-casein is recognized as the original beta-casein protein because it existed before a mutation caused the appearance of A1 beta-casein in European herds a few thousand years ago" http://www.betacasein.org





Cow's Milk: FDA Guideline

- All must be pasteurized
- Vitamins fortification should be performed prior to homogenization
- Authorized to increase Vitamin D amount July 18, 2016
 - Up to 84 IU/100g of vitamin D3 to milk (42IU in 2012), 84 IU/100g of vitamin D2 to plant-based beverages intended as milk alternatives and 89 IU/100g of vitamin D2 to plant-based yogurt alternatives
- FDA does not regulate homogenization;
 U.S. Public Health Service provides guidelines.

Grade "A" Pasteurized Milk Ordinance

(Includes provisions from the Grade "A" Condensed and Dry Milk Products and Condensed and Dry Whey–Supplement I to the Grade "A" PMO)

2009 Revision



U.S. Department of Health and Human Services

Public Health Service

Food and Drug Administration

Source: http://www.fda.gov/downloads/Food/GuidanceRegulation/UCM209789.pdf

Plant Based Protein: Soymilk

HORIBA Scientific

- June, 2016
- Coca-Cola has entered into an agreement to acquire Unilever's AdeS soy-based beverage business

http://www.coca-colacompany.com/coca-cola-unbottled/the-coca-cola-company-and-coca-cola-femsa-to-acquireades-soy-based-beverage-business-from-unilever

- May, 2016
- Plant-based food category tops 3.5B. 9% growth for the past 2 years versus 4% for other. Plant-based milk drives the sales, 14% growth.

http://www.soyfoods.org/blog/plant-powered-sales-top-3-5-billion

- July, 2016
- DANONE acquired WhiteWave Foods Co. for 10.4 billion, stock jumped 19%

http://www.wsj.com/articles/danone-boosts-u-s-business-with-whitewave-deal-1467870422



Plant-based Protein: Soymilk Processing

Traditionally: water extraction of soybeans



• Commercially:

Clean Dehull Pressure cook	Grinding (milky slurry)	Centrifuge	Fortification /flavoring	Pasteurized/ homogenized	Bottle aseptically (no air)	
-------------------------------	-------------------------------	------------	-----------------------------	-----------------------------	-----------------------------------	--

"Beany" flavor – arise through lipoxygenase activity of soybean oil.

- Heat it to deactivate enzyme (e.g. lipoxygenases, trypsin inhibitors)
- Use defatted soy flour, soybean protein concentrate, or isolated soybean protein
- Mask it with sugar or other flavors (e.g. coffee, chocolate)
- Use GMO those without lipoxidase activity

Source: http://www.madehow.com/Volume-5/Soy-Milk.html#ixzz4ED6SuLY7

© 2016 HORIBA, Ltd. All rights reserved




Plant-based Protein: Importance of PSD in soymilk

- Track grinding processing
- Particle Size determines taste and texture
- Mouth is sensitive to particles >30um
- Particle size also determines the stability and quality of the emulsion
- High pressure heat treatment denatures soy proteins (11S, 7S at different temperature). This leads to rearrangement and denaturation of proteins (40nm> PSD<1um)*
- Milk becomes more viscous with high pressure heat treatment (homogeneous, smooth, and creamy texture)** when compare to traditional method

Sources: *Malaki Nik, A., Tosh, S., Poysa, V., Woodrow, L., & Corredig, M. (2008). Physicochemical characterization of soymilk after step-wise centrifugation. Food Research International, 41: 286-294. **Zuo Feng a,b, Peng Xingyun a, Shi Xiaodi a, Guo Shuntang a (2016). Effects of high-temperature pressure cooking and traditional cooking on soymilk: Protein particles formation and sensory quality

Plan-based Protein Soymilk



© 2016 HORIBA, Ltd. All rights reserved.

Explore the future

HORIBA





Sugar: Under attack

- Sugar, linked to obesity and diabetes*
- FDA (declare added sugar), American Heart Association (<9tsp/day), World Health Organization (<5% of daily calories). Reality =22-26 tsp
- Sugar reduction strategies easy said than done.

"When the company announced last spring that it planned to remove aspartame from Diet Pepsi, it cited declining sales and health concerns stemming from scientific studies linking artificial sweeteners to obesity and cancer in lab rats"....10months later, they're bringing it back to "give customers a choice" **

 Sugar serves as a flavor enhancer, preservative (think jams), bulking agent, retains moisture, browning process (visual color), controls spread of the biscuit dough (powdered sugar)


<u>Source: *http://www.ift.org/food-technology/past-issues/2016/may/features/science-of-taste.aspx</u> **http://money.cnn.com/2016/06/27/news/companies/pepsi-diet-aspartame/





Sugar: Stevia

- 200-450 times sweeter than granulated white sugar* depending on the species
- Spray-dried or vacuum dried
- Particle size affects the flow properties of stevia and defines how much a "spoon full of stevia" actually weighs
- Increase in flowability when PSD > 200um** agglomeration technology is used
- "Dusty" and harder to mix when particles are <10um



Source: *http://www.google.com/patents/EP2498625A1?cl=en **Comparison of critical particle diameters and its effect on flowability in stevia and sucrose, International conference on science and technique based on applied and fundamental research T. Dozan, M. Benkovic, I. Bauman ICoSTAF, 2014

© 2016 HORIBA, Ltd. All rights reserved.





Sugar: Powdered Sugar

- (Hammer) Milled* from white granulated crystals (Application notes AN141, AN175)
 - 2X defined as 82%<200mesh (74um)
 - 4X- defined as 92%<200mesh
 - 6X defined as 93.5% <200mesh
 - 8X defined as 96% <200mesh
 - 10X defined as 98%<200mesh
 - 12X defined as 98% <325mesh (45um)
 - Silk Sugar defined as 97% <20.5um

• Production goals:


- The efficiency of milling from pilot size to full production
- Narrow particle size distributions
 - Uniformity minimizes separation
 - Dissolution/mixing
- **<u>Flowability</u>** (anti-caking agent 3%)
- Dry Dispersion:
 - Sampling (>100um)
 - Energy (Pressure Size Titration test)
 - Slope Dv90>Dv50>Dv10

Source: *http://www.hmicronpowder.com/industries/food/sugar





Sugar: Particle size distribution



Explore the future



Sugar: Particle size distribution

				Diameter (um)	ASTM Mesh	Frequency %	Cum %
				20	635	69.071	69.071
				25	500	10.644	79.715
				32	450	8.164	87.879
			1	38	400	3.865	91.744
T T T T T T T T T T T T T T T T T T T	T T			45	325	2.768	94.512
				53	270	1.872	96.383
	-	F	F	63	230	1.417	97.801
			-				
				75	200	0.983	98.783
				90	170	0.649	99.432
				106	140	0.331	99.763
				125	120	0.186	99.949
				150	100	0.051	100
					1000	-60 (%) -140 -20 -20 -20 -20 -20 -20 -20 -20 -20 -2	
0.100 1.000 10.00 100.0	1.000 10.00 100.0	10.00 100.0			1000	5000	
	Diamates (com)	Diameter (µm)					



Questions?

Email us: Labinfo@horiba.com Call us: 1800-446-7422 www.horiba.com

© 2016 HORIBA, Ltd. All rights reserved.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

30

HORIBA





Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA