LIGHT SCATTERING THEORY

Laser Diffraction (Static Light Scattering)







When a Light beam Strikes a Particle

- Some of the light is:
 - Diffracted
 - Reflected
 - Refracted
 - Absorbed and Reradiated



- Small particles require knowledge of optical properties:
 - Real Index (degree of refraction)
 - Imaginary Index (absorption of light within particle)
 - Light must be collected over large range of angles
 - Index values less significant for large particles

Diffraction Pattern







Explore the future © 2007 HORIBA, Ltd. All rights reserved.

Diffraction Patterns







Single vs. Double Slit Patterns





Single slit diffraction pattern Double slit diffraction pattern



Fraunhofer Diffraction: Young's Double Slit

Distance between slits is INVERSELY proportional to angle (θ)



Explore the future © 2007 HORIBA, Ltd. All rights reserved

Fraunhofer Diffraction: Pinhole

The diameter of the pinhole is INVERSELY proportional to angle (θ)



FRAUNHOFER DIFFRACTION - PARTICLE

Angle of scatter from a particle is INVERSELY proportional to angle (θ)



© 2007 HORIBA, Ltd. All rights reserved

Edge Scattering Phenomenon



Light Scatter occurs whether from a slit, a pinhole or a particle. It occurs at the edge of an object. A SLIT and PARTICLE of the same size produce the same diffraction pattern





HORIBA

Explore the future © 2007 HORIBA, Ltd. All rights reserved

Fraunhofer Diffraction: Particles



Diffraction Pattern: Large vs. Small Particles



- Low angle scatter
- Large signal



Narrow Pattern - High intensity



Wide Pattern - Low intensity SMALL PARTICLE:

- High Angle Scatter
- Small Signal



Angle of Scatter vs. Particle Size

Fraunhofer Diffraction describes forward scatter technology. Scatter angles are relatively small, less than 30°



- ANGLE OF SCATTER is INVERSELY proportional to particle size. SMALL particles scatter at larger angles than large particles.
- The <u>AMOUNT OF LIGHT</u> scattered is DIRECTLY proportional to particle size. LARGE particles scatter MORE light than small particles.

Light Scattering

$$I(\Theta) = \frac{I_0}{2k^2 a^2} \left\{ \left[S_1(\Theta) \right]^2 + \left[S_2(\Theta) \right]^2 \right\}$$

- $I(\Theta)$ is the total scattered intensity as a function of angle Θ ;
- *I*₀ is the intensity of the incident light;
- k is the wavenumber = $2\pi/\lambda$;
- λ is the wavelength of I_0 of the illuminating source in air;
- *a* is the distance from scatterer to detector;

 $S1(\Theta)$ and $S2(\Theta)$ are dimensionless, complex functions defined in general scattering theory, describing the change of amplitude in respectively the perpendicular and the parallel polarized light as a function of angle Θ with respect to the forward direction. Computer algorithms have been developed in order to allow computation of these functions and, thus, of $I(\Theta)$.

Plot of Airy Function



HORIBA

Explore the future © 2007 HORIBA, Ltd. All rights reserved.

Comparison of Models



Light energy scattering patterns for an arbitrary detector configuration against particle size (μm) and scattering angle (°) for equal volumes of particles (Fraunhofer theory)

Light energy scattering patterns for an arbitrary detector configuration against particle size (µm) and scattering angle (°) for equal volumes of particles (Mie theory, latex particles RI 1.60 - 0.0i, in water RI 1.33)

Fraunhofer (left) vs. Mie (right)



Fraunhofer Approximation

$$(S_1)^2 = (S_2)^2 = \alpha^4 \left[\frac{J_1(\alpha \sin \Theta)}{\alpha \sin \Theta} \right]^2$$
$$I(\Theta) = \frac{I_0}{k^2 a^2} \alpha^4 \left[\frac{J_1(\alpha \sin \Theta)}{\alpha \sin \Theta} \right]^2$$

dimensionless size parameter $\alpha = \pi x / \lambda$;

J1 is the Bessel function of the first kind of order unity.

Assumptions:

a) all particles are much larger than the light wavelength (only scattering at the contour of the particle is considered; this also means that the same scattering pattern is obtained as for thin two-dimensional circular disks)

b) only scattering in the near-forward direction is considered (Q is small).

Limitation: (diameter at least about 40 times the wavelength of the light, or a >>1)* If λ =680nm (.68 µm), then 40 x .68 = 27 µm If the particle size is larger than about 50 µm, then the Fraunhofer approximation gives good results.

Fraunhofer Approximation



Explore the future © 2007 HORIBA, Ltd. All rights reserved.



Explore the future © 2007 HORIBA, Ltd. All rights reserved.

Mie vs. Fraunhofer



Figure A.3 -- Comparison of scattering patterns of non-absorbing particles according to Fraunhofer and Mie calculations ($N_p = 1,59 - 0,0; n_{water} = 1,33;$ wavelength = 633 nm)

Extinction Efficiencies



Figure A.2 -- Extinction efficiencies in relation to particle size and refractive index (Mie prediction). (n_p and k_p as indicated; $n_{water} = 1,33$; wavelength = 633 nm) (Fraunhofer assumes an extinction efficiency of 2 for all particle sizes)

Good agreement for transparent particles >50 μm """""" opaque particles > 2 μm

Scattering intensity pattern for single particles in relation to size



Mie calculation; wavelength 633 nm; Np = 1,59 - 0,0 i; nm = 1,33

© 2007 HORIBA, Ltd. All rights reserved

Light intensity scattering patterns for equal particle volumes in relation to size



Mie calculation; wavelength 633 nm; Np = 1,59 - 0,0; nm = 1,33

Influence of particle size on angular light intensity scattering patterns



Explore the future © 2007 HORIBA, Ltd. All rights reserved.

Influence of imaginary parts of RI (absorbancies)



© 2007 HORIBA, Ltd. All rights reserved.

Influence of distribution width on angular light intensity scattering patterns



Explore the future

© 2007 HORIBA, Ltd. All rights reserved.



Result Calculation







Iterations





Iterations

- How many times to calculate new (better) results before reporting final result?
- For LA-930 standard = 30, sharp = 150
- For LA-950 standard = 15, sharp = 1000
- Fewer iterations = broader peak
- More iterations = narrow peak
- Too many iterations = fit to noise, over-resolved



Viewing Raw Data



Explore the future © 2007 HORIBA, Ltd. All rights reserved.

LA-950 Optics





Why 2 Wavelengths?



Explore the future © 2007 HORIBA, Ltd. All rights reserved.



Practical Application of Theory: Mie vs. Fraunhofer for Glass Beads



Explore the future © 2007 HORIBA, Ltd. All rights reserved.

Practical Application of Theory: Mie vs. Fraunhofer for CMP Slurry



Explore the future © 2007 HORIBA, Ltd. All rights reserved.

Conclusions

- It helps to know some of the theory to best use a laser diffraction particle size analyzer
 - Small particles wide angles
 - Large particles low angles
- Look at Intensity curves, R parameter & Chi square calculations
- Use Mie theory at all times (default whenever choosing an RI kernel other than Fraunhofer)