

# Fine Particle Characterization in Chemical Mechanical Planarization (CMP) Slurries with Fluorescence Correlation Spectroscopy



June 8th, 2021

- 1. Brief CMP description
  - 1. Particle size variation
  - 2. Simplified abrasion model
- 2. Motivation and scope
  - 1. Process variation
  - 2. Roadmap Fine Particle Focus
    - 1. IRDS
    - 2. SEMI Standard
- 3. 2 methods using fluorescence
  - 1. FCS
  - 2. NTA



### **Basic Description of CMP**



Dynamic process...focus on particle metrology



- Slurry: Liquid, additives, particles
  - Chemically modify surface
  - Mechanical abrasion of modified surface
- Compliant pad
- Wafer

### **Simplified Mechanical Abrasion Model**



- ~D20 and less
- <u>Difficult to measure</u>
- These can be useless filler if too high
- Essentially depleting the active abrasives near MPS
- Can remain on the wafer as post CMP defects
- High surface area may react with other additives uncontrollably



### **Motivation and Scope**

#### Metrology improvements to reduce variation from particles



Reactive: respond to process parameter excursions for slurry in MFG process
Proactive: New metrology and methods
Particle suppliers: improve PSD metrology → control PSD
Slurry R&D → Optimize abrasive performance → robust CMP performance
Commercial slurry → Avoid process excursions due to abrasive particle raw material variation

## Background – IRDS 2018 & 2020

#### 2020 Challenge #5 is a perennial for IRDS Yield Enhancement...

Table YE-2b Yield Enhancement Difficult Challenges--SECC

Near-Term Difficult Challenges: 2020–2027	Description
Challenge #1 Control of electrically active and other particles in ultrapure water and	Insufficient metrology and marginal treatment technology capability
liquid chemicals for advanced semiconductor manufacturing (logic key driver)	
Challenge #2 Control of electrically active and other particles in UPW of EUV mask	Insufficient metrology and marginal treatment technology capability
production	
Challenge #3 Critical organics in UPW	Insufficient metrology for online speciation of organics
Challenge #4 Metals in H <sub>2</sub> O <sub>2</sub>	Insufficient metrology capability. Characterization needed for specific critical metals for
	key process steps and segments.
Challenge #5 Slurry particle characterization and other impurity monitoring	Insufficient metrology capability

https://irds.ieee.org/images/files/pdf/2020/2020IRDS\_YE\_Tables.xlsx

First Year of IC Production		2018	2019	2020	2021
Logic 1/2 pitch, nm (contacted)	18	18	14	14	12
GENERAL					
Fluid purity impact on device yield / performance. Requirements are specific to the chemistry					
and the Process step that it is used for.					
Contaminant based process control ???					
CHEMICALS - Metrology					
Development of 40 nm particle counter [A]					
Development of 20 nm particle counter [A]					
Development of 10 nm particle counter [A]					
Particle characterization to identify source of contamination [A], (for components and					
systems)					
Micelle particle size and concentration characterization [B] for surfactant containing					
chemistries.					
Bubble discrimination of reactive chemistries [A]					
Improved metrology for concentration measurements [B]					
Anion & cation measurement in process chemicals, e.g. cleaning chemicals [A]					
Organic measurement (TOC and speciation) in process chemicals, e.g. H2O2, IPA [A]					
Contaminant characterization in CMP slurries, e.g. zeta potential, large particle size [A]*					
Contaminant characterization in CMP slurries. (Agglomertaes, Foreign Material) [B]					
Characterization of CMP slurries, e.g. Particle size distribution [B]					

Roadmap perspective:

Perpetual need for improvements to particle metrology...



### **Background – SEMI Standards**

Semi standards do not list FCS or NTA on the list of techniques for particle size...

*However,* they do acknowledge the needs for improved particle metrology as semiconductor technology advances.

#### SEMI C98-1219

#### GUIDE FOR CHEMICAL MECHANICAL PLANARIZATION (CMP) PARTICLE SIZE DISTRIBUTION (PSD) MEASUREMENT AND REPORTING USED IN SEMICONDUCTOR MANUFACTURING

#### "Response to Semiconductor Technology Advances

- As node sizes continue to decrease, new requirements for CMP PSD are becoming increasingly important.
- As CMP slurries move towards smaller particle sizes more nano-size sensitive PSD metrologies are needed.
- More PSD information is needed at smaller particle sizes (i.e., more particle size bins in the nanometer size range).
- Metrologies capable of differentiating the nature of the measured particles are needed...."

#### Next gen PSD requires: "Sensitivity in the nm range....Smaller bins in the nm range"



### **The Problem**

What happens when mean particle size approaches the limit of resolution?



Improve Process Control and meet Roadmap Challenges: Improve the resolution of fine particles



### **Fluorescent Enhancement Techniques**





### Outline

- Overview of Fluorescence Correlation Spectroscopy (FCS)
- Use of Multiple Dyes in FCS Analysis and Application of Fluorescence Lifetime Analysis (FLA)
- FCS Characterization of Different Dyes Adsorbed Silica Abrasives
- FCS Characterization of the Small Particle Tail in Silica Particle Size Distributions Using Fractionation via Differential Sedimentation Analysis of Supernatant Fractions
- Conclusions and Directions for Future Studies



### **Fluorescence Correlation Spectroscopy**

<u>Consider the following Thought Experiment</u>: Image the fluorescence from a dilute solution of dye molecules.





diameter of sample cell





### Single-Molecule Spectroscopic Analysis Fluorescence Correlation Spectroscopy (FCS)







#### FCS Measurements: Intensity Fluctuation Autocorrelation Function, $G(\tau)$



#### Evaluate the Diffusion Coefficient, $D_{T}$ , and Number of Fluorophores, N



 $\tau_D$  is related to the *translational diffusion coefficient*,  $D_{\rm T}$ :  $\tau_D = r_0^2 / 4D_{\rm T}$ 

where:  $r_0^2$  is the radius of focal volume (~0.5 µm)

 $D_{\rm T}$  is related to the *hydrodynamic diameter*,  $d_{\rm h}$ :

 $D_T = kT / 3\pi \eta d_h$  (Stokes-Einstein equation)

where: k is the Boltzmann constant. T is the absolute temperature.  $\eta$  is the solvent viscosity.

G(0) is related to the average number of fluorophores in the focal volume, N:

G(0) = 1 / N



#### Representative FCS Analysis of Dispersed Silica Abrasive Using Adsorbed Rhodamine 110



Jacobson, L.M.; Turner, D.K.; Wayman, A.E.; Rawat, A.K.; Carver, C.T.; Moinpour, M.; Remsen, E.E. *ECS J. Solid State Sci. Tech.* **2015**, *4*, P5053-P5057. Schorr, D.K.; Smith, M.A.; Rawat, A.K.; Carver, C.T.; Mansour, M.; Remsen, E.E. *ECS Trans.* **2016** *72*, 43-51.

#### Fluorophores Employed in FCS Studies of Silica Abrasive Dispersions







Alexa fluor 488 (A488)  $\lambda_{ex} = 490$  nm,  $\lambda_{em} = 525$  nm Rhodamine 110 (R110)  $\lambda_{ex} = 496$  nm,  $\lambda_{em} = 520$  nm Rhodamine 6G (R6G)  $\lambda_{ex} = 530$  nm,  $\lambda_{em} = 565$  nm



#### Overlay of Normalized FCS ACFs for 10 nM A488, R110 and R6G





τ (s)

#### Overlay of Normalized ACFs for A488 and Adsorbed A488 on PL-2





#### Overlay of Normalized ACFs for R110 and Adsorbed R110 on PL-2



#### Overlay of Normalized ACFs for R6G and Adsorbed R6G on PL-2





#### Overlay of Normalized ACFs for Adsorbed A488, R110 and R6G on PL-2





#### Overlay of Normalized ACFs for Supernatants from Centrifuged R110 + S3 Mixtures as Function of Spin Time at 14,500 rpm



#### Overlay of Normalized ACFs for Supernatants from Centrifuged R110 + S3 Mixtures as Function of Rotor Speed for 60 min



### Conclusions and Directions for Future Studies

- A488, R110 and R6G adsorbed to silica without significant static or dynamic quenching
- The adsorption of the dyes to silica increases in the following order: A488 < R110 < R6G</li>
- Differential sedimentation of silica allows isolation and analysis by FCS of the small particle tail of the particle size distribution (PSD)
- Use R6G adsorption in conjunction with differential sedimentation for highest resolution of the small particle tail of the PSD
- Use FCS to quantify the adsorption isotherm of the dye interaction with silica



### Outline

- Overview Nanoparticle Tracking Analysis (NTA, MANTA)
  - Ultramicroscope  $\rightarrow$  NTA
- NTA with 3 lasers  $\rightarrow$  How the Viewsizer is differentiated from classic NTA
  - Improved PSD
  - Tagging specific particles
- Concept: Enhancing the detection of fluorescent-tagged fine particles
- Choice of suitable particle and dye
- Series of experiments
- Conclusions and Directions for Future Studies

### Ultra Microscope Concepts $\rightarrow$ <u>NTA:</u> Nanoparticle Tracking Analysis

The projected scatter is detected as a 2D image captured by a microscope (or high resolution camera). Invented in 1912 for colloids in a fluid  $\rightarrow$  Nobel Prize in 1925



The incident beam is scattered by particles within the interrogated volume of the optical system. *Detection* by scattering intensity, *Sizing* by Brownian Diffuision

#### **Viewsizer 3000: Multispectral Nanoparticle Tracking Analysis**

NTA – Nanoparticle Tracking Analysis MANTA - Multi Spectral Nanoparticle Tracking Analysis

Sample cuvette





High resolution Camera & video capture

#### **Viewsizer 3000: Multispectral Nano Particle Tracking Analysis**

3 lasers produce a wide range of particle detection



28

### **Viewsizer 3000 Video Capture**



Colors added to help distinguish particle sizes





#### **Viewsizer 3000: Multispectral Nano Particle Tracking Analysis**

- Particles are tracked, sized, counted
- counts/ml for each definable bin



Discreet counts of a discreet volume; not estimations of a fitted curve. More video capture  $\rightarrow$  more individual counts  $\rightarrow$  true PSD



## **Multispectral NTA**

Simultaneous illumination with 3 lasers give an accurate size through the PSD compared to similar methods; fines and coarse are measured accurately.







### **Fluorescence for Selective Tagging**



Scatter and fluorescence measured orthogonally.

With 3 laser excitation sources, we have many options of fluorophores to selectively identify mixed particles systems. Long pass and short pass filters allow for isolation of fluorescent particles.

### **Detection Limits – Scattering vs Fluorescent Emission**



in water (n=1.337), laser 445 nm wavelength, polarized, objective NA=0.28, 80°÷100° integration

		Raylei	gh Scattering
Ţ	$=I_0\left(rac{1}{2} ight)$	$rac{+\cos^2 heta}{2R^2} ight)$	$\left(rac{2\pi}{\lambda} ight)^4 \left(rac{n^2-1}{n^2+2} ight)^2 \left(rac{d}{2} ight)^6$
	I Theta R Iamda n	= = = =	scattered intensity viewing angle viewing distance wavelength of light particle refractive index

particle diameter

Rayleigh scattering diminishes quickly ~ d<sup>6</sup>

d

Calculated limit for detecting scatter in silica is ~ 40 nm.

Fluorescent particle emission signal diminishes less rapidly,  $\sim d^2$ 



## **Concept for Improving Detection Limits**



combine the scattered signal with fluorescent emission



## Test the Concept: A Challenging Colloidal Silica

Aggregate ratio = 40nm / 15nm =2.7

D2 = 40 nm



With a ~ 40nm mean particle size and a wide range of aggregate sizes, Fuso PL-1 is a good challenge for the resolution of the Viewsizer 3000.



https://fusokk.co.jp/eng/wp/wp-content/uploads/2019/05/5d5d380df5aea451f62d038624cb8cb5.pdf

241374

Semiconducto

# **Test the Concept: An Efficient Dye**

Cationic dye can interact electrostatically with the negative surface charge colloidal silica

• Very efficient adsorption in the literature

Surface modification of colloidal silica particles using cationic surfactant and the resulting adsorption of dyes Asad M. Khan

Journal of Molecular Liquids, Volume 274, 15 January 2019

Selective adsorption of organic dyes by porous hydrophilic silica aerogels from aqueous system *Wei Wei* **Water Science & Technology | 78.2 | 2018** 

Adsorptive removal of methylene blue from aqueous solution using coal fly ash-derived mesoporous silica material *Ning Yuan* **Adsorption Science & Technology 2019, Vol. 37(3–4) 333–348** 

Rhodamine 6G (R6G)  $\lambda_{ex}$  = 530 nm,  $\lambda_{em}$  = 565 nm

Based on Prof. Remsen's results and literature supporting cationic dye adsorption on SiO2, and good fluorescence alignment with our green laser We chose a high purity colloidal silica **and Rhodamine 6G** 



# **Test the Concept: Series of Experiments**



Diameter, nm

# **Optimized Conditions**



# **Series of Experiments**



Comparing the calculated limit of scattering detection for the Viewsizer and the enhanced signal of the dyed particles, The improved detection is pronounced at sizes smaller than ~ 65nm.

### **Illustration of the Experiments**



Standard process Fines poorly detected



Dyed and diluted, we can optimize the optical settings to improve detection of fines

Dyed particles increase counts Need dilution to count all particles



# **Summary of Experiments**

Addition of Rhodamine 6G to Fuso PL-1:

- Consistently increases the counts
- Detection of finer particles upon optimal dilution



Opportunities for follow-up

- More efficient dyes
  - blue and red
- Lot-to-lot (good vs. bad)
- Dyes for other particle types
- Process for complex slurry formulations (interfering additives)



