

How to Achieve One-Coat, High Hiding Power Paints by Optimizing Titanium Dioxide Pigments

HORIBA
Scientific

How to Achieve One-Coat,
High Hiding Power Paints
by Optimizing
Titanium Dioxide Pigments




Glenn Evers

Webinar • October 1, 2020

Glenn R. Evers, IS2 LLC
(M: 302-379-1265 eversgr@gmail.com)

Basic Paint : Water, Resin & TiO_2



TiO₂ Optimization in Low VOC Paints, Promises Better:


- **TiO₂ Grades / Extenders**
- **Dispersion**
- **Zeta Potential**
- **Particle Size**
 - Instrument Pay Back in 3 months!
- **Gold Standard**

TiO₂ Quality Certificate of Analysis

(If you don't ask, you don't get)


1. **Particle Size Distribution (PSD)**
2. **GSD (Geometric Standard Distribution)**
3. **% Greater than 0.5 um**
4. **Compare to Gold Standard (GS)**
5. Nibs or Scatts and distribution
6. Slurry Solids (76%)/pH
7. Rheology
8. Grit (wt.%)

Effect of % TiO₂ (alumina treated) on Hiding



% Min. TiO₂	Oil Abs. (lb/ 100 lb)	20° Enamel Gloss	Average Hiding Low PVC	Average Hiding High PVC	Relative Chalk Resistance	Use
95	12	70	100	90	Good	General – High gloss
91	13	67	95	100	Very Good	Universal Interior/ Exterior
89	17	60	90	95	Excellent	Exterior Durable
80	32	45	85	115	Good	Flat


High Opacity & Whiteness by Spacing




Green Glass Loses Color and Becomes Whiter by Grinding & Sieving to 2mm

850 μm \rightarrow 177 μm \rightarrow μ

Beer Bottle Glass





Basic Paint : Water, Resin & TiO_2



Ideal Dispersion is Better than Random

(Ref: Diebold (2005))



NEED Dispersion Dispersion Dispersion

- **BEST Dispersible & Durable R706 type coating alumina/silica better than SuperDurable – eg Roof Coatings**
- **BEST Gloss**
- **Effects Tint Strength, Contrast Ratio, Hiding Power, Whiteness**

CHALKING

1. Light Degrades Resin Binder
2. Works with water, acid rain, environmental deposits and UV
3. Worst case Oxygen Free in the film



a) Particle size CONTROLLED

Sulfate = Hydrolysis

Chloride = Oxidation/ Pressure

b) Grind: Steam /Pigment control

c) High density TiO_2 spacing

TiO₂ for Low VOC Paint

TiO ₂ Grade #	Slurry Viscosity @ 80 wt. % (low Cp)	Particle % > 0.5 um (low %)	Median Particle Size & GSD	Draw Down Particle Size & # Particles (low #)	Resin & HEUR Comp.	Gloss (high)	Optical: Con. Ratio, Color Comp. & Tint S (high)
#3 GS	++	+	++	+	++	++	+
#3	-	-	0	--	0	-	-

Lab Equipment Need For TiO_2 Optimization

- **Zeta Potential / pH Meter**
- **Rheometer**
- **Particle Size Analyzer**

Zeta for TiO₂ :


- In Paint: Avoid Dangerous pH Zones
- In Pigment Grind: Shift IEP w/ Dispersants
- Disperse Near IEP

Table 3—Characteristics of Different Surface-Treated TiO₂ Pigments^a

Sample #	Isoelectric Point (IEP)	Surface Area (m ² /g)	Description of TiO ₂ Pigments (Surface Concentration wt%)
1	9.0	15	4.5 Al ₂ O ₃ TiO ₂
2	6.3	15	3.6 Al ₂ O ₃ – 6.5 SiO ₂ TiO ₂
3	4.51	12.2	1.25 ZrO ₂ 2.0 Al ₂ O ₃ 6.0 SiO ₂ TiO ₂
4	6.36	14.3	2.0 Al ₂ O ₃ 0.5 ZrO ₂ –2.0 SiO ₂ TiO ₂
5	8.49	14.1	2.0 Al ₂ O ₃ 0.5 ZrO ₂ –1.0 Al ₂ O ₃ TiO ₂
6	8.28	15.0	1.0 ZrO ₂ –1.0 Al ₂ O ₃ TiO ₂

(a) Surface composition of inorganic oxides is described according to the order of surface treatment, see text.

Universal Grade TiO_2 Slurry Zeta Potential as a Function of pH



Lab Equipment Need For TiO₂ Optimization

- Zeta Potential / pH Meter
- Rheometer
- Particle Size Analyzer

Horiba PSD Plots: iPhone Video Tracking TiO₂ Spacer Composite Agglomeration

Effect of mixing TiO₂ slurry with CaCO₃ coated TiO₂




Demonstrate

- Real Time Particle Growth
- GSD Peak Broadening
- Final Finger Print

GSD 1.7 @ 0 min.

GSD 2.2 @ 11 min.


GSD 4.5 @ 45 min.



Two Types of Flocculation

- 1. Bridging:** Polymer Molecules Connect TiO_2 Particles.
- 2. Depletion:** TiO_2 get squeezed out of the thickener solution
 - a) Lower Gloss & CR
 - b) Poorer Film Strength
 - c) Poorer Adhesion

Rheology Modifier / TiO2 Dispersant Phase Diagram



Great GSD needs a hydrophobic micro-codispersant (Helps Rheology Chemicals)

JCT Research, Vol. 3, No.3, July 2006






Figure 2—Adsorption isotherms of hydrophilic and hydrophobic dispersants on interior and exterior TiO₂. Exterior TiO₂ adsorbs more dispersant, probably because of a higher surface area. Surfaces are close to being saturated at the 0.6–0.8% dispersant level.


Effect of pH on TiO₂ Median Particle Size

(Better is < 0.4um, AMP-95 pH adjusted)



Effect of pH on TiO₂ Particle Size GSD

(Lower numbers are better, AMP- 95 pH adjusted)



Gold Standard: Narrow Particle Size
Over all pH's

High Opacity VOC Paints Requirements

- Minimize Thickeners & Surfactants
- High Solids Filler Slurries
 - ✓ Target : 70 wt. % solids (min.)
- TiO₂ Slurry Properties
 - 1) Small Particle Size vs. pH
 - 2) Narrow – GSD vs. pH
 - 3) Zeta Potential vs. pH
 - 4) High Solids / Low Viscosity: 380-480 Cp @ 79% Solids (Brookfield : #3 spindle, 100 rpm)
 - 5) Compatibility Tests & Phase Diagrams