

Kuba Tatarkiewicz PhD, R&D Global Director

# How to present and compare data obtained by particle tracking analysis and other related methods

3-6-2019

# Particle tracking analysis

- **Measuring sizes of individual particles**

*hydrodynamic diameter  $d_h > d_0$  due to diluent drag*

- **Counting tracked particles in a given volume**

*investigated volume depends on size and refractive index*

- **Testing small volume of a sample**

- statistical process of estimating particle size distribution
- volume tested about 100,000x smaller than sample volume


# Concentration measurements

- **Definition of  $C_N$ : counts per volume [part/mL]**  
*other defs: mass  $C_M$  [mg/mL] or volume  $C_V$  [ $\mu$ L/mL]*
- **Typically  $C_N$  is given for a range of sizes**  
*e.g. between 100 nm and 1  $\mu$ m diameters*
- **Hence use of plots with size bins**
  - bin widths not necessarily equal
  - some software do not support unequal bins (notably Excel)


# Statistical considerations

- **Error in counting proportional to  $\sqrt{N}$**   
*this applies to total counts and individual bin counts*
- **Number of bins depends on expected errors**  
 *$N=10,000$  and 100 bins, average error 10 counts/bin or 10%*
- **To obtain “nice”, smooth distributions:**
  - decrease number of bins and use wider bins
  - use narrow bins to calculate parameters of a distribution


# Example plots of size distributions



Counts



Density of PSD  
*logarithmic bins*



Cumulative

# Binning

- **Definition of a bin:**
  - $d \in [b_i, b_{i+1}[$  where  $i=1, \dots, N$
  - equivalent definition:  $b_i \leq d < b_{i+1}$  (non-overlapping bins)
  - typically  $b_1 = 0$  and  $b_N = \text{max size to be measured}$
- Strange bins definition encountered:

		Average particles per mL	CV, %
29	particle diameter bin		
30	( 50 to 59) nm		
31	( 60 to 69) nm		
32	( 70 to 79) nm		
33	( 80 to 89) nm		
34	( 90 to 99) nm		
35	( 100 to 109) nm		

Define the particle diameter bin range, for example:  $100 \text{ nm} \leq d \leq 109 \text{ nm}$ . The preferred bin width is 10 nm. The total range should cover from  $> 100 \text{ nm}$  to  $< 2 \mu\text{m}$ . The range may be narrower, depending the range of the instrument being used. An example range is given here, but may be modified as appropriate for the measurement/instrument.

# What is really plotted

- **Concentration is just a rational number, also per bin**
- **Histogram uses area to represent a distribution**
- **With unequal bins natural measurable is:**

*density of particle size distribution (PSD)*

*units: number of particles per bin width and  
per investigated volume [part/nm/mL]*

*total concentration  $\equiv$  area of a histogram*

# Dimensional analysis

- **When plotting concentration vs. size,**  
*DO NOT connect points\**
- **Total concentration is a sum of numbers, not an integral**
- **When plotting density of PSD, one can connect middle of bins** ➔ *linear approximation*
- **Total concentration  $C_N$  is then**  
*an area under the curve*



*\*There's no data between points – by definition of a bin*

# Practical procedure

- **Create a list of measured sizes**

*typically diameters of nanoparticles are given in [nm]*


- **Bin those sizes in specific binning**

*typically logarithmic binning*

$$(b_{i+2} - b_{i+1}) / (b_{i+1} - b_i) = \text{const}$$

- **Calculate density of PSD [part/nm/mL];  $V_0$  needed**
- **Plot as a histogram of density of PSD**
- **Calculate parameters of a distribution (AV, SD,  $D_{50}$  ...)**

# Example of real data in logarithmic binning



Concentration from 50 nm to 700 nm = area of histogram of density of PSD

# Real data vs. standards

- **Most standards of size are mono-disperse**

*there are no standards for concentration...*

- **Typically real life samples are poly-disperse**

*e.g. multiple sizes or continuous distributions*

- **Natural samples like sea water or blood:**

- highly poly-disperse with dominating small particles
- Junge distribution  $N(d) \sim d^{-x}$  where  $x=3.5\div 4$

# Fitted data problem (NanoSight FTLA)

- **Fitting unknown distribution of sizes**
- **Assumed binomial(s) distribution**
- **Height (or area) of different peaks  
is NOT conserved during fitting  
(not invariant of any fitting procedure)**
- **Looks great but lacks numerical accuracy**
- **Artificial peaks can be created**

cf. van der Pol et al. J Thrombosis and Haemostasis, 12, 1182 (2014)

# Parametric description of distributions

- **Concentration**

$$B_i = (b_{i+1} - b_i) \quad N_{total} = \sum_{i=1}^N \frac{n_i}{B_i} B_i = \sum_{i=1}^N n_i$$

- **Average size**

$$d_{average} = \frac{1}{N_{total}} \sum_{i=1}^N n_i * \left( b_i + \frac{B_i}{2} \right)$$

- **Standard deviation**


$$SD = \sqrt{\frac{1}{N_{total}} \sum_{i=1}^N n_i * \left[ d_{average} - \left( b_i + \frac{B_i}{2} \right) \right]^2}$$

- **D<sub>50</sub>**

$$D_k = \frac{1}{N_{total}} \sum_{i=1}^k n_i$$

*definition  $D_k = 0.5$  ⇨  $k$  ⇨  $b_k$  equal  $D_{50}$*

# D50 & mode depend on binning!



# Lies, damned lies, and statistics


- **Parametric description using AV and SD  
is not accurate or unique!**

*Anscombe's quartet – same AV and SD, various shapes*

- **Even higher moments do not give enough info**
- **Basic experimental question:**

*How similar are two measured distributions?*

# Kolmogorov-Smirnov statistics



Comparing parameters:

$d_{av} = 256 \text{ nm}$ ,  $SD = 145 \text{ nm}$ ,  $CV = 0.57$



$d_{av} = 163 \text{ nm}$ ,  $SD = 68 \text{ nm}$ ,  $CV = 0.42$

Non-parametric test:

*Kolmogorov-Smirnov statistics*

$D_{A,B}$	$\alpha$	$D_{A,B,\alpha}$	Reject?
0.2335	0.05	0.0338	<b>yes</b>

# Anderson-Darling statistics



*Area between two cumulatives is a good measure of distance between two or more distributions with unknown shape.*

# Distance between distributions

- **Normalize area between cumulatives**

*extreme case: a) particles at 10 nm, b) particles at 1000 nm*

- **Distance is a number from the range [0,1]**

*same distributions have distance 0*

*extreme case – distance 1*



- **If areas are calculated between distributions in same normalization, it can be a good measure**

# Thank you

Omoshiro-okashiku  
Joy and Fun



감사합니다      Cảm ơn

ありがとうございました

Dziękuję      धन्यवाद      Grazie

Merci      谢谢      நன்றி

ขอขอบคุณครับ      Obrigado

Σας ευχαριστούμε

شُكْرًا

Tack ska ni ha

Большое спасибо

Danke

Gracias